

ENDODONTICS

ELPLORERS/EXCAVATORS

Explorers: These instruments used in a procedures to preserve the health of the dental pulp and periapical tissues, and the treatment of the diseased pulp to enable the tooth (or part of the tooth) to be retained in function. Further, endodontic treatment may be carried out on a healthy pulp where the tooth has to be devitalized as part of the overall restorative treatment plan.

Excavators: The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber.

Endodontic-Explorers

Endodontic-Excavators



ENDODONTICS

ELPLORERS/EXCAVATORS

Explorers: These instruments used in a procedures to preserve the health of the dental pulp and periapical tissues, and the treatment of the diseased pulp to enable the tooth (or part of the tooth) to be retained in function. Further, endodontic treatment may be carried out on a healthy pulp where the tooth has to be devitalized as part of the overall restorative treatment plan.


Excavators: The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber.



ENDODONTICS

PLUGGERS

Used to plug and condense amalgam into cavity preparation. Single ended hand instruments used to plug and condense gutta percha in a root canal.



Endodontic Pluggers ISO Size



ISO Size 30
24mm

RCP30
1400

ISO Size 40
24mm

RCP40
1401

ISO size 50
24mm

RCP50
1402

ISO Size 60
24mm

RCP60
1403

Endodontic Luks



1/3
21mm

RCP1/3
1404

5/7
21mm

RCP5/7
1405

ENDODONTICS

PLUGGERS

Used to plug and condense amalgam into cavity preparation. Single ended hand instruments used to plug and condense gutta percha in a root canal.



ENDODONTICS

PLUGGERS

Used to compact filling material during vertical condensation.

Anterior



ENDODONTICS

PLUGGERS

Used to compact filling material during vertical condensation.

Posterior



ENDODONTICS

NICKLE TITANIUM ENDODONTICS SPERADERS

Used to compact filling material during lateral condensation. These instruments made of nickel titanium.



D11

D11T

D11TS

255

25

40

40S

MA57

3NT

RCSD11NT

RCSD11TNT

RCSD11TSNT

RCS255NT

RCS25NT

RCS40NT

RCS40SNT

RCSMA57NT

RCS3NT

1480

1481

1482

1483

1484

1485

1486

1487

1488

ENDODONTICS

SPREADERS

Hand instrument with a tapered, pointed end, used to condense and compress materials such as gutta percha in root canals.

ISO Sized Spreaders



- ISO
Size 30
24mm

RCS30
1540
- ISO
Size 40
24mm

RCS40
1541
- ISO
Size 50
28mm

RCS50
1542
- ISO
Size 60
28mm

RCS60
1543

Spreaders



- GP1
20mm

RCGP1
1544
- 1S
Wakai
22mm

RCSW1S
1545
- 2S
Wakai
22mm

RCSW2S
1546
- D11
Standard
23mm

RCSD11
1547

ENDODONTICS

SPREADERS

Hand instrument with a tapered, pointed end, used to condense and compress materials such as gutta percha in root canals.



ENDODONTICS

MICRO CURETTES/MICRO SURGICAL EXCAVATORS

Excavators: The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber.

Micro-Curettes



ENDODONTICS

APICAL PLUGGERS

Used to compact filling material during vertical condensation.

Pluggers-Straight



ENDODONTICS

APICAL PLUGGERS

Used to compact filling material during vertical condensation.



ENDODONTICS

APICAL BURNISHERS / RETRO FILLING

Used to condense, smooth and polish amalgam.

Retro Filling

The retro filling instruments are sized to work with ultrasonic preparation tips.

Burnishers-Straight



Left Apical
PFLA
1700

Right Apical
PFRA
1701

1 Abou-Rass Straight
AB-ARS1
1702

2 Abou-Rass Straight
AB-ARS2
1703

3 Abou-Rass Straight
AB-ARS3
1704

4 Abou-Rass Straight
AB-ARS4
1705

5 Abou-Rass Straight
AB-ARS5
1706

ENDODONTICS

APICAL BURNISHERS / RETRO FILLING

Used to condense, smooth and polish amalgam.

Burnishers-Angled



1
Abou-Rass
Angled

AB-ARS1

1707

2
Abou-Rass
Angled

AB-ARS2

1708

3
Abou-Rass
Angled

AB-ARS3

1709

4
Abou-Rass
Angled

AB-ARS4

1710

5
Abou-Rass
Angled

AB-ARS5


1711

ENDODONTICS

APICAL PLUGGERS / BURNISHERS

Pluggers: are used to condense filling materials into cavity preparation. Ends are serrated.

Burnisher: Used to condense, smooth and polish amalgam.



Small Left
Apical

APBSL

1712

Small Right
Apical

APBSR

1713

Medium Left
Apical

APBML

1714

Medium Right
Apical

APBMR

1715

Large Left
Apical

ABLL

1716

ENDODONTICS

APICAL PLUGGERS / BURNISHERS

Pluggers: are used to condense filling materials into cavity preparation. Ends are serrated.

Burnisher: Used to condense, smooth and polish amalgam.



Large Right
Apical

APBLR
1717

Universal
Small Apical

APBUS
1718

Universal
Medium Apical

APBUM
1719

Universal
Large Apical

APBUL
1720

ENDODONTICS

FORCEPS / PUNCHES / RUBBER DAM ACCESSORIES

Splinting Forceps: Immobilization of fractured bone ends by wiring, pinning or splinting.

R.D. Clamp Forceps: Instruments with self-locking device used to apply and remove rubber dam clamps.

R.D. Frame: Flexible frame used to hold the rubber dam in a taut state, away from the field of operation.

R.D. Punch: Instruments used to punch various sized holes in a sheet of rubber dam.

Martensite stainless steel, Clamps shown 50% actual.

Post And Silver Point Removal Forceps to remove posts from root canals




Rubber Dam Clamp Forceps

shown 40% size



Rubber Dam Frames



Rubber Dam Punches







































ENDODONTICS

RUBBER DAM CLAMPS






















Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls. these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

Molars - General

					
	8A		W8A		8AD
	RDC8A		RDCW8A		RDC8AD
	1830		1831		1832
					
	10		11		14
	RDC10		RDC11		RDC14
	1833		1834		1835
					
	14A		W14A		18
	RDC14A		RDCW14A		RDC18
	1836		1837		1838
					
	24		25		26N
	RDC24		RDC25		RDC26N
	1839		1840		1841

Upper Molars

			
	56		56S
	RDC56		RDC56S
	1842		1843
			
	4		W4
	RDC4		RDCW4
	1844		1845
			
	5		8
	RDC5		RDC8
	1846		1847
			
	27		
	RDC27		
	1848		

ENDODONTICS

RUBBER DAM CLAMPS

Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls. these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

Upper Molars	Lower Molars			Pre-molars
  30 RDC30 1849	  3 RDC3 1853	  W3 RDCW3 1854	  7 RDC7 1855	  0 RDC0 1862
  31 RDC31 1850	  W7 RDCW7 1856	  7A RDC7A 1857	  12A RDC12A 1858	  1 RDC1 1863
  201 RDC201 1851	  31A RDC31A 1859	  26 RDC26 1860	  28 RDC28 1861	  1A RDC1A 1864
  205 RDC205 1852				  2 RDC2 1865

ENDODONTICS

RUBBER DAM CLAMPS

Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls. these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

Pre-molars



Anterior

