

HOLE PUNCHING PLIERS



HOLE PUNCH PLIERS
Creates Half-Moon Cutouts in Clear Aligners
2426



TEAR DROP PUNCH PLIERS
Designed to Create a Notch at the
Gingival Margin of Clear Aligners
2427



HORIZONTAL PUNCH PLIERS
Horizontal Rectangular Shaped Punch
2428



VERTICAL PUNCH PLIERS
Vertical Tectangular Shaped Punch Pliers
2429

ORTHODONTICS

PIN & LIGATURE CUTTING PLIERS



**PIN & LIGATURE CUTTER
LONG HANDLE**
Cuts Ligature Wires upto 0.015"/0.38mm
2430



**PIN & LIGATURE CUTTER
MINI STYLE**
Cuts Ligature Wires upto 0.015"/0.38mm
2431



**PIN & LIGATURE CUTTER
NARROW TIP**
Cuts Ligature Wires upto 0.015"/0.38mm
2432



**PIN LIGATURE CUTTER
NARROW TIP 15, ANGLED**
Cuts Ligature Wires upto 0.015"/0.38mm
2433



**PIN & LIGATURE CUTTER
STANDARD**
Cuts Ligature Wires upto 0.015"/0.38mm
2434



**PIN & LIGATURE CUTTER
STANDARD 15 ANGLED**
Cuts Wires upto 0.015"/0.38mm
2435

HARD WIRE CUTTING PLIERS



ORTHO HARD WIRE CUTTER
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm
Precisely Angled Straight Tip
2436



ORTHO HARD WIRE CUTTER, 15 ANGLD
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm
2437



**DISTAL END CUTTER
MINI STYLE CUT-HOLD**
Safely Cuts & Holds Wires upto 0.025"/0.64mm
2438



**DISTAL END CUTTER
STANDARD CUT-HOLD**
Safely Cuts & Holds Wires upto 0.025"/0.64mm
2439



**DISTAL END CUTTER
LONG HANDLE, CUT-HOLD**
Safely Cuts & Holds Wires upto 0.025"/0.64mm
2440



**FLUSH CUT-HOLD DISTAL END CUTTER
STANDARD HANDLE**
Featuring a Tungsten Carbide Tip Insert
2441



**FLUSH CUT-HOLD DISTAL CUTTER
LONG HANDLE**
Featuring a Tungsten Carbide Tip Insert
2442

UTILITY PLIERS



**HOW PLIERS
STRAIGHT**
Precisely Serrated Jaws for a Firm Grip
2450



**DEBONDING PLIERS
STRAIGHT**
Precisely Aligned, Straight Tip
2451



**FORCE MODULE
SEPARATING PLIERS**
Smooth Reversible Action Control
2452



**COON STYLE
LIGATURE PLIERS**
Long Slim Tips with Reversible Action Control
2453



ADHESIVE REMOVING PLIERS
Double Sided Tungsten Carbide Tip Insert
Comes with Free Pack of Occlusal Pads
2454



**WEINGART PLIERS
LONG HANDLE**
Designed for Wires up to 0.030"/0.76mm
2555



**HOW PLIERS
CURVED**
Precisely Serrated Jaws for Firm Grip
2457



**POSTERIOR BAND REMOVING
PLIERS, LONG TIP**
Comes with a free Pack of Occlusal Pads
2458



**BRACKET REMOVING PLIERS,
STRAIGHT**
Slim edge Jaws for an easy Grip on Bracket
2459

ORTHODONTICS

UTILITY PLIERS



**BRACKET REMOVING PLIERS
ANGLED**
Slim Edge Jaws for an Easy Grip on Bracket
2460



**POSTERIOR BAND REMOVING PLIERS
STANDARD TIP**
Comes with a free pack of Occlusal pads
2461



**DEBONDING PLIERS
ANGLED**
60 Angled Tip
2462



**WEINGART PLIERS 90
ANGLED**
Designed for wires
upto 0.030"/0.76mm
2463



**WEINGART PLIERS
SLIM**
Designed for wires upto 0.030"/0.76mm
2464



**WEINGART UTILITY
PLIERS**
Designed for wires upto 0.030"/0.76mm
2465



**WEINGART PLIERS, UNIVERSAL
45 ANGLED TIP**
Designed for wires upto 0.030"/0.76mm
2466



**CROWN & BAND CONTOURING
PLIERS**
Featuring a concave and convex tip design
2467



**CROWN & BAND CRIMPING
PLIERS**
Slim design tips for precise crimping
2468



BALL HOOK CRIMPING PLIERS
Designed to crimp hooks, posts and
auxillary stops on arch wire
2469



**POSTERIOR BAND REMOVING
PLIERS, SHORT TIP**
Comes with a free pack of Occlusal pads
2470

WIRE FORMING PLIERS



BIRD BEAK PLIERS
Featuring a Cone & Pyramid shaped beak
Designed for wires upto 0.028"/0.71mm
2475



**LIGHT WIRE PLIERS
LONG BEAK**
Precisely built for wires upto 0.020"/0.51mm
2476



TWEED LOOP FORMING PLIERS
Designed for wire upto 0.020"/0.51mm
2477



LIGHT WIRE PLIERS WITH 3 GROOVES
Precisely built for wires upto 0.020"/0.51mm
2478



**LIGHT WIRE PLIERS WITH
CUTTER & 3 GROOVES**
Precisely built for wires upto 0.020"/0.51mm
2479



**LIGHT WIRE PLIERS
WITH CUTTER**
Precisely built for wires upto 0.020"/0.51mm
2480



JARABAK PLIERS
Featuring 3 grooves for a firm wire grip
2481



JARABAK PLIERS WITH CUTTER
Featuring 3 grooves & cutter at base
2482



LINGUAL ARCH FORMING PLIERS
Designed to create double back and triple back
bends in all wires upto 0.036"/0.91mm
2483

ORTHODONTICS

WIRE FORMING PLIERS



NANCE LOOP FORMING PLIERS
Four step beak to create 3,4,5 and 6mm loop
designed for all wires upto 0.022"/0.56mm
2484



De La Rosa PLIERS
Featuring a chamfered concave and convex tip
built for all wires upto 0.025"/0.64mm
2485



De La Rosa Pliers WITH 3 GROOVES
Featuring a chamfered concave and convex tip
built for all wires upto 0.025"/0.64mm
2486



HOLLOW CHOP CONTOURING PLIERS
Long concave & Convex tip
for intricate wire contouring
2487



OPTICAL PLIERS
Featuring a slim beak with concave end
designed for all wire upto 0.025"/0.64mm
2488



V-STOP PLIERS
Create V-shaped bends on
wires upto 0.022"/0.56mm
2489



YOUNG PLIERS
3 sectional beak with 3 grooves
designed for wires upto 0.028"/0.71mm
2490



TWEED ARCH FORMING PLIERS 1.75mm
Designed for wires upto 0.025"/0.66mm
2491



TWEED ARCH FORMING PLIERS 1.0mm
Designed for wires upto 0.025"/0.66mm
2492

ORTHODONTICS

WIRE FORMING PLIERS



TORQUING PLIERS FEMALE/OUTER
Places torques in all wires upto 0.030"/0.76mm alongside our inner torquing pliers
2493



TORQUING PLIERS MALE/INNER
Places torques in all wires upto 0.030"/0.76mm alongside our outer torquing pliers
2494



STEP PLIERS 1MM
Designed for all wires upto 0.028"/0.71mm
2495



STEP PLIERS 0.75MM
Designed for all wires upto 0.028"/0.71mm
2496



STEP PLIERS 0.5MM
Designed for all wires upto 0.028"/0.71mm
2497



NiTi DISTAL CINCH BACK PLIERS
Designed for SS & NiTi wires upto 0.025"/0.64mm
2498



NiTi DISTAL CINCH BACK PLIERS, LINGUAL
Designed for SS&NiTi wires upto 0.018"/0.46mm
2499



3 PRONG PLIERS STANDARD
Precisely built for all wires upto 0.030"/0.76mm
2500



3 PRONG PLIERS HEAVY
Precisely built for all wires upto 0.030"/0.76mm
2501



3 PRONG PLIERS NiTi
Precisely built for all wire upto 0.020"/0.5mm
2502



Ligature
Director/Band
Pusher
2506



Small
Plugger/
Picker
2507



Dull Scaler/
Band
Pusher
2508



Ligature
Director
Scaler
2509



Scaler/
Band
Pusher
2510



Lingual
Ligature
Director,
U Shape
2511

HAND INSTRUMENTS



Dull
Scaler

2 5 1 2



Dull
Scaler

2 5 1 3



Large
Plugger/
Picker

2 5 1 4



Picker/
Ligature
Director

2 5 1 5



Scaler/
Picker

2 5 1 6



Dull
Scaler

2 5 1 7




Ligature
Director/
Band Pusher

2 5 1 8


ORTHODONTICS

HAND INSTRUMENTS



Band Pusher with Short Tip

2 5 1 9



Band Pusher

2 5 2 0



Fillion Explorer/
Ligature Director

2 5 2 1



Band Pusher/
Scaler

2 5 2 2



Band Pusher/
Scaler

2 5 2 3



Ligature Director

2 5 2 4




Ligature Director

2 5 2 5

ORTHODONTICS

HAND INSTRUMENTS



Lingual Torquing Key

2 5 2 6

Straight .020" Applicator

2 5 2 7

Single End Ligature Director

2 5 2 8

Lingual Ligature Director, .018"

2 5 2 9

Lingual Ligature Director, .025"

2 5 3 0

Ligature Director/ Plugger

2 5 3 1

Ligature Director Tucker

2 5 3 2

Plugger/ Ligature Picker

2 5 3 3



Posterior Bracket Holder Tweezer, Long

2566



Direct Bond Tweezer with Aligner & Positioner

2567




Buccal Tube Tweezer

2568



Slim Buccal Tube Bonder Tweezer

2569



Direct Bond Bracket Tweezer, Lingual

2570



Direct Bond Tweezer, Aligner

2571



All Purposes Scissors
2572



Crown Scissors
2573



Crown Scissors, Straight
2574



Crown Scissors, Curved
2575

HEIGHT GAUGES



BONE BRACKET POSITIONING GAUGE
2580



BRACKET HEIGHT GAUGE FOR
0.022"/0.56mm Brackets
2581



BRACKET HEIGHT GAUGE FOR
0.018"/0.46mm Brackets
2582



BRACKET HEIGHT GAUGE, MOVEABLE HEAD
2583



MATHIEU NARROW TIP
Enhanced double spring function
2585



MATHIEU SMAHA FINE TIP
Enhanced double spring function
2586



MATHIEU NARROW TIP, TC
Enhanced double spring function
2587



MATHIEU STANDARD TIP
Enhanced double spring function
2588



MATHIEU STANDARD TIP, TC
Enhanced double spring function
2589