

## HOLE PUNCHING PLIERS



**HOLE PUNCH PLIERS**  
Creates Half-Moon Cutouts in Clear Aligners  
2426



**TEAR DROP PUNCH PLIERS**  
Designed to Create a Notch at the  
Gingival Margin of Clear Aligners  
2427



**HORIZONTAL PUNCH PLIERS**  
Horizontal Rectangular Shaped Punch  
2428



**VERTICAL PUNCH PLIERS**  
Vertical Rectangular Shaped Punch Pliers  
2429



## PIN & LIGATURE CUTTING PLIERS



**PIN & LIGATURE CUTTER  
LONG HANDLE**  
Cuts Ligature Wires upto 0.015"/0.38mm  
2430



**PIN & LIGATURE CUTTER  
MINI STYLE**  
Cuts Ligature Wires upto 0.015"/0.38mm  
2431



**PIN & LIGATURE CUTTER  
NARROW TIP**  
Cuts Ligature Wires upto 0.015"/0.38mm  
2432



**PIN LIGATURE CUTTER  
NARROW TIP 15, ANGLED**  
Cuts Ligature Wires upto 0.015"/0.38mm  
2433



**PIN & LIGATURE CUTTER  
STANDARD**  
Cuts Ligature Wires upto 0.015"/0.38mm  
2434



**PIN & LIGATURE CUTTER  
STANDARD 15 ANGLED**  
Cuts Wires upto 0.015"/0.38mm  
2435



# ORTHODONTICS

## HARD WIRE CUTTING PLIERS



**ORTHO HARD WIRE CUTTER**  
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm  
Precisely Angled Straight Tip  
2436



**ORTHO HARD WIRE CUTTER, 15 ANGLED**  
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm  
2437



**DISTAL END CUTTER MINI STYLE CUT-HOLD**  
Safely Cuts & Holds Wires upto 0.025"/0.64mm  
2438



**DISTAL END CUTTER STANDARD CUT-HOLD**  
Safely Cuts & Holds Wires upto 0.025"/0.64mm  
2439



**DISTAL END CUTTER LONG HANDLE, CUT-HOLD**  
Safely Cuts & Holds Wires upto 0.025"/0.64mm  
2440



**FLUSH CUT-HOLD DISTAL END CUTTER STANDARD HANDLE**  
Featuring a Tungsten Carbide Tip Insert  
2441



**FLUSH CUT-HOLD DISTAL CUTTER LONG HANDLE**  
Featuring a Tungsten Carbide Tip Insert  
2442



# ORTHODONTICS

## UTILITY PLIERS



**HOW PLIERS STRAIGHT**  
Precisely Serrated Jaws for a Firm Grip  
2450



**DEBONDING PLIERS STRAIGHT**  
Precisely Aligned, Straight Tip  
2451



**FORCE MODULE SEPARATING PLIERS**  
Smooth Reversible Action Control  
2452



**COON STYLE LIGATURE PLIERS**  
Long Slim Tips with Reversible Action Control  
2453



**ADHESIVE REMOVING PLIERS**  
Double Sided Tungsten Carbide Tip Insert  
Comes with Free Pack of Occlusal Pads  
2454



**WEINGART PLIERS LONG HANDLE**  
Designed for Wires up to 0.030"/0.76mm  
2555



**HOW PLIERS CURVED**  
Precisely Serrated Jaws for Firm Grip  
2457



**POSTERIOR BAND REMOVING PLIERS, LONG TIP**  
Comes with a free Pack of Occlusal Pads  
2458



**BRACKET REMOVING PLIERS, STRAIGHT**  
Slim edge Jaws for an easy Grip on Bracket  
2459



## UTILITY PLIERS



**BRACKET REMOVING PLIERS  
ANGLED**  
Slim Edge Jaws for an Easy Grip on Bracket  
2460



**POSTERIOR BAND REMOVING PLIERS  
STANDARD TIP**  
Comes with a free pack of Occusal pads  
2461



**DEBONDING PLIERS  
ANGLED**  
60 Angled Tip  
2462



**WEINGART PLIERS 90  
ANGLED**  
Designed for wires  
upto 0.030"/0.76mm  
2463



**WEINGART PLIERS  
SLIM**  
Designed for wires upto 0.030"/0.76mm  
2464



**WEINGART UTILITY  
PLIERS**  
Designed for wires upto 0.030"/0.76mm  
2465



**WEINGART PLIERS, UNIVERSAL  
45 ANGLED TIP**  
Designed for wires upto 0.030"/0.76mm  
2466



**CROWN & BAND CONTOURING  
PLIERS**  
Featuring a concave and convex tip design  
2467



**CROWN & BAND CRIMPING  
PLIERS**  
Slim design tips for precise crimping  
2468



**BALL HOOK CRIMPING PLIERS**  
Designed to crimp hooks, posts and  
auxillary stops on arch wire  
2469



**POSTERIOR BAND REMOVING  
PLIERS, SHORT TIP**  
Comes with a free pack of Occlusal pads  
2470



## WIRE FORMING PLIERS



**BIRD BEAK PLIERS**  
Featuring a Cone & Pyramid shaped beak  
Designed for wires upto 0.028"/0.71mm  
2475



**LIGHT WIRE PLIERS  
LONG BEAK**  
Precisely built for wires upto 0.020"/0.51mm  
2476



**TWEED LOOP FORMING PLIERS**  
Designed for wire upto 0.020"/0.51mm  
2477



**LIGHT WIRE PLIERS WITH 3 GROOVES**  
Precisely built for wires upto 0.020"/0.51mm  
2478



**LIGHT WIRE PLIERS WITH  
CUTTER & 3 GROOVES**  
Precisely built for wires upto 0.020"/0.51mm  
2479



**LIGHT WIRE PLIERS  
WITH CUTTER**  
Precisely built for wires upto 0.020"/0.51mm  
2480



**JARABAK PLIERS**  
Featuring 3 grooves for a firm wire grip  
2481



**JARABAK PLIERS WITH CUTTER**  
Featuring 3 grooves & cutter at base  
2482



**LINGUAL ARCH FORMING PLIERS**  
Designed to create double back and triple back  
bends in all wires upto 0.036"/0.91mm  
2483



## WIRE FORMING PLIERS



**NANCE LOOP FORMING PLIERS**  
Four step beak to create 3,4,5 and 6mm loop  
designed for all wires upto 0.022"/0.56mm  
2484



**De La Rosa PLIERS**  
Featuring a chamfered concave and convex tip  
built for all wires upto 0.025"/0.64mm  
2485



**De La Rosa Pliers WITH 3 GROOVES**  
Featuring a chamfered concave and convex tip  
built for all wires upto 0.025"/0.64mm  
2486



**HOLLOW CHOP CONTOURING PLIERS**  
Long concave & Convex tip  
for intricate wire contouring  
2487




**OPTICAL PLIERS**  
Featuring a slim beak with concave end  
designed for all wire upto 0.025"/0.64mm  
2488



**V-STOP PLIERS**  
Create V-shaped bends on  
wires upto 0.022"/0.56mm  
2489



**YOUNG PLIERS**  
3 sectional beak with 3 grooves  
designed for wires upto 0.028"/0.71mm  
2490



**TWEED ARCH FORMING PLIERS 1.75mm**  
Designed for wires upto 0.025"/0.66mm  
2491



**TWEED ARCH FORMING PLIERS 1.0mm**  
Designed for wires upto 0.025"/0.66mm  
2492



## WIRE FORMING PLIERS



**TORQUING PLIERS FEMALE/OUTER**  
Places torques in all wires upto 0.030"/0.76mm  
alongside our inner torquing pliers  
**2493**



**TORQUING PLIERS MALE/INNER**  
Places torques in all wires upto 0.030"/0.76mm  
alongside our outer torquing pliers  
**2494**



**STEP PLIERS 1MM**  
Designed for all wires upto 0.028"/0.71mm  
**2495**



**STEP PLIERS 0.75MM**  
Designed for all wires upto 0.028"/0.71mm  
**2496**



**STEP PLIERS 0.5MM**  
Designed for all wires upto 0.028"/0.71mm  
**2497**



**NiTi DISTAL CINCH BACK PLIERS**  
Designed for SS & NiTi wires upto 0.025"/0.64mm  
**2498**



**NiTi DISTAL CINCH BACK PLIERS, LINGUAL**  
Designed for SS&NiTi wires  
upto 0.018"/0.46mm  
**2499**



**3 PRONG PLIERS STANDARD**  
Precisely built for all wires  
upto 0.030"/0.76mm  
**2500**



**3 PRONG PLIERS HEAVY**  
Precisely built for all wires  
upto 0.030"/0.76mm  
**2501**



**3 PRONG PLIERS NiTi**  
Precisely built for all wire  
upto 0.020"/0.5mm  
**2502**





Ligature  
Director/Band  
Pusher  
2506



Small  
Plugger/  
Picker  
2507



Dull Scaler/  
Band  
Pusher  
2508



Ligature  
Director  
Scaler  
2509



Scaler/  
Band  
Pusher  
2510



Lingual  
Ligature  
Director,  
U Shape  
2511



# ORTHODONTICS

## HAND INSTRUMENTS



Dull Scaler  
2512



Dull Scaler  
2513



Large Plugger/  
Picker  
2514



Picker/  
Ligature  
Director  
2515



Scaler/  
Picker  
2516



Dull Scaler  
2517




Ligature  
Director/  
Band Pusher  
2518




# ORTHODONTICS

## HAND INSTRUMENTS



Band Pusher with Short Tip

2 5 1 9



Band Pusher

2 5 2 0



Fillion Explorer/  
Ligature Director

2 5 2 1



Band Pusher/  
Scaler

2 5 2 2



Band Pusher/  
Scaler

2 5 2 3



Ligature Director

2 5 2 4




Ligature Director

2 5 2 5



# ORTHODONTICS

## HAND INSTRUMENTS



Lingual Torquing Key

2 5 2 6

Straight .020" Applicator

2 5 2 7

Single End Ligature Director

2 5 2 8

Lingual Ligature Director, .018"

2 5 2 9

Lingual Ligature Director, .025"

2 5 3 0

Ligature Director/Plugger

2 5 3 1

Ligature Director Tucker

2 5 3 2

Plugger/Ligature Picker

2 5 3 3



# ORTHODONTICS

## TWEEZERS-SCISSORS



Posterior Bracket Holder Tweezer, Long

2566



Direct Bond Tweezer with Aligner & Positioner

2567




Buccal Tube Tweezer

2568



Slim Buccal Tube Bonder Tweezer

2569



Direct Bond Bracket Tweezer, Lingual

2570



Direct Bond Tweezer, Aligner

2571



All Purposes Scissors  
2572



Crown Scissors  
2573



Crown Scissors, Straight  
2574



Crown Scissors, Curved  
2575



## HEIGHT GAUGES



BONE BRACKET POSITIONING GAUGE  
2580



BRACKET HEIGHT GAUGE FOR  
0.022"/0.56mm Brackets  
2581



BRACKET HEIGHT GAUGE FOR  
0.018"/0.46mm Brackets  
2582



BRACKET HEIGHT GAUGE, MOVEABLE HEAD  
2583



## MATHIEU NEEDLE HOLDERS



MATHIEU NARROW TIP  
Enhanced double spring function  
2585



MATHIEU SMAHA FINE TIP  
Enhanced double spring function  
2586



MATHIEU NARROW TIP, TC  
Enhanced double spring function  
2587



MATHIEU STANDARD TIP  
Enhanced double spring function  
2588



MATHIEU STANDARD TIP, TC  
Enhanced double spring function  
2589