

DIAGNOSTIC MIRRORS

Cone Socket Handles



#M1 Cone Socket Mirror Handle **M1CS** 901



#M2 Cone Socket Mirror Handle **M2CS** 902



#M3 Cone Socket Mirror Handle **M3CS** 903



#M4 Cone Socket Mirror Handle **M4CS** 904

Simple Stem Handles



#M5 Simple Stem Mirror Handle **M5SS** 905



#M6 Simple Stem Mirror Handle **M6SS** 906



#M7 Simple Stem Mirror Handle **M7SS** 907



#M8 Simple Stem Mirror Handle **M8SS** 908



Simple Stem Mouth Mirror
#4-7/8" (22mm)

SSMM#4

910



Simple Stem Mouth Mirror
#5-15/16" (24mm)

SSMM#5

911



Simple Stem Mouth Mirror
Double Sided
#4-7/8" (22mm)

SSMMDS#4

913



Simple Stem Mouth Mirror
Double Sided
#5-15/16" (24mm)

SSMMDS#5

914



Cone Socket Mouth Mirror
#4-7/8" (22mm)

CSMM#4

915



Cone Socket Mouth Mirror
#5-15/16" (24mm)

CSMM#5

916

DIAGNOSTIC

MICRO SURGICAL MIRRORS



1
Small
Rectangular
MICMS1
1013



2
Small
Rectangular
MICMS2
1014



3
Small
Rectangular
MICMS3
1015




Small Round
MICMSR
1016



Large Round
MICMLR
1017

DIAGNOSTIC MICRO SURGICAL MIRRORS



DIAGNOSTIC

DRESSING PLIERS

Used to grasp or transfer materials in and out of the oral cavity. All dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment. Stainless steel made.



1
Perry
DP1
1071



2
College
DP2
1072



3
DP3
1073



17
DP17
1074



18
DP18
1075

DIAGNOSTIC

DRESSING PLIERS

Used to grasp or transfer materials in and out of the oral cavity. All dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment. Stainless steel made.



18L
Locking

DP18L

1076



19

DP19

1077



Utility
Pick-Up
Dressing Plier

DPU17

1078



317
Smooth
Beak

DP317

1079



318
Foi/Meriam
Serrated Beak


DP318

1080


DIAGNOSTIC HANDLE OPTIONS

Use for caries and calculus detection, exploration of pocket characteristics, furcation and restorations.

DOUBLE ENDED



SINGLE ENDED



DIAGNOSTIC EXPLORERS

tooth surface irregularities, defective margins on restorations, decalcified areas and carious lesions. For the quality assurance, we use the Germany made material. Moreover, explores are available in a variety of design types.

Single-End



DIAGNOSTIC

EXPLORERS



8
EXS8
1127

9
EXS9
1128

23
Shepherd's
Hook
EXS23
1129

54
EXS54
1130


TU17
Wilkins/Tufts
EXSTU17
1131

17
EXS17
1132

DIAGNOSTIC EXPLORERS

A dental explorer is an instrument in dentistry commonly used in the dental armamentarium. Use for caries and calculus detection, exploration of pocket characteristics, furcation and restorations.


Double-End



DIAGNOSTIC EXPLORERS

A dental explorer is an instrument in dentistry commonly used in the dental armamentarium. Use for caries and calculus detection, exploration of pocket characteristics, furcation and restorations.

Double-End



DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



UNC12
University of
North Carolina
CPUNC12
1207

UNC15
University of
North Carolina
CPUNC15
1208

PQ-O
Michigan
CPQO
1209

PQ-W
Williams
CPQW
1210

Moffitt/
Maryland
CPMM
1211

Williams
NovatechNT2
CPWNT2
1212

DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



DIAGNOSTIC COLOR PROBES

Double-Ended Screening Probes




DIAGNOSTIC COLOR PROBES



DIAGNOSTIC EX PROBES

What is an Ex-Probe? Ex-Probe is an excellent combination instrument for dental diagnostics. It has an tip on one end and a probe on the other.



10-9-8-7-5-3-2-1

23/QW
Williams

XP23/QW
1290

10-9-8-7-5-3-2-1

23/QOW
Thin
Williams

XP23/TW
1291

11-8-6-3

23/CP-8

XP23/CP8
1292

11-8-6-3

23/CP-11

XP23/CP11
1293

11.5-8.5-5.5-3.5

23/CP-11.5B
Screening
Probe

XP23/CP11.5B
1294

12-9-6-3

23/CP-12

XP23/CP12
1295

12-9-6-3

23/NT3CP12
Novatech

XP23/CP12NT3
1296


12-11-10-9-
8-7-6-5-4-3-2-1

23/UNC12
University of
North Carolina

XP23/UNC12
1297

DIAGNOSTIC EX PROBES

What is an Ex-Probe? Ex-Probe is an excellent combination instrument for dental diagnostics. It has an tip on one end and a probe tip on the other.



15-14-13-12-11-10-9-8-7-6-5-4-3-2-1

23/UNC15
University of North Carolina

23/CP12UNC
1298

15-14-13-12-11-10-9-8-7-6-5-4-3-2-1

TU17/
UNC15

XPTU17/UNC15
1299

10-9-8-7-5-3-2-1

23/W
Williams

XP23/W
1300

10-9-8-7-5-3-2-1

6/OW
Thin
Williams

XP6/OW
1301

10-9-8-7-5-3-2-1

6/W
Williams

XP6/W
1302

10-9-8-7-5-3-2-1

17/OW
Thin
Williams

XP17/OW
1303

15-14-13-12-11-10-9-8-7-6-5-4-3-2-1

3A/UNC15

XP3A/UNC15
1304

10-9-8-7-5-3-2-1
FLAT

23/GF
Goldman Fox

XP23/GF
1305

8-6-3

23/O
Michigan

XP23/O
1306