

Implantology Instruments



Implantology Instruments Index

Hard wire Cutting pliers				O3
Hole Punching Pliers				O1
Ortho Hand Instruments				O9-O12
Ortho height gauges				O14
Ortho Needle Holders				O15
Ortho Tweezers and Scissors				O13
Pin & Ligature Cutting pliers				O2
Utility Pliers				O5
Wire Forming Pliers				O6-O8

IMPLANTOLOGY

IMPLANT PERIOSTEALS

Hand instrument, usually double ended, used to separate the periosteum from the underlying bone.



1
LaBanca
Periosteal
IMPLAB1
1920

2
LaBanca
Modified
IMPLABM2
1921

Buser
IMPBUS
1922


P9/PA
IMP9/PA
1923

25
IMP25
1924

IMPLANTOLOGY

IMPLANT PERIODONTAL CHISELS

Hand instrument used to remove hard tissues such as enamel or bone by chipping, cleaving or paring. The beveled blade may be in line with the handle or at an angle to it.



1/2
Buser Modified
Periodontal
Chisel
3mm/4mm

IPCBM1/2
1970

2/3
Buser Modified
Periodontal
Chisel
4mm/5mm

IPCBM2/3
1971

1/3
Buser Modified
Periodontal
Chisel
3mm/5mm

IPCBM1/3
1972

Jovanovic
13KTG
Periodontal
Chisel
1mm/4mm

IPCJOV1/4
1973

Jovanovic
13KTG
Periodontal
Chisel
2mm/5mm

IPCJOV2/5
1974

Jovanovic
13KTG
Periodontal
Chisel
3mm/6mm

IPCJOV3/6
1975

P24G/13KL
Palacci Modified

IPC24G/13KL
1976

IMPLANTOLOGY

BONE SCOOPS & PLUGGERS

Bone Scoops: Hand instrument similar to a bone chisel but having a hollowed-out blade similar to coupland chisel.

Bone Plugger: Used to plug and condense grafted bone.

BONE SCOOPS



Jovaanovic
None Scoop
5mm

IBSJOV5mm
2020

Bone Graft
Scoop

IBSBG
2021

Labanca
Bone Curette
7mm

ICLBC
2022

BONE PLUGGERS



LaBanca
3x1.5mm

IBPLAB
2023

4
Jovanovic
3mm

IBPJ4
2024

Bone
Packer DE
5-7mm

IBPP5-7
2025

Bone Packer
Serrated
3-5mm

IBPS3-5
2026

Bone Grafting Pluggers/Condensers



Bone Grafting Plugger (1.0-1.7mm)
(8-11-13-15-18)
IBG1.0/1.7
2081

IBG1.0/1.7

2081



Bone Grafting Plugger (1.7-2.0mm)
(8-11-13-15-18)

IBG1.7/2.0

2082



Bone Grafting Plugger (2.4-2.7mm)
(8-11-13-15-18)

IBG2.4/2.7

2083



Grisdale Bone Condensor (3.3/4.0mm)
(8-11-13-15-18)

IBG3.3/4.0

2084



Grisdale Bone Condensor (3.3/4.0mm)
(8-11-13-15-18), Serrated

IBGS3.3/4.0

2085



TPI1-Tunneling Procedure Instrument

ITPI1
2090



TPI2-Tunneling Procedure Instrument

ITPI2
2091



TPI3-Tunneling Procedure Instrument

ITPI3
2092



TPI4-Tunneling Procedure Instrument

ITPI4
2093



ITPI5-Tunneling Procedure Instrument

ITPI5
2094



TPI6-Tunneling Procedure Instrument

ITPI6
2095



TPI7-Tunneling Procedure Instrument

ITPI7
2096

TUNNELING KNIVES



Urban Tunneling Knife

IUTK
2097



Tunneling Knife 1

ITK1
2098



Tunneling Knife 2 (Short Handle)

ITK2S
2099



Tunneling Knife 2

ITK2
2100



Giles De Quincey Knife 1

IGDQK1
2101



Giles De Quincey Knife 2

IGDQK2
2102



Steigmann Soft Tissue Splitter

ISSTS
2103



Steigmann Soft Tissue Splitter, Modified

ISSTSM
2104

IMPLANTOLOGY

OSTEOTOME CHISELS

Osteotome Chisels: Surgical instruments used to cut or divide bone.



7-10-13-15-18mm

Osteotome Chisel
Straight, Single Edge
4mm

IOCSSE4

2 1 1 0



7-10-13-15-18mm

Osteotome Chisel
Curved, Single Edge
4mm

IOCCSE4

2 1 1 1



7-10-13-15-18mm

Osteotome Chisel
Straight, Single Edge
7.5mm

IOCSSE7.5

2 1 1 2



7-10-13-15-18mm

Osteotome Chisel
Curved, Single Edge
7.5mm

IOCCSE7.5


2 1 1 3



Osteotome Chisel
Jovanovic
6mm

IOCJOV6

2 1 1 4




8-10-13-18mm

Osteotome Chisel
Single Edge
3mm

IOCSE3

2 1 1 5




8-10-13-18mm

Osteotome Chisel
8-10-13-18mm
Double edge, 3mm

IOCDE3

2 1 1 6




8-10-13-18mm

Osteotome Chisel
8-10-13-18mm
Single Edge, 6mm

IOCSE6

2 1 1 7



8-10-13-18mm

Osteotome Chisel
8-10-13-18mm
Double Edge, 6mm

IOCDE6

2 1 1 8




8-10-13-18mm

Osteotome Chisel
8-10-13-18mm
Single Edge, 9mm

IOCSE9

2 1 1 9



8-10-13-18mm

Osteotome Chisel
8-10-13-18mm
Double Edge, 9mm

IOCDE9

2 1 2 0

Martensite stainless steel tissue punches.



Tissue Punch STR
4mm
ITPS4
2 130



Tissue Punch STR
5mm
ITPS5
2 131



Tissue Punch STR
6mm
ITPS6
2 132



Tissue Punch CVD
2mm
ITPC2
2 133



Tissue Punch CVD
3mm
ITPC3
2 134



Tissue Punch CVD
4mm
ITPC4
2 135



Tissue Punch CVD
5mm
ITPC5
2 136



Tissue Punch CVD
6mm
ITPC6
2 137

Sinus Lift Instruments



Sinus Lift Cobra Shape 3.5mm (2-4-5-6-8-10mm)

ISC3.5
2152



Sinus Lift Cobra Shape In/Out (2-4-5-6-8-10mm)

ISCI/O
2153



Sinus Lift Mushroom Shape (2.8/3.8mm)
(2-4-5-6-8-10)

IMSL2.8/3.8
2154




Sinus Lift Mushroom Shape (4.8/5.8mm)
(2-4-5-6-8-10)

IMSL4.8/5.8
2155

IMPLANTOLOGY

SINUS LIFT INSTRUMENTS




IMPLANTOLOGY

SINUS LIFT INSTRUMENTS


Designed to separate the sinus lining from the bone without fraying or tearing.



Designed to separate the sinus lining from the bone without fraying or tearing.




Designed to separate the sinus lining from the bone without fraying or tearing.






2.0mm
ITCCSO2.0
2255

Tapered Concave Tip-Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm
ITCCSO3.5
2256

Tapered Concave Tip-Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm
ITCCSO4.3
2257

Tapered Concave Tip-Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm
ITCCSO5.0
2258

Tapered Concave Tip-Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm
ITCCSO5.5
2259

Tapered Concave Tip-Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm
ITCCAO2.0
2260

Tapered Concave Tip-Angulated - 2.0mm (4-6-8-10-13-16-18-20-23-26)



3.5mm
ITCCAO3.5
2261

Tapered Concave Tip-Angulated - 3.5mm (4-6-8-10-13-16-18-20-23-26)




4.3mm
ITCCAO4.3
2262

Tapered Concave Tip-Angulated - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm
ITCCAO5.0
2263

Tapered Concave Tip-Angulated - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm
ITCCAO5.5
2264

Tapered Concave Tip-Angulated - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm
ITCVSO2.0
2270

Tapered Convex- Tip-Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm
ITCVSO3.5
2271

Tapered Convex- Tip-Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm
ITCVSO4.3
2272

Tapered Convex- Tip-Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm
ITCVSO5.0
2273

Tapered Convex- Tip-Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm
ITCVSO5.5
2274

Tapered Convex- Tip-Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm
ITCVAO2.0
2275

Tapered Convex-Angulated Tip - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm
ITCVAO3.5
2276

Tapered Convex-Angulated Tip - 3.5mm (4-6-8-10-13-16-18-20-23-26)




4.3mm
ITCVAO4.3
2277

Tapered Convex-Angulated Tip - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm
ITCVAO5.0
2278

Tapered Convex-Angulated Tip - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm
ITCVA5.5
2279

Tapered Convex-Angulated Tip - 5.5mm (4-6-8-10-13-16-18-20-23-26)



2.0mm
IPSO2.0
2281

PUSHER - Tip-Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm
IPSO3.5
2282

PUSHER Tip-Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm
IPSO4.3
2283

PUSHER- Tip-Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm
IPSO5.0
2284

PUSHER- Tip-Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)



5.5mm
IPSO5.5
2285

PUSHER- Tip-Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm
IPA02.0
2286

PUSHER- Angulated Tip - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm
IPA03.5
2287

PUSHER- Angulated Tip - 3.5mm (4-6-8-10-13-16-18-20-23-26)




4.3mm
IPA04.3
2288

PUSHER- Angulated Tip - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm
IPA05.0
2289

PUSHER- Angulated Tip - 5.0mm (4-6-8-10-13-16-18-20-23-26)



5.5mm
IPA05.5
2290

PUSHER- Angulated Tip - 5.5mm (4-6-8-10-13-16-18-20-23-26)

Bone Spreading
Osteotome Kit



Interchangeable
Osteotome Handle

TAPERED CONVEX OSTETOMES


STRAIGHT

2.7 mm 3.2 mm 3.7 mm 4.2 mm 5.0 mm



ANGLED

2.7 mm 3.2 mm 3.7 mm 4.2 mm 5.0 mm




ICVOS2.7 2300 ICVOS3.2 2301 ICVOS3.7 2302 ICVOS4.2 2303 ICVOS5.0 2304

ICVOA2.7 2305 ICVOA3.2 2306 ICVOA3.7 2307 ICVOA4.2 2308 ICVOA5.0 2309


TAPERED CONCAVE OSTEOTOMES

STRAIGHT



ICCOS2.7 2310 ICCOS3.2 2311 ICCOS3.7 2312 ICCOS4.2 2313 ICCOS5.0 2314


ANGLED



ICCOA2.7 2315 ICCOA3.2 2316 ICCOA3.7 2317 ICCOA4.2 2318 ICCOA5.0 2319


BONE PUSHER OSTEOTOMES

STRAIGHT



IPOS2.7 2320 IPOS3.2 2321 IPOS3.7 2322 IPOS4.2 2323 IPOS5.0 2324

ANGLED



IPOA2.7 2325 IPOA3.2 2326 IPOA3.7 2327 IPOA4.2 2328 IPOA5.0 2329

IMPLANTOLOGY

CALIPERS & BONE GRAFTING SYRINGES



CASTROVIEJO
STR 17CM

ICVCS17
2382



Castroviejo
Calipers
Straight
0-40mm
18cm/7"

ICVCS18
2384



Castroviejo
Calipers
Angled
0-40mm
18cm/7"

ICVCCA18
2386



Bone Grafting
Syringe
7mm

IBGS7
2396

- 2389-1.0mm **IBGS1**
2389
- 2390-2.0mm **IBGS2**
2390
- 2391-2.5mm **IBGS2.5**
2391
- 2392-3.0mm **IBGS3**
2392
- 2393-4.0mm **IBGS4**
2393
- 2394-5.0mm **IBGS5**
2394
- 2395-6.0mm **IBGS6**
2395

NOTE: Also available in
1mm, 2mm, 2.5mm, 3mm, 4m, 5mm and 6mm



Weiss Modified
Castroviejo
Calipers
0-40mm
17cm/6 3/4"

IWMCV17
2383



CASTROVIEJO
STR 8.5cm
ICVC-ST8.5cm

ICVCS8.5
2385



CASTROVIEJO
CVD 17cm
ICVC-CV17cm

ICVCA17
2387



Cortical
Bone Clamp

ICBC
2388



Bone Grafting 1 Ring Style Syringe
Straight, 3.75mm

IBGS1RSS3.75
2400



Bone Grafting 1 Ring Style Syringe
Straight, 50.0mm

IBGS1RSS5.0
2402



Bone Grafting 1 Ring Style Syringe
Curved, 3.75mm

IBGS1RSC3.75
2401



Bone Grafting 1 Ring style Syringe
Curved, 5.0mm

IBGS1RSC5.0
2403



Ridge Mapping, Bone Width 0.25mm
Both Side Reading, 14cm

RMBC14
2410



Caliper Double sided Reading
Scissors Style

ICDSS
2411



Sliding Implant Caliper, 94mm
1.5mm Tip diameter

ISIC94
2412



Stainless Steel Ruler
15cm/6"

SSR
2413



Bone collector Chise, Straight

IBCC
2415



Bone Collector Blade, Straight

IBCBS
2416



Bone Collector Blade

IBCBS
2418

IBCCC
2417

Bone Collector Chisel, Curved

IMPLANTOLOGY


ENDO FORCEPS, BONE CRUSHERS & MALLETS

Micro Surgery Endo Forceps, 95mm

IMSEF95
2419

Micro Surgery Endo forceps, 110mm

IMSEF110
2420



Mead-2
Mallet

IM2M
2421



Bone Crusher

IBC
2424



130g
Mallet

I130M
2422



235g
Mallet

I235M
2423



Bone Morseriser

IBM
2425