

Orthodontic Instruments



HOLE PUNCHING PLIERS



HOLE PUNCH PLIERS
Creates Half-Moon Cutouts in Clear Aligners

2426



TEAR DROP PUNCH PLIERS
Designed to Create a Notch at the
Gingival Margin of Clear Aligners

2427



HORIZONTAL PUNCH PLIERS
Horizontal Rectangular Shaped Punch

2428



VERTICAL PUNCH PLIERS
Vertical Tectangular Shaped Punch Pliers

2429

PIN & LIGATURE CUTTING PLIERS



**PIN & LIGATURE CUTTER
LONG HANDLE**
Cuts Ligature Wires upto 0.015"/0.38mm

2430



**PIN & LIGATURE CUTTER
MINI STYLE**
Cuts Ligature Wires upto 0.015"/0.38mm

2431



**PIN & LIGATURE CUTTER
NARROW TIP**
Cuts Ligature Wires upto 0.015"/0.38mm

2432



**PIN LIGATURE CUTTER
NARROW TIP 15, ANGLED**
Cuts Ligature Wires upto 0.015"/0.38mm

2433



**PIN & LIGATURE CUTTER
STANDARD**
Cuts Ligature Wires upto 0.015"/0.38mm

2434



**PIN & LIGATURE CUTTER
STANDARD 15 ANGLED**
Cuts Wires upto 0.015"/0.38mm

2435

HARD WIRE CUTTING PLIERS



ORTHO HARD WIRE CUTTER
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm
Precisely Angled Straight Tip

2436



ORTHO HARD WIRE CUTTER, 15 ANGLD
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm

2437



**DISTAL END CUTTER
MINI STYLE CUT-HOLD**
Safely Cuts & Holds Wires upto 0.025"/0.64mm

2438



**DISTAL END CUTTER
STANDARD CUT-HOLD**
Safely Cuts & Holds Wires upto 0.025"/0.64mm

2439



**DISTAL END CUTTER
LONG HANDLE, CUT-HOLD**
Safely Cuts & Holds Wires upto 0.025"/0.64mm

2440



**FLUSH CUT-HOLD DISTAL END CUTTER
STANDARD HANDLE**
Featuring a Tungsten Carbide Tip Insert

2441



**FLUSH CUT-HOLD DISTAL CUTTER
LONG HANDLE**
Featuring a Tungsten Carbide Tip Insert

2442

UTILITY PLIERS



**HOW PLIERS
STRAIGHT**
Precisely Serrated Jaws for a Firm Grip

2450



**DEBONDING PLIERS
STRAIGHT**
Precisely Aligned, Straight Tip

2451



**FORCE MODULE
SEPARATING PLIERS**
Smooth Reversible Action Control

2452



**COON STYLE
LIGATURE PLIERS**
Long Slim Tips with Reversible Action Control

2453



ADHESIVE REMOVING PLIERS
Double Sided Tungsten Carbide Tip Insert
Comes with Free Pack of Occlusal Pads

2454



**WEINGART PLIERS
LONG HANDLE**
Designed for Wires upto 0.030"/0.76mm

2455



**HOW PLIERS
CURVED**
Precisely Serrated Jaws for Firm Grip

2457



**POSTERIOR BAND REMOVING
PLIERS, LONG TIP**
Comes with a free Pack of Occlusal Pads

2458



**BRACKET REMOVING PLIERS,
STRAIGHT**
Slim edge Jaws for an easy Grip on Bracket

2459

UTILITY PLIERS



**BRACKET REMOVING PLIERS
ANGLED**
Slim Edge Jaws for an Easy Grip on Bracket

2460



**POSTERIOR BAND REMOVING PLIERS
STANDARD TIP**
Comes with a free pack of Occlusal pads

2461



**DEBONDING PLIERS
ANGLED**
60 Angled Tip

2462



**WEINGART PLIERS 90
ANGLED**
Designed for wires
upto 0.030"/0.76mm

2463



**WEINGART PLIERS
SLIM**
Designed for wires upto 0.030"/0.76mm

2464



**WEINGART UTILITY
PLIERS**
Designed for wires upto 0.030"/0.76mm

2465



**WEINGART PLIERS, UNIVERSAL
45 ANGLED TIP**
Designed for wires upto 0.030"/0.76mm

2466



**CROWN & BAND CONTOURING
PLIERS**
Featuring a concave and convex tip design

2467



**CROWN & BAND CRIMPING
PLIERS**
Slim design tips for precise crimping

2468



BALL HOOK CRIMPING PLIERS
Designed to crimp hooks, posts and
auxillary stops on arch wire

2469



**POSTERIOR BAND REMOVING
PLIERS, SHORT TIP**
Comes with a free pack of Occlusal pads

2470

WIRE FORMING PLIERS



BIRD BEAK PLIERS
Featuring a Cone & Pyramid shaped beak
Designed for wires upto 0.028"/0.71mm

2475



**LIGHT WIRE PLIERS
LONG BEAK**
Precisely built for wires upto 0.020"/0.51mm

2476



TWEED LOOP FORMING PLIERS
Designed for wire upto 0.020"/0.51mm

2477



LIGHT WIRE PLIERS WITH 3 GROOVES
Precisely built for wires upto 0.020"/0.51mm

2478



**LIGHT WIRE PLIERS WITH
CUTTER & 3 GROOVES**
Precisely built for wires upto 0.020"/0.51mm

2479



**LIGHT WIRE PLIERS
WITH CUTTER**
Precisely built for wires upto 0.020"/0.51mm

2480



JARABAK PLIERS
Featuring 3 grooves for a firm wire grip

2481



JARABAK PLIERS WITH CUTTER
Featuring 3 grooves & cutter at base

2482



LINGUAL ARCH FORMING PLIERS
Designed to create double back and triple back
bends in all wires upto 0.036"/0.91mm

2483

WIRE FORMING PLIERS



NANCE LOOP FORMING PLIERS
Four step beak to create 3,4,5 and 6mm loop
designed for all wires upto 0.022"/0.56mm

2484



De La Rosa PLIERS
Featuring a chamfered concave and convex tip
built for all wires upto 0.025"/0.64mm

2485



De La Rosa Pliers WITH 3 GROOVES
Featuring a chamfered concave and convex tip
built for all wires upto 0.025"/0.64mm

2486



HOLLOW CHOP CONTOURING PLIERS
Long concave & Convex tip
for intricate wire contouring

2487



OPTICAL PLIERS
Featuring a slim beak with concave end
designed for all wire upto 0.025"/0.64mm

2488



V-STOP PLIERS
Create V-shaped bends on
wires upto 0.022"/0.56mm

2489



YOUNG PLIERS
3 sectional beak with 3 grooves
designed for wires upto 0.028"/0.71mm

2490



TWEED ARCH FORMING PLIERS 1.75mm
Designed for wires upto 0.025"/0.66mm

2491



TWEED ARCH FORMING PLIERS 1.0mm
Designed for wires upto 0.025"/0.66mm

2492

WIRE FORMING PLIERS



TORQUING PLIERS FEMALE/OUTER
Places torques in all wires upto 0.030"/0.76mm alongside our inner torquing pliers

2493



TORQUING PLIERS MALE/INNER
Places torques in all wires upto 0.030"/0.76mm alongside our outer torquing pliers

2494



STEP PLIERS 1MM
Designed for all wires upto 0.028"/0.71mm

2495



STEP PLIERS 0.75MM
Designed for all wires upto 0.028"/0.71mm

2496



STEP PLIERS 0.5MM
Designed for all wires upto 0.028"/0.71mm

2497



NiTi DISTAL CINCH BACK PLIERS
Designed for SS & NiTi wires upto 0.025"/0.64mm

2498



NiTi DISTAL CINCH BACK PLIERS, LINGUAL
Designed for SS&NiTi wires upto 0.018"/0.46mm

2499



3 PRONG PLIERS STANDARD
Precisely built for all wires upto 0.030"/0.76mm

2500



3 PRONG PLIERS HEAVY
Precisely built for all wires upto 0.030"/0.76mm

2501



3 PRONG PLIERS NiTi
Precisely built for all wire upto 0.020"/0.5mm

2502

HAND INSTRUMENTS



Ligature
Director/Band
Pusher
2506



Small
Plugger/
Picker
2507



Dull Scaler
Band Pusher
2508



Ligature
Director
Scaler
2509



Scaler/
Band
Pusher
2510



Lingual
Ligature
Director,
U Shape
2511

ORTHODONTICS

HAND INSTRUMENTS



Dull
Scaler
2512



Dull
Scaler
2513



Large
Plugger/Picker
2514



Picker/
Ligature
Director
2515



Scaler/
Picker
2516




Dull
Scaler
2517



Ligature
Director/
Band Pusher
2518


ORTHODONTICS

HAND INSTRUMENTS



Band Pusher
with Short
Tip

2519



Band Pusher

2520



Fillion Explorer/
Ligature
Director

2521



Band Pusher/
Scaler

2522



Band Pusher/
Scaler

2523



Ligature
Director

2524




Ligature
Director

2525

ORTHODONTICS

HAND INSTRUMENTS



Lingual
Torquing
Key

2526

Straight
.020"
Applicator

2527

Single End
Ligature
Director

2528

Lingual
Ligature
Director,
.018"

2529

Lingual
Ligature
Director,
.025"

2530

Ligature
Director/
Plugger

2531

Ligature
Director
Tucker


2532

Plugger/
Ligature
Picker

2533

ORTHODONTICS

TWEEZERS-SCISSORS



Posterior Bracket Holder Tweezer, Long

2566



Direct Bond Bracket Tweezer with Aligner & Positioner

2567




Buccal Tube Tweezer

2568



Slim Buccal Tube Bonder Tweezer

2569



Direct Bond Bracket Tweezer, Lingual

2570



Direct Bond Tweezer, Aligner

2571



All Purposes Scissors

2572



Crown Scissors

2573



Crown Scissors, Straight

2574



Crown Scissors, Curved

2575

HEIGHT GAUGES



BONE BRACKET POSITIONING GAUGE

2580



BRACKET HEIGHT GAUGE FOR
0.022"/0.56mm Brackets

2581



BRACKET HEIGHT GAUGE FOR
0.018"/0.46mm Brackets

2582



BRACKET HEIGHT GAUGE, MOVEABLE HEAD

2583

MATHIEU NEEDLE HOLDERS



MATHIEU NARROW TIP
Enhanced double spring function

2585



MATHIEU SMAHA FINE TIP
Enhanced double spring function

2586



MATHIEU NARROW TIP, TC
Enhanced double spring function

2587



MATHIEU STANDARD TIP
Enhanced double spring function

2588



MATHIEU STANDARD TIP, TC
Enhanced double spring function

2589