

EZ-LUXATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.



Spade Elevator, Straight
SELSS
 4180




Spade Elevator, Mesially Angled
SELSMA
 4181



Spade Elevator, Distally Angled
SELSDA
 4182



Spade Elevator, Left
SELSL
 4183



Spade Elevator, Right
SELSR
 4184



2.5mm, Small, Straight
SELS2.5
 4185



2.5mm, Small, Curved
SELSC2.5
 4186



2.5mm, Small, Mesially Angled
SELSMA2.5
 4187



2.5mm, Small, Distally Angled
SELSDA2.5
 4188

EZ-LUXATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.



3.75mm, Large, Straight

SELLS3.75
4190



3.75mm, Large, Curved

SELLC3.75
4191



3.75mm, Large, Mesially Angled

SELLMA3.75
4192



3.75mm, Large, Distally Angled

SELLDA3.75
4193



4.5mm, Large, Straight

SELLS4.5
4194



4.5mm, Large, Curved

SELLC4.5
4195



4.5mm, Large, Mesially Angled

SELLMA4.5
4196



4.5mm, Large, Distally Angled

SELLDA4.5
4197



Small, Tight, Mesial Angle, 2.5mm

SELSTMA2.5
4198



Tight Mesial Spade, 2.5mm

SELTMS2.5
4199

BERNARD & LINDO



Lindo Levin, Straight
3mm
LLS3
4204

Lindo Levin, Straight
4mm
LLS4
4205



11 Straight Bernard
SBERS11
BS11
4207



6 Left Bernard
SBERL6
BR6
4208



7 Right Bernard
SBERR7
BL7
4209



Lindo Levin, Curved
3mm
LLC3
4206



Bernard Spade, Offset IN
SBERO-IN
BSOIN
4210



Bernard Spade, Offset OUT
SBERO-OUT
BSOOT
4211

SURGICAL

HEIDBRINK

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Heidbrink



1
Heidbrink
Straight
E1
4213



2
Heidbrink
Straight
E2
4214



46
E46
4217



46R
E46R
4218



40
Heidbrink
Straight
E40
4215



41
Heidbrink
Straight
E41
4216



77
E77
4219



77R
E77R
4220

CRYER

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Suggested Pair.



21 Cryer Mini

22 Cryer Mini

EM21
4221

EM22
4222

Cryer:

Cryer's right and left elevators with triangular sharp-pointed blades and hollow handles.



44
Cryer Small

E44
4223



45
Cryer Small

E45
4224



25
Cryer Medium

E25
4225



26
Cryer Medium

E26
4226



A30
Cryer

EA30
4227



A31
Cryer

EA31
4228



31
Cryer Large

EA31
4229



32
Cryer Large

EA32
4230



LUXATING ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Luxating Elevator: The luxating elevators have thin, sharp blades that are used to cut periodontal ligaments. These instruments are used by rocking and cutting verses and elevator that lifts and pries. The importance of luxating elevator is to reduce pressure on adjacent tooth.

LUXATOR ELEVATORS



Luxating Elevator
Straight, 2mm

SLES2
4231



Luxating Elevator
Curved, 2mm

SLEC2
4232



Luxating Elevator
Straight, 2.5mm

SLES2.5
4233



Luxating Elevator
Curved, 2.5mm

SLEC2.5
4234



Luxating Elevator
Straight, 3mm

SLES3
4235



Luxating Elevator
Curved, 3mm

SLEC3
4236



Luxating Elevator
Straight, 4mm

SLES4
4237



Luxating Elevator
Curved, 4mm

SLEC4
4238



Luxating Elevator
Straight, 5mm

SLES5
4239



Luxating Elevator
Curved, 5mm

SLEC5
4240

SELDIN ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Seldin



APEXO ELEVATORS

Apexo



71 Miller
Apexo

E71
4252

72 Miller
Apexo

E72
4253

73 Miller
Apexo

E73
4254

74 Miller
Apexo

E74
4255



4
Schmeckebeer
Apexo

E4
4256

5
Schmeckebeer
Apexo

E5
4257

Cogswell & Crane Pick



8
Crane Pick

E8
4261

A
Cogswell

EA
4262

B
Cogswell

EB
4263



301
Apexo

E301
4258

302
Apexo


E302
4259

303
Apexo

E303
4260

COUPLAND & MACMILLAN

Coupland Gouge



 1C Coupland Gouge EGIC 4264	 2C Coupland Gouge EG2C 4265	 3C Coupland Gouge EG3C 4266
 A1C Coupland Gouge EGA1C 4267	 A2C Coupland Gouge EGA2C 4268	 A3C Coupland Gouge EGA3C 4269

Macmillan



 11M MacMillan Gouge EG11M 42611	 12M MacMillan Gouge EG12M 42621	 13M MacMillan Gouge EG13M 42631
--	--	--

SURGICAL

MISCELLANEOUS ELEVATORS



19

E19
4270

20

E20
4271

23

E23
4272

Friedman



31F
Friedman
[Craig]
Gouge Serrated

E31F
4278

32F
Friedman
[Craig]
Gouge Serrated

E32F
4279

Curtis



51 Curtis

EC51
4280

52 Curtis

EC52
4281



38
Serrated

E38
4273

39
Serrated

E38
4274

Bein



559/1
Bein, Large

EBE559/1
4275

559/2
Bein, Medium

EBE559/2
4276

559/3
Bein, Small

EBE559/3
4277

Hourigan



Suggested
Pair

H5
Hourigan, Serrated

H5S
4282

H6
Hourigan, Serrated

H6S
4283

SURGICAL

MISCELLANEOUS ELEVATORS

		Woodward		Stout			
		Suggested Pair					
							
6X Potts	7X Potts	15 Woodward	16 Woodward	11 Stout	11A Stout	80	81
E6X 4284	E7X 4285	EWW15 4286	EWW16 4287	E11 4288	E11A 4289	E80 4290	E81 4291

						
81A	81E	92 Serrated	93 Serrated	98	190	191
E81A 4292	E81E 4293	E92 4294	E93 4295	E98 4296	E190 4297	E191 4298

SURGICAL

T-BAR HANDLE ELEVATORS




Warwick-James Straight **EWJS** Handle #504
 Warwick-James Left **EWJL** Handle #504
 Warwick-James Right **EWJR** Handle #504

SURGICAL

ROOT TIP PICKS ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.




SURGICAL

ROOT TIP PICKS ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

West



Howard

Slightly heavier design than Heidbrink picks.




SURGICAL

ROOT TIP PICKS ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Davis Teasers



PERIOTOME ELEVATORS

Thin sharp blades used to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament. May be used wherever an extraction is indicated, but are particularly important when considering the placement of a dental implant.



1 Posterior
SPT1
4347



3 Posterior-Right
SPTR3
4348



4 Posterior-Left
SPTL4
4349



2 Anterior
SAT2
4350



5 Anterior-Angled
SATA5
4351



6 Anterior-Straight
SATS6
4352



Wise
SPW
4353

SURGICAL

PERIOSTEALS


These instruments are used for retracting periosteum after removing gingival tissue.



SURGICAL

RERIOSTEALS

These instruments are used for retracting periosteum after removing gingival tissue.



SURGICAL

RERIOSTEALS

These instruments are used for retracting periosteum after removing gingival tissue.



24G

SPS24G

4 4 1 7

Tin Coating
24G

SPS24GT

4 4 1 8



Benque

SPSBEN

4 4 1 9



Buser

SPSBUS

4 4 2 0

Tin Coating
Buser

SPSBUST

4 4 2 1



KN-1

SPSKN1

4 4 2 2



KN-152

SPSKN152

4 4 2 3



Hourigan-2

SPSHOURI

4 4 2 4

SURGICAL

RERIOSTEALS

These instruments are used for retracting periosteum after removing gingival tissue.



HF-20


SPSHF20

4 4 2 5

Tin Coating
HF-20

SPSHF20T

4 4 2 6



Molt-5

SPSM5

4 4 2 7



Selden-23

SPSSEL23


4 4 2 8



Freer

SPSFR

4 4 2 9

OHL-10
Periosteal
Elevator

SPSOHL10

4 4 3 0

Woodson-1
Modify

SPSWSONM

4 4 3 1

Allen Periosteal
Elevator

SPPAEL


4432

SURGICAL

SURGICAL CURETTES (SPOON)

These instruments are used for curettage and cyst removal and debridement of tooth socket.

Lucas-Curettes




SURGICAL

SURGICAL CURETTES (SPOON)

These instruments are used for curettage and cyst removal and debridement of tooth socket.

Miller Currettes



Molt-Currettes



**BONE CHISELS
SINGLE BEVEL**



2 Gardner
4mm
SBCSE2G4
4541



5 Gardner
4.3mm
SBCSE5G4.3
4542



6 Gardner
4mm
SBCSE6G4
4543



9 Gardner
4mm
SBCSE9G4
4544



1 Chandler
5mm
SBCSE1CH5
4545



2 Chandler
5mm
SBCSE2CH5
4546



4 Chandler
5mm
SBCSE4CH5
4547



DOUBLE BEVEL

1 Gardner
3mm
SBCSE1G3
4540



7 Chandler
4.7mm
SBCDB7CH4.7
4548



8 Chandler
6.1mm
SBCDB8CH6.1
4549



9 Chandler
7.1mm
SBCDB9CH7.1
4550



1W Wakefield
2.7mm
SBCDB1W2.7
4553



2W Wakefield
3mm
SBCDB2W3
4554

4.2mm
SBCDB18M4.2
4555

SURGICAL

BONE FILES

Bone Files: For final smoothing of bone by pushing out.

Miller Collburn-Cross Cut

The small paddle shaped end is half the size of the larger paddle shaped end.



1X
Miller-Collburn
Cross-Cut

SBF1X

4 6 0 0

Curved working end combined with a straight working end.



2X
Miller-Collburn
Cross-Cut

SBF2X

4 6 0 1

Curved working end combined with a small paddle shaped working end



3X
Miller-Collburn
Cross-Cut

SBF3X

4 6 0 2

Miller-Straight Cut

Same design as the 1X Cross-Cut but with Straight-Cut.



21
Miller

SBF21M

4 6 0 3

Same design as the 2X Cross-Cut but with Straight-Cut.



45
Miller

SBF45M

4 6 0 4

SURGICAL

BONE FILES

Bone Files: For final smoothing of bone by pushing out.

Miller-Straight Cut



52
Miller

SBF52M

4 6 0 5



64
Miller

SBF64M

4 6 0 6

Wahl-Straight Cut



1
Wahl

SBF1W

4 6 0 7



2
Wahl

SBF2W

4 6 0 8

Howard-Straight Cut



12
Howard

SBF12H

4 6 0 9



67
Howard

SBF67H

4 6 1 0

SURGICAL


SCISSORS

Used for trimming tissue or cutting sutures.
All Scissors shown 50% size.

Wagner



Kelly



SURGICAL

SCISSORS

Used for trimming tissue or cutting sutures.
All Scissors shown 50% size.



SURGICAL

SCISSORS

Used for trimming tissue or cutting sutures.
All Scissors shown 50% size.



31 Castroviejo Curved
10cm/4"

SS31CASC
4740



17 Iris Straight/Delicate
11.5cm/4½"

SS17ISD
4741



18 Iris Curved/Delicate
11.5cm/4½"

SS18ICD
4742



Metzenbaum Baby
12.5cm/5"

SSMB
4743

SURGICAL

SCISSORS

Used for trimming tissue or cutting sutures.
All Scissors shown 50% size.


General Surgical Scissors



21
Straight/
Pointed
14.5cm/5³/₄"

SS21SP

4 7 4 4



22
Straight/
Blunt-Sharp
14.5cm/5³/₄"

SS22SBS

4 7 4 5



23
Curved/
Pointed
14.5cm/5³/₄"

SS24CP

4 7 4 6



Lister
Bandage
14.5cm/5³/₄"

SSLB

4 7 4 7

Suture Scissors



13
Suture
15cm/6"

SS13S

4 7 4 8



13S
Suture
15cm/4³/₄"

SS13SS

4 7 4 9



Crown &
Gold
Straight
15cm/4³/₄"

SSC&GS

4 7 5 0



Crown &
Gold
Curved
15cm/4³/₄"

SSC&GC

4 7 5 1



Crown &
Gold Curved
Pedodontic
15cm/3.5"

SSC&GP

4 7 5 2

Crown & Gold Remover

SURGICAL

TUNGSTON CARBIDE TIPPED SCISSORS

Super-Cut Scissors: The razor-sharp cutting edge cuts effortlessly and softly also through close material. The serrated cutting edge on the second blade prevents a slipping from tissue and sutures. Super-Cut is not practical for wires, cartilages and other hard material. As identification, one handle ring is gold plated.

All Scissors shown 50% size.

Super cut scissors



10
Double-Curved
13.5cm/5 1/4"

SS10DC

4800



14
Double-Curved
11.5cm/5 1/4"

SS14DC

4801



16
Goldman-Fox
12.5cm/5"

SS16GFTC

4802



17
Iris
Straight/Delicate
11.5cm/4 1/2"

SS17ISDTC

4803



18
Iris
Curved/Delicate
11.5cm/4 1/2"

SS18ICDTC

4804



Metzenbaum
Straight/Blunt
14.5cm/5 1/2"

SSMSBTC

4805



Metzenbaum
Curved/Blunt
14.5cm/5 1/2"

SSMCBTC

4806



Metzenbaum
Straight/
Pointed
14.5cm/5 1/2"

SSMSPTC

4807



Metzenbaum
Curved/
Pointed
14.5cm/5 1/2"

SSMCPTC

4808



Metzenbaum
Delicate-
Curved/
Blunt
14.5cm/5 1/2"

SSMDCBTC

4809

SURGICAL

SCISSORS

All Scissors shown 50% size.



SURGICAL

TISSUE PLIERS AND FORCEPS

These pliers are used to grasp and stabilize soft tissue flaps during suturing.

All Tissue pliers and forceps shown 50% size.

Semkin-Taylor

<p>31 Semkin-Taylor Straight 12.5cm/5"</p> <p>TP31STS 4870</p>	<p>32 Semkin-Taylor Curved 12.5cm/5"</p> <p>TP32STC 4871</p>	<p>33 Semkin-Taylor 1x2 Straight 12.5cm/5"</p> <p>TP33STS 4872</p>	<p>34 Semkin-Taylor 1x2 Curved 12.5cm/5"</p> <p>TP34STC 4873</p>
---	---	---	---

<p>16 3x4 15cm/6"</p> <p>TP16 4874</p>	<p>50-L Fox Locking 3x3 15cm/6"</p> <p>TP50L 4875</p>	<p>Kramer-Nevins Tissue Pliers 1x2 Curved</p> <p>TPKN 4876</p>
---	--	---

Suture Pliers

<p>20 Corn Suture Pliers 15cm/6"</p> <p>TP20C 4877</p>

<p>41 Adson Plain 12cm/4 3/4"</p> <p>TPAP41 4878</p>	<p>42 Adson 1x2 12cm/4 3/4"</p> <p>TPA42 4879</p>	<p>43 Adson-Brown 7x7 12cm/4 3/4"</p> <p>TPAB43 4880</p>	<p>Russian 15cm/6"</p> <p>TPRN15 4881</p>	<p>Russian 20cm/8"</p> <p>TPRN20 4882</p>	<p>Gerald Plain 18cm/7" Straight/Serrated</p> <p>TPGP18 4883</p>	<p>Gerald 1x2 15cm/7" Straight</p> <p>TPG15 4884</p>
---	--	---	--	--	---	---

Allison

<p>Allison Baby 4x5 12.5cm/5"</p> <p>FALIB12.5 4885</p>	<p>Allison Straight 5x6 15cm/6"</p> <p>FALIS15 4886</p>
--	--

SURGICAL

TISSUE PLIERS / FORCEPS

All Tissue Pliers and Forceps shown 50% size.

Adson



41

Adson
12cm/4³/₄"

TPAD12TC
4 8 8 7



42

Adson 1x2
12cm/4³/₄"

TPAD1X2TC
4 8 8 8



Adson
Brown
12cm/4³/₄"

TPADB12TC
4 8 8 9

General



General
15cm/6"

TPG15TC
4 8 9 0



General
1x2
15cm/6"

TPG1X2TC
4 8 9 1

Gerald



Gerald
18cm/7"

TPGRD18TC
4 8 9 2



Gerald 1x2
18cm/7"

TPGRD1X2TC
4 8 9 3



Gushing
18cm/7"

TPGUSHTC
4 8 9 4

Cushing



Cushing
Angled
18cm/7"

SLB
4 8 9 5



Cushing
Bayone
19cm/7¹/₂"

SLB
4 8 9 6



Potts Smith
18cm/7"
(Shown)

SLB
4 8 9 7

Potts Smith
14.5cm/5-1/4"

SLB
4 8 9 8



De Bakay
15cm/6"
(Shown)

SLB
4 8 9 9

De Bakay
19cm/7¹/₂"

SLB
4 9 0 0

SURGICAL

NEEDLE HOLDERS

Martensite stainless steel instruments are used for grasp and guide the needle while suturing.

All Needle Holders shown 50% size.



Hegar-Baumgartner (Gardner 5)
14cm/5 1/2"
NHHB
4950

Webster Smooth
13cm/5"
NHWS13
4951

Mayo-Hegar (Gardner 6) Grooved
16cm/6 1/4"
NHMH
4952

Crile-Wood Grooved
15cm/6"
NHCW
4953

Olsen-Hegar Needle Holder/Scissors
17cm/6 3/4"
NHOHHS
4954

Boynton (Hu-Friedy)
12cm/4 3/4"
NHBOY
4955

Castroviejo
14cm/5 1/2"
NHCV
4956

Collier
12.5cm/5"
NHC
4957

Derf Grooved
12.5cm/5"
NHDG
4958

Mathieu (Machatt)
14cm/5 1/2"
NHMM
4959

SURGICAL

TC-NEEDLE HOLDERS

Used for grasp and guide the needle when suturing.
All Needle Holders shown 50% size.



Derf
11.5cm/4 1/2"

NHDTTC
4960



Hegar-
Baumgartner
11.5cm/4 1/2"

NHHBTC
4961

Halsey



Hasley
Smooth
13cm/5"

NHHSTC
4962



Halsey
13cm/5"

NHHTC
4963

FLAT HANDLE

Castroviejo TC Needle Holders



Castroviejo
Straight
Flat Handle
13cm/5.15"

NHCVSFH13TC

4970



Castroviejo
Curved
Flat Handle
13cm/5.15"

NHCVFHC13TC

4971



Castroviejo
Straight
Flat Handle
14cm/5.5"

NHCVSFH14TC

4972



Castroviejo
Curved
Flat Handle
14cm/5.5"

NHCVFHC14TC

4973



Castroviejo
Straight
Flat Handle NH
15cm/5.9"

NHCVFHS15TC

4974



Castroviejo
Curved
Flat Handle NH
15cm/5.9"

NHCVFHC15TC

4975



Castroviejo
Straight
Flat Handle
16cm/6.3"

NHCVFHS16TC

4976



Castroviejo
Curved
Flat Handle
16cm/6.3"

NHCVFHC16TC

4977



Castroviejo
Straight
Flat Handle
18cm/7"

NHCVFHS18TC

4978



Castroviejo
Curved
Flat Handle
18cm/7"

NHCVFHC18TC

4979

SURGICAL

ROUND HANDLE Castroviejo TC Needle Holders



Castroviejo
Straight
Round Handle
13cm/5"
NHCVSRH13TC
4982



Castroviejo
Curved
Round handle
13cm/5"
NHCVCRH13TC
4983



Castroviejo
Straight
Round Handle
14cm/5.5"
NHCVSRH14TC
4984



Castroviejo
Curved
Round Handle
14cm/5.5"
NHCVCRH14TC
4985

4987



Castroviejo
NH, ST
Round Handle
15cm/5.9"
NHCVSRH15TC
4986



Castroviejo
NH, CV
Round Handle
15cm/5.9"
NHCVCRH15TC
4987



Castroviejo
Straight
Round Handle
16cm/1.6"
NHCVSRH16TC
4988



Castroviejo
Curved
Round Handle
16cm/1.6"
NHCVCRH16TC
4989



Castroviejo
Straight
Round Handle
18cm/7"
NHCVSRH18TC
4990




Castroviejo
Curved
Round Handle
18cm/7"
NHCVCRH18TC
4991

SURGICAL

NEEDLE HOLDERS

All Needle Holders shown 50% size.



Arruga
Serrated
16cm/6 1/2"
NHARRUS16TC
4995



Webster
Smooth
11.5cm/4 1/2"
NHWEBS11.5TC
4996



Crile-Wood
15cm/6"
NHCW15TC
4997



Mayo-Hegar
16cm/6 1/4"
NHMH16TC
5001

Olsen-Hegar (Needle Holder)



Olsen-Hegar
Needle Holder
Scissors
17cm/6 3/4"
NHSH17TC
4998



Olsen-Hegar
Needle Holder
Scissors
14cm/5 3/4"
NHSH14TC
4999



DeBakey
18cm/7"
Delicate
NHDEB18TC
5000



Mayo-Hegar
18cm/7"
NHMH18TC
5002

SURGICAL

NEEDLE HOLDERS

Used for grasp and guide the needle while suturing.
All Needle Holders shown 50% size.



Matheu-Kocher
14cm/5"

NHMK14TC

5005



Matheu
14cm/5 1/2"

NHM14TC

5006



HMT
Thin Jaw
14cm/5 1/2"

NHHMT14TC

5006



Fine Swedish
Pattern
15cm/6"

NHFSP15TC

5007



Micro-Vascular
Pattern
15cm/6"

NHMVP15TC

5008



Micro-Vascular
Pattern
18cm/7"

NHMVP18TC

5009



Mini-Ryder
Pattern
15cm/6"

NHMRP15TC

5010



Mini-Ryder
Pattern
18cm/7"

NHMRP18TC

5011

SURGICAL

HEMOSTATIC FORCEPS

A multipurpose instrument used to clamp off blood vessels, remove small root tips and grasp loose objects. All Hemostats

Hartman



Hartman
Straight
10cm/4"

HHS10
5 0 3 0



Hartman
Curved
10cm/4"


HHC10
5 0 3 1



5
1x2 Teeth
Straight
12.5cm/5"

H5S1X2
5 0 3 4

Adson




Adson
Straight
19cm/7"

HAS19
5 0 3 2



Adson
Curved
19cm/7-1/2"

HAC19
5 0 3 3



Rochester-pean
Curved
16cm/6"

HRPC16
5 0 3 5

SURGICAL

HEMOSTATIC FORCEPS

A multipurpose instrument used to clamp off blood vessels, remove small root tips and grasp loose objects. All Hemostats

Halsted-Mosquito



Kelly



SURGICAL

CHEEK RETRACTORS/TISSUE RETRACTORS



Middeldorf
Retractor, 17x14mm
SMID17x14
5050



Middeldorf
Retractor, 22x20mm
SMID22x20
5051



9 Ramus
Retractor, 16x80mm
SR9
5052



8 Surgical
Retractors, 11x70mm
SR8
5053



1 Surgical
Retractor, 7x25mm
SR17x25
5059



1 Surgical
Retractor, 10x42mm
SR110x42
5060



1 Surgical
Retractors, 12x55mm
SR1-12x55
5061



1 Austin
S1A
5066



2 Austin
S2A
5067



1 Henahan
SR1HENN
5068



2 Henahan
SR2HENN
5069




3 Henahan
SR3HENN
5070

SURGICAL

CHEEK RETRACTORS/TISSUE RETRACTORS

Retractors: For holding back the lips and cheeks, and tongue.


Mouth Gag: To prop open the mouth and prevent it from closing.



SURGICAL

CHEEK RETRACTORS/TISSUE RETRACTORS

Retractors: For holding back the lips and cheeks, and tongue.



Skin Hooks: For hooking skin during surgery.
Retractors: For holding back the lips and cheeks.



SURGICAL

SCALPEL HANDLES

Scalpel Handle: For holding a mess blade.



SURGICAL

ASPIRATORS

These tools are used for precision section at surgical sites.

Coupland Handle
25cm/9.9"

SASPCH
5249



Actual tip size (Internal diameter)

5250	SASPCT1	<input type="radio"/>	1 Coupland Tip	1.5mm
5251	SASPCT2	<input type="radio"/>	2 Coupland Tip	2.0mm
5252	SASPCT3	<input type="radio"/>	3 Coupland Tip	2.5mm
5253	SASPCT4	<input type="radio"/>	4 Coupland Tip	3.0mm

Actual tip size (Internal diameter)

5254	SASPCT1SH	<input type="radio"/>	1 Coupland Tip/Side Holes	1.5mm
5255	SASPCT2SH	<input type="radio"/>	2 Coupland Tip/Side Holes	2.0mm
5256	SASPCT3SH	<input type="radio"/>	3 Coupland Tip/Side Holes	2.5mm
5257	SASPCT4SH	<input type="radio"/>	4 Coupland Tip/Side Holes	3.0mm

Shown 50% size.

BYRD



Byrd
Self-Cleaning
20cm/8"

Actual tip size (Internal diameter)

5259	SASPB1	<input type="radio"/>	B1 Byrd	1.5mm
5260	SASPB2	<input type="radio"/>	B2 Byrd	2.0mm
5261	SASPB3	<input type="radio"/>	B3 Byrd	2.5mm

Cogswell 3
Suction Tube
23cm/9"

SASPCG3
5262



Shown 50% size.

Frazier
Suction Tube

Actual tip size (Internal diameter)

5263	SASPF6	<input type="radio"/>	6 Frazier	1.5mm
5264	SASPF8	<input type="radio"/>	8 Frazier	2.0mm
5265	SASPF10	<input type="radio"/>	10 Frazier	2.5mm
5266	SASPF12	<input type="radio"/>	12 Frazier	3.0mm

Used in oral surgery



1mm

2mm

Shown 50% size. 3mm

4mm

SASP1

SASP2

SASP3

SASP4

5270

5271

5272

5273

Yankauer
26.5cm/10.5"

SASPYK

5275

Shown 50% size.



Andrew-Pychon
23cm/9"

SASPAP

5276

Shown 50% size.



Locking Pliers
16.5cm/6 1/2"

LP

5280



SURGICAL

SYRINGES

Syringes: An instrument for injecting aesthetics for a cartridge.
 3-Way Syringes - Supplies air, water, or a mist through separate interval channels created by combining pressurized air with water flow.



Intraligamentary
Syringe
1.8ml

SY11.8ml
5320



Finger Loop Style

Syringe
1.8ml
SY1.8ml
5322



Finger Rest Style

Syringe
1.8ml
SY111.8ml
5323



Syringe
1.8ml
SYS1.8
5324

Cook-Waite Style



Paroject

SYPARS
5328



Pen Style with Blue Handle

SYPARP
5329

SURGICAL

MISCELLANEOUS

Miscellaneous tools for laboratory.



Stainless Steel 6" Rule

SS6R
5386



Tofflemire
Junior

TJ
5383




Tofflemire
Senior

TS
5384



Tofflemire
Universal

TU
5385



Boley
Gauge

BG
5387

SURGICAL

ATRAUMATIC EXTRACTING FORCEPS SET

Our Unique tip design makes it mold to the individual shape of your patient's tooth. The smooth handle with openings, allows for more thorough and easier cleaning. The newly styled handle allows a more uniform transmission while applying pressure, pulling and for rotary movements. More comfortable design can be used by right or left-handed practitioners! Our handles are lighter and are easier to use.



F1

UPPER
Universal

EF1

5430



F2

LOWER
Universal

EF2

5431



F3

Upper
Anterior

EF3

5432



F5

Lower
Anterior

EF4N

5433



F6

Upper
3rd Molar

EF6

5434



F7C

Lower Molar, modified
Cow-Horn Bead

EF7C

5435



46XL

Lower Root
Extra Long

E46XL

5436



51XL

Upper Root

E51XL

5437



222

Lower
Wisdoms

E222

5438



301

Lower
Roots

E301

5439



300

Upper Roots

E30

5440

SURGICAL

EXTRACTION FORCEPS AMERICAN PATTERN

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



5
Lower Third
Molars
E5
5441

15
Lower
Molars
E15
5442

18L
LOWER MOLARS
LEFT SIDE
E18L
5447

18R
Upper Molars
Right Side
E18R
5448

16S
Lower Molars
E16S
5443

17
Lower Molars
E17
5444

6
Splitting
Forceps
E6
5449

1
Upper Centrals
Laterals &
Cuspids
E1
5450

16
Lower Molars
Cow-Horn Beak
E16
5445

1AF
Upper
Anteriors
E1AE
5446

SURGICAL

EXTRACTION FORCEPS AMERICAN PATTERN

Forceps: For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



99C
Upper
Premolars
& Incisors

E99C


5451



10S
Upper Molars
Straight-Handle

E10S

5454



53L
Upper Molars
Left Side

E53L

5455



53R
Upper Molars
Right Side

E53R

5456



23
Lower Molars
Cow-Horn Beak

E23

5457



10H
UPPER MOLARS
HOOK HANDLE

E10H

5458



69
Fragments
& Small Roots
Upper & Lower

E69

5459



65
Upper Incisors
& Roots

E65

5460



Universal
for Upper
Teeth & Roots

EU

5461



32
Upper Molars
& Premolars
Strong Beak

E32

5462

SURGICAL

EXTRACTION FORCEPS AMERICAN PATTERN

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



151A
Lower
Premolars
Incisors & Roots
E151A
5463



151XAS Children
Lower Premolars
Incisors &
Roots Canines
E151XAS
5464



151AS
Split Lower
Premolars
Incisors &
Roots Canines
E151AS
5465



151AF Apical
Beak Lower
Premolars
Incisors &
Roots Canines
E151AE
5466



101
Universal
for Premolars
E101
5467



150
Upper Premolars
Incisors &
Roots Canines
E150
5468



150AF
Apical Beak
Upper Premolars
Incisors &
Roots Canines
150AF
5469



150AS
Split Beak
Upper Premolars
Incisors &
Roots Canines
E150AS
5470



150XAS
Children Split
Upper Premolars
Incisors &
Roots Canines
E150XAS
5471



150A
Upper PreMolars
Incisors &
Roots Canines
E150A
5472



150S
Upper Premolars
Incisors &
Roots for
Children
E150S
5473



151
Upper Premolars
Incisors &
Roots
E151
5474

SURGICAL

EXTRACTION FORCEPS AMERICAN PATTERN

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



151S
Lower Premolars
Incisors &
Roots for Children

E151S

5 4 7 5



210H
Upper Wisdoms
with Hook Handle

E210H

5 4 7 6



210S Upper
Wisdoms
with Straight
Handle

E210S

5 4 7 7



222
Lower Wisdoms

E222

5 4 7 8



222AF
Apical Beak
Lower Wisdoms

E222AF

5 4 7 9



286
Upper Roots
Premolars
& Incisors

E286

5 4 8 0

SURGICAL

Extracting Forceps European Style

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



30
Upper Roots
& Incisors

E30

5481



22
Lower Molars
Hawks Bill

E22

5482



21
Lower Molars

E21

5483



20
Lower Wisdoms

E20

5484



8
Lower
Premolars

E8


5485



7
Upper
Premolars

E7


5486



3
Upper Incisors
& Canines

E3

5487



1
Upper Central
& Canines

E1


5488



19
Upper Wisdoms

E19

5489



18A
Upper Molars

E18A

5490



18
Upper Molars
Left

E18

5491



17
Upper Molars
Right

E17

5492


SURGICAL

Extracting Forceps European Style

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



67 AF
Apical Beak
Upper Wisdoms
E67AF
5493



68
Lower Roots
E68
5494




73
Lower Molars
(Hawks Bill)
E73
5495



74
Lower Incisors
& Roots
E74
5496



31
Lower Roots
E32
5497



32
Lower Molars
E31
5498



33
Lower Roots
E33
5499



33M
Lower Roots
With Long Beak
E33M
5500



67
Upper Wisdoms
E67
5501



51A
Upper Roots
Narrow Beak
E51A
5502



45
Lower Roots
E45
5503



44
Lower Roots
E44
5504

SURGICAL

Extracting Forceps European Style

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



95
Upper Molars
Left

E95

5505



121
Upper Wisdoms

E121


5506



167 Universal
Patter for
Lower Teeth
Molars & Premolars

E167

5507



168
Upper
Canines

E168

5508



74N
Lower Incisors
& Roots
Canines

E74N

5509



74NF
Upper Incisors
& Roots

E74NF


5510



79
Lower Wisdoms

E79

5511



79A
Lower Wisdoms

E79A

5512



86C
Lower Molars

E86C

5513



89
Upper Molars
Right

E89

5514



90
Upper Molars
Left

E90

5515



94
Upper Molars
Right

E94

5516

SURGICAL

Extracting Forceps European Style

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



F-1
Upper
Universal
EF-1
5517



F-2
Lower
Universal
EE-2
5518



F-3
Upper Anterior
EF-3
5519



F-4N
Lower Anterior
EF-4N
5520



F-4
Lower Anterior
EF-4
5521



F-5
Lower Anterior
Parallel Beaks
EF-5
5522



F-6
Upper
3rd Molars
EF-6
5523



F-7
Lower
3rd Molars
EF-7
5524



F-7C
Lower Molar
Modify
Cow-Horn Beak
EF-7C
5525



F-8
Upper Universal
EF-8
5526



F-9
Lower Universal
EF-9
5527



F-10
Lower Molar
EF-10
5528


SURGICAL

Extracting Forceps European Style

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



F11
Upper Molar
Left
EF-11
5529



F12
Upper Molar
Right
EF-12
5530



F13
Upper Molar
Right
EF-13
5531



F14
Upper Molar
Left
EF-14
5532



169
Lower Cnines
E169
5533



233
Lower Roots
Parallel Beaks
E233
5534



235
Upper Roots
Parallel Beaks
E235
5535




MD-1
Universal
Upper Wisdom
Premolars &
Roots
EMD-1
5536



MD-4
Lower Molars
EMD-4
5537



22L
Lower Wisdoms
Left
E22L
5538



22R
Lower Wisdoms
Right
E22R
5539

SURGICAL

Extracting Forceps European Style

Forceps: For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

Root Splinter And Extraction Tweezer: For removing the residual part of tooth root.

SPONGE & DRESSING/HEMOSTATES: For grasping off blood vessels and removing small pieces of root tips.



F15
Upper
Premolars
Universal

EF-15

5540



F16
Upper
Root

EF-16

5541



F17
Lower Root

EF-17

5542



F18
Upper Roots
Fragment Narrow

EF-18

5543



F 18TN
Upper Roots
Fragment Narrow

EF-18TN

5544



F 19
Lower Roots
Fragment Narrow

EF-19

5545



F 19NW
Lower Roots
Fragment Narrow

EF-19NW

5546



46XL
Lower Root
Extra Long Beaks

EF46XL

5547



F 20
Upper
Premolar
1ST & 2RD Molars

EF-20

5548



49
Upper Roots

EF49

5549



51XL
Upper Roots

EF51XL

5550

Atraumatic Forceps Children's / Pedodontic Pattern



F1
Upper Universal
Pedo

5551



F2
Lower Universal
Pedo

5552



F3
Upper Anterior
Pedo

5553



F4
Lower Anterior
Pedo

5554



F5
Lower Anterior
Pedo

5555



F6
Upper 3rd
Pedo

5556



F7
Lower 3rd
Pedo

5557



F16
Lower Root
Pedo

5558



F17
Lower Root
Pedo

5559



B23
Lower Molar
Pedo

5560



150-SK
Upper Universal
Pedo

E150-SK

5561



151-SK
Lower Universal
Pedo

E151-SK

5562

RONGEURS

FRIEDMAN



Friedman
14cm/5-1/2"

RFM

6000

BLUMENTHAL



Blumenthal
Standard 30°
14.5cm/5-3/4"

RBIS

6003

CLEVELAND



1A Mead
16.5cm/6-1/2"

RIA

6006



3H Dean
17cm/6-3/4"

R3HD

6008



4 Cleveland
14cm/5-1/2"

R4C

6010



Mini Friedman
14cm/5-1/2"

RMFM

6001



Blumenthal 45°
14.5cm/5-3/4"

RBL45

6004



4A Cleveland
17cm (6-3/4")

R4AC

6007



5 Cleveland
Side Cutting
17cm (6-3/4")

R5C

6012



5S Cleveland
Side Cutting
15cm (6")

R5SC


6013



Mini Friedman 90°
14.5cm/5-1/2"

RMFM90

6002



Blumenthal 90°
14.5cm/5-3/4"

RBM90

6005



Kramer-Nevis
14.5cm (5-3/4")

RKN

6014

RONGEURS



Kerrison Rongeur 3mm
for widening the sinus window

RK3mm

6020



Kerrison Rongeur 4mm
for widening the sinus window

RK4mm

6022

DOUBLE ACTION HINGED RONGEURS



18M- (Double Action)
Side-Cutting
17mm (6-7/8")

R18MDA17

6025



Beyer- (Double Action)
18cm (7")

RBEYDA18

6030



Boehler (Double Action)
15cm (6")

RBOEDA15

6035



Boehler (Double Action)
16cm (6-1/4")

RBOEDA16

6040