


## Diagnostic Instruments



## Diagnostic Instruments Index

Endo Explorers				E1	
Endo Excavators				E2	
Endo Pluggers				E3-E6	
Endo Spreaders				E8-E9	
Endo Nickel Titanium spreaders				E7	
Endo Micro Curettes				E10	
Endo apical Pluggers				E11-E12	
Endo Apical Burnishers				E13-E14	
Endo Apical Pluggers & Burnishers Combo				E15-E16	
Splitting Forceps				E17	
Rubber Dam Clamp Frame				E17	
Ruber Dam Punches				E17	
Ruber dam Clamps				E18-E20	
Rubber Dam Holder				E20	

# DIAGNOSTIC MIRRORS

## Cone Socket Handles



#M1 Cone Socket Mirror Handle **M1CS** 901



#M2 Cone Socket Mirror Handle **M2CS** 902



#M3 Cone Socket Mirror Handle **M3CS** 903



#M4 Cone Socket Mirror Handle **M4CS** 904

## Simple Stem Handles



#M5 Simple Stem Mirror Handle **M5SS** 905



#M6 Simple Stem Mirror Handle **M6SS** 906



#M7 Simple Stem Mirror Handle **M7SS** 907



#M8 Simple Stem Mirror Handle **M8SS** 908



Simple Stem Mouth Mirror  
#4-7/8" (22mm)

**SSMM#4**

910



Simple Stem Mouth Mirror  
#5-15/16" (24mm)

**SSMM#5**

911



Simple Stem Mouth Mirror  
Double Sided  
#4-7/8" (22mm)

**SSMMDS#4**

913



Simple Stem Mouth Mirror  
Double Sided  
#5-15/16" (24mm)

**SSMMDS#5**

914



Cone Socket Mouth Mirror  
#4-7/8" (22mm)

**CSMM#4**

915



Cone Socket Mouth Mirror  
#5-15/16" (24mm)

**CSMM#5**

916

# DIAGNOSTIC

## MICRO SURGICAL MIRRORS



1  
Small  
Rectangular  
**MICMS1**  
1013



2  
Small  
Rectangular  
**MICMS2**  
1014




3  
Small  
Rectangular  
**MICMS3**  
1015



Small Round  
**MICMSR**  
1016



Large Round  
**MICMLR**  
1017

**DIAGNOSTIC****MICRO SURGICAL MIRRORS**

# DIAGNOSTIC

## DRESSING PLIERS

Used to grasp or transfer materials in and out of the oral cavity. All dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment. Stainless steel made.



1  
Perry  
**DP1**  
1071



2  
College  
**DP2**  
1072



3  
**DP3**  
1073



17  
**DP17**  
1074



18  
**DP18**  
1075

# DIAGNOSTIC

## DRESSING PLIERS

Used to grasp or transfer materials in and out of the oral cavity. All dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment. Stainless steel made.



18L  
Locking

**DP18L**

1076



19

**DP19**

1077



Utility  
Pick-Up  
Dressing Plier

**DPU17**

1078



317  
Smooth  
Beak

**DP317**

1079



318  
Foi/Meriam  
Serrated Beak


**DP318**

1080


# DIAGNOSTIC HANDLE OPTIONS

Use for caries and calculus detection, exploration of pocket characteristics, furcation and restorations.

## DOUBLE ENDED



## SINGLE ENDED



# DIAGNOSTIC EXPLORERS

tooth surface irregularities, defective margins on restorations, decalcified areas and carious lesions. For the quality assurance, we use the Germany made material. Moreover, explores are available in a variety of design types.

Single-End



# DIAGNOSTIC

# EXPLORERS



8  
**EXS8**  
1127

9  
**EXS9**  
1128

23  
Shepherd's  
Hook  
**EXS23**  
1129

54  
**EXS54**  
1130


TU17  
Wilkins/Tufts  
**EXSTU17**  
1131

17  
**EXS17**  
1132

# DIAGNOSTIC EXPLORERS

A dental explorer is an instrument in dentistry commonly used in the dental armamentarium. Use for caries and calculus detection, exploration of pocket characteristics, furcation and restorations.


Double-End



# DIAGNOSTIC EXPLORERS

A dental explorer is an instrument in dentistry commonly used in the dental armamentarium. Use for caries and calculus detection, exploration of pocket characteristics, furcation and restorations.

Double-End



# DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



# DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



UNC12  
University of  
North Carolina  
**CPUNC12**  
1207

UNC15  
University of  
North Carolina  
**CPUNC15**  
1208

PQ-O  
Michigan  
**CPQO**  
1209

PQ-W  
Williams  
**CPQW**  
1210

Moffitt/  
Maryland  
**CPMM**  
1211

Williams  
NovatechNT2  
**CPWNT2**  
1212

# DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



# DIAGNOSTIC COLOR PROBES

The calibrated periodontal probe is an instruments of assessment. It has a rod-shaped working end that is marked in millimeter units. The calibrated probe is actually similar to the ruler. It is used to... **1.** determine the consistency of the gingival tissue. **2.** Survey sulcus and pockets depths. **3.** Measure clinical attachment levels **4.** Measure the width of attached gingival. **5.** Verify whether bleeding and/or purulent exudate is in presence. These probes are made of spring stainless steel



# DIAGNOSTIC COLOR PROBES

## Double-Ended Screening Probes




# DIAGNOSTIC

# COLOR PROBES



# DIAGNOSTIC EX PROBES

**What is an Ex-Probe?** Ex-Probe is an excellent combination instrument for dental diagnostics. It has an tip on one end and a probe on the other.



10-9-8-7-5-3-2-1

23/QW  
Williams**XP23/QW**  
1290

10-9-8-7-5-3-2-1

23/QOW  
Thin  
Williams**XP23/TW**  
1291

11-8-6-3

23/CP-8

**XP23/CP8**  
1292

11-8-6-3

23/CP-11

**XP23/CP11**  
1293

11.5-8.5-5.5-3.5

23/CP-11.5B  
Screening  
Probe**XP23/CP11.5B**  
1294

12-9-6-3

23/CP-12


**XP23/CP12**  
1295

12-9-6-3

23/NT3CP12  
Novatech**XP23/CP12NT3**  
129612-11-10-9-  
8-7-6-5-4-3-2-123/UNC12  
University of  
North Carolina**XP23/UNC12**  
1297

# DIAGNOSTIC EX PROBES

**What is an Ex-Probe?** Ex-Probe is an excellent combination instrument for dental diagnostics. It has an tip on one end and a probe tip on the other.



15-14-13-12-11-10-9-8-7-6-5-4-3-2-1

23/UNC15  
University of North Carolina

**23/CP12UNC**  
1298

15-14-13-12-11-10-9-8-7-6-5-4-3-2-1

TU17/  
UNC15

**XPTU17/UNC15**  
1299

10-9-8-7-5-3-2-1

23/W  
Williams

**XP23/W**  
1300

10-9-8-7-5-3-2-1

6/OW  
Thin  
Williams

**XP6/OW**  
1301

10-9-8-7-5-3-2-1

6/W  
Williams

**XP6/W**  
1302

10-9-8-7-5-3-2-1

17/OW  
Thin  
Williams

**XP17/OW**  
1303

15-14-13-12-11-10-9-8-7-6-5-4-3-2-1

3A/UNC15

**XP3A/UNC15**  
1304

10-9-8-7-5-3-2-1  
FLAT

23/GF  
Goldman Fox

**XP23/GF**  
1305

8-6-3

23/O  
Michigan

**XP23/O**  
1306

## Endodontic Instruments



## Endodontic Instruments Index

Bone crushers				I21
Bone Grafting Instruments				I4
Bone Grafting Syringes				I19-I20
Bone Pluggers				I3
Bone Scoops				I3
Bone Scrapers & Calipers				I20
Bone Spreading Osteotome Kit				I17-I18
Bone Spreading Osteotomes				I15-I16
Calipers				I19
Implant Periosteal Chisels				I2
Implant Periosteals				I1
Osteotome Chisels				I7
Sinus Lift Instruments				I10-I13
Sinus Lift/Bone Graft				I9
Tissue Punches				I8
Tunneling Knives				I6
Tunneling Procedure Instruments				I5

# ENDODONTICS

## ELPLORERS/EXCAVATORS

**Explorers:** These instruments used in a procedures to preserve the health of the dental pulp and periapical tissues, and the treatment of the diseased pulp to enable the tooth (or part of the tooth) to be retained in function. Further, endodontic treatment may be carried out on a healthy pulp where the tooth has to be devitalized as part of the overall restorative treatment plan.

**Excavators:** The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber.

### Endodontic-Explorers



DG16  
**EXDG16**  
1350

DG16/23  
**EXDG16/23**  
1351

DG16/17  
**EXDG16/17**  
1352

1  
Abou-Rass  
Explorer  
**EX1AR**  
1353

### Endodontic-Excavators



1.5  
2  
Glick  
**EXCGL2**  
1354

1.0  
31L  
**EXC31L**  
1355

# ENDODONTICS

## ELPLORERS/EXCAVATORS

**Explorers:** These instruments used in a procedures to preserve the health of the dental pulp and periapical tissues, and the treatment of the diseased pulp to enable the tooth (or part of the tooth) to be retained in function. Further, endodontic treatment may be carried out on a healthy pulp where the tooth has to be devitalized as part of the overall restorative treatment plan.

**Excavators:** The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber.



32L

EXC32L

1356



33L

EXC33L

1357



31W

Oval Spoon

EXC31W

1358

# ENDODONTICS

## PLUGGERS

Used to plug and condense amalgam into cavity preparation. Single ended hand instruments used to plug and condense gutta percha in a root canal.

Endodontic Pluggers ISO Size



ISO  
Size 30  
24mm

**RCP30**  
1400

ISO  
Size 40  
24mm

**RCP40**  
1401

ISO  
size 50  
24mm

**RCP50**  
1402

ISO  
Size 60  
24mm

**RCP60**  
1403

Endodontic Luks



1/3  
21mm

**RCP1/3**  
1404

5/7  
21mm

**RCP5/7**  
1405

# ENDODONTICS

## PLUGGERS

Used to plug and condense amalgam into cavity preparation. Single ended hand instruments used to plug and condense gutta percha in a root canal.



# ENDODONTICS

## PLUGGERS

Used to compact filling material during vertical condensation.

Anterior



8A 29mm	8 1/2 A 29mm	9A 29mm	9 1/2 A 29mm	10A 29mm	10 1/2 A 29mm	11A 29mm	11 1/2 A 29mm	12A 29mm
------------	-----------------	------------	-----------------	-------------	------------------	-------------	------------------	-------------

<b>ARCP8</b> 14 12	<b>ARCP8-1/2</b> 14 13	<b>ARCP9</b> 14 14	<b>ARCP9-1/2</b> 14 15	<b>ARCP10</b> 14 16	<b>ARCP10-1/2</b> 14 17	<b>ARCP11</b> 14 18	<b>ARCP11-1/2</b> 14 19	<b>ARCP12</b> 14 20
-----------------------	---------------------------	-----------------------	---------------------------	------------------------	----------------------------	------------------------	----------------------------	------------------------

# ENDODONTICS

## PLUGGERS

Used to compact filling material during vertical condensation.


### Posterior



# ENDODONTICS

## NICKEL TITANIUM ENDODONTIC SPERADERS

Used to compact filling material during lateral condensation. These instruments made of nickel titanium.



<b>RCSD11NT</b>	<b>RCSD11TNT</b>	<b>RCSD11TSNT</b>	<b>RCS255NT</b>	<b>RCS25NT</b>	<b>RCS40NT</b>	<b>RCS40SNT</b>	<b>RCSMA57NT</b>	<b>RCS3NT</b>
1480	1481	1482	1483	1484	1485	1486	1487	1488

# ENDODONTICS

## SPREADERS

Hand instrument with a tapered, pointed end, used to condense and compress materials such as gutta percha in root canals.

ISO Sized Spreaders



Spreaders



# ENDODONTICS

## SPREADERS

Hand instrument with a tapered, pointed end, used to condense and compress materials such as gutta percha in root canals.



# ENDODONTICS

## MICRO CURETTES/MICRO SURGICAL EXCAVATORS

**Excavators:** The shank is long to reach canals. To curettage inside of tooth to base of pulp chamber.

Micro-Curettes



1  
Abou-Rass  
Straight Apical

**MCARAS1**

1600

2  
Abou-Rass  
Straight Apical

**MCARAS2**

1601

3  
Abou-Rass  
Straight Apical

**MCARAS3**

1602

4  
Abou-Rass  
Straight Apical

**MCARAS4**

1603

4  
Velart  
Apical

**MCAVA4**

1604

1  
Abou-Rass

**MECARA1**

1605

2  
Abou-Rass

**MECARA2**

1606

Abou-Rass

**MECARA3**

1607

4  
Abou-Rass

**MECARAD4**

1608

# ENDODONTICS

## APICAL PLUGGERS

Used to compact filling material during vertical condensation.


### Pluggers-Straight



# ENDODONTICS

## APICAL PLUGGERS

Used to compact filling material during vertical condensation.



# ENDODONTICS

## APICAL BURNISHERS / RETRO FILLING

Used to condense, smooth and polish amalgam.

### Retro Filling

The retro filling instruments are sized to work with ultrasonic preparation tips.

### Burnishers-Straight



Left Apical  
**PFLA**  
1700

Right Apical  
**PFRA**  
1701

1 Abou-Rass Straight  
**AB-ARS1**  
1702

2 Abou-Rass Straight  
**AB-ARS2**  
1703

3 Abou-Rass Straight  
**AB-ARS3**  
1704

4 Abou-Rass Straight  
**AB-ARS4**  
1705

5 Abou-Rass Straight  
**AB-ARS5**  
1706

# ENDODONTICS

## APICAL BURNISHERS / RETRO FILLING

Used to condense, smooth and polish amalgam.

### Burnishers-Angled



1  
Abou-Rass  
Angled

**AB-ARS1**  
1707

2  
Abou-Rass  
Angled

**AB-ARS2**  
1708

3  
Abou-Rass  
Angled

**AB-ARS3**  
1709

4  
Abou-Rass  
Angled

**AB-ARS4**  
1710

5  
Abou-Rass  
Angled


**AB-ARS5**  
1711

# ENDODONTICS

## APICAL PLUGGERS / BURNISHERS

**Pluggers:** are used to condense filling materials into cavity preparation. Ends are serrated.

**Burnisher:** Used to condense, smooth and polish amalgam.



Small Left Apical  
**APBSL**  
 1712

Small Right Apical  
**APBSR**  
 1713

Medium Left Apical  
**APBML**  
 1714

Medium Right Apical  
**APBMR**  
 1715

Large Left Apical  
**ABLL**  
 1716

# ENDODONTICS

## APICAL PLUGGERS / BURNISHERS

**Pluggers:** are used to condense filling materials into cavity preparation. Ends are serrated.

**Burnisher:** Used to condense, smooth and polish amalgam.



# ENDODONTICS

## SPLINTING FORCEPS RUBBER DAM PUNCHES & R.D. ACCESSORIES

**Splinting Forceps:** Immobilization of fractured bone ends by wiring, pinning or splinting.

**R.D. Clamp Forceps:** Instruments with self-locking device used to apply and remove rubber dam clamps.

**R.D. Frame:** Flexible frame used to hold the rubber dam in a taut state, away from the field of operation.

**R.D. Punch:** Instruments used to punch various sized holes in a sheet of rubber dam.

Martensite stainless steel, Clamps shown 50% actual.

Post And Silver Point Removal Forceps to remove posts from root canals

Four types of Steiglitz forceps are shown, each with a diagram of its tip. The forceps are: Steiglitz Straight (14cm/5 1/2"), Steiglitz (Archer) 45° (13.5cm/5 1/4"), Steiglitz 90° (12cm/5 3/4"), and Peer 45° (12cm/4 3/4").

<p>Steiglitz Straight 14cm/5 1/2"</p> <p><b>RFSTR</b> 1770</p>	<p>Steiglitz (Archer)45° 13.5cm/5 1/4"</p> <p><b>RF45</b> 1771</p>	<p>Steiglitz 90° 12cm/5 3/4"</p> <p><b>RF90</b> 1772</p>	<p>Peer 45° 12cm/4 3/4"</p> <p><b>RFP</b> 1773</p>
--	--	--	--

Rubber Dam Clamp Forceps

shown 40% size

Two types of Rubber Dam Clamp Forceps are shown. The first is RDF4 (17.5cm/7") and the second is RDF (17cm/6 1/2").

<p>4 17.5cm/7"</p> <p><b>RDF4</b> 1774</p>	<p>Rubber Dam Forceps 17cm/6 1/2"</p> <p><b>RDF</b> 1775</p>
--	--

Rubber Dam Frames

Two types of Rubber Dam Frames are shown: RDCF5 (13cm/5") and RDAF6 (15cm/6").

<p>Rubber Dam Frame-Child 13cm/5"</p> <p><b>RDCF5</b> 1776</p>	<p>Rubber Dam Frame-Adult 15cm/6"</p> <p><b>RDAF6</b> 1777</p>
--	--

Rubber Dam Punches

Two types of Rubber Dam Punches are shown: RDPA (17cm/6 3/4") and RDP (16cm/6 1/4").

<p>Ainsworth 17cm/6 3/4"</p> <p><b>RDPA</b> 1778</p>	<p>Rubber Dam Punch 16cm/6 1/4"</p> <p><b>RDP</b> 1779</p>
--	--





# ENDODONTICS

## RUBBER DAM CLAMPS



Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls, these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

### Molars - General

					
	8A		W8A		8AD
	<b>RDC8A</b>		<b>RDCW8A</b>		<b>RDC8AD</b>
	1830		1831		1832
					
	10		11		14
	<b>RDC10</b>		<b>RDC11</b>		<b>RDC14</b>
	1833		1834		1835
					
	14A		W14A		18
	<b>RDC14A</b>		<b>RDCW14A</b>		<b>RDC18</b>
	1836		1837		1838
					
	24		25		26N
	<b>RDC24</b>		<b>RDC25</b>		<b>RDC26N</b>
	1839		1840		1841

### Upper Molars

			
	56		56S
	<b>RDC56</b>		<b>RDC56S</b>
	1842		1843
			
	4		W4
	<b>RDC4</b>		<b>RDCW4</b>
	1844		1845
			
	5		8
	<b>RDC5</b>		<b>RDC8</b>
	1846		1847
			
	27		
	<b>RDC27</b>		
	1848		

# ENDODONTICS

## RUBBER DAM CLAMPS

Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls. these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

### Upper Molars

### Lower Molars

### Pre-molars

	30
<b>RDC30</b>	1849
	31
<b>RDC31</b>	1850
	201
<b>RDC201</b>	1851
	205
<b>RDC205</b>	1852

	3		W3		7
<b>RDC3</b>	1853	<b>RDCW3</b>	1854	<b>RDC7</b>	1855
	7		7A		12A
<b>RDC7</b>	1855	<b>RDC7A</b>	1857	<b>RDC12A</b>	1858
	7		7A		12A
<b>RDC7</b>	1855	<b>RDC7A</b>	1857	<b>RDC12A</b>	1858
	31A		26		28
<b>RDC31A</b>	1859	<b>RDC26</b>	1860	<b>RDC28</b>	1861

	0
<b>RDC0</b>	1862
	1
<b>RDC1</b>	1863
	1A
<b>RDC1A</b>	1864
	2
<b>RDC2</b>	1865


# ENDODONTICS

## RUBBER DAM CLAMPS


Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls. these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

Pre-molars



Anterior



Rubber Dam Clamp Organizing Board  
1877

# Implantology Instruments



## Implantology Instruments Index

Hard wire Cutting pliers				O3
Hole Punching Pliers				O1
Ortho Hand Instruments				O9-O12
Ortho height gauges				O14
Ortho Needle Holders				O15
Ortho Tweezers and Scissors				O13
Pin & Ligature Cutting pliers				O2
Utility Pliers				O5
Wire Forming Pliers				O6-O8

# IMPLANTOLOGY

## IMPLANT PERIOSTEALS

Hand instrument, usually double ended, used to separate the periosteum from the underlying bone.



# IMPLANTOLOGY

## IMPLANT PERIODONTAL CHISELS

Hand instrument used to remove hard tissues such as enamel or bone by chipping, cleaving or paring. The beveled blade may be in line with the handle or at an angle to it.



1/2  
Buser Modified  
Periodontal  
Chisel  
3mm/4mm

**IPCBM1/2**  
1970

2/3  
Buser Modified  
Periodontal  
Chisel  
4mm/5mm

**IPCBM2/3**  
1971

1/3  
Buser Modified  
Periodontal  
Chisel  
3mm/5mm

**IPCBM1/3**  
1972

Jovanovic  
13KTG  
Periodontal  
Chisel  
1mm/4mm

**IPCJOV1/4**  
1973

Jovanovic  
13KTG  
Periodontal  
Chisel  
2mm/5mm

**IPCJOV2/5**  
1974

Jovanovic  
13KTG  
Periodontal  
Chisel  
3mm/6mm

**IPCJOV3/6**  
1975

P24G/13KL  
Palacci Modified

**IPC24G/13KL**  
1976

# IMPLANTOLOGY

## BONE SCOOPS & PLUGGERS

**Bone Scoops:** Hand instrument similar to a bone chisel but having a hollowed-out blade similar to coupland chisel.

**Bone Plugger:** Used to plug and condense grafted bone.

### BONE SCOOPS



Jovaanovic  
None Scoop  
5mm

**IBSJOV5mm**  
2020

Bone Graft  
Scoop

**IBSBG**  
2021

Labanca  
Bone Curette  
7mm

**ICLBC**  
2022

### BONE PLUGGERS



LaBanca  
3x1.5mm

**IBPLAB**  
2023

4  
Jovanovic  
3mm

**IBPJ4**  
2024

Bone  
Packer DE  
5-7mm

**IBPP5-7**  
2025

Bone Packer  
Serrated  
3-5mm

**IBPS3-5**  
2026

## Bone Grafting Pluggers/Condensors



Bone Grafting Plugger (1.0-1.7mm)  
(8-11-13-15-18)  
IBG1.0/1.7  
2081

IBG1.0/1.7  
2081



Bone Grafting Plugger (1.7-2.0mm)  
(8-11-13-15-18)

IBG1.7/2.0  
2082



Bone Grafting Plugger (2.4-2.7mm)  
(8-11-13-15-18)

IBG2.4/2.7  
2083



Grisdale Bone Condensor (3.3/4.0mm)  
(8-11-13-15-18)

IBG3.3/4.0  
2084



Grisdale Bone Condensor (3.3/4.0mm)  
(8-11-13-15-18), Serrated

IBGS3.3/4.0  
2085



TPI1-Tunneling Procedure Instrument

ITPI1  
2090



TPI2-Tunneling Procedure Instrument

ITPI2  
2091



TPI3-Tunneling Procedure Instrument

ITPI3  
2092



TPI4-Tunneling Procedure Instrument

ITPI4  
2093



ITPI5-Tunneling Procedure Instrument

ITPI5  
2094



TPI6-Tunneling Procedure Instrument

ITPI6  
2095



TPI7-Tunneling Procedure Instrument

ITPI7  
2096

## TUNNELING KNIVES



Urban Tunneling Knife

IUTK  
2097



Tunneling Knife 1

ITK1  
2098



Tunneling Knife 2 (Short Handle)

ITK2S  
2099



Tunneling Knife 2

ITK2  
2100



Giles De Quincey Knife 1

IGDQK1  
2101



Giles De Quincey Knife 2

IGDQK2  
2102



Steigmann Soft Tissue Splitter

ISSTS  
2103



Steigmann Soft Tissue Splitter, Modified

ISSTSM  
2104

# IMPLANTOLOGY

## OSTEOTOME CHISELS

**Osteotome Chisels:** Surgical instruments used to cut or divide bone.



7-10-13-15-18mm

Osteotome Chisel  
Straight, Single Edge  
4mm

**IOCSSE4**

2 1 1 0



7-10-13-15-18mm

Osteotome Chisel  
Curved, Single Edge  
4mm

**IOCCSE4**

2 1 1 1



7-10-13-15-18mm

Osteotome Chisel  
Straight, Single Edge  
7.5mm

**IOCSSE7.5**

2 1 1 2



7-10-13-15-18mm

Osteotome Chisel  
Curved, Single Edge  
7.5mm

**IOCCSE7.5**


2 1 1 3



Osteotome Chisel  
Jovanovic  
6mm

**IOCJOV6**

2 1 1 4




8-10-13-18mm

Osteotome Chisel  
Single Edge  
3mm

**IOCSE3**

2 1 1 5




8-10-13-18mm

Osteotome Chisel  
8-10-13-18mm  
Double edge, 3mm

**IOCDE3**

2 1 1 6




8-10-13-18mm

Osteotome Chisel  
8-10-13-18mm  
Single Edge, 6mm

**IOCSE6**

2 1 1 7




8-10-13-18mm

Osteotome Chisel  
8-10-13-18mm  
Double Edge, 6mm

**IOCDE6**

2 1 1 8




8-10-13-18mm

Osteotome Chisel  
8-10-13-18mm  
Single Edge, 9mm

**IOCSE9**

2 1 1 9



8-10-13-18mm

Osteotome Chisel  
8-10-13-18mm  
Double Edge, 9mm

**IOCDE9**

2 1 2 0

# IMPLANTOLOGY

## TISSUE PUNCHES

Martensite stainless steel tissue punches.



Tissue Punch STR  
4mm  
**ITPS4**  
2 130



Tissue Punch STR  
5mm  
**ITPS5**  
2 131



Tissue Punch STR  
6mm  
**ITPS6**  
2 132



Tissue Punch CVD  
2mm  
**ITPC2**  
2 133



Tissue Punch CVD  
3mm  
**ITPC3**  
2 134



Tissue Punch CVD  
4mm  
**ITPC4**  
2 135



Tissue Punch CVD  
5mm  
**ITPC5**  
2 136



Tissue Punch CVD  
6mm  
**ITPC6**  
2 137

### Sinus Lift Instruments



Sinus Lift Cobra Shape 3.5mm (2-4-5-6-8-10mm)

ISC3.5  
2152



Sinus Lift Cobra Shape In/Out (2-4-5-6-8-10mm)

ISCI/O  
2153



Sinus Lift Mushroom Shape (2.8/3.8mm)  
(2-4-5-6-8-10)

IMSL2.8/3.8  
2154




Sinus Lift Mushroom Shape (4.8/5.8mm)  
(2-4-5-6-8-10)

IMSL4.8/5.8  
2155

# IMPLANTOLOGY

## SINUS LIFT INSTRUMENTS




# IMPLANTOLOGY

## SINUS LIFT INSTRUMENTS


Designed to separate the sinus lining from the bone without fraying or tearing.



Designed to separate the sinus lining from the bone without fraying or tearing.




Designed to separate the sinus lining from the bone without fraying or tearing.





2.0mm  
**ITCCSO2.0**  
2255

Tapered Concave Tip-Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm  
**ITCCSO3.5**  
2256

Tapered Concave Tip-Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm  
**ITCCSO4.3**  
2257

Tapered Concave Tip-Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm  
**ITCCSO5.0**  
2258

Tapered Concave Tip-Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm  
**ITCCSO5.5**  
2259

Tapered Concave Tip-Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm  
**ITCCA02.0**  
2260

Tapered Concave Tip-Angulated - 2.0mm (4-6-8-10-13-16-18-20-23-26)



3.5mm  
**ITCCA03.5**  
2261

Tapered Concave Tip-Angulated - 3.5mm (4-6-8-10-13-16-18-20-23-26)




4.3mm  
**ITCCA04.3**  
2262

Tapered Concave Tip-Angulated - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm  
**ITCCA05.0**  
2263

Tapered Concave Tip-Angulated - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm  
**ITCCA05.5**  
2264

Tapered Concave Tip-Angulated - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm  
**ITCVSO2.0**  
2270

Tapered Convex- Tip-Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm  
**ITCVSO3.5**  
2271

Tapered Convex- Tip-Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm  
**ITCVSO4.3**  
2272

Tapered Convex- Tip-Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)




5.0mm  
**ITCVSO5.0**  
2273

Tapered Convex- Tip-Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm  
**ITCVSO5.5**  
2274

Tapered Convex- Tip-Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)



2.0mm  
**ITCVAO2.0**  
2275

Tapered Convex-Angulated Tip - 2.0mm (4-6-8-10-13-16-18-20-23-26)



3.5mm  
**ITCVAO3.5**  
2276

Tapered Convex-Angulated Tip - 3.5mm (4-6-8-10-13-16-18-20-23-26)




4.3mm  
**ITCVAO4.3**  
2277

Tapered Convex-Angulated Tip - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm  
**ITCVAO5.0**  
2278

Tapered Convex-Angulated Tip - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm  
**ITCVA5.5**  
2279

Tapered Convex-Angulated Tip - 5.5mm (4-6-8-10-13-16-18-20-23-26)



2.0mm  
**IPSO2.0**  
2281

PUSHER - Tip-Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)



3.5mm  
**IPSO3.5**  
2282

PUSHER Tip-Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm  
**IPSO4.3**  
2283

PUSHER- Tip-Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm  
**IPSO5.0**  
2284

PUSHER- Tip-Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)




5.5mm  
**IPSO5.5**  
2285

PUSHER- Tip-Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)




2.0mm  
**IPA02.0**  
2286

PUSHER- Angulated Tip - 2.0mm (4-6-8-10-13-16-18-20-23-26)




3.5mm  
**IPA03.5**  
2287

PUSHER- Angulated Tip - 3.5mm (4-6-8-10-13-16-18-20-23-26)



4.3mm  
**IPA04.3**  
2288

PUSHER- Angulated Tip - 4.3mm (4-6-8-10-13-16-18-20-23-26)



5.0mm  
**IPA05.0**  
2289

PUSHER- Angulated Tip - 5.0mm (4-6-8-10-13-16-18-20-23-26)



5.5mm  
**IPA05.5**  
2290

PUSHER- Angulated Tip - 5.5mm (4-6-8-10-13-16-18-20-23-26)

Bone Spreading  
Ostetome Kit




Interchangeable  
Osteotome Handle

### TAPERED CONVEX OSTETOMES

#### STRAIGHT

2.7 mm    3.2 mm    3.7 mm    4.2 mm    5.0 mm



#### ANGLED

2.7 mm    3.2 mm    3.7 mm    4.2 mm    5.0 mm




ICVOS2.7 2300    ICVOS3.2 2301    ICVOS3.7 2302    ICVOS4.2 2303    ICVOS5.0 2304

ICVOA2.7 2305    ICVOA3.2 2306    ICVOA3.7 2307    ICVOA4.2 2308    ICVOA5.0 2309


### TAPERED CONCAVE OSTEOTOMES

#### STRAIGHT



ICCOS2.7 2310    ICCOS3.2 2311    ICCOS3.7 2312    ICCOS4.2 2313    ICCOS5.0 2314


#### ANGLED



ICCOA2.7 2315    ICCOA3.2 2316    ICCOA3.7 2317    ICCOA4.2 2318    ICCOA5.0 2319


### BONE PUSHER OSTEOTOMES

#### STRAIGHT



IPOS2.7 2320    IPOS3.2 2321    IPOS3.7 2322    IPOS4.2 2323    IPOS5.0 2324

#### ANGLED



IPOA2.7 2325    IPOA3.2 2326    IPOA3.7 2327    IPOA4.2 2328    IPOA5.0 2329

# IMPLANTOLOGY

## CALIPERS & BONE GRAFTING SYRINGES



CASTROVIEJO  
STR 17CM

**ICVCS17**  
2382



Castroviejo  
Calipers  
Straight  
0-40mm  
18cm/7"

**ICVCS18**  
2384



Castroviejo  
Calipers  
Angled  
0-40mm  
18cm/7"

**ICVCCA18**  
2386



Bone Grafting  
Syringe  
7mm

**IBGS7**  
2396

- 2389-1.0mm **IBGS1**  
2389
- 2390-2.0mm **IBGS2**  
2390
- 2391-2.5mm **IBGS2.5**  
2391
- 2392-3.0mm **IBGS3**  
2392
- 2393-4.0mm **IBGS4**  
2393
- 2394-5.0mm **IBGS5**  
2394
- 2395-6.0mm **IBGS6**  
2395

NOTE: Also available in  
1mm, 2mm, 2.5mm, 3mm, 4m, 5mm and 6mm



Weiss Modified  
Castroviejo  
Calipers  
0-40mm  
17cm/6<sup>3</sup>/<sub>4</sub>"

**IWMCV17**  
2383



CASTROVIEJO  
STR 8.5cm  
ICVC-ST8.5cm

**ICVCS8.5**  
2385



CASTROVIEJO  
CVD 17cm  
ICVC-CV17cm

**ICVCA17**  
2387



Cortical  
Bone Clamp

**ICBC**  
2388



Bone Grafting 1 Ring Style Syringe  
Straight, 3.75mm

**IBGS1RSS3.75**  
2400



Bone Grafting 1 Ring Style Syringe  
Straight, 50.0mm

**IBGS1RSS5.0**  
2402



Bone Grafting 1 Ring Style Syringe  
Curved, 3.75mm

**IBGS1RSC3.75**  
2401



Bone Grafting 1 Ring style Syringe  
Curved, 5.0mm

**IBGS1RSC5.0**  
2403



Ridge Mapping, Bone Width 0.25mm  
Both Side Reading, 14cm

**RMBC14**  
2410



Caliper Double sided Reading  
Scissors Style

**ICDSS**  
2411



Sliding Implant Caliper, 94mm  
1.5mm Tip diameter

**ISIC94**  
2412




Stainless Steel Ruler  
15cm/6"

**SSR**  
2413



Bone collector Chise, Straight

**IBCC**  
2415



Bone Collector Blade, Straight

**IBCBS**  
2416



Bone Collector Blade

**IBCBS**  
2418

**IBCCC**  
2417

Bone Collector Chisel, Curved

# IMPLANTOLOGY


## ENDO FORCEPS, BONE CRUSHERS & MALLETS

Micro Surgery Endo Forceps, 95mm

IMSEF95  
2419

Micro Surgery Endo forceps, 110mm

IMSEF110  
2420




Mead-2  
Mallet

IM2M  
2421



Bone Crusher

IBC  
2424



130g  
Mallet

I130M  
2422



235g  
Mallet

I235M  
2423



Bone Morseriser

IBM  
2425

## Orthodontic Instruments



## Orthodontic Instruments Index

Hard wire Cutting pliers					O3
Hole Punching Pliers					O1
Ortho Hand Instruments					O9-O12
Ortho height gauges					O14
Ortho Needle Holders					O15
Ortho Tweezers and Scissors					O13
Pin & Ligature Cutting pliers					O2
Utility Pliers					O5
Wire Forming Pliers					O6-O8


## HOLE PUNCHING PLIERS



**HOLE PUNCH PLIERS**  
Creates Half-Moon Cutouts in Clear Aligners

2426



**TEAR DROP PUNCH PLIERS**  
Designed to Create a Notch at the  
Gingival Margin of Clear Aligners

2427



**HORIZONTAL PUNCH PLIERS**  
Horizontal Rectangular Shaped Punch

2428



**VERTICAL PUNCH PLIERS**  
Vertical Rectangular Shaped Punch Pliers

2429

## PIN & LIGATURE CUTTING PLIERS



**PIN & LIGATURE CUTTER  
LONG HANDLE**  
Cuts Ligature Wires upto 0.015"/0.38mm

2430



**PIN & LIGATURE CUTTER  
MINI STYLE**  
Cuts Ligature Wires upto 0.015"/0.38mm

2431



**PIN & LIGATURE CUTTER  
NARROW TIP**  
Cuts Ligature Wires upto 0.015"/0.38mm

2432



**PIN LIGATURE CUTTER  
NARROW TIP 15, ANGLED**  
Cuts Ligature Wires upto 0.015"/0.38mm

2433



**PIN & LIGATURE CUTTER  
STANDARD**  
Cuts Ligature Wires upto 0.015"/0.38mm

2434



**PIN & LIGATURE CUTTER  
STANDARD 15 ANGLED**  
Cuts Wires upto 0.015"/0.38mm

2435

## HARD WIRE CUTTING PLIERS



**ORTHO HARD WIRE CUTTER**  
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm  
Precisely Angled Straight Tip

2436



**ORTHO HARD WIRE CUTTER, 15 ANGLD**  
Cuts Hard Orthodontic Wires upto 0.028"/0.71mm

2437



**DISTAL END CUTTER MINI STYLE CUT-HOLD**  
Safely Cuts & Holds Wires upto 0.025"/0.64mm

2438



**DISTAL END CUTTER STANDARD CUT-HOLD**  
Safely Cuts & Holds Wires upto 0.025"/0.64mm

2439



**DISTAL END CUTTER LONG HANDLE, CUT-HOLD**  
Safely Cuts & Holds Wires upto 0.025"/0.64mm

2440



**FLUSH CUT-HOLD DISTAL END CUTTER STANDARD HANDLE**  
Featuring a Tungsten Carbide Tip Insert

2441



**FLUSH CUT-HOLD DISTAL CUTTER LONG HANDLE**  
Featuring a Tungsten Carbide Tip Insert

2442

## UTILITY PLIERS



**HOW PLIERS  
STRAIGHT**  
Precisely Serrated Jaws for a Firm Grip

2450



**DEBONDING PLIERS  
STRAIGHT**  
Precisely Aligned, Straight Tip

2451



**FORCE MODULE  
SEPARATING PLIERS**  
Smooth Reversible Action Control

2452



**COON STYLE  
LIGATURE PLIERS**  
Long Slim Tips with Reversible Action Control

2453



**ADHESIVE REMOVING PLIERS**  
Double Sided Tungsten Carbide Tip Insert  
Comes with Free Pack of Occlusal Pads

2454



**WEINGART PLIERS  
LONG HANDLE**  
Designed for Wires upto 0.030"/0.76mm

2455



**HOW PLIERS  
CURVED**  
Precisely Serrated Jaws for Firm Grip

2457



**POSTERIOR BAND REMOVING  
PLIERS, LONG TIP**  
Comes with a free Pack of Occlusal Pads

2458



**BRACKET REMOVING PLIERS,  
STRAIGHT**  
Slim edge Jaws for an easy Grip on Bracket

2459

## UTILITY PLIERS



**BRACKET REMOVING PLIERS  
ANGLED**  
Slim Edge Jaws for an Easy Grip on Bracket

2460



**POSTERIOR BAND REMOVING PLIERS  
STANDARD TIP**  
Comes with a free pack of Occlusal pads

2461



**DEBONDING PLIERS  
ANGLED**  
60 Angled Tip

2462



**WEINGART PLIERS 90  
ANGLED**  
Designed for wires  
upto 0.030"/0.76mm

2463



**WEINGART PLIERS  
SLIM**  
Designed for wires upto 0.030"/0.76mm

2464



**WEINGART UTILITY  
PLIERS**  
Designed for wires upto 0.030"/0.76mm

2465



**WEINGART PLIERS, UNIVERSAL  
45 ANGLED TIP**  
Designed for wires upto 0.030"/0.76mm

2466



**CROWN & BAND CONTOURING  
PLIERS**  
Featuring a concave and convex tip design

2467



**CROWN & BAND CRIMPING  
PLIERS**  
Slim design tips for precise crimping

2468



**BALL HOOK CRIMPING PLIERS**  
Designed to crimp hooks, posts and  
auxillary stops on arch wire

2469



**POSTERIOR BAND REMOVING  
PLIERS, SHORT TIP**  
Comes with a free pack of Occlusal pads

2470

## WIRE FORMING PLIERS



**BIRD BEAK PLIERS**  
Featuring a Cone & Pyramid shaped beak  
Designed for wires upto 0.028"/0.71mm

2475



**LIGHT WIRE PLIERS  
LONG BEAK**  
Precisely built for wires upto 0.020"/0.51mm

2476



**TWEED LOOP FORMING PLIERS**  
Designed for wire upto 0.020"/0.51mm

2477



**LIGHT WIRE PLIERS WITH 3 GROOVES**  
Precisely built for wires upto 0.020"/0.51mm

2478



**LIGHT WIRE PLIERS WITH  
CUTTER & 3 GROOVES**  
Precisely built for wires upto 0.020"/0.51mm

2479



**LIGHT WIRE PLIERS  
WITH CUTTER**  
Precisely built for wires upto 0.020"/0.51mm

2480



**JARABAK PLIERS**  
Featuring 3 grooves for a firm wire grip

2481



**JARABAK PLIERS WITH CUTTER**  
Featuring 3 grooves & cutter at base

2482



**LINGUAL ARCH FORMING PLIERS**  
Designed to create double back and triple back  
bends in all wires upto 0.036"/0.91mm

2483

## WIRE FORMING PLIERS



**NANCE LOOP FORMING PLIERS**  
Four step beak to create 3,4,5 and 6mm loop  
designed for all wires upto 0.022"/0.56mm

2484



**De La Rosa PLIERS**  
Featuring a chamfered concave and convex tip  
built for all wires upto 0.025"/0.64mm

2485



**De La Rosa Pliers WITH 3 GROOVES**  
Featuring a chamfered concave and convex tip  
built for all wires upto 0.025"/0.64mm

2486



**HOLLOW CHOP CONTOURING PLIERS**  
Long concave & Convex tip  
for intricate wire contouring

2487



**OPTICAL PLIERS**  
Featuring a slim beak with concave end  
designed for all wire upto 0.025"/0.64mm

2488



**V-STOP PLIERS**  
Create V-shaped bends on  
wires upto 0.022"/0.56mm

2489



**YOUNG PLIERS**  
3 sectional beak with 3 grooves  
designed for wires upto 0.028"/0.71mm

2490



**TWEED ARCH FORMING PLIERS 1.75mm**  
Designed for wires upto 0.025"/0.66mm

2491



**TWEED ARCH FORMING PLIERS 1.0mm**  
Designed for wires upto 0.025"/0.66mm

2492

## WIRE FORMING PLIERS



**TORQUING PLIERS FEMALE/OUTER**  
Places torques in all wires upto 0.030"/0.76mm alongside our inner torquing pliers

2493



**TORQUING PLIERS MALE/INNER**  
Places torques in all wires upto 0.030"/0.76mm alongside our outer torquing pliers

2494



**STEP PLIERS 1MM**  
Designed for all wires upto 0.028"/0.71mm

2495



**STEP PLIERS 0.75MM**  
Designed for all wires upto 0.028"/0.71mm

2496



**STEP PLIERS 0.5MM**  
Designed for all wires upto 0.028"/0.71mm

2497



**NiTi DISTAL CINCH BACK PLIERS**  
Designed for SS & NiTi wires upto 0.025"/0.64mm

2498



**NiTi DISTAL CINCH BACK PLIERS, LINGUAL**  
Designed for SS&NiTi wires upto 0.018"/0.46mm

2499



**3 PRONG PLIERS STANDARD**  
Precisely built for all wires upto 0.030"/0.76mm

2500



**3 PRONG PLIERS HEAVY**  
Precisely built for all wires upto 0.030"/0.76mm

2501



**3 PRONG PLIERS NiTi**  
Precisely built for all wire upto 0.020"/0.5mm

2502

## HAND INSTRUMENTS



Ligature  
Director/Band  
Pusher  
2506



Small  
Plugger/  
Picker  
2507



Dull Scaler  
Band Pusher  
2508



Ligature  
Director  
Scaler  
2509



Scaler/  
Band  
Pusher  
2510



Lingual  
Ligature  
Director,  
U Shape  
2511

# ORTHODONTICS

## HAND INSTRUMENTS



Dull  
Scaler  
2512



Dull  
Scaler  
2513



Large  
Plugger/Picker  
2514



Picker/  
Ligature  
Director  
2515



Scaler/  
Picker  
2516




Dull  
Scaler  
2517



Ligature  
Director/  
Band Pusher  
2518


# ORTHODONTICS

## HAND INSTRUMENTS



Band Pusher  
with Short  
Tip

2519



Band Pusher

2520



Fillion Explorer/  
Ligature  
Director

2521



Band Pusher/  
Scaler

2522



Band Pusher/  
Scaler

2523



Ligature  
Director

2524




Ligature  
Director

2525

# ORTHODONTICS

## HAND INSTRUMENTS



Lingual  
Torquing  
Key

2526

Straight  
.020"  
Applicator

2527

Single End  
Ligature  
Director

2528

Lingual  
Ligature  
Director,  
.018"

2529

Lingual  
Ligature  
Director,  
.025"

2530

Ligature  
Director/  
Plugger

2531

Ligature  
Director  
Tucker

2532

Plugger/  
Ligature  
Picker

2533

# ORTHODONTICS

## TWEEZERS-SCISSORS



Posterior Bracket Holder Tweezer, Long

2566



Direct Bond Bracket Tweezer with Aligner & Positioner

2567



Buccal Tube Tweezer

2568



Slim Buccal Tube Bonder Tweezer

2569



Direct Bond Bracket Tweezer, Lingual

2570



Direct Bond Tweezer, Aligner

2571



All Purposes Scissors

2572



Crown Scissors

2573



Crown Scissors, Straight

2574



Crown Scissors, Curved

2575

## HEIGHT GAUGES



BONE BRACKET POSITIONING GAUGE

2580



BRACKET HEIGHT GAUGE FOR  
0.022"/0.56mm Brackets

2581



BRACKET HEIGHT GAUGE FOR  
0.018"/0.46mm Brackets

2582



BRACKET HEIGHT GAUGE, MOVEABLE HEAD

2583

## MATHIEU NEEDLE HOLDERS



MATHIEU NARROW TIP  
Enhanced double spring function

2585



MATHIEU SMAHA FINE TIP  
Enhanced double spring function

2586



MATHIEU NARROW TIP, TC  
Enhanced double spring function

2587



MATHIEU STANDARD TIP  
Enhanced double spring function


2588



MATHIEU STANDARD TIP, TC  
Enhanced double spring function

2589

## Periodontic Instruments



## Periodontal Instruments Index

File Scalers				P19
Gracey A5 Curettes				p17
Gracey M5 Curettes				p18
Gracey Prophy Curettes				P16
Gracey Standard Curettes				p14-p15
Hoe Scalers				P20
Periodontal Chisels				P25-P27
Periodontal Files and Chisels				P24
Periodontal Handle Options				P1
Periodontal Knives				P21-P22
Periodontal Surgical Curettes				P23
Scaler Curette Sets-Goldman-Fox				P9
Sckel Scalers				P2-P8
Universal Curettes				p10-p13

# PERIODONTAL

## HANDLE OPTIONS

### DOUBLE ENDED

#1		9.5mm (3/8")
	PATENTED	
#2		9.5mm (3/8")
	PATENTED	
#3		9.5mm (3/8")
	PATENTED	
#4		9.5mm (3/8")

### SINGLE ENDED

#1		9.5mm (3/8")
#2		9.5mm (3/8")
#3		9.5mm (3/8")
#4		9.5mm (3/8")

# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.



H6/H7  
Hygienist  
**SH6/7**  
2 6 3 0



6/7  
Goldman  
**SH6/7**  
2 6 3 1



H5/U15  
Hygienist-Towner  
**SH5/U15**  
2 6 3 2



U15/30  
Towner-Jacquette  
**SU15/30**  
2 6 3 3



H5/33  
Hygienist  
Jacquette  
**SH5/33**  
2 6 3 4



U15/33  
Towner-Jacquette  
**SU15/33**  
2 6 3 5



U15  
Towner  
**SU15**  
2636

# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.

### Nevi Scalars



Nevi1  
Anterior  
**SCNEV19**  
2 6 3 7



Nevi 1/H5  
Anterior  
**SNEV1/H59**  
2 6 3 8



Nevi2  
Posterior  
**SCNEV129**  
2 6 3 9



Nevi3  
Posterior  
**SCNEV139**  
2 6 4 0



Nevi4  
Posterior  
**SCNEV149**  
2 6 4 1

# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.

### Jacquette Scalers



30/33  
Jacquette  
**SJ30/33**  
2 6 4 2



31/32  
Jacquette  
**SJ31/32**  
2 6 4 3



34/35  
Jacquette  
**SJ34/35**  
2 6 4 4



H3/H4  
Jacquette  
**SH3/4**  
2 6 4 5



14/15  
Posterior  
Jacquette  
**SJ14/156**  
2 6 4 6

# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.



204

**S204**  
2647



204S

**S204S**  
2648



204SD

**S204SD**  
2649



204IUFW

**SIUFW204**  
2650



1  
Springfield  
Anterior

**SSP1**  
2651



2/3  
Springfield

**SSP2/3**  
2652

# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.

### Posterior Scalers



### Anterior Scalers



# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.

### Anterior Scalers



# PERIODONTAL

## SICKLE SCALERS

**Sickle scalers:** Sickle scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70 to 80 degree.



- |                 |                   |                              |                               |                                 |                    |                  |
|-----------------|-------------------|------------------------------|-------------------------------|---------------------------------|--------------------|------------------|
| 11/12<br>Scaler | 11A/12A<br>Scaler | 11/12S<br>Modified<br>McCall | 13S/14S<br>Pointed<br>Curette | MTC 13/14<br>Pointed<br>Curette | 17S/18S<br>Curette | 19/20<br>Curette |
| <b>PSS11/12</b> | <b>PSS11/12A</b>  | <b>PSS11/12S</b>             | <b>PSC13/14S</b>              | <b>PSCTC13/14</b>               | <b>PSC17/18</b>    | <b>PSC19/20</b>  |
| 2705            | 2706              | 2707                         | 2708                          | 2709                            | 2710               | 2711             |

# PERIODONTAL

## SCALAR / CURETTE SETS

One side is sickle scaler and the other is a universal curette.



### Goldman-Fox Set



GF-1  
Universal Scaler

**PSGF1**  
2760



GF-2  
Curette

**PSGF2**  
2761



GF-3  
Curette

**PSGF3**  
2762



GF-4  
Curette

**PSGF**  
2763



GF-5 HOE  
BUCCAL/lingual

**PSGF5**  
2764



GF-6  
Hoe

**PSGF6**  
2765

# PERIODONTAL

## UNIVERSAL CURETTES

Blades are sharpened on both sides. Blade curved at 90° degree to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

The Langer currettes combine the shank design of the gracey currettes with the 90° cutting edge of the universal currettes, edged on both sides. As the Langer curette 1-2 is applicable for the treatment of lower bicuspids and molars, and the Langer curette 3-4 for upper bicuspids and molars, as well as the Langer curette 5-6 for upper and lower front teeth, the model 17-18 is especially adapted for all surfaces in the area of upper an lower jaw, specifically at small stoma.



- |                          |                          |                        |                         |
|--------------------------|--------------------------|------------------------|-------------------------|
| 1/2<br>Langer<br>(shown) | 1/2<br>Langer<br>After 5 | 1/2<br>Langer<br>Rigid | 1/2<br>Langer<br>Mini 5 |
| <b>PUCL1/2</b>           | <b>PUCL1/2A5</b>         | <b>PUCL1/2R</b>        | <b>PUCSL1/2M5</b>       |
| 2820                     | 2821                     | 2822                   | 2823                    |



- |                          |                          |                        |                         |
|--------------------------|--------------------------|------------------------|-------------------------|
| 3/4<br>Langer<br>(shown) | 3/4<br>Langer<br>After 5 | 3/4<br>Langer<br>Rigid | 3/4<br>Langer<br>Mini 5 |
| <b>PUCL3/4</b>           | <b>PUCL3/4A5</b>         | <b>PUCL3/4R</b>        | <b>PUCL3/4M5</b>        |
| 2824                     | 2825                     | 2826                   | 2827                    |



- |                          |                        |                          |                         |
|--------------------------|------------------------|--------------------------|-------------------------|
| 5/6<br>Langer<br>(shown) | 5/6<br>Langer<br>Rigid | 5/6<br>Langer<br>After 5 | 5/6<br>Langer<br>Mini 5 |
| <b>PUCL5/6</b>           | <b>PUCL5/6R</b>        | <b>PUCL5/6A5</b>         | <b>PUCL5/6M5</b>        |
| 2828                     | 2829                   | 2830                     | 2831                    |



- |                            |                          |
|----------------------------|--------------------------|
| 17/18<br>Langer<br>(shown) | 17/18<br>Langer<br>Rigid |
| <b>PUCL17/18</b>           | <b>PUCL17/18R</b>        |
| 2832                       | 2833                     |

# PERIODONTAL

## UNIVERSAL CURETTES

Blades are sharpened on both sides. Blade curved at 90° degree to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

Anterior



IU13/14  
Indiana  
University

**PUCIU13/14**

2834

Posterior



IU17/18  
Indiana  
University

**PUCIU17/18**

2835

Curette-Chisel Combo



1/2  
Bunting

**PCCCB1/2**

2836



5/6  
Bunting

**PUCBN5/6**

2837



YG7/8  
Younger-Good

**PUCYG7-8**

2838

# PERIODONTAL

## UNIVERSAL CURETTES

Blades are sharpened on both sides. Blade curved at 90° degree to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

### Barnhart Universal Curettes



- 1/2  
Barnhart

**PUCBH1/2**  
2839
- 1/2  
Barnhart  
Rigid

**PUCBH1/2R**  
2840



- 5/6  
Barnhart

**PUCBH5/6**  
2841
- 5/6  
Langer M-5

**PUCBH5/6R**  
2842

### Columbia Universal Curettes



- 2R/2L  
Columbia

**PUC2R/2L**  
2843



- 4R/4L  
Columbia

**PUC4R/4L**  
2844



- 13/14  
Columbia

**PUC13/14**  
2845
- 13/14  
Columbia  
Rigid

**PUC13/14R**  
2846
- 13/14  
Columbia  
Curved Tip

**PUC13/14C**  
2847

# PERIODONTAL

## UNIVERSAL CURETTES

Blades are sharpened on both sides. Blade curved at 90° degree to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.



10/11 Loma Linda (Shown)	10/11 Loma Linda Rigid	10/11 Loma Linda Slim
<b>PUCLL10/11</b> 2853	<b>PUCLL10/11R</b> 2854	<b>PUCLL10/11S</b> 2855



413/414
<b>PUC413/414</b> 2856



1/2 Mallery
<b>PUCMAL1/2</b> 2857



31/32 Goldman
<b>PUCG31/32</b> 2858



ME-1 [10-11]
<b>PUCME1</b> 2859

# PERIODONTAL


## STANDARD GRACEY CURETTES

Curettes are used supragingivally and subgingivally on enamel and root surfaces. There are two major sub types of curettes; universal curettes and area specific curettes. Universal curettes can be used to remove medium-size and small size calculus deposits. Area specific curettes are used to remove light calculus and endotoxins from root surfaces. Gracey curettes is an area-specific curette. It has offset blade that provides a perfect working angulation for the teeth surfaces.



### Standard Gracey Curettes:

The exact application areas for the different Gracey figures you will find in the follow chart:



- For all surfaces of the front teeth roots.
- For all surfaces of the front teeth Bicuspids roots.
- For all facial and oral surfaces of the Bicuspids and molars.
- For all mesial surfaces of the Bicuspids and molars.
- For all mesial surfaces of the Bicuspids and molars.
- For all distal surfaces of the bicuspids and molars.
- For all distal surfaces of the bicuspids and molars.


1/2 (Shown) <b>PCG1/2</b> 2910	3/4 (Shown) <b>PCG3/4</b> 2912	5/6 (Shown) <b>PCG5/6</b> 2914	7/8 (Shown) <b>PCG7/8</b> 2917	9/10 (Shown) <b>PCG9/10</b> 2920
1/2 Rigid <b>PCG1/2R</b> 2911	3/4 Rigid <b>PCG3/4R</b> 2913	5/6 Rigid <b>PCG5/6R</b> 2915	7/8 Rigid <b>PCG7/8R</b> 2918	9/10 Rigid <b>PCG9/10R</b> 2921
		5/6 Extra Rigid <b>PCG5/6ER</b> 2916	7/8 Extra Rigid <b>PCG7/8ER</b> 2919	

# PERIODONTAL

## STANDARD GRACEY CURETTES



GRACEY CURETTE



11/12 (Shown)	13/14 (Shown)	15/16 (Shown)	17/18 (Shown)
<b>PCG11/12</b> 2922	<b>PCG13/14</b> 2925	<b>PCG15/16</b> 2928	<b>PCG17/18</b> 2930
11/12 Rigid	13/14 Rigid	15/16 Rigid	17/18 Rigid
<b>PCG11/12R</b> 2923	<b>PCG13/14R</b> 2926	<b>PCG15/16R</b> 2929	<b>PCG17/18R</b> 2931
11/12 Extra Rigid	13/14 Extra Rigid		
<b>PCG11/12ER</b> 2924	<b>PCG13/14ER</b> 2927		



Curettes are used supragingivally and subgingivally on enamel and root surfaces. There are two major sub types of curettes; universal curettes and area specific curettes. Universal curettes can be used to remove medium-size and small size calculus deposits. Area specific curettes are used to remove light calculus and endotoxins from root surfaces. Gracey curettes is an area-specific curette. It has offset blade that provides a perfect working angulation for the teeth surfaces.

# PERIODONTAL

## Gracey Mesial-Distal and Prophy Curettes

**Rigid Gracey:** The rigid curette has a stronger and less flexible shank than the standard one. So it can be used for heavy calculus removal without using additional scalers.

### Gracey Mesial-Distal Curettes

 11/14 Mesial Distal (Shown) <b>PGMDC11/14</b> 2932	 12/13 Mesial Distal (Shown) <b>PGMDC12/13</b> 2934
 11/14 Mesial Distal Rigid <b>PGMDC11/14R</b> 2933	 12/13 Mesial Distal Rigid <b>PGMDC12/13R</b> 2935


### Gracey Prophy Curettes

 P1/2 <b>PGPC1/2</b> 2936	 P3/4 <b>PGPC5/6</b> 2937	 P5/6 <b>PGPC3/4</b> 2938	 P7/8 <b>PGPC7/8</b> 2939
--	---	--	--

# PERIODONTAL

## A-5 GRACEY CURETTES

- A-5 Curettes:**
1. Terminal shank elongated by 3mm for access into deep periodontal pockets and root surfaces of 5mm or more and better clearance around the crown.
  2. Blade thinned by 10% to ease gingival insertion and reduce tissue distention.
  3. Area-specific designed patterned after the standard Graceys offers exact adaption and control.
  4. Available in Rigid and Extra Rigid designs for greater strength when needed.




1/2	3/4	5/6	7/8	11/12	13/14	15/16	11/14 Mesial/Distal	12/13 Mesial/Distal
<b>A5G1/2</b> 2940	<b>A5G3/4</b> 2942	<b>A5G5/6</b> 2944	<b>A5G7/8</b> 2947	<b>A5G11/12</b> 2950	<b>A5G13/14</b> 2953	<b>A5G15/16</b> 2956	<b>A5G11/14</b> 2958	<b>A5G12/13</b> 2960
Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid
<b>A5G1/2R</b> 2941	<b>A5G3/4R</b> 2943	<b>A5G5/6R</b> 2945	<b>A5G7/8R</b> 2948	<b>A5G11/12R</b> 2951	<b>A5G13/14R</b> 2954	<b>A5G15/16R</b> 2957	<b>A5G11/14R</b> 2959	<b>A5G12/13R</b> 2961
		Extra Rigid	Extra Rigid	Extra Rigid	Extra Rigid			
		<b>A5G5/6ER</b> 2946	<b>A5G7/8ER</b> 2949	<b>A5G11/12ER</b> 2952	<b>A5G13/14ER</b> 2955			

# PERIODONTAL

## M-5 GRACEY CURETTES

- M-5 Curettes:**
1. Terminal shank elongated by 3mm for access into deep periodontal pockets and root surfaces of 5mm or more and better clearance around the crown.
  2. Blade thinned by 10% to ease gingival insertion and reduce tissue distention.
  3. Area-specific designed patterned after the standard Graceys offers exact adaption and control.
  4. Available in Rigid and Extra Rigid designs for greater strength when needed.



1/2	3/4	5/6	7/8	11/12	13/14	15/16	11/14 Mesial/Distal	12/13 Mesial/Distal
<b>M5G1/2</b> 2962	<b>M5G3/4</b> 2964	<b>M5G5/6</b> 2966	<b>M5G7/8</b> 2968	<b>M5G11/12</b> 2970	<b>M5G13/14</b> 2972	<b>M5G15/16</b> 2974	<b>M5G11/14</b> 2976	<b>M5G12/13</b> 2977
Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid		
<b>M5G1/2R</b> 2963	<b>M5G3/4R</b> 2965	<b>M5G5/6R</b> 2967	<b>M5G7/8R</b> 2969	<b>M5G11/12R</b> 2971	<b>M5G13/14R</b> 2973	<b>M5G15/16R</b> 2975		

# PERIODONTAL

## FILE SCALERS

**Hoe Scalers:** Use to dislodge heavy supramarginal calculus deposits use a pull stroke.

**File Scalers:** Use to crush and remove heavy calculus deposits. Use a pull stroke.

### File Scalers



# PERIODONTAL

## HOE SCALERS

**Hoe Scalars:** Use to dislodge heavy supramarginal calculus deposits use a pull stroke.  
**File Scalars:** Use to crush and remove heavy calculus deposits. Use a pull stroke.

Hoe Scalars



XS-15  
Towner  
**PHSXS15**  
3035

7K  
**PHS7K**  
3036

6/7  
Orban  
Mesial/Distal  
**PHSO6/7**  
3037

8/9  
Orban  
Buccal/Lingual  
**PHSO8/9**  
3038

U15/13R  
**PHSU15/13R**  
3039

# PERIODONTAL

## PERIODONTAL KNIVES

For incision that remove or recontour soft tissue. Also useful to excise inter-proximal tissue.



15K/16K  
Kirkland  
**PK15K/16K**  
3080



7  
Goldman-Fox  
**PKGF7**  
3081



8  
Goldman-Fox  
**PKGF8**  
3082



9  
Goldman-Fox  
**PKGF9**  
3083



11  
Goldman-Fox  
**PKGF11**  
3084



3/4  
Buck  
**PKB3/4**  
3085



5/6  
Buck  
**PKB5/6**  
3086

# PERIODONTAL

## PERIODONTAL KNIVES

For incision that remove or recontour soft tissue. Also useful to excise inter-proximal tissue.



1/2  
Solt

**PKSOLT1/2**  
3087



1/2  
Sanders

**PKS1/2**  
3088



1/2  
Orban

**PKO1/2**  
3089



USC  
Towner  
19/20

**PKUSCT19/20**  
3090



1/2  
Orban  
Knife

**PKOK1/2**  
3091

30° thinner than 1/2 Orban Knife



CK-3  
Crane  
Kaplan

**PKCK3**  
3092



Dr Allen  
End-Cutting  
Instrument  
Knife

**PKDA**  
3093

# PERIODONTAL

## PERIODONTAL SURGICAL CURETTES

### Kramer

### Sugrman



PR1/2  
Prichard  
**PSCPR1/2**  
3 140

KRA-1  
Kramer  
Mesial  
**PSCKRA1**  
3 141

KRA-2  
Kramer  
Distal  
**PSCKRA2**  
3 142

KRA-3  
Kramer  
**PSCKRA3**  
3 143

3/4  
Kramer  
Nevins  
**PSCKN3/4**  
3 144

8K/9K  
Kirkland  
**PSCK8/9K**  
3 145

5S/6S  
Sugrman  
Curette  
**PSCS5/6S**  
3 146

9S/10S  
Sugrman  
Curette  
**PSCS9/10S**  
3 147

11S/12S  
Sugrman  
Curette  
**PSCS11/12S**  
3 148

# PERIODONTAL

## PERIODONTAL FILES & CHISELS

**Periodontal Files:** File surfaces on both sides allow for push/pull application. These periodontal files are used interproximally.

**Periodontal Chisels:** These stainless steel chisels are used for the moving and reshaping bone.


Periodontal Files

Kirkland Periodontal Chisels



# PERIODONTAL

## PERIODONTAL CHISELS



1/2  
Kramer  
Nevis  
**PCKN1/2**  
3206

36/37  
Rhodes  
Rack-Action  
**PC36/37**  
3207

44/45  
Remington  
Hoe  
**PC44/45**  
3208

TG  
**PCTG**  
3209

TG-0  
Modified TG  
with Ochdenbein  
Cutouts  
**PCTG-0**  
3210

3/4  
Solt  
**PCSOLT3/4**  
3211

# PERIODONTAL

## PERIODONTAL CHISELS

### Fedi Chisels



### Ochsenbein Chisels



# PERIODONTAL

## PERIODONTAL CHISELS

### Kramer-Nevins Chisels



### Buser Chisels



### Ochsenbein Chisels



3 Kramer-Nevins  
7.5mm  
PKN3

**PKN3**

3221

55 Kramer-Nevins  
5.5mm  
PKN55

**PKN55**

3222

75 Kramer-Nevins  
7.5mm  
PKN75

**PKN75**

3223

4 Buser  
Round Tip-4mm  
PBS4

**PBS4**

3224

6 Buser  
Pointed Tip-6mm  
PBS6

**PBS6**

3225

1 Ochsenbein  
PCO1

**PCO1**

3226

2 Ochsenbein  
PCO2

**PCO2**

3227

## Restorative Instruments



## Restorative Instruments Index

Amalgam Carriers				R16-R18
Articulating Paper Forceps				R15
Burnishers				R28-R29
Carvers				R23-R25
Chisels & Margin Trimmers				R8
Copposite/Plastic Filling Instruments				R4-R6
Crown Splittler & Remoer Instruments				R32-R33
Discoid-Cleoids Carvers				R26-R27
Excavators				R11-R14
Gingival Cord Packers				R31
Ginival Retractors				R30
GoldFogle Composite/Plastic Filling Instruments				R3
Goldstein Composite Instruments				R1
Hatchets				R7
Hoes & Angle Formers				R9
Ipression Trays				R35
Margin Trimmers				R10
Placement Instruments				R15
Pluggers & Condensors				R19-R22
Spatulas				R34

# RESTORATIVES

## Goldstein Flexible-Thin Composite Instruments

Flexible, Thin, Highly Polished, Non Stick Stainless Steel Blades.




1  
RCIGFT1  
RCIGFT1  
3230




2  
RCIGFT2  
RCIGFT2  
3231




6  
RCIGFT6  
RCIGFT6  
3235



3  
RCIGFT3  
RCIGFT3  
3232



Micro-Mini  
RTNCIPCS  
RTNCIPCS  
3236



4  
RCIGFT4  
RCIGFT4  
3233



1(Mini)  
RCIGFTM1  
RCIGFTM1  
3237



5  
RCIGFT5  
RCIGFT5  
3234



3(Mini)  
RCIGFTM3  
RCIGFTM3  
3238



4(Mini)  
RCIGFTM4  
RCIGFTM4  
3239

# RESTORATIVES

## Composite/Plastic Filling Instruments

Highly Polished Stainless Steel Instruments Used for Composite Placement and Contouring



AB1  
Boghosian  
RPF1AB1  
3245



AB2  
Boghosian  
RPF1AB2  
3246



Interproximal Carver  
RPF1IC  
3250




W3  
RPF1W3  
3247



3Tufts  
RPF13T  
3251



4F  
Felt/Tufts  
RPF14F  
3248



A6 (156)  
RPF1A6  
3252



UF  
RPF1UF  
3249



4/5 Gregg  
RPF14/5G  
3253



PFI  
RPF1  
3254


## GOLDFOGEL COMPOSITE/PLASTIC FILLING INSTRUMENTS

### ANTERIOR




A  
Cosmetic Contouring

RCCIA  
3255



B  
Cosmetic Contouring

RCCIB  
3256




E  
Cosmetic Contouring

RCCIE  
3259




C  
Cosmetic Contouring

RCCIC  
3257




F  
Cosmetic Contouring

RCCIF  
3260



D  
Cosmetic Contouring

RCCID  
3258



G  
Cosmetic Contouring

RCCIG  
3261



H  
Cosmetic Contouring

RCCIH  
3262



I  
Cosmetic Contouring

RCCII  
3263

# RESTORATIVES

## COMPOSITE / PLASTIC FILLING INSTRUMENTS

### ANTERIOR



11  
RPF111  
3265

179  
RPF1179  
3266

G  
RPF1G  
3267

6SH  
RPF16SH  
3268

1 Woodson  
RPF1WDS1  
3269

8A  
RPF18A  
3264

PF1  
Heidemann  
Spatula  
Small  
RPF1HS1  
3270

PF2  
Heidemann  
Spatua  
Medium  
RPF1HM2  
3271

PF3  
Heidemann  
Large  
RPF1HL3  
3272

48  
RPF148  
3273

W1  
RPF1W1  
3274

118  
RPF1118  
3275

# RESTORATIVES

## COMPOSITE / PLASTIC FILLING INSTRUMENTS

Martensite stainless steel instruments used for composite placement and contouring.

### POSTERIOR



Novadent 2  
RPFIND2  
3276



Interproximal  
Off-Angle  
RCVIPCOA  
3277



49 Baldwin  
RPF149  
3278



UP1  
RPF1UP1  
3279



21  
Burnisher  
RBB21  
3280



21B  
Burnisher  
RBB21B  
3281



PKT-3R  
Round Cone  
RPKT3R  
3282

# RESTORATIVES

## COMPOSITE / PLASTIC FILLING INSTRUMENTS

Martensite stainless steel instruments used for composite placement and contouring.

### UNIVERSAL



1  
Loesche  
RPFIL1  
3286



2  
Loesche  
RPFIL2  
3287



3T  
Tufts  
RPFIL3T  
3288



1  
RPFIL1  
RPFIL1  
3289



1  
Woodson  
RPFWD1  
3290



2  
Woodson  
RPFWDS2  
3291



3  
Woodson  
RPFWDS3  
3292



S4  
Woodson  
RPFWDS4  
3293

BEQ1

RPFIBEQ1  
3285

# RESTORATIVES

## HATCHETS

Martensite stainless steel instruments used for composite placement and contouring.

**8/9H**  
[10-7-14]  
NOVADENT USA

**Suggested Pair.**

**14/14**  
Off Angle  
[15-8-14]  
**RH14/14**  
3401

**15/15**  
Off Angle  
[15-8-14]  
**RH15/15**  
3402

**14/14-0**  
Off Angle  
[15-10-16]  
**RH14/14-0**  
3403

**15/15-0**  
Off Angle  
[15-10-16]  
**RH15/15-0**  
3404

**44S**  
Off Angle  
**RH44S**  
3405

**45S**  
Off Angle  
**RH45S**  
3406

**13/14**  
[20-9-14]  
**RH13/14**  
3407

**15/16**  
[15-8-14]  
**RH15/16**  
3408

**17/18**  
[10-6-14]  
**RH17/18**  
3409

**51/52**  
[15-8-12]  
**RH51/52**  
3410

**53/54**  
[10-6-12]  
**RH53/54**  
3411

**11/12**  
[15-8-8]  
**RH11/12**  
3412

# RESTORATIVES

## CHISELS & MARGIN TRIMMERS

Margin Trimmer: These instruments used to produce proper bevel on enamel margins. These tools are similar to a hatchet except the blade curved and the cutting edge angled.

Ceramic crowns and Onlays-Chisels/Margin trimmers

Chisels			Margin Trimmers			
Used to refine the cavity preparation						
1/2 Wedelstaedt [20-15-3]	3/4 Wedelstaedt [11.5-15-3]	5/6 Wedelstaedt [15-15-3]	7/10 Straight [15][20]	8/9 Binangle [20-9-8]	11-12 Binangle [15-8-8]	40/41 Binangle [18-10-16]
<b>CP1/2</b> 3 4 4 0	<b>CP3/4</b> 3 4 4 1	<b>CP5/6</b> 3 4 4 2	<b>CP7/10</b> 3 4 4 3	<b>CP8/9</b> 3 4 4 4	<b>CP11/12</b> 3 4 4 5	<b>CP40/41</b> 3 4 4 6

# RESTORATIVES


## HOES & ANGLE FORMERS

**Angle Former:** For defining line angles, obtaining retentive form in dentin and placing bevels on enamel margins.

### Angle Formers

### Hoes

Forming line angles on anterior preparations.



# RESTORATIVES

## MARGIN TRIMMERS

**Margin Trimmers:** Used to produce proper bevel on enamel margins. Similar to a hatchet except the blade is curved and the cutting edge angled.

Margin Trimmers




# RESTORATIVES

## EXCAVATORS

These instruments made of marten stainless steel and used for carious dentin.

### Spoons




# RESTORATIVES

## EXCAVATORS

These instruments made of marten stainless steel and used for carious dentin.

### Oval Spoons




# RESTORATIVES

## EXCAVATORS

These instruments made of martensite stainless steel and used for carious dentin.  
These instruments are English patterns.

Blades

English Patterns




36/37 [15-8-14]	40-41 [11-7-14]	44/45 [20-9-14]	61/62 [20-9-12]	63/64 [15-8-12]	65/66 [10-6-12]	125/126	127/128
<b>EXC36/37</b> 3 5 6 3	<b>EXC40/41</b> 3 5 6 4	<b>EXC44/45</b> 3 5 6 5	<b>EXC61/62</b> 3 5 6 6	<b>EXC63/64</b> 3 5 6 7	<b>EXC65/66</b> 3 5 6 8	<b>EXC125/6</b> 3 5 6 9	<b>EXC127/8</b> 3 5 7 0

# RESTORATIVES

## EXCAVATORS

These instruments made of martensite stainless steel and used for carious dentin.  
 These instruments are English patterns.



# RESTORATIVES

## PLACEMENT INSTRUMENTS (Aesthetics)

These instruments are used to place base or liner within cavity preparations.



Calcium  
Hydroxide  
Placer

**PIC**  
3630

6061  
Mini Spatula/  
Placer

**SP6061**  
3631

10  
Novatech  
Placer

**PINT10**  
3632

11  
Novatech  
Placer

**PINT11**  
3633

1  
Composite  
Brush Handle

**HCB1**  
3634

Miller  
Articulating  
Paper Forcep

**APF2**  
3635

Miller  
Articulating  
Paper Forcep

**APF4**  
3636

Miller Mod  
Articulating  
Paper Forcep

**APF3**  
3637

## AMALGAM CARRIERS

Single Ended  
Steel & Teflon Tips

Amalgam Carrier 1.5mm, Teflon Tip

ACSTF1.5

3665



Amalgam Carrier 2.0mm, Teflon Tip

ACSTF2.0

3666



Amalgam Carrier 2.5mm, Teflon Tip

ACSTF2.5

3667



Amalgam Carrier 3.2mm, Teflon Tip

ACSTF3.2

3668



Amalgam Carrier 1.5mm, SS Tip

ACSS1.5

3669



Amalgam Carrier 2.0mm, SS Tip

ACSS2.0

3670



Amalgam Carrier 2.5mm, SS Tip

ACSS2.5

3671



Amalgam Carrier 3.2mm, SS Tip

ACSS3.2

3672

# RESTORATIVES

## AMALGAM CARRIERS



Amalgam Carrier Mini 1.5mm/1.5mm, Teflon Tip

ACMTF1.5/1.5  
3673



Amalgam Carrier Regular 2/2mm, Teflon Tip

ACRTF2/2  
3674



Amalgam Carrier Large 2.5/2.5mm, Teflon Tip

ACLTF2.5/2.5  
3675



Amalgam Carrier Jumbo 3.2/3.2mm, Teflon Tip

ACJTF3.2/3.2  
3676



Amalgam Carrier Mini 1.5/1.5mm, S. Tip

ACMSS1.5/1.5  
3677



Amalgam Carrier Regular 2/2mm, S. Tip

ACRSS2/2  
3678



Amalgam Carrier Large 2.5/2.5mm, S. Tip

ACLSS2.5/2.5  
3679



Amalgam Carrier Jumbo 3.2/3.2mm, S. Tip

ACJSS3.2/3.2  
3680

# RESTORATIVES

## Amalgam Carriers



Amalgam Carrier Mini 1.5/Regular 2.0mm, Teflon Tip

ACDTF1.5/2.0  
3681



Amalgam Carrier Large 2.5mm/Jumbo 3.2mm, Teflon Tip

ACDTF2.5/3.2  
3682



Amalgam Carrier Mini 1.5mm/Regular 2.0mm, SS. Tip

ACDSS1.5/2.0  
3683




Amalgam Carrier Large 2.5mm/Jumbo 3.2mm, SS Tip

ACDSS2.5/3.2  
3684



Amalgam Carrier, Chrome Plated with SS tip 2mm

ACCPSS2  
3685



Amalgam Carrier, Chrome Plated 2mm with Teflon Tip

ACCTF2  
3686



Amalgam Carrier, Chrome Plated 3.2mm with SS Tip

ACCPSS3.2  
3687



Amalgam Carrier, Chrome Plated 3.2mm with Teflon Tip

ACCTF3.2  
3688




Amalgam Well

AW  
3695

# RESTORATIVES

## PLUGGERS / CONDENSERS

The tips of instruments are very precisely machined and used to pack filling materials into cavity preparations.




½ Andrews	½ Black [15-7-12][25-7-12]	½ Serrated	H1 Hollenback	H2 Hollenback	H3 Hollenback	H4 Hollenback	0/1 Marquette Condenser	0/1 Marquette
<b>PLGA1/2</b> 3750	<b>PLG1/2NS</b> 3751	<b>PLG1/2</b> 3752	<b>PLGH1</b> 3753	<b>PLGH2</b> 3754	<b>PLGH3</b> 3755	<b>PLGH4</b> 3756	<b>PLG0/1</b> 3757	<b>PLG0/1NS</b> 3758

# RESTORATIVES

## PLUGGERS / CONDENSERS

The tips of instruments are very precisely machined and used to pack filling materials into cavity preparations.




2 Mortonson	2 Smith Serrated	SPO Serrated	1M Markley	2M Markley	4M Markley	5M Markley	6M Markley
<b>PLGMO2</b> 3759	<b>PLGS2</b> 3760	<b>PLGSP0</b> 3761	<b>PLG1M</b> 3762	<b>PLG2M</b> 3763	<b>PLG4M</b> 3764	<b>PLG5M</b> 3765	<b>PLG6M</b> 3766

# RESTORATIVES

## PLUGGERS / CONDENSERS

The tips of instruments are very precisely machined and used to pack filling materials into cavity preparations.



# RESTORATIVES

## PLUGGERS / CONDENSERS

The tips of instruments are very precisely machined and used to pack filling materials into cavity preparations.



# RESTORATIVES

## CARVERS (Amalgam)

These martensite stainless steel instruments are used to carve anatomical features and trim excess materials.



Inter-  
proximal

**CVIPC**  
3830



Inter-  
proximal  
Off-Angle

**CVIPCOA**  
3831



1/2-3  
Hollenback

**CVHL1/2**  
3832



3S  
Hollenback

**CVHL3S**  
3833



3  
Hollenback

**CVHL3**  
3834




8  
Willand  
Minn. IPC1

**CVW186**  
3835

# RESTORATIVES

## CARVERS (Amalgam)

These martensite stainless steel instruments are used to carve anatomical features and trim excess materials.



1  
Andrew

2  
Andrew

3  
Andrew

4  
Clappison

2/3  
Frahm

104

1 Hollenback  
2 Modified  
Cleoid Small

1 Hollenback  
3 Modified  
Cleoid Small

**CVA1**  
3836

**CVA2**  
3837

**CVA3**  
3838

**CVCL4**  
3839

**CVFR2/3**  
3840

**CV104**  
3841

**CVH1/MC2**  
3842

**CVH1/MC3**  
3843

# RESTORATIVES

## CARVERS (Amalgam)

These martensite stainless steel instruments are used to carve anatomical features and trim excess materials.




# RESTORATIVES

## CARVERS/DISCOID-CLEOIDS

These martensite stainless steel instruments are used to carve anatomical features and trim excess materials.

### Carvers



3  
Tanner

**CV3T**  
3851

4  
Tanner

**CV4T**  
3852

5  
Tanner

**CV5T**  
3853

6  
Tanner

**CV6T**  
3854

1  
Wall

**CVWL1**  
3855

3  
Wall

**CVWL3**  
3856

1  
Ward

**CVWR1**  
3857

2  
Ward


**CVWR2**  
3858

# RESTORATIVES

## CARVERS/DISCOID-CLEOIDS

These martensite stainless steel instruments are used to carve anatomical features and trim excess materials.

### Discoid-Cleoids



# RESTORATIVES

## BURNISHERS

Highly precision tools are used to condense, smooth, carve and polish amalgam.



# RESTORATIVES

## BURNISHERS

Highly precision tools are used to condense, smooth, carve and polish amalgam.

Acorn

Ladmore



2.7

1.7



3.0

2.7



1.9

1.3



1.9

1.3



21  
**BB21**  
3916

21B  
**BB21B**  
3917

2  
Ladmore  
**BBL2**  
3918

PKT-3R  
Rounded  
Cone  
**PKT3R**  
3919

Romerowski  
**BBROM**  
3920

Ladmore  
**BBL3**  
3921

29/34  
**RBB29/34**  
3922

# RESTORATIVES

## GINGIVAL RETRACTORS

These instruments are used to protect tissue during cavity preparations and placement of composite restorations.



Kincheloe

**GRK1**

3970



GP 10  
Golkman-Fox

**TRGF10**

3971



1

Meinershagen

**GRM1**

3972



2

Meinershagen

**GRM2**

3973



3

Meinershagen

**GRM3**

3974



4

Meinershagen

**GRM4**

3975

# RESTORATIVES

## GINGIVAL CORD PACKERS

Tips are made of martensite stainless steel of these instruments for atraumatic and accurate cord placement.



CSI-1  
Serrated

**GCPCSI1**  
4030



113  
Serrated

**GCP113**  
4031



BN1

**GCRBN1**  
4032



S6

**GCPS6**  
4033



Balshi

**GCPBAL**  
4034



1  
Yardley

**GCPYD1**  
4035



7  
Guyer  
Serrated

**GCPG7**  
4036

# RESTORATIVES

## Crown Splitter & Remover Instruments



Crown Spreader Forceps

RCSF  
4040



Morrel Crown Remover  
W/3 Tips and Key

RMCR  
4044



Dental Crown Gripper, TC

RDCGTC  
4042



Automatic Gun Type Crown Remover  
with 3 Tips and key

RAGTCR  
4046



Dental Crown Splitter

RDCS  
4048

# RESTORATIVES

## SPATULAS



Christensen Crown Remover  
90°

RCCR90  
4050



Christensen Crown Remover  
Straight

RCCRS  
4052



Crown Spreader

RCS  
4054



Goldstein  
Crown Remover  
Straight

RGSCRS  
4056



Goldstein  
Crown Remover  
45° Angle

RGSCR45  
4058



Goldstein  
Crown Remover  
90°

RGSCR90  
4060



Goldstein  
Occlusal Surface

RDSOCS  
4062

# RESTORATIVES

## SPATULAS

Very accurately punched blades are flexible and are used to mix and load cement and other materials into crowns or inlay/onlay preparations.



24 Flexible 1<sup>3</sup>/<sub>4</sub>"  
44mm

**CS24**

4 0 9 0



324 Rigid2"  
51mm

**CS324**

4 0 9 1



A6 Rigid1"  
25mm

**CSA6**


4 0 9 2



5 Novatech  
Long/Fluted

**CSNT5**

4 0 9 3



6 Novatech  
Long/Blade

**CSNT6**

4 0 9 4



7 Novatech  
Short/Blade

**CSNT7**


4 0 9 5



8 Novatech  
Long

**CSNT8**

4 0 9 6



9 Novatech  
Short

**CSNT9**

4 0 9 7

# RESTORATIVES

## IMPRESSION TRAYS

Used to mix and load cement and other materials into crowns or inlays/onlays preparations

Set of 10 ea  
**RITS10**  
4100

Set of 8 ea  
**RITS8**  
4101



Upper XS  
Upper S  
Upper M  
Upper L  
Upper XL

Set of 6 ea  
**RITS6**  
4102

Set of 10EA  
Upper-XS,S,M,L,XL  
Lower-XS,S,M,L,XL

Set of 8EA  
Upper-XS,S,M,L  
Lower-XS,S,M,L

Set of 6EA  
Upper-XS,S,M,L  
Lower-XS,S,M,L



Lower XS  
Lower S  
Lower M  
Lower L  
Lower XL

Set of 10 ea  
**RITPS10**  
4103

Set of 8 ea  
**RITPS8**  
4104



Upper XS  
Upper S  
Upper M  
Upper L  
Upper XL

Set of 6 ea  
**RITPS6**  
4105


Set of 10EA  
Upper-XS,S,M,L,XL  
Lower-XS,S,M,L,XL


Set of 8EA  
Upper-XS,S,M,L  
Lower-XS,S,M,L

Set of 6EA  
Upper-XS,S,M,L  
Lower-XS,S,M,L



Lower XS  
Lower S  
Lower M  
Lower L  
Lower XL

		U-1(XS)	U-2(S)	U-3(M)	U-4(L)	U-5(XL)
	A	52	54	58	59	63
	B	65	67	72	75	84

		L-1(XS)	L-2(S)	L-3(M)	L-4(L)	L-5(XL)
	A	50	53	57	66	60
	B	65	68	75	78	84

Set of 8 ea  
**RITS8**  
4110


Set of 8EA  
Upper-S,M,XL  
Lower-S,M,XL



Upper S  
Upper M  
Upper L  
Upper XL



Lower S  
Lower M  
Lower L  
Lower XL

		U-1(S)	U-2(M)	U-3(L)	U-4(XL)
	A	58	63	64	68
	B	64	70	77	78

Set of 8 ea  
**RITPS8**  
4112


Set of 8EA  
Upper-S,M,XL  
Lower-S,M,XL



Upper S  
Upper M  
Upper L  
Upper XL



Lower S  
Lower M  
Lower L  
Lower XL

		L-1(S)	L-2(M)	L-3(L)	L-4(XL)
	A	56	60	62	63
	B	65	72	75	83



Upper  
65x39 mm  
**RITU**  
4114



Lower  
64x39 mm  
**RITL**  
4116



Upper  
65x39 mm  
**RITPU**  
4118



Lower  
64x39 mm  
**RITPL**  
4120



Upper Right  
Lower Left  
**RITPUR**  
4122



Upper Left  
Lower Right  
**RITPLR**  
4123



Adjustable  
Perforated  
**RITPA**  
4124



Adjustable  
Solid  
**RITSA**  
4125



Bite Checa Tray  
**RITBCT**  
4130

## Surgical Instruments



# Surgical Instruments Index

Aspirators			S46
Bone Chisels			S23
Bone Files			S24-S25
Cheak &Tissue Retractors			S43-S44
Elevator Luxating			S6
Elevators Apexo			S8
Elevators Bernard			S3
Elevators Cogswell & Crane Pick			S8
Elevators Coupland & Macmillan			S9
Elevators cryer			S5
Elevators EZ Luxators			S1-S2
Elevators Heidbrink			S4
Elevators Lindo			S3
Elevators Miscelanoius			S10-S11
Elevators Periotomes			S16
Elevators Root Tip Picks			S13-S15
Elevators Seldin			S7
Elevators T-Bar Handle			S12
Extracting Forceps -American Styles			S50-S53
Extracting Forceps Atraumatic			S49
Extracting Forceps- European Style			S55-S59
Extracting Forceps-Children			S60
Hemostatic Forceps			S41
Matrix Retainers			S48
Measuring Instruments			S48
Mouth Gauges			S43
Mouth Retractors			S42
Needle Holders			S34
Needle Holders -TC Tipped			S35-S39
Periosteals			S17-S20
Rongeurs			S61-S62
Scalpel Handles			S45
Scissors			S26-S29
Scissors-TC Tipped			S30-S31
Surgical Curettes			S21-S22
Syringes			S47
Tissue Pliers & Forceps			S32
Tissue Pliers and Forceps-TC Tipped			S33

## EZ-LUXATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.



Spade Elevator, Straight

**SELSS**  
4180




Spade Elevator, Mesially Angled

**SELSMA**  
4181




Spade Elevator, Distally Angled

**SELSDA**  
4182



Spade Elevator, Left

**SELSL**  
4183



Spade Elevator, Right

**SELSR**  
4184



2.5mm, Small, Straight

**SELSS2.5**  
4185



2.5mm, Small, Curved

**SELS2.5**  
4186



2.5mm, Small, Mesially Angled

**SELSMA2.5**  
4187



2.5mm, Small, Distally Angled

**SELSDA2.5**  
4188

## EZ-LUXATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.



3.75mm, Large, Straight

**SELLS3.75**  
4190



3.75mm, Large, Curved

**SELLC3.75**  
4191



3.75mm, Large, Mesially Angled

**SELLMA3.75**  
4192



3.75mm, Large, Distally Angled

**SELLDA3.75**  
4193



4.5mm, Large, Straight

**SELLS4.5**  
4194



4.5mm, Large, Curved

**SELLC4.5**  
4195



4.5mm, Large, Mesially Angled

**SELLMA4.5**  
4196



4.5mm, Large, Distally Angled

**SELLDA4.5**  
4197



Small, Tight, Mesial Angle, 2.5mm

**SELSTMA2.5**  
4198



Tight Mesial Spade, 2.5mm

**SELTMS2.5**  
4199

## BERNARD & LINDO



Lindo Levin, Straight  
3mm  
**LLS3**  
4204

Lindo Levin, Straight  
4mm  
**LLS4**  
4205



11 Straight Bernard  
SBERS11  
**BS11**  
4207



6 Left Bernard  
SBERL6  
**BR6**  
4208



7 Right Bernard  
SBERR7  
**BL7**  
4209



Lindo Levin, Curved  
3mm  
**LLC3**  
4206



Bernard Spade, Offset IN  
SBERO-IN  
**BSOIN**  
4210



Bernard Spade, Offset OUT  
SBERO-OUT  
**BSOOT**  
4211

## HEIDBRINK

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Heidbrink



1  
Heidbrink  
Straight

**E1**  
4213



2  
Heidbrink  
Straight

**E2**  
4214



46

**E46**  
4217



46R

**E46R**  
4218



40  
Heidbrink  
Straight

**E40**  
4215



41  
Heidbrink  
Straight

**E41**  
4216



77

**E77**  
4219



77R

**E77R**  
4220

## CRYER

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Suggested Pair.



21 Cryer Mini

EM21  
4221

22 Cryer Mini

EM22  
4222

Cryer:

Cryer's right and left elevators with triangular sharp-pointed blades and hollow handles.



44 Cryer Small

E44  
4223



45 Cryer Small

E45  
4224



25 Cryer Medium

E25  
4225



26 Cryer Medium

E26  
4226



A30 Cryer

EA30  
4227



A31 Cryer

EA31  
4228



31 Cryer Large

EA31  
4229



32 Cryer Large

EA32  
4230



## LUXATING ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

**Luxating Elevator:** The luxating elevators have thin, sharp blades that are used to cut periodontal ligaments. These instruments are used by rocking and cutting verses and elevator that lifts and pries. The importance of luxating elevator is to reduce pressure on adjacent tooth.

### LUXATOR ELEVATORS



Luxating Elevator  
Straight, 2mm

**SLES2**  
4231



Luxating Elevator  
Curved, 2mm

**SLEC2**  
4232



Luxating Elevator  
Straight, 2.5mm

**SLES2.5**  
4233



Luxating Elevator  
Curved, 2.5mm

**SLEC2.5**  
4234



Luxating Elevator  
Straight, 3mm

**SLES3**  
4235



Luxating Elevator  
Curved, 3mm

**SLEC3**  
4236



Luxating Elevator  
Straight, 4mm

**SLES4**  
4237



Luxating Elevator  
Curved, 4mm

**SLEC4**  
4238



Luxating Elevator  
Straight, 5mm

**SLES5**  
4239



Luxating Elevator  
Curved, 5mm

**SLEC5**  
4240

## SELDIN ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Seldin



## APEXO ELEVATORS

### Apexo



71 Miller  
Apexo

**E71**  
4252

72 Miller  
Apexo

**E72**  
4253

73 Miller  
Apexo

**E73**  
4254

74 Miller  
Apexo

**E74**  
4255



4  
Schmeckebeer  
Apexo

**E4**  
4256

5  
Schmeckebeer  
Apexo

**E5**  
4257

### Cogswell & Crane Pick



8  
Crane Pick

**E8**  
4261

A  
Cogswell

**EA**  
4262

B  
Cogswell

**EB**  
4263



301  
Apexo

**E301**  
4258

302  
Apexo

**E302**  
4259

303  
Apexo

**E303**  
4260

## COUPLAND & MACMILLAN

Coupland Gouge

 1C Coupland Gouge <b>EG1C</b> 4264	 2C Coupland Gouge <b>EG2C</b> 4265	 3C Coupland Gouge <b>EG3C</b> 4266
 A1C Coupland Gouge <b>EGA1C</b> 4267	 A2C Coupland Gouge <b>EGA2C</b> 4268	 A3C Coupland Gouge <b>EGA3C</b> 4269

Macmillan

 11M MacMillan Gouge <b>EG11M</b> 42611	 12M MacMillan Gouge <b>EG12M</b> 42621	 13M MacMillan Gouge <b>EG13M</b> 42631
---	--	--

**SURGICAL**

**MISCELLANEOUS ELEVATORS**



19

**E19**  
4270

20

**E20**  
4271

23

**E23**  
4272

**Friedman**



31F  
Friedman  
[Craig]  
Gouge Serrated

**E31F**  
4278

32F  
Friedman  
[Craig]  
Gouge Serrated

**E32F**  
4279

**Curtis**



51 Curtis

**EC51**  
4280

52 Curtis

**EC52**  
4281



38  
Serrated

**E38**  
4273

39  
Serrated

**E38**  
4274

**Bein**



559/1  
Bein, Large

**EBE559/1**  
4275

559/2  
Bein, Medium

**EBE559/2**  
4276

559/3  
Bein, Small

**EBE559/3**  
4277

**Hourigan**



Suggested  
Pair

H5  
Hourigan, Serrated

**H5S**  
4282

H6  
Hourigan, Serrated

**H6S**  
4283

# SURGICAL


## MISCELLANEOUS ELEVATORS

		Woodward				Stout			
		Suggested Pair							
									
6X Potts	7X Potts	15 Woodward	16 Woodward	11 Stout	11A Stout	80	81		
<b>E6X</b> 4284	<b>E7X</b> 4285	<b>EWW15</b> 4286	<b>EWW16</b> 4287	<b>E11</b> 4288	<b>E11A</b> 4289	<b>E80</b> 4290	<b>E81</b> 4291		

						
81A	81E	92 Serrated	93 Serrated	98	190	191
<b>E81A</b> 4292	<b>E81E</b> 4293	<b>E92</b> 4294	<b>E93</b> 4295	<b>E98</b> 4296	<b>E190</b> 4297	<b>E191</b> 4298

# SURGICAL

## T-BAR HANDLE ELEVATORS



7X Potts <b>E7XTBH</b> 4301	3 Hourigan <b>E3HTBH</b> 4302	4 Hourigan <b>E4HTBH</b> 4303	1LX <b>E1LXCBH</b> 4307	1RX <b>E1RXC BH</b> 4308	11LX <b>E11LXCBH</b> 4309	11RX <b>E11RXC BH</b> 4310
--------------------------------------	-------------------------------------	-------------------------------------	-------------------------------	--------------------------------	---------------------------------	----------------------------------

Warwick James Straight <b>WWJS</b> 4304	Warwick James Left <b>WWJL</b> 4305	Warwick James Right <b>WWJW</b> 4306	12LX <b>E12LXCBH</b> 4311	12RX <b>E12RXC BH</b> 4312	14LX <b>E14LXCBH</b> 4 3 1 3	14RX <b>E14RXC BH</b> 4 3 1 4
--	--	---	---------------------------------	----------------------------------	------------------------------------	-------------------------------------

6X  
**E6XTBH**  
4300

Warwick-James  
Straight  
**EWJS**  
Handle #504

Warwick-James  
Left  
**EWJL**  
Handle #504


Warwick-James  
Right  
**EWJR**  
Handle #504



# SURGICAL

## ROOT TIP PICKS ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.




# SURGICAL

## ROOT TIP PICKS ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

West



Howard

Slightly heavier design than Heidbrink picks.




# SURGICAL

## ROOT TIP PICKS ELEVATORS

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Davis Teasers



## PERIOTOME ELEVATORS

Thin sharp blades used to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament. May be used wherever an extraction is indicated, but are particularly important when considering the placement of a dental implant.



1 Posterior  
SPT1  
4347



3 Posterior-Right  
SPTR3  
4348



4 Posterior-Left  
SPTL4  
4349



2 Anterior  
SAT2  
4350



5 Anterior-Angled  
SATA5  
4351



6 Anterior-Straight  
SATS6  
4352




Wise  
SPW  
4353

# SURGICAL

## PERIOSTEALS


These instruments are used for retracting periosteum after removing gingival tissue.



# SURGICAL

## RERIOSTEALS

These instruments are used for retracting periosteum after removing gingival tissue.



# SURGICAL

## RERIOSTEALS

These instruments are used for retracting periosteum after removing gingival tissue.



24G

SPS24G

4 4 1 7

Tin Coating  
24G

SPS24GT

4 4 1 8



Benque

SPSBEN

4 4 1 9



Buser

SPSBUS

4 4 2 0

Tin Coating  
Buser

SPSBUST

4 4 2 1



KN-1

SPSKN1

4 4 2 2



KN-152

SPSKN152

4 4 2 3



Hourigan-2

SPSHOURI

4 4 2 4

# SURGICAL

## RERIOSTEALS

These instruments are used for retracting periosteum after removing gingival tissue.



HF-20

SPSHF20

4425

Tin Coating  
HF-20

SPSHF20T

4426



Molt-5

SPSM5

4427



Selden-23

SPSSEL23

4428



Freer

SPSFR

4429

OHL-10  
Periosteal  
Elevator

SPSOHL10

4430

Woodson-1  
Modify

SPSWSONM

4431

Allen Periosteal  
Elevator

SPPAEL

4432

# SURGICAL

## SURGICAL CURETTES (SPOON)

These instruments are used for curettage and cyst removal and debridement of tooth socket.

Lucas-Curettes



# SURGICAL

## SURGICAL CURETTES (SPOON)

These instruments are used for curettage and cyst removal and debridement of tooth socket.

### Miller Currettes



### Molt-Currettes



**BONE CHISELS  
SINGLE BEVEL**



2 Gardner  
4mm

**SBCSE2G4**  
4541

5 Gardner  
4.3mm

**SBCSE5G4.3**  
4542

6 Gardner  
4mm

**SBCSE6G4**  
4543

9 Gardner  
4mm

**SBCSE9G4**  
4544

1 Chandler  
5mm

**SBCSE1CH5**  
4545

2 Chandler  
5mm

**SBCSE2CH5**  
4546

4 Chandler  
5mm

**SBCSE4CH5**  
4547

**DOUBLE BEVEL**



1 Gardner  
3mm

**SBCSE1G3**  
4540

7 Chandler  
4.7mm

**SBCDB7CH4.7**  
4548

8 Chandler  
6.1mm

**SBCDB8CH6.1**  
4549

9 Chandler  
7.1mm

**SBCDB9CH7.1**  
4550

1W Wakefield  
2.7mm

**SBCDB1W2.7**  
4553

2W Wakefield  
3mm

**SBCDB2W3**  
4554

4.2mm

**SBCDB18M4.2**  
4555

# SURGICAL

## BONE FILES

Bone Files: For final smoothing of bone by pushing out.

### Miller Collburn-Cross Cut

The small paddle shaped end is half the size of the larger paddle shaped end.



1X  
Miller-Collburn  
Cross-Cut

**SBF1X**  
4 6 0 0

Curved working end combined with a straight working end.



2X  
Miller-Collburn  
Cross-Cut

**SBF2X**  
4 6 0 1

Curved working end combined with a small paddle shaped working end



3X  
Miller-Collburn  
Cross-Cut

**SBF3X**  
4 6 0 2

### Miller-Straight Cut

Same design as the 1X Cross-Cut but with Straight-Cut.



21  
Miller

**SBF21M**  
4 6 0 3

Same design as the 2X Cross-Cut but with Straight-Cut.



45  
Miller

**SBF45M**  
4 6 0 4

# SURGICAL

## BONE FILES

Bone Files: For final smoothing of bone by pushing out.

### Miller-Straight Cut



52  
Miller

**SBF52M**

4605



64  
Miller

**SBF64M**

4606

### Wahl-Straight Cut



1  
Wahl

**SBF1W**

4607



2  
Wahl

**SBF2W**

4608

### Howard-Straight Cut



12  
Howard

**SBF12H**

4609



67  
Howard

**SBF67H**

4610

# SURGICAL


## SCISSORS

Used for trimming tissue or cutting sutures.  
All Scissors shown 50% size.

Wagner



Kelly



# SURGICAL

## SCISSORS

Used for trimming tissue or cutting sutures.  
All Scissors shown 50% size.



# SURGICAL

## SCISSORS

Used for trimming tissue or cutting sutures.  
All Scissors shown 50% size.



31 Castroviejo Curved  
10cm/4"

**SS31CASC**  
4740



17 Iris Straight/Delicate  
11.5cm/4½"

**SS17ISD**  
4741



18 Iris Curved/Delicate  
11.5cm/4½"

**SS18ICD**  
4742



Metzenbaum Baby  
12.5cm/5"

**SSMB**  
4743

# SURGICAL

## SCISSORS

Used for trimming tissue or cutting sutures.  
All Scissors shown 50% size.


### General Surgical Scissors



21  
Straight/  
Pointed  
14.5cm/5<sup>3</sup>/<sub>4</sub>"

**SS21SP**

4 7 4 4



22  
Straight/  
Blunt-Sharp  
14.5cm/5<sup>3</sup>/<sub>4</sub>"

**SS22SBS**

4 7 4 5



23  
Curved/  
Pointed  
14.5cm/5<sup>3</sup>/<sub>4</sub>"

**SS24CP**

4 7 4 6



Lister  
Bandage  
14.5cm/5<sup>3</sup>/<sub>4</sub>"

**SSLB**

4 7 4 7

### Suture Scissors



13  
Suture  
15cm/6"

**SS13S**

4 7 4 8



13S  
Suture  
15cm/4<sup>3</sup>/<sub>4</sub>"

**SS13SS**

4 7 4 9



Crown &  
Gold  
Straight  
15cm/4<sup>3</sup>/<sub>4</sub>"

**SSC&GS**

4 7 5 0



Crown &  
Gold  
Curved  
15cm/4<sup>3</sup>/<sub>4</sub>"

**SSC&GC**

4 7 5 1



Crown &  
Gold Curved  
Pedodontic  
15cm/3.5"

**SSC&GP**

4 7 5 2

### Crown & Gold Remover

# SURGICAL

## TUNGSTON CARBIDE TIPPED SCISSORS

**Super-Cut Scissors:** The razor-sharp cutting edge cuts effortlessly and softly also through close material. The serrated cutting edge on the second blade prevents a slipping from tissue and sutures. Super-Cut is not practical for wires, cartilages and other hard material. As identification, one handle ring is gold plated.

All Scissors shown 50% size.

Super cut scissors



10  
Double-Curved  
13.5cm/5 1/4"

SS10DC

4800



14  
Double-Curved  
11.5cm/5 1/4"

SS14DC

4801



16  
Goldman-Fox  
12.5cm/5"

SS16GFTC

4802



17  
Iris  
Straight/Delicate  
11.5cm/4 1/2"

SS17ISDTC

4803



18  
Iris  
Curved/Delicate  
11.5cm/4 1/2"

SS18ICDTC

4804



Metzenbaum  
Straight/Blunt  
14.5cm/5 1/2"

SSMSBTC

4805



Metzenbaum  
Curved/Blunt  
14.5cm/5 1/2"

SSMCBTC

4806



Metzenbaum  
Straight/  
Pointed  
14.5cm/5 1/2"

SSMSPTC

4807



Metzenbaum  
Curved/  
Pointed  
14.5cm/5 1/2"

SSMCPTC

4808



Metzenbaum  
Delicate-  
Curved/  
Blunt  
14.5cm/5 1/2"

SSMDCBTC

4809

# SURGICAL

## SCISSORS

All Scissors shown 50% size.



# SURGICAL

## TISSUE PLIERS AND FORCEPS

These pliers are used to grasp and stabilize soft tissue flaps during suturing.

All Tissue pliers and forceps shown 50% size.

### Semkin-Taylor

<p>31 Semkin-Taylor Straight 12.5cm/5"</p> <p><b>TP31STS</b> 4870</p>	<p>32 Semkin-Taylor Curved 12.5cm/5"</p> <p><b>TP32STC</b> 4871</p>	<p>33 Semkin-Taylor 1x2 Straight 12.5cm/5"</p> <p><b>TP33STS</b> 4872</p>	<p>34 Semkin-Taylor 1x2 Curved 12.5cm/5"</p> <p><b>TP34STC</b> 4873</p>
---	---	---	---

<p>16 3x4 15cm/6"</p> <p><b>TP16</b> 4874</p>	<p>50-L Fox Locking 3x3 15cm/6"</p> <p><b>TP50L</b> 4875</p>	<p>Kramer-Nevins Tissue Pliers 1x2 Curved</p> <p><b>TPKN</b> 4876</p>
---	--	---

### Suture Pliers

<p>20 Corn Suture Pliers 15cm/6"</p> <p><b>TP20C</b> 4877</p>
---

<p>41 Adson Plain 12cm/4<sup>3</sup>/<sub>4</sub>"</p> <p><b>TPAP41</b> 4878</p>	<p>42 Adson 1x2 12cm/4<sup>3</sup>/<sub>4</sub>"</p> <p><b>TPA42</b> 4879</p>	<p>43 Adson-Brown 7x7 12cm/4<sup>3</sup>/<sub>4</sub>"</p> <p><b>TPAB43</b> 4880</p>	<p>Russian 15cm/6"</p> <p><b>TPRN15</b> 4881</p>	<p>Russian 20cm/8"</p> <p><b>TPRN20</b> 4882</p>	<p>Gerald Plain 18cm/7" Straight/Serrated</p> <p><b>TPGP18</b> 4883</p>	<p>Gerald 1x2 15cm/7" Straight</p> <p><b>TPG15</b> 4884</p>
--	---	--	--	--	---	---

### Allison

<p>Allison Baby 4x5 12.5cm/5"</p> <p><b>FALIB12.5</b> 4885</p>	<p>Allison Straight 5x6 15cm/6"</p> <p><b>FALIS15</b> 4886</p>
--	--

# SURGICAL

## TISSUE PLIERS / FORCEPS

All Tissue Pliers and Forceps shown 50% size.

Adson



41

Adson  
12cm/4<sup>3</sup>/<sub>4</sub>"

**TPAD12TC**  
4 8 8 7



42

Adson 1x2  
12cm/4<sup>3</sup>/<sub>4</sub>"

**TPAD1X2TC**  
4 8 8 8



Adson  
Brown  
12cm/4<sup>3</sup>/<sub>4</sub>"

**TPADB12TC**  
4 8 8 9

General



General  
15cm/6"

**TPG15TC**  
4 8 9 0



General  
1x2  
15cm/6"

**TPG1X2TC**  
4 8 9 1

Gerald



Gerald  
18cm/7"

**TPGRD18TC**  
4 8 9 2



Gerald 1x2  
18cm/7"

**TPGRD1X2TC**  
4 8 9 3



Gushing  
18cm/7"

**TPGUSHTC**  
4 8 9 4

Cushing



Cushing  
Angled  
18cm/7"

**SLB**  
4 8 9 5



Cushing  
Bayone  
19cm/7<sup>1</sup>/<sub>2</sub>"

**SLB**  
4 8 9 6



Potts Smith  
18cm/7"  
(Shown)

**SLB**  
4 8 9 7

Potts Smith  
14.5cm/5-1/4"

**SLB**  
4 8 9 8



De Bakay  
15cm/6"  
(Shown)

**SLB**  
4 8 9 9

De Bakay  
19cm/7<sup>1</sup>/<sub>2</sub>"


**SLB**  
4 9 0 0

# SURGICAL

## NEEDLE HOLDERS

Martensite stainless steel instruments are used for grasp and guide the needle while suturing.

All Needle Holders shown 50% size.



Hegar-Baumgartner (Gardner 5)  
14cm/5 1/2"  
**NHHB**  
4950

Webster Smooth  
13cm/5"  
**NHWS13**  
4951

Mayo-Hegar (Gardner 6) Grooved  
16cm/6 1/4"  
**NHMH**  
4952

Crile-Wood Grooved  
15cm/6"  
**NHCW**  
4953

Olsen-Hegar Needle Holder/Scissors  
17cm/6 3/4"  
**NHOHHS**  
4954

Boynton (Hu-Friedy)  
12cm/4 1/4"  
**NHBOY**  
4955

Castroviejo  
14cm/5 1/2"  
**NHCV**  
4956

Collier  
12.5cm/5"  
**NHC**  
4957

Derf Grooved  
12.5cm/5"  
**NHDG**  
4958

Mathieu (Machatt)  
14cm/5 1/2"  
**NHMM**  
4959

# SURGICAL

## TC-NEEDLE HOLDERS

Used for grasp and guide the needle when suturing.  
All Needle Holders shown 50% size.



Derf  
11.5cm/4 1/2"

**NHDTTC**  
4960



Hegar-  
Baumgartner  
11.5cm/4 1/2"

**NHHBTC**  
4961

Halsey



Hasley  
Smooth  
13cm/5"

**NHHSTC**  
4962



Halsey  
13cm/5"

**NHHTC**  
4963

## FLAT HANDLE

### Castroviejo TC Needle Holders



Castroviejo  
Straight  
Flat Handle  
13cm/5.15"

**NHCVSFH13TC**

4970



Castroviejo  
Curved  
Flat Handle  
13cm/5.15"

**NHCVFHC13TC**

4971



Castroviejo  
Straight  
Flat Handle  
14cm/5.5"

**NHCVSFH14TC**

4972



Castroviejo  
Curved  
Flat Handle  
14cm/5.5"

**NHCVFHC14TC**

4973



Castroviejo  
Straight  
Flat Handle NH  
15cm/5.9"

**NHCVFHS15TC**

4974



Castroviejo  
Curved  
Flat Handle NH  
15cm/5.9"

**NHCVFHC15TC**

4975



Castroviejo  
Straight  
Flat Handle  
16cm/6.3"

**NHCVFHS16TC**

4976



Castroviejo  
Curved  
Flat Handle  
16cm/6.3"

**NHCVFHC16TC**

4977



Castroviejo  
Straight  
Flat Handle  
18cm/7"

**NHCVFHS18TC**

4978



Castroviejo  
Curved  
Flat Handle  
18cm/7"

**NHCVFHC18TC**

4979

# SURGICAL

## ROUND HANDLE Castroviejo TC Needle Holders



Castroviejo  
Straight  
Round Handle  
13cm/5"  
**NHCVSRH13TC**  
4982



Castroviejo  
Curved  
Round handle  
13cm/5"  
**NHCVCRH13TC**  
4983



Castroviejo  
Straight  
Round Handle  
14cm/5.5"  
**NHCVSRH14TC**  
4984



Castroviejo  
Curved  
Round Handle  
14cm/5.5"  
**NHCVCRH14TC**  
4985

4987



Castroviejo  
NH, ST  
Round Handle  
15cm/5.9"  
**NHCVSRH15TC**  
4986



Castroviejo  
NH, CV  
Round Handle  
15cm/5.9"  
**NHCVCRH15TC**  
4987



Castroviejo  
Straight  
Round Handle  
16cm/1.6"  
**NHCVSRH16TC**  
4988



Castroviejo  
Curved  
Round Handle  
16cm/1.6"  
**NHCVCRH16TC**  
4989



Castroviejo  
Straight  
Round Handle  
18cm/7"  
**NHCVSRH18TC**  
4990




Castroviejo  
Curved  
Round Handle  
18cm/7"  
**NHCVCRH18TC**  
4991

# SURGICAL

## NEEDLE HOLDERS

All Needle Holders shown 50% size.



Arruga  
Serrated  
16cm/6 1/2"  
**NHARRUS16TC**  
4995



Webster  
Smooth  
11.5cm/4 1/2"  
**NHWEBS11.5TC**  
4996



Crile-Wood  
15cm/6"  
**NHCW15TC**  
4997



Mayo-Hegar  
16cm/6 1/4"  
**NHMH16TC**  
5001

Olsen-Hegar (Needle Holder)



Olsen-Hegar  
Needle Holder  
Scissors  
17cm/6 3/4"  
**NHSH17TC**  
4998



Olsen-Hegar  
Needle Holder  
Scissors  
14cm/5 3/4"  
**NHSH14TC**  
4999



DeBakey  
18cm/7"  
Delicate  
**NHDEB18TC**  
5000



Mayo-Hegar  
18cm/7"  
**NHMH18TC**  
5002

# SURGICAL

## NEEDLE HOLDERS

Used for grasp and guide the needle while suturing.  
All Needle Holders shown 50% size.



Matheu-Kocher  
14cm/5"

**NHMK14TC**

5005



Matheu  
14cm/5 1/2"

**NHM14TC**

5006



HMT  
Thin Jaw  
14cm/5 1/2"

**NHHMT14TC**

5006



Fine Swedish  
Pattern  
15cm/6"

**NHFSP15TC**

5007



Micro-Vascular  
Pattern  
15cm/6"

**NHMVP15TC**

5008



Micro-Vascular  
Pattern  
18cm/7"

**NHMVP18TC**

5009



Mini-Ryder  
Pattern  
15cm/6"

**NHMRP15TC**

5010



Mini-Ryder  
Pattern  
18cm/7"

**NHMRP18TC**

5011

# SURGICAL

## HEMOSTATIC FORCEPS

A multipurpose instrument used to clamp off blood vessels, remove small root tips and grasp loose objects. All Hemostats

Hartman



Hartman  
Straight  
10cm/4"

**HHS10**  
5 0 3 0



Hartman  
Curved  
10cm/4"


**HHC10**  
5 0 3 1



5  
1x2 Teeth  
Straight  
12.5cm/5"

**H5S1X2**  
5 0 3 4

Adson




Adson  
Straight  
19cm/7"

**HAS19**  
5 0 3 2



Adson  
Curved  
19cm/7-1/2"

**HAC19**  
5 0 3 3



Rochester-pean  
Curved  
16cm/6"

**HRPC16**  
5 0 3 5

# SURGICAL

## HEMOSTATIC FORCEPS

A multipurpose instrument used to clamp off blood vessels, remove small root tips and grasp loose objects. All Hemostats

Halsted-Mosquito



Kelly



# SURGICAL

## CHEEK RETRACTORS/TISSUE RETRACTORS



Middeldorf  
Retractor, 17x14mm  
SMID17x14  
5050



Middeldorf  
Retractor, 22x20mm  
SMID22x20  
5051



9 Ramus  
Retractor, 16x80mm  
SR9  
5052



8 Surgical  
Retractors, 11x70mm  
SR8  
5053



1 Surgical  
Retractor, 7x25mm  
SR17x25  
5059



1 Surgical  
Retractor, 10x42mm  
SR110x42  
5060



1 Surgical  
Retractors, 12x55mm  
SR1-12x55  
5061



1 Austin  
S1A  
5066



2 Austin  
S2A  
5067



1 Henahan  
SR1HENN  
5068



2 Henahan  
SR2HENN  
5069




3 Henahan  
SR3HENN  
5070

# SURGICAL

## CHEEK RETRACTORS/TISSUE RETRACTORS

**Retractors:** For holding back the lips and cheeks, and tongue.


**Mouth Gag:** To prop open the mouth and prevent it from closing.



# SURGICAL

## CHEEK RETRACTORS/TISSUE RETRACTORS

**Retractors:** For holding back the lips and cheeks, and tongue.



**Skin Hooks:** For hooking skin during surgery.  
**Retractors:** For holding back the lips and cheeks.



# SURGICAL

## SCALPEL HANDLES

**Scalpel Handle:** For holding a mess blade.



# SURGICAL

## ASPIRATORS

These tools are used for precision section at surgical sites.

Coupland Handle  
25cm/9.9"

**SASPCH**  
5249



Actual tip size (Internal diameter)

5250	<b>SASPCT1</b>	<input type="radio"/>	1 Coupland Tip	1.5mm
5251	<b>SASPCT2</b>	<input type="radio"/>	2 Coupland Tip	2.0mm
5252	<b>SASPCT3</b>	<input type="radio"/>	3 Coupland Tip	2.5mm
5253	<b>SASPCT4</b>	<input type="radio"/>	4 Coupland Tip	3.0mm

Actual tip size (Internal diameter)

5254	<b>SASPCT1SH</b>	<input type="radio"/>	1 Coupland Tip/Side Holes	1.5mm
5255	<b>SASPCT2SH</b>	<input type="radio"/>	2 Coupland Tip/Side Holes	2.0mm
5256	<b>SASPCT3SH</b>	<input type="radio"/>	3 Coupland Tip/Side Holes	2.5mm
5257	<b>SASPCT4SH</b>	<input type="radio"/>	4 Coupland Tip/Side Holes	3.0mm

Shown 50% size.

BYRD



Byrd Self-Cleaning  
20cm/8"

Actual tip size (Internal diameter)

5259	<b>SASPB1</b>	<input type="radio"/>	B1 Byrd	1.5mm
5260	<b>SASPB2</b>	<input type="radio"/>	B2 Byrd	2.0mm
5261	<b>SASPB3</b>	<input type="radio"/>	B3 Byrd	2.5mm

Cogswell 3 Suction Tube  
23cm/9"

**SASPCG3**  
5262




Shown 50% size.

Actual tip size (Internal diameter)

5263	<b>SASPF6</b>	<input type="radio"/>	6 Frazier	1.5mm
5264	<b>SASPF8</b>	<input type="radio"/>	8 Frazier	2.0mm
5265	<b>SASPF10</b>	<input type="radio"/>	10 Frazier	2.5mm
5266	<b>SASPF12</b>	<input type="radio"/>	12 Frazier	3.0mm

Frazier Suction Tube



Used in oral surgery



**SASP1**  
5270

**SASP2**  
5271

**SASP3**  
5272

**SASP4**  
5273

Yankauer  
26.5cm/10.5"

**SASPYK**  
5275

Shown 50% size.



Andrew-Pychon  
23cm/9"

**SASPAP**  
5276

Shown 50% size.



Locking Pliers  
16.5cm/6 1/2"

**LP**

5280



# SURGICAL

## SYRINGES

**Syringes:** An instrument for injecting aesthetics for a cartridge.  
 3-Way Syringes - Supplies air, water, or a mist through separate interval channels created by combining pressurized air with water flow.



Intraligamentary  
Syringe  
1.8ml

**SY11.8ml**  
5320



Finger Loop Style

Syringe  
1.8ml  
**SY1.8ml**  
5322



Finger Rest Style

Syringe  
1.8ml  
**SY111.8ml**  
5323



Syringe  
1.8ml  
**SYS1.8**  
5324

Cook-Waite Style



Paroject

**SYPARS**  
5328



Pen Style with Blue Handle

**SYPARP**  
5329

# SURGICAL

## MISCELLANEOUS

Miscellaneous tools for laboratory.



Stainless Steel 6" Rule

**SS6R**  
5386



Tofflemire  
Junior

**TJ**  
5383




Tofflemire  
Senior

**TS**  
5384



Tofflemire  
Universal

**TU**  
5385




Shown at 50% size.

Boley  
Gauge  
**BG**  
5387

# SURGICAL

## ATRAUMATIC EXTRACTING FORCEPS SET

Our Unique tip design makes it mold to the individual shape of your patient's tooth. The smooth handle with openings, allows for more thorough and easier cleaning. The newly styled handle allows a more uniform transmission while applying pressure, pulling and for rotary movements. More comfortable design can be used by right or left-handed practitioners! Our handles are lighter and are easier to use.



F1

UPPER  
Universal

EF1

5430



F2

LOWER  
Universal

EF2

5431



F3

Upper  
Anterior

EF3

5432



F5

Lower  
Anterior

EF4N

5433



F6

Upper  
3rd Molar

EF6

5434



F7C

Lower Molar, modified  
Cow-Horn Bead

EF7C

5435



46XL

Lower Root  
Extra Long

E46XL

5436



51XL

Upper Root

E51XL

5437



222

Lower  
Wisdoms

E222

5438



301

Lower  
Roots

E301

5439



300

Upper Roots

E30

5440

# SURGICAL

## EXTRACTION FORCEPS AMERICAN PATTERN

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.

The image displays various dental extraction forceps and root splinter tweezers. Each tool is accompanied by a diagram of the tooth it is designed to extract, with 'R-L' indicating right and left sides. The tools are arranged in a grid-like fashion around a central, larger image of a root splinter tweezer.

- 5 Lower Third Molars** (E5, 5441)
- 15 Lower Molars** (E15, 5442)
- 18L LOWER MOLARS LEFT SIDE** (E18L, 5447)
- 18R Upper Molars Right Side** (E18R, 5448)
- 16S Lower Molars** (E16S, 5443)
- 17 Lower Molars** (E17, 5444)
- 6 Splitting Forceps** (E6, 5449)
- 1 Upper Centrals Laterals & Cuspids** (E1, 5450)
- 16 Lower Molars Cow-Horn Beak** (E16, 5445)
- 1AF Upper Anteriors** (E1AE, 5446)

# SURGICAL

## EXTRACTION FORCEPS AMERICAN PATTERN

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



99C  
Upper  
Premolars  
& Incisors

E99C

5451



10S  
Upper Molars  
Straight-Handle

E10S

5454



53L  
Upper Molars  
Left Side

E53L

5455



53R  
Upper Molars  
Right Side

E53R

5456



23  
Lower Molars  
Cow-Horn Beak

E23

5457



10H  
UPPER MOLARS  
HOOK HANDLE

E10H

5458



69  
Fragments  
& Small Roots  
Upper & Lower

E69

5459



65  
Upper Incisors  
& Roots

E65


5460



Universal  
for Upper  
Teeth & Roots

EU

5461



32  
Upper Molars  
& Premolars  
Strong Beak

E32

5462

# SURGICAL

## EXTRACTION FORCEPS AMERICAN PATTERN

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



151A  
Lower  
Premolars  
Incisors & Roots

E151A

5463



151XAS Children  
Lower Premolars  
Incisors &  
Roots Canines

E151XAS

5464



151AS  
Split Lower  
Premolars  
Incisors &  
Roots Canines

E151AS

5465



151AF Apical  
Beak Lower  
Premolars  
Incisors &  
Roots Canines

E151AE

5466



101  
Universal  
for Premolars

E101

5467



150  
Upper Premolars  
Incisors &  
Roots Canines

E150

5468



150AF  
Apical Beak  
Upper Premolars  
Incisors &  
Roots Canines

E150AF

5469



150AS  
Split Beak  
Upper Premolars  
Incisors &  
Roots Canines

E150AS

5470



150XAS  
Children Split  
Upper Premolars  
Incisors &  
Roots Canines

E150XAS

5471



150A  
Upper PreMolars  
Incisors &  
Roots Canines

E150A

5472



150S  
Upper Premolars  
Incisors &  
Roots for  
Children

E150S

5473



151  
Upper Premolars  
Incisors &  
Roots

E151

5474

# SURGICAL

## EXTRACTION FORCEPS AMERICAN PATTERN

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



151S  
Lower Premolars  
Incisors &  
Roots for Children

**E151S**  
5 4 7 5



210H  
Upper Wisdoms  
with Hook Handle

**E210H**  
5 4 7 6




210S Upper  
Wisdoms  
with Straight  
Handle

**E210S**  
5 4 7 7



222  
Lower Wisdoms

**E222**  
5 4 7 8



222AF  
Apical Beak  
Lower Wisdoms

**E222AF**  
5 4 7 9



286  
Upper Roots  
Premolars  
& Incisors

**E286**  
5 4 8 0

# SURGICAL

## Extracting Forceps European Style

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



30  
Upper Roots  
& Incisors

E30

5481



22  
Lower Molars  
Hawks Bill

E22

5482



21  
Lower Molars

E21

5483



20  
Lower Wisdoms

E20

5484



8  
Lower  
Premolars

E8


5485



7  
Upper  
Premolars

E7

5486



3  
Upper Incisors  
& Canines

E3

5487



1  
Upper Central  
& Canines

E1


5488



19  
Upper Wisdoms

E19

5489



18A  
Upper Molars

E18A

5490



18  
Upper Molars  
Left

E18

5491



17  
Upper Molars  
Right

E17

5492


# SURGICAL

## Extracting Forceps European Style

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



67 AF  
Apical Beak  
Upper Wisdoms  
**E67AF**  
5493



68  
Lower Roots  
**E68**  
5494



73  
Lower Molars  
(Hawks Bill)  
**E73**  
5495



74  
Lower Incisors  
& Roots  
**E74**  
5496



31  
Lower Roots  
**E32**  
5497



32  
Lower Molars  
**E31**  
5498



33  
Lower Roots  
**E33**  
5499



33M  
Lower Roots  
With Long Beak  
**E33M**  
5500



67  
Upper Wisdoms  
**E67**  
5501



51A  
Upper Roots  
Narrow Beak  
**E51A**  
5502



45  
Lower Roots  
**E45**  
5503



44  
Lower Roots  
**E44**  
5504

# SURGICAL

## Extracting Forceps European Style

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



95  
Upper Molars  
Left

E95

5505



121  
Upper Wisdoms

E121

5506



167 Universal  
Patter for  
Lower Teeth  
Molars & Premolars

E167

5507



168  
Upper  
Canines

E168

5508



74N  
Lower Incisors  
& Roots  
Canines

E74N

5509



74NF  
Upper Incisors  
& Roots

E74NF

5510



79  
Lower Wisdoms

E79

5511



79A  
Lower Wisdoms

E79A

5512



86C  
Lower Molars

E86C

5513



89  
Upper Molars  
Right

E89

5514



90  
Upper Molars  
Left

E90

5515



94  
Upper Molars  
Right

E94

5516

# SURGICAL

## Extracting Forceps European Style

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.


**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



F-1  
Upper  
Universal  
EF-1  
5517



F-2  
Lower  
Universal  
EE-2  
5518



F-3  
Upper Anterior  
EF-3  
5519



F-4N  
Lower Anterior  
EF-4N  
5520



F-4  
Lower Anterior  
EF-4  
5521



F-5  
Lower Anterior  
Parallel Beaks  
EF-5  
5522



F-6  
Upper  
3rd Molars  
EF-6  
5523



F-7  
Lower  
3rd Molars  
EF-7  
5524



F-7C  
Lower Molar  
Modify  
Cow-Horn Beak  
EF-7C  
5525



F-8  
Upper Universal  
EF-8  
5526



F-9  
Lower Universal  
EF-9  
5527



F-10  
Lower Molar  
EF-10  
5528

# SURGICAL

## Extracting Forceps European Style

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



F11  
Upper Molar  
Left  
EF-11  
5529



F12  
Upper Molar  
Right  
EF-12  
5530



F13  
Upper Molar  
Right  
EF-13  
5531



F14  
Upper Molar  
Left  
EF-14  
5532



169  
Lower Cnines  
E169  
5533



233  
Lower Roots  
Parallel Beaks  
E233  
5534



235  
Upper Roots  
Parallel Beaks  
E235  
5535



MD-1  
Universal  
Upper Wisdom  
Premolars &  
Roots  
EMD-1  
5536



MD-4  
Lower Molars  
EMD-4  
5537



22L  
Lower Wisdoms  
Left  
E22L  
5538



22R  
Lower Wisdoms  
Right  
E22R  
5539

## Extracting Forceps European Style

**Forceps:** For extraction teeth form alveolar bone. American Pattern/English Pattern/Unique Design/Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**SPONGE & DRESSING/HEMOSTATES:** For grasping off blood vessels and removing small pieces of root tips.



F15  
Upper  
Premolars  
Universal

EF-15

5540



F16  
Upper  
Root

EF-16

5541



F17  
Lower Root

EF-17

5542



F18  
Upper Roots  
Fragment Narrow

EF-18

5543



F 18TN  
Upper Roots  
Fragment Narrow

EF-18TN

5544



F 19  
Lower Roots  
Fragment Narrow

EF-19

5545



F 19NW  
Lower Roots  
Fragment Narrow

EF-19NW

5546



46XL  
Lower Root  
Extra Long Beaks

EF46XL

5547



F 20  
Upper  
Premolar  
1ST & 2RD Molars

EF-20

5548



49  
Upper Roots

EF49

5549



51XL  
Upper Roots

EF51XL

5550

Atraumatic Forceps Children's / Pedodontic Pattern



F1  
Upper Universal  
Pedo

5551



F2  
Lower Universal  
Pedo

5552



F3  
Upper Anterior  
Pedo

5553



F4  
Lower Anterior  
Pedo

5554



F5  
Lower Anterior  
Pedo

5555



F6  
Upper 3rd  
Pedo

5556



F7  
Lower 3rd  
Pedo

5557



F16  
Lower Root  
Pedo

5558



F17  
Lower Root  
Pedo

5559



B23  
Lower Molar  
Pedo

5560



150-SK  
Upper Universal  
Pedo

E150-SK

5561



151-SK  
Lower Universal  
Pedo

E151-SK

5562

## RONGEURS

## FRIEDMAN



Friedman  
14cm/5-1/2"

RFM

6000

## BLUMENTHAL



Blumenthal  
Standard 30°  
14.5cm/5-3/4"

RBIS

6003

## CLEVELAND



1A Mead  
16.5cm/6-1/2"

RIA

6006



3H Dean  
17cm/6-3/4"

R3HD

6008



4 Cleveland  
14cm/5-1/2"

R4C

6010



Mini Friedman  
14cm/5-1/2"

RMFM

6001



Blumenthal 45°  
14.5cm/5-3/4"

RBL45

6004



4A Cleveland  
17cm (6-3/4")

R4AC

6007



5 Cleveland  
Side Cutting  
17cm (6-3/4")

R5C

6012



5S Cleveland  
Side Cutting  
15cm (6")

R5SC


6013



Mini Friedman 90°  
14.5cm/5-1/2"

RMFM90

6002



Blumenthal 90°  
14.5cm/5-3/4"

RBM90

6005



Kramer-Nevis  
14.5cm (5-3/4")

RKN

6014

## RONGEURS



Kerrison Rongeur 3mm  
for widening the sinus window

RK3mm

6020



Kerrison Rongeur 4mm  
for widening the sinus window

RK4mm

6022

### DOUBLE ACTION HINGED RONGEURS



18M- (Double Action)  
Side-Cutting  
17mm (6-7/8")

R18MDA17

6025



Beyer- (Double Action)  
18cm (7")

RBEYDA18

6030



Boehler (Double Action)  
15cm (6")

RBOEDA15

6035



Boehler (Double Action)  
16cm (6-1/4")

RBOEDA16

6040

## Black Line Instruments




## Black Line Instruments

Black Line Cheek Retractors			BL29
Black Line Comp/Plastic Instruments			BL13-BL19
Black Line Diagnostic Explorers			BL-1
Black Line Endo Excavators			BL-2
Black Line Endo Explorers			BL-2
Black Line Endo Pluggers			BL4
Black Line Endo Spreaders			BL-3
Black Line Implantology Pluggers & Condensers			BL8
Black Line Luxating Elevators			BL30
Black Line Luxating Elevators (Thin Version)			BL31
Black Line PDL Periosteals			BL27-BL28
Black Line Periodontal Bone chisels			BL10-BL11
Black Line Periodontal Knives			BL12
Black Line Periosteals			BL23-BL25
Black Line Periotome Elevators			BL34
Black Line Root Tip Pickes			BL32-BL33
Black Line Rubber Dam Clamps			BL5-BL7
Black Line Sinus Lift Instrumens			BL9
Black Line Surgical Curettes			BL26

# BLACK LINE

## DIAGNOSTIC



2  
DEXD2  
**DEXD2**  
BL1000

3CH  
DEXD3CH  
**DEXD3CH**  
BL1001

5  
DEXD5  
**DEXD5**  
BL1002

TU17/23  
Wilkins/Tufts  
DEXTU17/23  
**DEXDTU17/23**  
BL1003


11/12  
Old Dominion University  
DEXD11/12ODU  
**DEXD11/12ODU**  
BL1004

3-6-9-12  
23/CP-12  
DEP23/CP12  
**DEXD23/CP12**  
BL1005

# BLACK LINE

## Explorers and Excavators

### Endo Explorers



DG16

**EDG16**

BL1020

DG16/17

**EDG16/17**


BL1021

DG16/23

**EDG16/23**

BL1022

### Excavators



1.0mm

31L

**E31L**

BL1023

1.3mm

31W

**E31W**

BL1024

1.5mm

32L

**E32L**

BL1025

2.0mm

33L

**E33L**

BL1026

## SPREADERS

0.20



D11TS  
Thinnest 21mm

**ERCSD11TS**  
BL1030

0.20



MA57  
Heat Carrier 27mm

**ERC SMA57**  
BL1031

0.55mm



00  
25mm

**ERC S00**  
BL1033

0.25



GP1-20MM

**ERC SGP1**  
BL1034

0.23



D11T  
Thinner 22mm

**ERCSD11T**  
BL1032

0.40




D11  
Standard 23mm

**ERCSD11**  
BL1035

# BLACK LINE

## PLUGGERS



# BLACK LINE

Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some cases cotton wool rolls. these are different and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

## Molars - General

	8A		W8A		8AD
<b>RDC8A</b>		<b>RDCW8A</b>		<b>RDC8AD</b>	
BL1830		BL1831		BL1832	
	10		11		14
<b>RDC10</b>		<b>RDC11</b>		<b>RDC14</b>	
BL1833		BL1834		BL1835	
	14A		W14A		18
<b>RDC14A</b>		<b>RDCW14A</b>		<b>RDC18</b>	
BL1836		BL1837		BL1838	
	24		25		26N
<b>RDC24</b>		<b>RDC25</b>		<b>RDC26N</b>	
BL1839		BL1840		BL1841	

## Upper Molars

	56		56S
<b>RDC56</b>		<b>RDC56S</b>	
BL1842		BL1843	
	4		W4
<b>RDC4</b>		<b>RDCW4</b>	
BL1844		BL1845	
	5		8
<b>RDC5</b>		<b>RDC8</b>	
BL1846		BL1847	
	27		
<b>RDC27</b>			
BL1848			

## Rubber Dam Clamps


Martensite stainless steel, Clamps shown 75% actual.

Upper Molars	Lower Molars	Pre-molars
 <p>30</p> <p><b>ERDC30</b> BL1850</p>	 <p>3</p> <p><b>ERDC3</b> BL1854</p>	 <p>0</p> <p><b>ERDC0</b> BL1863</p>
 <p>31</p> <p><b>ERDC31</b> BL1851</p>	 <p>W3</p> <p><b>ERDCW3</b> BL1855</p>	 <p>7</p> <p><b>ERDC7</b> BL1856</p>
 <p>31A</p> <p><b>ERDC31A</b> BL1860</p>	 <p>7A</p> <p><b>ERDC7A</b> BL1858</p>	 <p>12A</p> <p><b>ERDC12A</b> BL1859</p>
 <p>201</p> <p><b>ERDC201</b> BL1852</p>	 <p>26</p> <p><b>ERDC26</b> BL1861</p>	 <p>1</p> <p><b>ERDC1</b> BL1864</p>
 <p>205</p> <p><b>ERDC205</b> BL1853</p>	 <p>28</p> <p><b>ERDC28</b> BL1862</p>	 <p>1A</p> <p><b>ERDC1A</b> BL1865</p>
		 <p>2</p> <p><b>ERDC2</b> BL1866</p>


# BLACK LINE

## Rubber Dam Clamps

Pre-molars



Anterior



## Pluggers and Condensers



Rompen Plugger, Mesial/Distal

**BLRPM/D**


BL1900



Rompen Plugger, Buccal/Lingual

**BLRPB/L**

BL1901



Jovanovic Plugger

**BLJP**

BL1902



Labanca Plugger

**BLLBP**

BL1903



Gridsdale Bone Condenser

**BLGBP**

BL1904

8-11-13-15-18mm

**SINUS LIFT INSTRUMENTS**



Sinus Curettes 1  
ISC1  
BL1915



Sinus Curettes 2  
ISC2  
BL1916



Kramer-Nevins  
Obtuse  
IKNO  
BL1917



Kramer-Nevins  
Acute  
IKNA  
BL1918



Sinus Lift 300  
ISL300  
BL1919

Sinus Curettes 1  
Sinus Curettes 2  
Kramer-Nevins  
Obtuse  
Kramer-Nevins  
Acute

Periodontal Bone Chisels



1/2 Buser Modified  
Periodontal Chisel  
4 mm/5 mm  
P19AUSC1127

1/2 Buser Modified  
Periodontal Chisel (4mm/5mm)  
**I1/2BMPBC**  
BL1930



2/3 Buser Mo  
Periodontal C  
5 mm/6 mm  
P19AUSC1128

2/3 Buser Modified  
Periodontal Chisel (5mm/6mm)  
**I2/3BMPBC**  
BL1931



1/3 Buser Mod  
Periodontal C  
4 mm/6 mm  
P19AUSC1129

1/3 Buser Modified  
Periodontal Chisel (4mm/6mm)  
**I1/3BMPBC**  
BL1932



P24G/13KL  
Palacci Modified  
**1PAL13KG**  
BL1933



3/4 Solt  
**ISOLT3/4**  
BL1934

Sinus Curettage  
Sinus Curette 2  
Kramer-Meyers  
Osteop  
Removal

Periodontal Bone Chisels



15/16  
Kirkland Knife  
**BLKK15/16**  
BL1940



Allen End-Cutting  
Intrасulcular Knife  
**BLKAECI**  
BL1941



1/2 Allen Orban Knife  
Round  
**BLKAOR1/2**  
BL1942



1/2 Allen Orban Knife  
**BLKAO1/2**  
BL1943



11 Goldman-Fox Knife  
**BL11GFK**  
BL1944



1/2 Orban Knife  
**BL1/2OK**  
BL1945

sinus curette 2  
sinus curette 2  
curette  
parodontal

Periodontal Knives



Tunneling Knife1  
**BLTK1**  
BL1950



Tunneling Knife 2  
**BLTK2**  
BL1951




Giles De Quincey  
Tunneling Knife 1  
**BLGDQTK1**  
BL1952



Gile De Quincey  
Tunneling Knife 2  
**BLGDQTK2**  
BL1953

GOLDSTEIN FLEXIBLE THIN COMPOSITE INSTRUMENTS



1  
**RTNCIGFT1**  
BL3230



6 (Flexible Thin)  
**RTNCIGFT6**  
BL3235




2  
**RTNCIGFT2**  
BL3231




Micro-Mini  
**RTNCIPCS**  
BL3236



3 (Extra Flexible)  
**RTNCIGFT3**  
BL3232




Mini 1  
**RTNCIGTM1**  
BL3237




4 (Extra-Flexible)  
**RTNCIGFT4**  
BL3233



Mini 3 (Extra-Thin)  
**RTNCIGTM3**  
BL3238



5 (Flexible-Thin)  
**RTNCIGF5**  
BL3234




Mini 4 (Extra-Thin)  
**RTNCIGTM4**  
BL3239



AB1  
Boghosian  
**RTNPF1AB1**  
BL3240

AB1  
Boghosian  
RTNPF1AB1




AB2  
Boghosian  
**RTNPF1AB2**  
BL3241

AB2  
Boghosian  
RTNPF1AB2



Interproximal  
Carver  
**RTNCVIPC**  
BL3242

Interproximal  
Carver  
RTNCVIPC



W3  
**RTNPF1W3**  
BL3243

W3  
RTNPF1W3



A6  
**RTNPF1A6**  
BL3244

A6  
RTNPF1A6



4/5 Gregg  
**RTNPF18A**  
BL3245

4/5  
RTNPF18A



Contouring Instrument  
**RTNCF1R/L**  
BL3246

Contouring  
Instrument  
RTNCF1R/L



**RTNPF18A**

8A  
RTNPF18A  
BL3247

8A  
RTNPF18A

BURNISH

### BURNISHERS



21B  
**RTNBB21B**  
BL3250



27/29  
**RTNBB2729**  
BL3251



2 Ladmore  
**RTNBBL2**  
BL3252




3 Ladmore  
**RTNBBL3**  
BL3253



Small/Medium  
Ball Burnisher  
**RTNBBS/M**  
BL3254

### FREEMAN BURNISHERS



Freedman (Duckhead)  
**RTNPCCI**  
BL3255



Freedman (Small)  
Contact Forming  
**RTNFCIS**  
BL3256



Freedman (Large)  
Contact Forming  
**RTNFCIL**  
BL3257

**GOLDFOGEL FREEHAND INSTRUMENTS**



A  
Cosmetic Contouring

**RTNCCIA**  
BL3260



B  
Cosmetic Contouring

**RTNCCIB**  
BL3261



F  
Cosmetic Contouring

**RTNCCIF**  
BL3262



C  
Cosmetic Contouring

**RTNCCIC**  
BL3263



G  
Marginal Ridge & Embrasure Shaping Instrument

**RTNCCIG**  
BL3264




D  
Cosmetic Contouring

**RTNCCID**  
BL3265



H  
Occlusal Anatomy Instrument

**RTNCCIH**  
BL3266



E  
Cosmetic Contouring

**RTNCCIE**  
BL3267



I  
Cosmetic Packing Instrument


**RTNCCII**  
BL3268

**ANTERIOR KIT**

**POSTERIOR KIT**



3  
Extra-Flex  
RTNCIGFT3  
BL3270




OT  
Tanner  
RTNPLGOT  
BL3271




Mini 1  
RTNCIGFTMI1  
BL3272




3  
Hollenback  
RTNPLGH3  
BL3273




Micro-Mini  
RTNCIPCS  
BL3274




5A  
RTNPLG5A  
BL3275



Medium Placing/Condensing  
RTNCIPCM  
BL3276



Small/Medium  
Contact Forming  
RTNCFISM  
BL3277



Large  
Placing/Condensing  
RTNCIPCL  
BL3278



Medium/Large  
Contact Forming  
RTNCFIM/L  
BL3279

**COMPOSITE INSTRUMENTS**

Composite Spatula  
Small, 1-1.5mm  
**CS-1/1.5**  
3300



Composite Spatula  
Large  
**CSL**  
3301



Former and Spatula  
**CF&S**  
3302



Composite Spatula 1.7mm  
Special Bend for Mesial and Distal Areas  
**CS1.7**  
3303



Composite Spatula  
Combo Flattened-Beavertail  
**CSCFB**  
3304



1.6mm  
**CS1.6**  
3305



Plastic Instrument/Composite Fissure Former Black  
**CICFFB**  
3306

## COMPOSITE INSTRUMENTS



Spatula 2.2mm/Burnisher 1.5mm

**S2.2B1.5**  
3310



Composite Pluggers

**CP1**  
3311



Composite Pluggers

**CP2**  
3312



Spatula 2.4mm/Burnisher 1.9mm

**S2.4/B1.9**  
3313



Composite Plugger

**CP3**  
3314

Plastic  
Instrument/Composite Plugger

**PICP**  
3315



Plastic Instrument, specail Plugger

**PISP**  
3316

# BLACK LINE



Probe 0.2mm

**BLPT23**

BL3582



Heidemann Spatula

**BLHST2**

BL3580



Probe 0.2mm

**BLPT11**

BL3581



Plugger (Wiland Shape)

0.4mm

**BLPWST11**

BL3585



Spatula T1

**BLST1**

BL3586



Spatula (Cobra Shape)

0.4mm

**BLSCST2**

BL3587



Spatula T4

**BLST4**

BL3588

PERIOSTEAL



Buser Periosteal Elevator

**BLBPE**

BL3650



Allen Periosteal Elevator

**BLAPE**

BL3651



Allen Periosteal Elevator, Anterior

**BLAPEA**

BL3652



9 Molt Periosteal Elevator

**BL9MPE**

BL3653



3 Prichard Periosteal Elevator

**BL3PPE**

BL3654



7mm

3.5mm

2/4 Molt Curette

**BLMC2/4**

BL3655

## PERIOSTEAL



Swiss Perio  
Modified P24GSP Periosteal

**BLSPMP24GSPP**

BL3670



Swiss Perio  
Modified P24GSIC Periosteal

**BLSPMP24GSICP**

BL3671



Swiss Periosteal  
Modified P24GSP Periosteal 2

**BLSPMP24GSP2**

BL3672



Swiss Periosteal  
Modified Schlee Periosteal

**BLSPMSP**

BL3673



Swiss Periosteal  
Modified PH2 Periosteal

**BLSPMPSH2**

BL3674



Modified Labanca Periosteal

**BLMLP**

BL3675

## PERIOSTEAL



Urban Mini Me Periosteal

**BLUMMP**

BL3680



Buser/Modified PH2MBH (Small)

**BLBMPH2MBHS**

BL3681



Buser/Modified PH2MBH (Large)

**BLBMPH2MBHL**

BL3682



152-Kramer-Nevins Periosteal

**BLKN152P**

BL3683



24G Periosteal

**BL24GP**

BL3684



14 Goldman-Fox Periosteal

**BL14GFP**

BL3685

# BLACK LINE

## SURGICAL CURETTES



86 Lucas  
Surgical Curette (2.8mm)

**BLLSC86**  
BL3690



84 Lucas  
Surgical Curette (2.0mm)

**BLLSC84**  
BL3691



85 Lucas  
Surgical Curette (2.5mm)

**BLLSC85**  
BL3692



88 Lucas  
Surgical Curette (4.7mm)

**BLLSC88**  
BL3693



87 Lucas  
Surgical Curette (3.5mm)

**BLLSC87**  
BL3694

# BLACK LINE

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.



Spade Elevator, Straight

**BLEZLSS**  
BL4180



Spade Elevator, Mesially Angled

**BLEZLSMA**  
BL4181



Spade Elevator, Distally Angled

**BLEZLSDA**  
BL4182



Spade Elevator, Left

**BLEZLSL**  
BL4183



Spade Elevator, Right

**BLEZLSR**  
BL4184



2.5mm, Small, Straight

**BLEZLSS2.5**  
BL4185



2.5mm, Small, Curved

**BLEZLSC2.5**  
BL4186



2.5mm, Small, Mesially Angled

**BLEZLSMA2.5**  
BL4187




2.5mm, Small, Distally Angled

**BLEZLSDA2.5**  
BL4188



# BLACK LINE

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.



# BLACK LINE

## RETRACTOTRS



# BLACK LINE

## LUXATING ELEVATORS



Luxator, Straight  
2mm

**BLLXS2**

BL4220



Luxator, Curved  
2mm

**BLLXC2**

BL4221



Luxator, Straight  
2.5mm

**BLLXS2.5**

BL4222



Luxator, Curved  
2.5mm

**BLLXC2.5**

BL4223



Luxator, Straight  
3mm

**BLLXS3**

BL4224



Luxator, Curved  
3mm

**BLLXC3**

BL4225



Luxator, Straight  
4mm

**BLLXS4**

BL4226



Luxator, Curved  
4mm

**BLLXC4**

BL4227



Luxator, Straight  
5mm

**BLLXS5**

BL4228



Luxator, Curved  
5mm

**BLLXC5**

BL4229

# BLACK LINE

## LUXATING ELEVATORS



Slim Luxating Elevator  
Long, Straight  
2mm  
**SLELS2**  
BL4250



Slim Luxating Elevator  
Long, Curved  
2mm  
**SLELC2**  
BL4251



Slim Luxating Elevator  
Long, Straight  
3mm  
**SLELS3**  
BL4252



Slim Luxating Elevator  
Long, Curved  
3mm  
**SLELC3**  
BL4253



Slim Luxating Elevator  
Long, Straight  
4mm  
**SLELS4**  
BL4255



Slim Luxating Elevator  
Long, Curved  
4mm  
**SLELC4**  
BL4256



Slim Luxating Elevator  
Long, Straight  
5mm  
**SLELS5**  
BL4257



Slim Luxating Elevator  
Long, Curved  
5mm  
**SLELC5**  
BL4258

# BLACK LINE

## HEIDBRINK

## APICAL



1  
Heidbrink  
**BLEHB1**  
BL4316



2  
Heidbrink  
**BLEHB2**  
BL4317



3  
Heidbrink  
**BLEHB3**  
BL4318



13/14  
Heidbrink  
**BLEHB13/14**  
BL4321



9R  
**BLE9R**  
BL4323



9L  
**BLE9L**  
BL4324

ROOT TIP PICKS

DAVIS

HOWARD



10  
Davis  
Double-End  
**BLED10**  
BL4343


11  
Davis  
Double-End  
**BLED11**  
BL4344

2  
Howard  
**BLEHO2**  
BL4332

3  
Howard  
**BLEHO3**  
BL4333

# BLACK LINE

## PERIOTOMES



Periotome,  
Posterior (DE)

**BLPPDE**

BL4347

Periotome  
Posterior, Right (SE)

**BLPPRSE**

BL4348

Periotome  
Posterior, Left (SE)

**BLPPLSE**

BL4349

Periotome  
Anterior (DE)

**BLPADE**

BL4350

Periotome  
Anterior, Angled (SE)

**BLPAASE**

BL4351

Periotome  
Anterior, Straight (SE)

**BLPASDE**

BL4352

Periotome  
Wise

**BLPW**

BL4353

Item No.	Item Descriptions	Item Code	Page
901	M1-Cone Socket Handle	M1CS	D-1
902	M2-Cone Socket Handle	M2CS	D-1
903	M3-Cone Socket Handle	M3CS	D-1
904	M4-Cone Socket Handle	M4CS	D-1
905	M5-Simple Stem Handle	M5SS	D-1
906	M6-Simple Stem Handle	M6SS	D-1
907	M7-Simple Stem Handle	M7SS	D1
908	M8-Simple Stem Handle	M8SS	D1
911	Simple stem Mouth Mirror #4 (7/8 inch-22mm)	CSMWHM2	D-1
912	Simple stem Mouth Mirror #5 (15/16 inch-24mm)	CSMWHM3	D-1
913	Simple Stem Double Sided, 30° Angulation Mouth Mirror #4 (7/8 inch-22mm)	CSMWHM4	D-1
914	Simple Stem Double Sided, 30° Angulation Mouth Mirror #5 (15/16 inch-24mm)	SSMWHM5	D-1
915	Cone Socket Mirror #4 (7/8 inch-22mm)	SSMWHM6	D-1
916	Cone Socket Mirror #5 (15/16 inch-24mm)	SS	
1013	Micro surgical Mirrors-1 small Rectangular	MICMS1	D-2
1014	Micro surgical Mirrors-2 small Rectangular	MICMS2	D-2
1015	Micro surgical Mirrors-3 small Rectangular	MICMS3	D-2
1016	Micro surgical Mirrors- small round	MICMSR	D-2
1017	Micro surgical Mirrors - Large Round	MICMLR	D-2
1018	Micro surgical Mirrors - 1 Vervart	MICMV1	D-3
1019	Micro surgical Mirrors - 2 Vervart	MICMV2	D-3
1020	Micro Mirror Round 3mm	MICMR3	D-3
1021	Micro Mirror Round 5mm	MICMR5	D-3
1022	Micro Mirror 2mmx7mm	MICM2X7	D-3
1023	Micro Mirror 3mmx9mm	MICM3X9	D-3
1071	Dressing Pliers - 1 Perry	DP1P	D-5
1072	Dressing Pliers - 2 College	DP2C	D-5
1073	Dressing Pliers - 3	DP3	D-5
1074	Dressing Pliers - 17	DP17	D-5
1075	Dressing Pliers - 18	DP18	D-5
1076	Dressing Pliers - 18L Locking	DP18L	D-6
1077	Dressing Pliers - 19	DP19	D-6
1078	Dressing Pliers - Utility Pick -Up Dressing Plier	DPU17	D-6
1079	Dressing Pliers - 317 Smooth Beak	DP317	D-6
1080	Dressing Pliers - 318 Foil/Meriam, Serrated Beaks	DP318	D-6
1121	Explorers - 3	EXS3	D-7
1122	Explorers - 3A	EXS3A	D-7
1123	Explorers - 6 6.5mm	EXS6	D-7
1124	Explorers - 6A 12mm	EXS6A	D-7
1125	Explorers - 6L 9mm	EXS6L	D-7
1126	Explorers - 6XL 12.5mm	EXSXL	D-7
1127	Explorers - 8	EXS8	D-8
1128	Explorers - 9	EXS9	D-8
1129	Explorers - 23 Shepherd's Hook	EXS23	D-8
1130	Explorers - 54	EXS54	D-8
1131	Explorers - TU17 Wilkins /Tufts	EXSTU17	D-8

1132	Explorers - 17	EXS17	D-8
1133	Explorers - 1	EXD1	D-9
1134	Explorers - 2	EXD2	D-9
1135	Explorers - 2A Pigtail	EXD2A	D-9
1136	Explorers - 2H	EXD2H	D-9
1137	Explorers - 3	EXD3	D-9
1138	Explorers - TU 17/23 Wilkins /Tufts	EXDTU17/23	D-9
1139	Explorers - 3ES	EXD3ES	D-9
1140	Explorers - 3CH Cowhrn	EXD3CH	D-9
1141	Explorers - 5	EXD5	D-10
1142	Explorers - 2R/2L	EXD2R/2L	D-10
1143	Explorers - 6/3	EXD6/3	D-10
1144	Explorers - 23/6	EXD23/6	D-10
1145	Explorers - EN1	EXDEN1	D-10
1146	Explorers - 3A/TU17	EXD3A/TU17	D-10
1147	Explorers - 11/12ODU (Old Dominion University)	EXD11/12ODU	D-10
1148	Explorers - 11/12 A-5	EXD11/12A5	D-10
1200	Color Probes - CP-2	CP2	D-11
1201	Color Probes - CP-8	CP8	D-11
1202	Color Probes - CP-10	CP10	D-11
1203	Color Probes - CP-11	CP11	D-11
1204	Color Probes - CP-12	CP12	D-11
1205	Color Probes - CP-18	CP18	D-11
1206	Color Probes - CP-11.5B Screening Probe	CP11.5B	D-11
1207	Color Probes - CP-UNC12 (University of North Carolina)	CPUNC12	D-12
1208	Color Probes - CP-UNC15 (University of North Carolina)	CPUNC15	D-12
1209	Color Probes - PQ-O Michigan	CPQ-O	D-12
1210	Color Probes - PQ-W Williams	CPQW	D-12
1211	Color Probes - Moffitt/Maryland	CPMM	D-12
1212	Color Probes - Williams Novatec NT2	CPWNT2	D-12
1213	Color Probes - CP-11 Novatec	CP11NT	D-13
1214	Color Probes - CP-12 Novatec	CP12NT	D-13
1215	Color Probes - UNC12 Novatec	CPUNC12NT	D-13
1216	Color Probes - UNC15 Novatec	CPUNC15NT	D-13
1217	Color Probes - Thin William' s Novatec	CPWNT	D-13
1218	Color Probes - Glickman 26G	CP26G	D-13
1219	Color Probes - Goldman Fox-Flat	CPGF	D-14
1220	Color Probes - Merritt B	CPMB	D-14
1221	Color Probes - Michigan "O"	CPM-O	D-14
1222	Color Probes - Thin Williams	CPWT	D-14
1223	Color Probes - Williams	CPW	D-14
1224	Color Probes - CP-11.5B Novatec Screening Probe	CP11.5BNT	D-14
1225	Color Probes - CP-12/ Screening Probe	CP12/11.5B	D-15
1226	Color Probes - CP-UNC15/ Screening Probe	CPUNC15/11.5B	D-15
1227	Color Probes - Thin Williams / Screening Probe	CPTW/11.5B	D-15
1228	Color Probes - 3N Nabers	CP3N	D-15
1229	Color Probes - Williams DE Off set	CPWOSDE	D-15

1230	Color Probes - CP-12/Thin Williams	CP12/WT	D-15
1231	Color Probes - UT3/4 Cattoni	CP12UT3/4	D-16
1232	Color Probes - Goldman-Fox Williams	CPGF/W	D-16
1233	Color Probes - Goldman-Fox/Thin Williams	CPGFWT	D-16
1234	Color Probes - 1N Nabers	CP1N	D-16
1235	Color Probes - 2N Nabers	CP2N	D-16
1236	Color Probes - Q-2N Nabers	CPQ2N	D-16
1290	EX-PROBES - 23/QW Williams	XP23/QW	D-17
1291	EX-PROBES - 23/QOW Thin Williams	XP23/WT	D-17
1292	EX-PROBES - 23/CP-8	XP23/CP8	D-17
1293	EX-PROBES - 23/CP-11	XP23/CP11	D-17
1294	EX-PROBES - 23/CP-11.5B Screening Probe	XP23/CP11.5B	D-17
1295	EX-PROBES - 23/CP-12	XP23/CP12	D-17
1296	EX-PROBES - NT3 CP-12 Novatech	XP23/CP12NT3	D-17
1297	EX-PROBES - 23/UNC12 (University of North Carolina)	XP23/UNC12	D-17
1298	EX-PROBES - 23/UNC15 (University of North Carolina)	23/UNC15	D-18
1299	EX-PROBES - TU17/UNC15	XPTU17/UNC15	D-18
1300	EX-PROBES - 23/W Williams	XP23/WT	D-18
1301	EX-PROBES - 6/OW Thin Williams	XP6/WT	D-18
1302	EX-PROBES - 6/W Williams	XP6/WT	D-18
1303	EX-PROBES - 17/OW Thin Williams	XP17/WT	D-18
1304	EX-PROBES - 3A/UNC15	XP3A/UNC15	D-18
1305	EX-PROBES - 23/GF Goldman Fox	XP23/GF	D-18
1306	EX-PROBES - 23/O Michigan	XP23/O	D-18
1350	ENDO EXPLORERS - DG16	EXDG16	E-1
1351	ENDO EXPLORERS - DG16/23	EXDG16/23	E-1
1352	ENDO EXPLORERS - DG16/17	EXDG16/17	E-1
1353	ENDO ENDO XPLORERS - 1 Abou Rass Explorer	EX1AR	E-1
1354	ENDO EXCAVATORS - 2 Glick	EXCGL2	E-1
1355	ENDO EXCAVATORS - 31L	EXC31L	E-1
1356	ENDO EXCAVATORS - 32L	EXC32L	E-2
1357	ENDO EXCAVATORS - 33L	EXC33L	E-2
1358	ENDO EXCAVATORS - 31W Oval Spoon	EXC31W	E-2
1400	RCP30 - ISO Size 30 24mm	RCP30	E-3
1401	RCP40 - ISO Size 40 24mm	RCP40	E-3
1402	RCP50 - ISO size 50 24mm	RCP50	E-3
1403	RCP60 - ISO Size 60 24mm	RCP60	E-3
1404	RCP1/3 - 1/3 21mm	RCP1/3	E-3
1405	RCP5/7 - 5/7 21mm	RCP5/7	E-3
1406	RCP9/11 - 9/11 21mm	RCP9/11	E-4
1407	RCPGL1 - 1 Glick, Blade/Pluggler	RCPGL1	E-4
1408	RCPL1 - L1 18mm	RCPL1	E-4
1409	RCPL2 - L2 18mm	RCPL2	E-4
1410	RCPL3 - L3 18mm	RCPL3	E-4
1411	RCPL4 - L4 18mm	RCPL4	E-4
1412	RCP8A - 29mm	RCP8A	E-5
1413	RCP8-1/2, 29mm	RCP8A-1/2	E-5

1414	RCP9A - 29mm	RCP9A	E-5
1415	RCP9-1/2A - 29mm	RCP9A-1/2	E-5
1416	ARCP10A - 29mm	RCP10A	E-5
1417	RCP10-1/2A, 29mm	RCP10A-1/2	E-5
1418	RCP11A - 29mm	RCP11A	E-5
1419	RCP11-1/2A, 29mm	RCP11-1/2	E-5
1420	RCP12A - 29mm	RCP12A	E-5
1421	RCP8P - 22mm	RCP8	E-6
1422	RCP8-1/2P, 22mm	RCP8-1/2	E-6
1423	RCP9P - 9 22mm	RCP9	E-6
1424	RCP9-1/2P - 9 1/2 22mm	RCP9-1/2	E-6
1425	RCP10P - 10 22mm	RCP10	E-6
1426	RCP10-1/2P - 10 1/2 22mm	RCP10-1/2	E-6
1427	RCP11P - 11 22mm	RCP11	E-6
1428	RCP11-1/2P - 11 1/2 22mm	RCP11-1/2	E-6
1429	RCP12P - 12 22mm	RCP12A	E-6
1480	RCSD11NT- D11 (NITI)	RCSD11NT	E-7
1481	RCSD11TNT- D11T (NITI)	RCSD11TNT	E-7
1482	RCSD11TSNT- D11TS (NITI)	RCSD11TSNT	E-7
1483	RCS25SNT - 255 (NITI)	RCS25SNT	E-7
1484	RCS25NT- 25 (NITI)	RCS25NT	E-7
1485	RCS40NT- 40 (NITI)	RCS40NT	E-7
1486	RCS40SNT- 40S (NITI)	RCS40SNT	E-7
1487	RCSMA57NT- MA57 (NITI)	RCSMA57NT	E-7
1488	RCS3NT - 3NT (NITI)	RCS3NT	E-7
1540	RCS30 - ISO Size 30 24mm	RCS30	E-8
1541	RCS40- ISO Size 40 24mm	RCS40	E-8
1542	RCS50- ISO Size50 28mm	RCS50	E-8
1543	RCS60- ISO Size 60 28mm	RCS60	E-8
1544	RCSGP1- GP1 20mm	RCSGP1	E-8
1545	RCS1S- 1S Wakai 22mm	RCSW1S	E-8
1546	RCS2S- 2S Wakai 22mm	RCSW2S	E-8
1547	RCD11- D11 Standard 23mm	RCSD11	E-8
1548	RCSD11S- D11S Thin 23mm	RCSD11S	E-9
1549	RCSD11T- D11T Thinner 22mm	RCSD11T	E-9
1550	RCSD11TS- D11TS Thinner 21mm	RCSD11TS	E-9
1551	RCS3- 3 28mm	RCS3	E-9
1552	RCSGP2- GP2 25mm	RCSGP2	E-9
1553	RCSOP - O 25mm	RCSO	E-9
1554	RCS00 - OO 25mm	RCSOO	E-9
1555	RCSMA57 - MA57 27mm	RCSM57	E-9
1556	RCSGP3 - GP3 30mm	RCSGP3	E-9
1600	MCARAS1-S1 Abou-Rass Straight Apical	MCARAS1	E-10
1601	MCARAS2-S2 Abou-Rass Straight Apical	MCARAS2	E-10
1602	MCARAS4-S4 Abou-Rass Straight Apical	MCARAS3	E-10
1603	MCARAS3-S3 Abou-Rass Straight Apical	MCARAS4	E-10
1604	MCVA4-4 Velart Apica	MCVA4	E-10

1605	MICRO SURGICAL EXCAVATORS-1 Abou-Rass	MEXCAR1	E-10
1606	MICRO SURGICAL EXCAVATORS-2 Abou-Rass	MEXCAR2	E-10
1607	MICRO SURGICAL EXCAVATORS-3 Abou-Rass	MEXCAR3	E-10
1608	MICRO SURGICAL EXCAVATORS-4 Abou-Rass	MEXCAR4	E-10
1650	APICAL PUGGERS - 0.5 Abou-Rass Straight	APARS0.5	E-11
1651	APICAL PUGGERS - 1 Abou-Rass Straight	APABS1.0	E-11
1652	APICAL PUGGERS - 1.5 Abou-Rass Straight	APARS1.5	E-11
1653	APICAL PUGGERS - 2 Abou-Rass Straight	APARS2.0	E-11
1654	APICAL PUGGERS - 3 Abou-Rass Straight	APARS3.0	E-11
1655	APICAL PUGGERS - 0.5 Abou-Rass Angled	APAAR0.5/0.6	E-12
1656	APICAL PUGGERS - 1 Abou-Rass Angled	APAAR0.8	E-12
1657	APICAL PUGGERS - 1.5 Abou-Rass Angled	APAAR0.1.5/1.5	E-12
1658	APICAL PUGGERS - 2 Abou-Rass Angled	APAAR1.5/2.0	E-12
1659	APICAL PUGGERS - 3 Abou-Rass Angled	APAAR2.0/2.0	E-12
1700	RETRO FILLING - Left Apical	PFLA	E-13
1701	RETRO FILLING - Right Apical	PFRA	E-13
1702	APICAL BURNISHERS, Straight - 1 Abou-Rass Straight	ABARS1	E-13
1703	APICAL BURNISHERS, Straight- 2 Abou-Rass Straight	ABARS2	E-13
1704	APICAL BURNISHERS, Straight- 3 Abou-Rass Straight	ABARS3	E-13
1705	APICAL BURNISHERS, Straight- 4 Abou-Rass Straight	ABARS4	E-13
1706	APICAL BURNISHERS Straight- 5 Abou-Rass Straight	ABARS5	E-13
1707	APICAL BURNISHERS, Angled - 1 Abou-Rass Angled	ABAAR10.5	E-14
1708	APICAL BURNISHERS, Angled - 2 Abou-Rass Angled	ABAAR21.0	E-14
1709	APICAL BURNISHERS, Angled - 3 Abou-Rass Angled	ABAAR31.25	E-14
1710	APICAL BURNISHERS, Angled - 4 Abou-Rass Angled	ABAAR41.75	E-14
1711	APICAL BURNISHERS, Angled - 5 Abou-Rass Angled	ABAAR52.0	E-14
1712	APICAL PLUGGERS/BURNISHERS - Small Left Apical	APBSL	E-15
1713	APICAL PLUGGERS/BURNISHERS - Small Right Apical	APBSR	E-15
1714	APICAL PLUGGERS/BURNISHERS - Medium Left Apical	APBML	E-15
1715	APICAL ANGLED PLUGGERS/BURNISHERS - Medium Right Apical	APBMR	E-15
1716	APICAL ANGLED PLUGGERS/BURNISHERS - Large Left Apical	ABLL	E-15
1717	APICAL ANGLED PLUGGERS/BURNISHERS - Large Right Apical	APBLR	E-16
1718	APICAL PLUGGERS/BURNISHERS - Universal Small Apical	APBUS	E-16
1719	APICAL PLUGGERS/BURNISHERS - Universal Medium Apical	APBUM	E-16
1720	APICAL PLUGGERS/BURNISHERS - Universal Large Apical	APBUL	E-16
1770	Steiglitz, Straight 14cm/5-1/2"	RFSTR	E-17
1771	Steiglitz, 45°, 13cm/5-1/4"	RF45	E-17
1772	Steiglitz, 90°, 12cm/4-3/2"	RF90	E-17
1773	Peet, 45°, 12cm/4-3/4 "	RFP	E-17
1774	Rubber Dam Clamp Forceps, (4), 17.5cm7"	RDF4	E-17
1775	Rubber Dam Clamp Forceps, 17cm/6-1/2"	RDF	E-17
1776	Rubber Dam Frame-Child, 13cm/5"	RDCF5	E-17
1777	Rubber Dam Frame-Adult 15cm/6"	RDAF6	E17
1778	RUBBER DAM PUNCH - Ainsworth 17cm/6 3/4 "	RDPA	E17
1779	Rubber Dam Punch1 16cm/6.75"	RDP	E-18
1830	RUBBER DAM CLAMP - 8A	RDC8A	E-18
1831	RUBBER DAM CLAMP - W8A	RDCW8A	E-18

1832	RUBBER DAM CLAMP - 8AD	RDC8AD	E-18
1833	RUBBER DAM CLAMP - 10	RDC10	E-18
1834	RUBBER DAM CLAMP - 11	RDC11	E-18
1835	RUBBER DAM CLAMP - 14	RDC14	E-18
1836	RUBBER DAM CLAMP - 14A	RDC14A	E-18
1837	RUBBER DAM CLAMP - W14A	RDCW14A	E-18
1838	RUBBER DAM CLAMP - 18	RDC18	E-18
1839	RUBBER DAM CLAMP - 24	RDC24	E-19
1840	RUBBER DAM CLAMP - 25	RDC25	E-19
1841	RUBBER DAM CLAMP - 26N	RDC26N	E-19
1842	RUBBER DAM CLAMP - 56	RDC56	E-19
1843	RUBBER DAM CLAMP - 56S	RDC56S	E-19
1844	RUBBER DAM CLAMP - 4	RDC4	E-19
1845	RUBBER DAM CLAMP - W4	RDCW4	E-19
1846	RUBBER DAM CLAMP - 5	RDC5	E-19
1847	RUBBER DAM CLAMP - 8	RDC6	E-19
1848	RUBBER DAM CLAMP - 27	RDC27	E-19
1849	RUBBER DAM CLAMP - 30	RDC30	E-19
1850	RUBBER DAM CLAMP - 31	RDC31	E-19
1851	RUBBER DAM CLAMP - 201	RDC201	E-19
1852	RUBBER DAM CLAMP - 205	RDC205	E-19
1853	RUBBER DAM CLAMP - 3	RDC3	E-19
1854	RUBBER DAM CLAMP - W3	RDCW3	E-19
1855	RUBBER DAM CLAMP - 7	RDC7	E-19
1856	RUBBER DAM CLAMP - W7	RDCW7	E-19
1857	RUBBER DAM CLAMP - 7A	RDC7A	E-19
1858	RUBBER DAM CLAMP - 12A	RDC12A	E-19
1859	RUBBER DAM CLAMP - 31A	RDC31A	E-19
1860	RUBBER DAM CLAMP - 26	RDC26N	E-19
1861	RUBBER DAM CLAMP - 28	RDC28	E-19
1862	RUBBER DAM CLAMP - 0	RDC0	E-19
1863	RUBBER DAM CLAMP - 1	RDC1	E-19
1864	RUBBER DAM CLAMP - 1A	RDC1A	E-19
1865	RUBBER DAM CLAMP - 2	RDC2	E-19
1866	RUBBER DAM CLAMP - W2	RDCW2	E-20
1867	RUBBER DAM CLAMP - 2A	RDC2A	E-20
1868	RUBBER DAM CLAMP - 2AS	RDC2AS	E-20
1869	RUBBER DAM CLAMP - 27N	RDC27N	E-20
1870	RUBBER DAM CLAMP - 29	RDC29	E-20
1871	RUBBER DAM CLAMP - 6	RDC6	E-20
1872	RUBBER DAM CLAMP - 9	RDC9	E-20
1873	RUBBER DAM CLAMP - 9S	RDC9S	E-20
1874	RUBBER DAM CLAMP - 212	RDC212	E-20
1875	RUBBER DAM CLAMP - 212SA	RDC212SA	E-20
1876	RUBBER DAM CLAMP - 00	RDC00	E-20
1877	Rubber Dam Clamp Organizing Board	RDCOB	E-20
1920	IMPLANT PERIOSTEALS - 1 LaBanca Periosteal	IMPLAB1	I-1

1921	IMPLANT PERIOSTEALS - 2 LaBanca Modified	IMPLABM2	I-1
1922	IMPLANT PERIOSTEALS - Buser	IMPBUS	I-1
1923	IMPLANT PERIOSTEALS - P9/PA	IMP9/PA	I-1
1924	IMPLANT PERIOSTEALS - 25	IMP25	I-1
1970	1/2 Buser Modified Periodontal Chisel 3mm/4mm	IPCBM1/2	I-2
1971	2/3 Buser Modified Periodontal Chisel 4mm/5mm	IPCBM2/3	I-2
1972	1/3 Buser Modified Periodontal Chisel 3mm/5mm	IPCBM1/3	I-2
1973	1mm/4mm Jovanovic 13KTG Periodontal Chisel	IPCJOV1/4	I-2
1974	2mm/5mm Jovanovic 13KTG Periodontal Chise	IPCJOV2/5	I-2
1975	3mm/6mm Jovanovic 13KTG Periodontal Chise	IPCJOV3/6	I-2
1976	P24G/13KL Palacci Modified	IPC24G/13KL	I-2
2020	Bone Scoop Jovanovic None Scoop 5mm	IBSJOV5	I-3
2021	BONE SCOOPS - Bone Graft Scoop	IBSBG	I-3
2022	BONE SCOOPS - Labanca Bone Curette 7mm	ICLBC	I-3
2023	Bone Pluggers LaBanca, 3x1.5mm	IBPLAB	I-3
2024	BONE PLUGGERS - 3mm, No.4 Jovanovic	IBPJ4	I-3
2025	BONE PLUGGERS- Bone Packer DE 5-7mm	IBPP5-7	I-3
2026	BONE PLUGGERS - Bone Packer DE 3-5mm	IBPS3-5	I-3
2052	Sinus Lift Cobra 3.5mm	ISC3.5	I-4
2053	Sinus Lift Cobra In/Out (2-4-5-6-8-10mm)	ISCI/O	I-4
2054	Sinus Lift Mushroom 2.8/3.8mm (2-4-5-6-8-10mm)	IMSL2.8/3.8	I-4
2055	Sinus Lift Mushroom 4.8/5.8mm (2-4-5-6-8-10mm)	IMSL4.8/5.8	I-4
2081	Bone Graft Pluggers 1.0mm/1.7mm, (8-11-13-15-18)	IBG1.0/1.7	I-5
2082	Bone Graft Pluggers 1.7mm/2.0mm, (8-11-13-15-18)	IBG1.7/2.0	I-5
2083	Bone Graft Pluggers 2.4mm/2.7mm, (8-11-13-15-18)	IBG2.4/2.7	I-5
2084	Bone Graft Pluggers 3.3mm/4.0mm - (8-11-13-15-18)	IBG3.3/4.0	I-5
2085	Bone Graft Pluggers 3.3mm/4.0mm, Serrated (8-11-13-15-18)	IBGS3.3/4.0	I-5
2090	TPI1-Tunneling Instrument -1	ITPI1	I-6
2091	TPI2-Tunneling Instrument -2	ITPI2	I-6
2092	TPI3-Tunneling Instrument -3	ITPI3	I-6
2093	TPI4-Tunneling Instrument -4 Curved	ITPI4	I-6
2094	TPI5-Tunneling Instrument -5	ITPI5	I-6
2095	TPI6-Tunneling Instrument -4 Spade	ITPI6	I-6
2096	TPI7-Tunneling Instrument -6	ITPI7	I-6
2097	Urban Tunneling Knife	IUTK	I-7
2098	Tunneling Knife 1	ITK1	I-7
2099	Tunneling Knife 2 (Short Handle)	ITK2S	I-7
2100	Tunneling Knife 2	ITK2	I-7
2101	Giles de Quincey Tunneling Instrument 1	IGDQK1	I-7
2102	Giles de Quincey Tunneling Instrument 2	IGDQK2	I-7
2103	Steigmann Soft Tissue Splitter	ISSTS	I-7
2104	Steigmann Soft Tissue Splitter, Modified	ISSTSM	I-7
2110	OSTEOTOME CHISEL- Single Edge Straight 4mm, (7-10-13-15-18mm)	IOCSSE4	I-8
2111	OSTEOTOME CHISEL- Single Edge Curved 4mm, (7-10-13-15-18mm)	IOCCSE4	I-8
2112	OSTEOTOME CHISEL- Single Edge Straight 7.5mm, (7-10-13-15-18mm)	IOCSSE7.5	I-8
2113	OSTEOTOME CHISEL- Single Edge Curved 7.5mm, (7-10-13-15-18mm)	IOCCSE7.5	I-8
2114	OSTEOTOME CHISEL- Jovanovic Osteotome 6mm	IOCJOV6	I-8

2115	OSTEOTOME CHISEL, Single Edge 3mm-(7-10-13-15-18mm)	IOCSE3	I-8
2116	OSTEOTOME CHISEL, Double Edge 3mm-(7-10-13-15-18mm)	IOCDE3	I-8
2117	OSTEOTOME CHISEL- Single Edge 6mm -(7-10-13-15-18mm)	IOCSE6	I-8
2118	OSTEOTOME CHISEL-Double Edge 6mm- (7-10-13-15-18mm)	IOCDE6	I-8
2119	OSTEOTOME CHISEL-Single Edge 9mm-(7-10-13-15-18mm)	IOCSE9	I-8
2120	OSTEOTOME CHISEL-Double Edge 9mm STR (7-10-13-15-18mm)	IOCDE9	I-8
2130	Tissue Punch STR 4mm	ITPS4	I-9
2131	Tissue Punch STR 5mm	ITPS5	I-9
2132	Tissue Punch STR 6mm	ITPS6	I-9
2133	Tissue Punch CVD at 90°, 2mm	ITPC2	I-9
2134	Tissue Punch CVD at 90°, 3mm	ITPC3	I-9
2135	Tissue Punch CVD at 90°, 4mm	ITPC4	I-9
2136	Tissue Punch CVD at 90°, 5mm	ITPC5	I-9
2137	Tissue Punch CVD at 90°, 6mm	ITPC6	I-9
2230	SINUS LIFT INSTRUMENT - 1 MILLS SINUS	IMS1	I-10
2231	SINUS LIFT INSTRUMENT - 2 MILLS SINUS	IMS2	I-10
2232	SINUS LIFT INSTRUMENT - 3 MILLS SINUS	IMS3	I-10
2233	SINUS LIFT INSTRUMENT - 1 PIKOS SINUS	IPS1	I-10
2234	SINUS LIFT INSTRUMENT - 2 PIKOS SINUS	IPS2	I-10
2235	SINUS LIFT INSTRUMENT - 3 PIKOS SINUS	IPS3	I-10
2236	SINUS LIFT INSTRUMENT - 4 PIKOS SINUS	IPS4	I-11
2237	SINUS LIFT INSTRUMENT - 5 PIKOS SINUS	IPS5	I-11
2238	SINUS LIFT INSTRUMENT - 6 PIKOS SINUS	IPS6	I-11
2239	SINUS LIFT INSTRUMENT - 1 SINUS	IS1	I-11
2240	SINUS LIFT INSTRUMENT - 2 SINUS	IS2	I-11
2241	SINUS LIFT INSTRUMENT - 4 SINUS	IS3	I-12
2242	SINUS LIFT INSTRUMENT - 5 SINUS	IS4	I-12
2243	SINUS LIFT INSTRUMENT - 6 SINUS	IS5	I-12
2244	SINUS LIFT INSTRUMENT - 1 SINUS FREER	ISFR1	I-12
2245	SINUS LIFT INSTRUMENT - SINUS FREER-CU	ISF-CU	I-13
2246	SINUS LIFT INSTRUMENT - MISCH SINUS-1	IMISCH/S-1	I-13
2247	SINUS LIFT INSTRUMENT - WALLACE SINUS-1	IWS-1	I-13
2248	SINUS LIFT INSTRUMENT - WALLACE SINUS-2	IWS-2	I-13
2255	Osteotomes- Tapered Concave Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)	ITCCSO2.0	I-14
2256	Osteotomes- Tapered Concave Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)	ITCCSO3.5	I-14
2257	Osteotomes- Tapered Concave Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)	ITCCSO4.3	I-14
2258	Osteotomes- Tapered Concave Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)	ITCCSO5.0	I-14
2259	Osteotomes- Tapered Concave Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)	ITCCSO5.5	I-14
2260	Osteotomes- Tapered Concave-Angulated -2.0mm (4-6-8-10-13-16-18-20-23-26)	ITCCAO2.0	I-14
2261	Osteotomes- Tapered Concave-Angulated - 3.5mm (4-6-8-10-13-16-18-20-23-26)	ITCCAO3.5	I-14
2262	Osteotomes- Tapered Concave-Angulated - 4.3mm (4-6-8-10-13-16-18-20-23-26)	ITCCAO4.3	I-14
2263	Osteotomes- Tapered Concave-Angulated - 5.0mm (4-6-8-10-13-16-18-20-23-26)	ITCCAO5.0	I-14
2264	Osteotomes- Tapered Concave-Angulated - 5.5mm (4-6-8-10-13-16-18-20-23-26)	ITCCAO5.5	I-14
2270	Osteotomes- Tapered Convex -Straight - 2.0mm (4-6-8-10-13-16-18-20-23-26)	ITCVSO2.0	I-15
2271	Osteotomes- Tapered Convex -Straight - 3.5mm (4-6-8-10-13-16-18-20-23-26)	ITCVSO3.5	I-15
2272	Osteotomes- Tapered Convex -Straight - 4.3mm (4-6-8-10-13-16-18-20-23-26)	ITCVSO4.3	I-15
2273	Osteotomes- Tapered Convex -Straight - 5.0mm (4-6-8-10-13-16-18-20-23-26)	ITCVSO5.0	I-15

2274	Osteotomes- Tapered Convex -Straight - 5.5mm (4-6-8-10-13-16-18-20-23-26)	ITCVSO5.5	I-15
2275	Osteotomes- Tapered Convex -Angulated - 2.0mm (4-6-8-10-13-16-18-20-23-26)	ITCVAO2.0	I-15
2276	Osteotomes- Tapered Convex -Angulated - 3.5mm (4-6-8-10-13-16-18-20-23-26)	ITCVAO3.5	I-15
2277	Osteotomes- Tapered Convex -Angulated - 4.3mm (4-6-8-10-13-16-18-20-23-26)	ITCVAO4.3	I-15
2278	Osteotomes- Tapered Convex -Angulated - 5.0mm (4-6-8-10-13-16-18-20-23-26)	ITCVAO5.0	I-15
2279	Osteotomes- Tapered convex -Angulated-5.5mm (4-6-8-10-13-16-18-20-23-26)	ITCVAO5.5	I-15
2281	Osteotomes- Bone Pusher, Straight (Slightly Convex) - 2.0mm (4-6-8-10-13-16-18-20-23-26)	IBPS2.0	I-16
2282	Osteotomes- Bone Pusher, Straight (Slightly Convex) - 3.5mm (4-6-8-10-13-16-18-20-23-26)	IBPS3.5	I-16
2283	Osteotomes- Bone Pusher, Straight (Slightly Convex) - 4.3mm (4-6-8-10-13-16-18-20-23-26)	IBPS4.3	I-16
2284	Osteotomes- Bone Pusher, Straight (Slightly Convex) - 5.0mm (4-6-8-10-13-16-18-20-23-26)	IBPS5.0	I-16
2285	Osteotomes- Bone Pusher, Straight (Slightly Convex) - 5.5mm (4-6-8-10-13-16-18-20-23-26)	IBPS5.5	I-16
2286	Osteotomes- Bone Pusher, Angulated (Slightly Convex) -2.0mm (4-6-8-10-13-16-18-20-23-26)	IBPA2.0	I-16
2287	Osteotomes- Bone Pusher, Angulated (Slightly Convex) - 3.5mm (4-6-8-10-13-16-18-20-23-26)	IBPA3.5	I-16
2288	Osteotomes- Bone Pusher, Angulated (Slightly Convex) - 4.3mm (4-6-8-10-13-16-18-20-23-26)	IBPA4.3	I-16
2289	Osteotomes- Bone Pusher, Angulated (Slightly Convex) - 5.0mm (4-6-8-10-13-16-18-20-23-26)	IBPA5.0	I-16
2290	Osteotomes- Bone Pusher, Angulated (Slightly Convex) - 5.5mm (4-6-8-10-13-16-18-20-23-26)	IBPA5.5	I-16
2299	Interchangeable Osteotome Handle	IBSKOIH	I-17
2300	Osteotomes- Tapered Convex Straight Tip - 2.7mm (7-10-13-15-18-20)	ICVOS2.7	I-17
2301	Osteotomes- Tapered Convex Straight Tip - 3.2mm (7-10-13-15-18-20)	ICVOS3.2	I-17
2302	Osteotomes- Tapered Convex Straight Tip - 3.7mm (7-10-13-15-18-20)	ICVOS3.7	I-17
2303	Osteotomes- Tapered Convex Straight Tip - 4.2mm (7-10-13-15-18-20)	ICVOS4.2	I-17
2304	Osteotomes- Tapered Convex Straight Tip - 5.0mm (7-10-13-15-18-20)	ICVOS5.0	I-17
2305	Osteotomes- Tapered Convex-Angulated Tip -2.7mm (7-10-13-15-18-20)	ICVOA2.7	I-17
2306	Osteotomes- Tapered Convex-Angulated Tip - 3.2mm (7-10-13-15-18-20)	ICVOA3.2	I-17
2307	Osteotomes- Tapered Convex-Angulated Tip - 3.7mm (7-10-13-15-18-20)	ICVOA3.7	I-17
2308	Osteotomes- Tapered Convex-Angulated Tip - 4.2mm (7-10-13-15-18-20)	ICVOA4.2	I-17
2309	Osteotomes- Tapered Convex-Angulated Tip - 5.0mm (7-10-13-15-18-20)	ICVOA5.0	I-17
2310	Osteotomes- Tapered Concave Tip-Straight - 2.7mm (7-10-13-15-18-20)	ICCOS2.7	I-18
2311	Osteotomes- Tapered Concave Tip-Straight - 3.2mm (7-10-13-15-18-20)	ICCOS3.2	I-18
2312	Osteotomes- Tapered Concave Tip-Straight - 3.7mm (7-10-13-15-18-20)	ICCOS3.7	I-18
2313	Osteotomes- Tapered Concave Tip-Straight - 4.2mm (7-10-13-15-18-20)	ICCOS4.2	I-18
2314	Osteotomes- Tapered Concave Tip-Straight - 5.0mm (7-10-13-15-18-20)	ICCOS5.0	I-18
2315	Osteotomes- Tapered Concave Tip-Angulated - 2.7mm (7-10-13-15-18-20)	ICCOA2.7	I-18
2316	Osteotomes- Tapered Concave Tip-Angulated - 3.2mm (7-10-13-15-18-20)	ICCOA3.2	I-18
2317	Osteotomes- Tapered Concave Tip-Angulated - 3.7mm (7-10-13-15-18-20)	ICCOA3.7	I-18
2318	Osteotomes- Tapered Concave Tip-Angulated - 4.2mm (7-10-13-15-18-20)	ICCOA4.2	I-18
2319	Osteotomes- Tapered Concave Tip-Angulated-5.0mm (7-10-13-15-18-20)	ICCOA5.0	I-18
2320	Bone Pusher-Slightly Convex Tip-Straight - 2.7mm (7-10-13-15-18-20)	IBPOS2.7	I-18
2321	Bone Pusher-Slightly Convex Tip-Straight - 3.2mm (7-10-13-15-18-20)	IBPOS3.2	I-18
2322	Bone Pusher-Slightly Convex Tip-Straight - 3.7mm (7-10-13-15-18-20)	IBPOS3.7	I-18
2323	Bone Pusher-Slightly Convex Tip-Straight - 4.2mm (7-10-13-15-18-20)	IBPOS4.2	I-18
2324	Bone Pusher-Slightly Convex Tip-Straight - 5.0mm (7-10-13-15-18-20)	IBPOS5.0	I-18
2325	Bone Pusher-Slightly Convex Tip-Angulated - 2.7mm (7-10-13-15-18-20)	IBPOA2.7	I-18
2326	Bone Pusher-Slightly Convex Tip-Angulated - 3.2mm (7-10-13-15-18-20)	IBPOA3.2	I-18
2327	Bone Pusher-Slightly Convex Tip-Angulated - 3.7mm (7-10-13-15-18-20)	IBPOA3.7	I-18
2328	Bone Pusher-Slightly Convex Tip-Angulated - 4.2mm (7-10-13-15-18-20)	IBPOA4.2	I-18
2329	Bone Pusher-Slightly Convex Tip-Angulated - 5.0mm (7-10-13-15-18-20)	IBPOA5.0	I-18

2382	Castroviejo Caliper, Straight, 17cm/6-3/4"	ICVCS17	I-19
2383	Weiss Modified Castroviejo Calipers 0-40mm 17cm/6 3/4 "	IWMCVC17	I-19
2384	Castroviejo Calipers-Straight 0-40mm 18cm/7"	ICVCS18	I-19
2385	Castroviejo Calipers-Straight 0-40mm 8.5cm/3.35"	ICVCS8.5	I-19
2386	Castroviejo Calipers Angled 0-40mm 18cm/7"	ICVCA18	I-19
2387	Castroviejo Calipers-Curved 17cm/6.75"	ICVCA17	I-19
2388	Cortical Bone Clamp	ICBC	I-19
2389	Bone Grafting Triple Ring Style Syringe 1mm	IBGS1	I-19
2390	Bone Grafting Triple Ring Style Syringe 2mm	IBGS2	I-19
2391	Bone Grafting Triple Ring Style Syringe 2.5mm	IBGS2.5	I-19
2392	Bone Grafting Triple Ring Style Syringe 3mm	IBGS3	I-19
2393	Bone Grafting Triple Ring Style Syringe 4mm	IBGS4	I-19
2394	Bone Grafting Triple Ring Style Syringe 5mm	IBGS5	I-19
2395	Bone Grafting Triple Ring Style Syringe 6mm	IBGS6	I-19
2396	Bone Grafting Triple Ring Style Syringe 7mm	IBGS7	I-19
2400	Bone Grafting 1 Ring Style Syringe 3.75mm Straight	IBGSIRSS3.75	I-20
2401	Bone Grafting 1 Ring Style Syringe 3.75mm Curved	IBGSIRSC3.76	I-20
2402	Bone Grafting 1 Ring Style Syringe 5.0mm Straight	IBGSIRSS5.0	I-20
2403	Bone Grafting 1 Ring Style Syringe 5.0mm Curved	IBGSIRSC5.0	I-20
2410	Ridge Mapping Bone Width Caliper 0-25mm, Bone side reading 14cm	RBMC14	I-20
2411	Ridge Mapping Caliper Scissors Style with double sided scale 15cm	ICDSS	I-20
2412	Sliding Implant Caliper 94mm, 1.5mm tip diameter	ISIC94	I-20
2413	Stainless Steel Ruler 15cm/6"	SSR	I-20
2415	Bone Scraper/Harvesting Chisel with Replaceable Blade, Straight	IBCC	I-20
2416	Bone Collector Blade, Straight	IBCBS	I-20
2417	Bone Scraper/Harvesting Chisel with Replaceable Blade, Curved	IBCCC	I-20
2418	Bone Collector Blade, Straight	IBCBS	I-20
2419	Micro Surgery Endo Forceps 95mm	IMSEF95	I-21
2420	Micro Surgery Endo Forceps 110mm	IMSEF110	I-21
2421	Mead-2 Mallet	IM2M	I-21
2422	130g Mallet	I130GM	I-21
2423	235g Mallet	I235GM	I-21
2424	Bone Crusher	IBC	I-21
2425	Bone Morselizer	IBM	I-21
2426	Hole Punch Pliers	OHPP	O-1
2427	Tear Drop Punch Pliers	OTDPP	O-1
2428	Horizontal Punch Pliers	OHPP	O-1
2429	Vertical Punch Pliers	OVPP	O-1
2430	Pin & Ligature Cutter (Long Handle)	OPLCL	O-2
2431	Pin & Ligature Cutter, Mini Style	OPLCMS	O-2
2432	Pin & Ligature Cutter, Narrow Tip	OPLCNT	O-2
2433	Pin & Ligature Cutter, Narrow Tip at 15° Angled	OPLCNT15	O-2
2434	Pin & Ligature Cutter, Standard	OPLCS	O-2
2435	Pin & Ligature Cutter, Standard at 15° Angled	OPLCS15	O-2
2436	Ortho Hard Wire Cutter	OOHWC	O-3
2437	Ortho hard Wire Cutter at 15° Angled	OOHWC15	O-3
2438	Distal End Cutter, Mini Style Cut&Hold	ODECMSC&H	O-3

2439	Distal End Cutter, Standard Cut&Hold	ODECSC&H	O-3
2440	Distal End Cutter, Cut & Hold, (Long Handle)	ODECC7HLH	O-3
2441	Flush Cut&Hold Distal End Cutter, (Standard Handle)	OFCHDECSH	O-3
2442	Flush Cut&Hold Distal End Cutter, (Long Handle)	OFCHDECLH	O-3
2450	Utility Pliers - How Pliers, Straight	OHPS	O-4
2451	Utility Pliers - Debonding Pliers, Straight	OHOP	O-4
2452	Utility Pliers - Force Module Separating Pliers	OFMSP	O-4
2453	Utility Pliers - Coon Stle Liigature Pliers	OCSLP	O-4
2454	Utility Pliers - Adhesive Removing Pliers	OARP	O-4
2455	Utility Pliers - Weingart Pliers, Long Handle	OWPLH	O-4
2457	Utility Pliers - How Pliers, Curved	OHPC	O-4
2458	Utility Pliers - Posterior Band Removing Pliers, Long Tip	OPBRP	O-4
2459	Utility Pliers - Bracket Removing Pliers, Straight	OBRPS	O-4
2460	Utility Pliers - Bracket Removing Pliers, Angled	OBRPC	O-5
2461	Utility Pliers - Posterior Band Removing Pliers, Standard Tip	OPBRPST	O-5
2462	Utility Pliers - Debonding Pliers, Angled	ODPA	O-5
2463	Utility Pliers - Weingart Pliers, at 90° Angled	<u>OWP90</u>	O-5
2464	Utility Pliers - Weingart Pliers , Slim	OWPS	O-5
2465	Utility Pliers - Weingart Utility Pliers	OWUP	O-5
2466	Utility Pliers - Weingart Pliers, Universal at 45° Angle Tip	OWPU45AT	O-5
2467	Utility Pliers - Crown & Band Contouring Pliers	OCBCP	O-5
2468	Utility Pliers - Crown & Band Crimping Pliers	OCBCRIMP	O-5
2469	Utility Pliers - Ball Hook Crimping Pliers	OBHCP	O-5
2470	Utility Pliers - Posterior Band Removing Pliers, Short Tip	OPBRPST	O-5
2475	Bird Beak Pliers	OBBP	O-6
2476	Light Wire Pliers, Long Beak	OLWPLB	O-6
2477	Tweed Loop Forming Pliers	OTLFP	O-6
2478	Light Wire Pliers with 3 Grooves	OLWW3G	O-6
2479	Light wire Pliers with Cutter and 3 Grooves	OLWPWC3G	O-6
2480	Light wire Pliers with Cutter	OLWPWC3G	O-6
2481	Jarabak Pliers	OJBP	O-6
2482	Jarabak Pliers with Cutter	OJBPWC	O-6
2483	Ligual Arch Forming Pliers	OLAFP	O-6
2484	Nance Loop Forming Pliers	ONLFP	O-7
2485	De La Rosa Pliers	ODLARP	O-7
2486	De La Rosa Pliers with 3 Grooves	ODLARPW3G	O-7
2487	Hollow Chop Contouring Pliers	OHCCP	O-7
2488	Optical Pliers	OOP	O-7
2489	V-Stop Pliers	OVSP	O-7
2490	Young Pliers	OYP	O-7
2491	Tweed Arch Forming Pliers 1.75mm	OTAFP1.75	O-7
2492	Tweed Arch Forming Pliers 1.0mm	OTAFP1.0	O-7
2493	Torquing Pliers Female/Outer Places	OTPFO	O-8
2494	Torquing Pliers Male/Inner	OTPMI	O-8
2495	Step Pliers 1mm	OSP1	O-8
2496	Step Pliers 0.75mm	OSP0.75	O-8
2497	step Pliers 0.5mm	OSP0.5	O-8

2498	NITI Distal Cinch Back Pliers	ONTDCBP	O-8
2499	NITI Distal Cinch Back Pliers, Lingual	ONTDCB	O-8
2500	3 Prong Pliers, Standard	O3PPS	O-8
2501	3 Prong Pliers, Heavy	O3PPH	O-8
2502	3 Prong Pliers, NITI	O3PPNT	O-8
2506	Ligature Director/Band Pusher	OLDBP	O-9
2507	Small Plugger/Picker	OSPP	O-9
2508	Dull Scaler/Band Pusher	ODSBP	O-9
2509	Ligature Director Scaler	OLDS	O-9
2510	Scaler Band Pusher	OSBP	O-9
2511	Ligal Ligature Director, U-Shape	OLLDUS	O-9
2512	Dull Scaler	ODS	O-10
2513	Dull Scaler -2	ODSC	O-10
2514	Large Plugger/Picker	OLPP	O-10
2515	Picker Ligature Director	OPLD	O-10
2516	Scaler Picker	OSP	O-10
2517	Dull Scaler, Single Ended	ODSSE	O-10
2518	Ligature Director/Band Pusher	OLDBP	O-10
2519	Band Pusher with Short Tip	OBPWST	O-11
2520	Band Pusher with Long Tip	OBP	O-11
2521	Fillion Explorer/Ligature	OFELD	O-11
2522	Band Pusher/Scaler	OBPS	O-11
2523	Band Pusher/Scaler long Tip	OBPS	O-11
2524	Ligature Director with Short Tip	OLD	O-11
2525	Ligature Director with Long Tip	OLDS	O-11
2526	Ligature Torquing Ker, Single ended	OLTK	O-12
2527	Straight Applicator	OS0.02	O-12
2528	Ligature Director, Single Ended	OSELD	O-12
2529	Lingual Ligature Director 0.018"	OLLD.018	O-12
2530	Lingual Ligature Director 0.025"	OLLD.025	O-12
2531	Ligature Director/Plugger	OLDP	O-12
2532	Ligature Director Tucker	OLDT	O-12
2533	Plugger/Ligature Picker	OPLP	O-12
2566	POSTERIOR Bracket Holder Tweezer, Long	OPBTL	O-13
2567	Direct Bond Bracket Tweezer with Aligner & Positioner	ODBBTAP	O-13
2568	Buccal Tube Tweezer	OBTT	O-13
2569	Slim Buccal Tube Bonder Tweezer	OSBTBT	O-13
2570	Direct Bond Bracket Tweezer, Lingual	ODBBTL	O-13
2571	Direct Bond Bracket Tweezer, Aligner	ODBTA	O-13
2572	All Purposes Scissors	OAPS	O-13
2573	Crown Scissors	OCS	O-13
2574	Crown Scissors, Straight	OCSS	O-13
2575	Crown Scissors, Curved	OCSC	O-13
2580	Bone Bracket Positioning Gauge	OBPG	O-14
2581	Bracket Height Gauge for 0.022"/0.56mm Brackets	OBHG.022	O-14
2582	Bracket Height Gauge for 0.018"/0.46mm Brackets	OBHG.018	O-14
2583	Bracket Height Gauge, Moveable Head	OBHGMH	O-14

2585	Mathieu Narrow Tip, Enhance Double Spring Function	OMNT	O-15
2586	Mathieu Smaha Fine Tip, Enhance Double Spring Function	OMSFT	O-15
2587	Mathieu Narrow TC Tip, Enhance Double Spring Function	OMNTTC	O-15
2588	Mathieu Standard Tip, Enhance Double Spring Function	OMST	O-15
2589	Mathieu Standard TC Tip, Enhance Double Spring Function	OMSTTC	O-15
2630	SICKLE SCALERS - H6/H7 Hygienist	SH6/7	P-2
2631	SICKLE SCALERS - 6/7 Goldman	SHG6/7	P-2
2632	SICKLE SCALERS - H5/U15 Hygienist -Towner	SH5/U15	P-2
2633	SICKLE SCALERS - U15/30 Towner-Jacquette	SU15/30	P-2
2634	SICKLE SCALERS - H5/33 Hygienist Jacquette	SH5/33	P-2
2635	SICKLE SCALERS - U15/33 Towner - Jacquette	SU15/33	P-2
2636	SICKLE SCALERS - U15 Towner	SU15	P-2
2637	SICKLE SCALERS - Nevi1 Anterior	SCNEV119	P-3
2638	SICKLE SCALERS - Nevi1/H5 Anterior	SNEV1/H59	P-3
2639	SICKLE SCALERS - Nevi2 Posterior	SCNEV126	P-3
2640	SICKLE SCALERS - Nevi3 Posterior	SCNEV139	P-3
2641	SICKLE SCALERS - Nevi4 Posterior	SCNEV149	P-3
2642	SICKLE SCALERS - 30/33 Jacquette	SJ30/33	P-4
2643	SICKLE SCALERS - 31/32 Jacquette	SJ31/32	P-4
2644	SICKLE SCALERS - 34/35 Jacquette	SJ34/35	P-4
2645	SICKLE SCALERS - H3/H4 Jacquette	SH3/4	P-4
2646	SICKLE SCALERS - 14/15 Posterior Jacquette	SJ14/156	P-4
2647	SICKLE SCALERS - 204	S204	P-5
2648	SICKLE SCALERS - 204S	S204S	P-5
2649	SICKLE SCALERS - 204SD	S204SD	P-5
2650	SICKLE SCALERS - 204IUFW	SIUFW204	P-5
2651	SICKLE SCALERS - 1 Springfield Anterior	SSP1	P-5
2652	SICKLE SCALERS - 2/3 Springfield	SSP2/3	P-5
2653	SICKLE SCALERS - 135	SN135	P-6
2654	SICKLE SCALERS - REM 3/4 Remington	SRM3/4	P-6
2655	SICKLE SCALERS - 129	S129	P-6
2656	SICKLE SCALERS - T2/T3 Taylor	ST2/3	P-6
2657	SICKLE SCALERS - 128 Nebraska	SN128	P-6
2658	SICKLE SCALERS - USC-128	SUSC128	P-6
2659	SICKLE SCALERS - WH-2 Whiteside	SWH2	P-7
2660	SICKLE SCALERS - YG-15 Younger -Good	SYG15	P-7
2661	SICKLE SCALERS - 137	SN137	P-7
2662	SICKLE SCALERS - 137 Mini Five	SN13M	P-7
2663	SICKLE SCALERS - 0/00 Morse	SM0/00	P-7
2664	SICKLE SCALERS - OD1/33 Odont Jacquette	SOD1/33	P-7
2694	SICKLE SCALERS - 152 Cumine	PSSCM152	P-8
2695	SICKLE SCALERS - CI -2/3	PSSCI-2/3	P-8
2696	SICKLE SCALERS - CI -1	PSSCI-1	P-8
2697	Ball SICKLE SCALERS - B-1	PSSB-1	P-8
2698	Ball SICKLE SCALERS - B-4	PSSB	P-8
2699	Ball SICKLE SCALERS - B-5/6	PSSB-5/6	P-8
2700	Ball SICKLE SCALERS - B-2/3	PSSB-2/3	P-8

2701	SICKLE SCALERS - CK-6	PSSCK-6	P-8
2702	SICKLE SCALERS - CK-6D	PSSCK-6D	P-8
2705	SICKLE SCALERS - 11/12 Scale	PSS11/12	P-9
2706	SICKLE SCALERS - 11A/12A Scaler	PSS11/12A	P-9
2707	SICKLE SCALERS - 11/12S Modified McCall	PSS11/12S	P-9
2708	SICKLE SCALERS - 13S/14S Pointed Curette	PSC13/14S	P-9
2709	SICKLE SCALERS - MTC 13/14 Pointed Curette	PSCTC13/14	P-9
2710	SICKLE SCALERS - 17S/18S Curette	PSC17/18	P-9
2711	SICKLE SCALERS - 19/20 Curette	PSC19/20	P-9
2760	SICKLE SCALERS Goldman-Fox - GF-1 Universal Scaler	PSGF1	P-10
2761	SICKLE SCALERS Goldman-Fox - GF-5 Hoe Buc cal / Lingual	PSGF5	P-10
2762	SICKLE SCALERS Goldman-Fox - GF-2 Curettre	PSGF2C	P-10
2763	SICKLE SCALERS Goldman-Fox - GF-6 Hoe	PSGF6	P-10
2764	SICKLE SCALERS Goldman-Fox - GF-3 Curette	PSGF3C	P-10
2765	SICKLE SCALERS Goldman-Fox - GF-21	PSGF21	P-10
2766	SICKLE SCALERS Goldman-Fox - GF-4 Curette	PSGF4C	P-11
2767	SICKLE SCALERS Goldman-Fox - GF-21S/33 Comination	PSGF21S/33C	P-11
2768	SICKLE SCALERS Cattoni Scalers - 103/306	PS103/106	P-11
2769	SICKLE SCALERS Cattoni Scalers - 107/108	PS107/108	P-11
2770	SICKLE SCALERS MJ Scaler (Montana Jack Scaler)	PS104/105HC	P-11
2771	SICKLE SCALERS Cattoni Scalers - 107/108S Modified	PS107/108SM	P-11
2820	UNIVERSAL CURETTE - 1/2 Langer ( shown)	PUCL1/2	P-12
2821	UNIVERSAL CURETTE - 1/2 Langer Rigid	PUCL1/2R	P-12
2822	UNIVERSAL CURETTE - 1/2 Langer A-5	PUCL1/2A5	P-12
2823	UNIVERSAL CURETTE - 1/2 Langer M-5	PUCSL1/2M5	P-12
2824	UNIVERSAL CURETTE - 3/4 Langer ( shown)	PUCL3/4	P-12
2825	UNIVERSAL CURETTE - 3/4 Langer Rigid	PUCL3/4R	P-12
2826	UNIVERSAL CURETTE - 3/4 Langer A-5	PUCL3/4A5	P-12
2827	UNIVERSAL CURETTE - 3/4 Langer M-5	PUCL3/4M5	P-12
2828	UNIVERSAL CURETTE - 5/6 Langer ( shown)	PUCL5/6	P-12
2829	UNIVERSAL CURETTE - 5/6 Langer Rigid	PUCL5/6R	P-12
2830	UNIVERSAL CURETTE - 5/6 Langer A-5	PUCL5/6A5	P-12
2831	UNIVERSAL CURETTE - 5/6 Langer M-5	PUCL5/6M5	P-12
2832	UNIVERSAL CURETTE - 17/18 Langer ( shown)	PUCL17/18	P-12
2833	UNIVERSAL CURETTE - 17/18 Langer Rigid	PUCL17/18R	P-12
2834	UNIVERSAL CURETTE - IU13/14 Indiana University	PUCIU13/14	P-13
2835	UNIVERSAL CURETTE - IU17/18 Indiana University	PUCIU17/18	P-13
2836	UNIVERSAL CURETTE - 1/2 Bunting	PCCCB1/2	P-13
2837	UNIVERSAL CURETTE - 5/6 Bunting	PUCBN5/6	P-13
2838	UNIVERSAL CURETTE - YG7/8 Younger -Good	PUCYG7-8	P-13
2839	UNIVERSAL CURETTE - 1/2 Barnhart (Shown)	PUCBH1/2	P-14
2840	UNIVERSAL CURETTE - 1/2 Barnhart Rigid	PUCBH1/2R	P-14
2841	UNIVERSAL CURETTE - 5/6 Barnhar t (Shown)	PUCBH5/6	P-14
2842	UNIVERSAL CURETTE - 5/6 Langer Mini Five	PUCBH5/6R	P-14
2843	UNIVERSAL CURETTE - 2R/2L Columbia	PUCC2R/2L	P-14
2844	UNIVERSAL CURETTE - 4R/4L Columbia	PUCC4R/4L	P-14
2845	UNIVERSAL CURETTE - 13/14 Columbia (Shown)	PUCC13/14	P-14

2846	UNIVERSAL CURETTE - 13/14 Columbia Rigid	PUCC13/14R	P-14
2847	UNIVERSAL CURETTE - 13/14 Columbia Curved Tip	PUCC13/14C	P-14
2848	UNIVERSAL CURETTE - 13S/14S Bates	PUCBT13/14S	P-15
2849	UNIVERSAL CURETTE - HU1/2 Hut chinson	PUCHU1/2	P-15
2850	UNIVERSAL CURETTE - 3/4 Rule [23-24]	PUCR3/4	P-15
2851	UNIVERSAL CURETTE - UC3/4 Rat cliff	PUCUCR3/4	P-15
2852	UNIVERSAL CURETTE - 411/412	PUC411/412	P-15
2853	UNIVERSAL CURETTE - 10/11 Loma Linda (Shown)	PUCLL10/11	P-16
2854	UNIVERSAL CURETTE - 10/11 Loma Linda Rigid	PUCLL10/11R	P-16
2855	UNIVERSAL CURETTE - 10/11 Loma Linda Slim	PUCLL10/11S	P-16
2856	UNIVERSAL CURETTE - 413/414	PUC413/414	P-16
2857	UNIVERSAL CURETTE - 1/2 Mallery	PUCMAL1/2	P-16
2858	UNIVERSAL CURETTE - 31/32 Goldman	PUCG31/32	P-16
2859	UNIVERSAL CURETTE - ME-1 [10-11]	PUCME1	P-16
2910	GRACEY CURETTE - 1/2	PCG1/2	P-17
2911	GRACEY CURETTE - 1/2 Rigid	PCG1/2R	P-17
2912	GRACEY CURETTE - 3/4	PCG3/4	P-17
2913	GRACEY CURETTE - 3/4 Rigid	PCG3/4R	P-17
2914	GRACEY CURETTE - 5/6 (Shown)	PCG5/6	P-17
2915	GRACEY CURETTE - 5/6 Rigid	PCG56R	P-17
2916	GRACEY CURETTE - 5/6 Extra Rigid	PCG5/6ER	P-17
2917	GRACEY CURETTE - 7/8 (Shown)	PCG7/8	P-17
2918	GRACEY CURETTE - 7/8 Rigid	PCG7/8R	P-17
2919	GRACEY CURETTE - 7/8 Extra Rigid	PCG7/8ER	P-17
2920	GRACEY CURETTE - 9/10 (Shown)	PCG9/10	P-17
2921	GRACEY CURETTE - 9/10 Rigid	PCG9/10R	P-17
2922	GRACEY CURETTE - 11/12 (Shown)	PCG11/12	P-18
2923	GRACEY CURETTE - 11/12 Rigid	PCG11/12R	P-18
2924	GRACEY CURETTE - 11/12 Extra Rigid	PCG11/12ER	P-18
2925	GRACEY CURETTE - 13/14 (Shown)	PCG13/14	P-18
2926	GRACEY CURETTE - 13/14 Rigid	PCG13/14R	P-18
2927	GRACEY CURETTE - 13/14 Extra Rigid	PCG13/14ER	P-18
2928	GRACEY CURETTE - 15/16 (Shown)	PCG15/16	P-18
2929	GRACEY CURETTE - 15/16 Rigid	PCG15/16R	P-18
2930	GRACEY CURETTE - 17/18 (Shown)	PCG17/18	P-18
2931	GRACEY CURETTE - 17/18 Rigid	PCG17/18R	P-18
2932	GRACEY CURETTE - 11/14 Mesial Distal (Shown)	PGMDC11/14	P-19
2933	GRACEY CURETTE - 11/14 Mesial Distal Rigid	PGMDC11/14R	P-19
2934	GRACEY CURETTE - 12/13 Mesial Distal (Shown)	PGMDC12/13	P-19
2935	GRACEY CURETTE - 12/13 Mesial Distal Rigid	PGMDC12/13R	P-19
2936	GRACEY CURETTE - P 1/2	PGPC1/2	P-19
2937	GRACEY CURETTE - P 3/4	PGPC5/6	P-19
2938	GRACEY CURETTE - P 5/6	PGPC3/4	P-19
2939	GRACEY CURETTE - GRACEY / CURETTES - P 7/8	PGPC7/8	P-19
2940	GRACEY CURETTE - A-5 GC1/2	A5G1/2	P-20
2941	GRACEY CURETTE - A-5 GC 1/2, Rigid	A5G1/2R	P-20
2942	GRACEY CURETTE - A-5 GC 3/4	A5G3/4	P-20

2943	GRACEY CURETTE - A-5 GC 3/4, Rigid	A5G3/4R	P-20
2944	GRACEY CURETTE - A-5 GC 5/6	A5G5/6	P-20
2945	GRACEY CURETTE - A-5 GC 5/6, Rigid	A5G5/6R	P-20
2946	GRACEY CURETTE - A-F GC 5/7, Extra Rigid	A5G5/6ER	P-20
2947	GRACEY CURETTE - A-5 GC 7/8	A5G7/8	P-20
2948	GRACEY CURETTE - A-F GC 7/8, Rigid	A5G7/8R	P-20
2949	GRACEY CURETTE - A-5 GC 7/8, Extra Rigid	A5G7/8ER	P-20
2950	GRACEY CURETTE - A-5 GC 11/12	A5G11/12	P-20
2951	GRACEY CURETTE - A-5 GC 11/2, Rigid	A5G11/12R	P-20
2952	GRACEY CURETTE - A-5 GC 11/12, Extra Rigid	A5G11/12ER	P-20
2953	GRACEY CURETTE - A-5 GC 13/14	A5G13/14	P-20
2954	GRACEY CURETTE - A-5 GC 13/14, Rigid	A5G13/14R	P-20
2955	GRACEY CURETTE - A-5 GC 13/14, Extra Rigid	ASG13/14ER	P-20
2956	GRACEY CURETTE - A-5 GC 15/16	A5G15/16	P-20
2957	GRACEY CURETTE - A-5 GC 15/16, Rigid	A5G15/16R	P-20
2958	GRACEY CURETTE - A-5 GC 11/14 Mesial /Distal	A5G11/14	P-20
2959	GRACEY CURETTE - A-5 GC 11/14 Mesial /Distal, Rigid	A5G11/14R	P-20
2960	GRACEY CURETTE - A-5 GC 12/13 Mesial /Distal	A5G12/13	P-20
2961	GRACEY CURETTE - A-5 GC 12/13 Mesial /Distal, Rigid	A5G12/13R	P-20
2962	GRACEY CURETTE - M-5 GC 1/2	M5G1/2	P-21
2963	GRACEY CURETTE - M-5 GC 1/2, Rigid	M5G1/2R	P-21
2964	GRACEY CURETTE -M-5 GC 3/4	M5G3/4	P-21
2965	GRACEY CURETTE - M-5 GC 3/4, Rigid	M5G3/4R	P-21
2966	GRACEY CURETTE - M-5 GC 5/6	M5G5/6	P-21
2967	GRACEY CURETTE - M-5 GC 5/6, Rigid	M5G5/GR	P-21
2968	GRACEY CURETTE - M-5 GC 7/8	M5G7/8	P-21
2969	GRACEY CURETTE -M-5 GC 7/8, Rigid	M5G7/8R	P-21
2970	GRACEY CURETTES -M-5 GC 11/12	M5G11/12	P-21
2971	GRACEY CURETTE - M-5 GC 11/12, Rigid	M5G11/12R	P-21
2972	GRACEY CURETTE - M-5 GC 13/14	M5G13/14	P-21
2973	GRACEY CURETTE - M-5 GC 13/14, Rigid	M5G13/14R	P-21
2974	GRACEY CURETTE - M-5 GC 15/16	M5G15/16	P-21
2975	GRACEY CURETTE - M-5 GC 15/16, Rigid	M5G15/16R	P-21
2976	GRACEY CURETTE - M-5 GC 11/14 Mesial /Distal	M5G11/14	P-21
2977	GRACEY CURETTE - M-5 12/13 Mesial /Distal	M5G12/13	P-21
3030	FILE SCALERS - 3/7 Orban Buccal/Lingual	PFS-OBL3/7	P-22
3031	FILE SCALERS - 6/7 Orban Buccal/Lingual	PSF-OBL6/7	P-22
3032	FILE SCALERS - 5/11 Hirschteld Mesial/Distal	PFS-OMD12/13	P-22
3033	FILE SCALERS - 12/13 Orban Mesial/Distal	PFS-HMD5/11	P-22
3034	FILE SCALERS - 9/10 Hirschteld Anterior	PFS-HA9/10	P-22
3035	HOE SCALERS - XS-15 Towner	PHSXS15	P-23
3036	HOE SCALERS - 7K	PHS7K	P-23
3037	HOE SCALERS - 6/7 Orban Mesial /Distal	PHSOS6/7	P-23
3038	HOE SCALERS - 8/9 Orban Buccal /Lingual	PHSO8/9	P-23
3039	HOE SCALERS - U 15/13R	PHSU15/13R	P-23
3080	PERIODONTAL KNIVES - 15K/16K Kirkland	PK15K/16K	P-24
3081	PERIODONTAL KNIVES - 7 Gold-Fox	PKGF7	P-24

3082	PERIODONTAL KNIVES - 8 Gold-Fox	PKGF8	P-24
3083	PERIODONTAL KNIVES - 9 Gold-Fox	PKGF9	P-24
3084	PERIODONTAL KNIVES - 11 Goldman-Fox	PKGF11	P-24
3085	PERIODONTAL KNIVES - 3/4 Buck	PKB3/4	P-24
3086	PERIODONTAL KNIVES - 5/6 Buck	PKB5/6	P-24
3087	PERIODONTAL KNIVES - 1/2 Solt	PKSOLT1/2	P-25
3088	PERIODONTAL KNIVES - 1/2 Sanders	PKS1/2	P-25
3089	PERIODONTAL KNIVES - 1/2 Orban	PKO1/2	P-25
3090	PERIODONTAL KNIVES - USC Towner 19/20	PKUSCT19/20	P-25
3091	PERIODONTAL KNIVES - 1/2 Orban Knife	PKOK1/2	P-25
3092	PERIODONTAL KNIVES - CK-3 Crane Kaplan	PKCK3	P-25
3093	PERIODONTAL KNIVES - Dr Allen End-Cutting Intralucular Knife	PKDA	P-25
3140	PSC - PR1/2 Prichard	PSCPR1/2	P-26
3141	PSC - KRA-1 Kramer Mesial	PSCKRA1	P-26
3142	PSC - KRA-2 Kramer Distal	PSCKRA2	P-26
3143	PSC - KRA-3 Kramer	PSCKRA3	P-26
3144	PSC - 3/4 Kramer Nevins	PSCKN3/4	P-26
3145	PSC - 8K/9K Kirkland	PSCK8/9K	P-26
3146	PSC - 5S/6S Sugarman Curette	PSCS5/6S	P-26
3147	PSC - 9S/10S Sugarman Curette	PSCS9/10S	P-26
3148	PSC - 11S/12S Sugarman Curette	PSCS11/12S	P-26
3200	PSC - 9/10 Schluger Curved File	PFSC9/10	P-27
3201	PSC - 11/12 Buck Straight File	PFBS11/12	P-27
3202	PSC - 1S/2S Sugarman File Mesial /Distal	PFSMD1S/2S	P-27
3203	PERIODONTAL FILE & CHISELS - 3S/4S Sugarman File Buccal /Lingual	PFSBL3S/4S	P-27
3204	PERIODONTAL FILE & CHISELS - 13K/TG	PC13K/TG	P-27
3205	PERIODONTAL FILE & CHISELS - 13K/13KL	PC13K/13KL	P-27
3206	PERIODONTAL CHISELS - 1/2 Kramer Nevins	PCKN1/2	P-28
3207	PERIODONTAL CHISELS - 36/37 Rhodes Rack -Action	PC36/37	P-28
3208	PERIODONTAL CHISELS - 44/45 Remington Hoe	PC44/45	P-28
3209	PERIODONTAL CHISELS - TG	PCTG	P-28
3210	PERIODONTAL CHISELS - TG-0 Modified TG with Ochsenbein Cutouts	PCTG-0	P-28
3211	PERIODONTAL CHISELS - 3/4 Solt	PCSOLT3/4	P-28
3214	PERIODONTAL CHISELS - 3 Ochsenbein	PCO3	P-29
3215	PERIODONTAL CHISELS - 4 Ochsenbein	PCO4	P-29
3216	PERIODONTAL CHISELS - 1 Fedi	PCF1	P-29
3217	PERIODONTAL CHISELS - 2 Fedi	PCF2	P-29
3218	PERIODONTAL CHISELS - 3 Fedi	PCF3	P-29
3219	PERIODONTAL CHISELS - 4 Fedi	PCF4	P-29
3221	3 Kramer-Nevins, 7.5mm (Periodontal Chisel)	PKN3	P-30
3222	55 Kramer-Nevins, 5.5mm (Periodontal Chisel)	PKN55	P-30
3223	75 Kramer-Nevins, short edge 7.5mm (Periodontal Chisel)	PKN75	P-30
3224	4 Buser, Round Tip-4mm (Periodontal Chisel)	PBS4	P-30
3225	6 Buser, Pointed Tip-6mm (Periodontal Chisel)	PBS6	P-30
3226	1 Ochsenbein (Periodontal Chisel)	PCO1	P-30
3227	2 Ochsenbein (Periodontal Chisel)	PCO2	P-30
3230	Goldstein Flexible, Thin Composite Instruments - #1	RCIGFT1	R-1

3231	Goldstein Flexible, Thin Composite Instrument - #2	RCIGFT2	R-1
3232	Goldstein Flexible, Thin Composite Instrument - #3 Extra Flexible	RCIGFT3	R-1
3233	Goldstein Flexible, Thin Composite Instrument - #4 Extra Flexible	RCIGFT4	R-1
3234	Goldstein Flexible, Thin Composite Instrument - #5 Flexi -Thin	RCIGFT5	R-1
3235	Goldstein Flexible, Thin Composite Instrument - #6 Flexi -Thin	RCIGFT6	R-1
3236	Goldstein Flexible, Thin Composite Instrument - Micro Mini	RCIGFTMM	R-1
3237	Goldstein Flexible, Thin Composite Instrument - Mini 1	RCIGFTM1	R-1
3238	Goldstein Flexible, Thin Composite Instrument - Mini 3	RCIGFTM3	R-1
3239	Goldstein Flexible, Thin Composite Instrument -Mini 4	RCIGFTM4	R-1
3245	Plastic Filling Instruments - AB1 Boghosian	RPFIAB1	R-2
3246	Plastic Filling Instruments - AB2 Boghosian	RPFIAB2	R-2
3247	Plastic Filling Instruments - W3	RPFIW3	R-2
3248	Composite Instruments - 4F Tufts	RCI4FT	R-2
3249	Plastic Filling Instruments - UF	RPFIUF	R-2
3250	Composite Instruments - Interproximal Carver, Long	RCVIPC	R-2
3251	Composite Instruments - 3 Tufts	RCI3T	R-2
3252	Plastic Filling Instruments - A6 (156)	RPFIA6	R-2
3253	Plastic Filling Instruments - 4/5 Gregg	RPFI4/5G	R-2
3254	Plastic Filling Instruments - PF1	RPFIPF1	R-2
3255	Goldfogel Cosmatic Contouring Instruments - A Cosmatic Contouring	RGFCCA	R-3
3256	Goldfogel Cosmatic Contouring Instruments - B Cosmatic Contouring	RGFCCB	R-3
3257	Goldfogel Cosmatic Contouring Instruments - C Cosmatic Contouring	RGFCCB	R-3
3258	Goldfogel Cosmatic Contouring Instruments - D Cosmatic Contouring	RGFCCD	R-3
3259	Goldfogel Cosmatic Contouring Instruments - E Cosmatic Contouring	RGFCCE	R-3
3260	Goldfogel Cosmatic Contouring Instruments - F Cosmatic Contouring	RGFCCF	R-3
3261	Goldfogel Cosmatic Contouring Instruments - G Cosmatic Contouring	RGFCCG	R-3
3262	Goldfogel Cosmatic Contouring Instruments - H Cosmatic Contouring	RGFCCH	R-3
3263	Goldfogel Cosmatic Contouring Instruments - I Cosmatic Contouring	RGFCCI	R-3
3264	Plastic Filling Instruments - 8A	RPFI8A	R-4
3265	Plastic Filling Instruments - 11	RPFI11	R-4
3266	FILLING INSTRUMENTS - 179	RPFI179	R-4
3267	FILLING INSTRUMENTS - G	RPFIG	R-4
3268	FILLING INSTRUMENTS - 6SH	RPFI6SH	R-4
3269	Plastic Filling Instruments - 1 Woodson	RPFI1W	R-4
3270	FILLING INSTRUMENTS - PF#1 Heidmann Spatula, Small	RPFI1HSS	R-4
3271	Plastic Filling Instruments - PF#2 Heidmann Spatula, Medium	RPFI2HSM	R-4
3272	Plastic Filling Instruments - PF#3 Heidmann Spatula, Large	RPFI3HSL	R-4
3273	Plastic Filling Instruments - 48	RPFI48	R-4
3274	Plastic Filling Instruments - W1	RPFIW1	R-4
3275	Plastic Filling Instruments - 118	RPFI118	R-4
3276	Plastic Filling Instruments - Novadent 2	RPFIND2	R-5
3277	Plastic Filling Instruments - Interproximal Off-Angle	RPFIIOA	R-5
3278	Plastic Filling Instruments - 49 Baldwin	RPFI49	R-5
3279	Plastic Filling Instruments - UP1	RPFIUP1	R-5
3280	Plastic Filling Instruments - 21 Burnisher	RPFI21	R-5
3281	Plastic Filling Instruments - 21B Burnisher	RPFI21B	R-5
3282	Plastic Filling Instruments - PKT Round Cone	RPFIPKT	R-5

3285	Plastic Filling Instruments - BEQ1 Quet in	RPFIBEQ1	R-6
3286	Plastic Filling Instruments - 1 Loesche	RPF11L	R-6
3287	Plastic Filling Instruments - 2 Loesche	RPF12L	R-6
3288	Plastic Filling Instruments - 3T	RPF13T	R-6
3289	Plastic Filling Instruments - 1RPF1		
3290	Plastic Filling Instruments - 1 Woodson	RPF12W	R-6
3291	Plastic Filling Instruments - 2 Woodson	RPF13W	R-6
3292	Plastic Filling Instruments - 3 Woodson	RPF1	R-6
3293	Plastic Filling Instruments - S4 Woodson		
3400	Hatchet 8/9H [10-7-14]	RH8/9H	R-7
3401	Hatchet 14/14 Off-Angle [15-8-14]	RH14/14OA	R-7
3402	Hatchet 15/15 Off-Angle [15-8-14]	RH15/15OA	R-7
3403	Hatchet 14/14-0 Off-Angle [15-10-16]	RH14/14-0OA	R-7
3404	Hatchet 15/15-0 Off-Angle [15-10-16]	RH15/15-0OA	R-7
3405	Hatchet 44S Off Angle	RH44SOA	R-7
3406	Hatchet 45S Off Angle	RH45SOA	R-7
3407	Hatchet 13/14 [20-9-14]	RH13/14	R-7
3408	Hatchet 15/16 [15-8-14]	RH15/16	R-7
3409	Hatchet 17/18 [10-6-14]	RH17/18	R-7
3410	Hatchet 51/52 [15-8-12]	RH51/52	R-7
3411	Hatchet 53/54 [10-6-12]	RH53/54	R-7
3412	Hatchet 11/12 [15-8-8]	RH11/12	R-7
3440	RESTORATIVE CHISEL - ½ Wedelstaedt [20-15-3]	RC1/2WD	R-8
3441	RESTORATIVE CHISEL - ¾ Wedelstaedt [11.5-15-3]	RC3/4WD	R-8
3442	RESTORATIVE CHISEL - 5/6 Wedelstaedt [15-15-3]	RC5/6WD	R-8
3443	RESTORATIVE CHISEL - 7/10 St raight [15] [20]	RCT7/10S	R-8
3444	RESTORATIVE CHISEL - 8/9 Binangle [20-9-8]	RCT8/9B	R-8
3445	RESTORATIVE CHISEL - 11/12 Binangle [15-8-8]	RCT11/12	R-8
3446	RESTORATIVE CHISEL - 40/41 Binangle [18-10-16]	RCT40/41	R-8
3447	RESTORATIVE HOES - 20 [14-6-8]	RH20	R-9
3448	RESTORATIVE HOES - 21 [10-4-8]	RH21	R-9
3449	RESTORATIVE HOES - 22 [10-4-14]	RH22	R-9
3450	RESTORATIVE HOES - 23 [6.5-2.5-9]	RH23	R-9
3460	Angle Former 30/31 [12-80-5.8]	RAF30/31	R-9
3462	Angle Former 32/33 [9-80-408]	RAF32/33	R-9
3464	Angle Former 34/35 [7-80-2.5-5.9]	RAF34/35	R-9
3500	Restorative Margin Trimmers - 26 [13-95-8-14] Distal	RMT26D	R-10
3501	Restorative Margin Trimmers - 27 [13-80-8-14] Mesial	RMT27M	R-10
3502	Restorative Margin Trimmers - 28 [10-95-7-14] Distal	RMT28D	R-10
3503	Restorative Margin Trimmers - 29 [10-80-7-14] Mesial	RMT29M	R-10
3504	Restorative Margin Trimmers - 77/78 {15-95-8-12} Distal	RMT77/78D	R-10
3505	Restorative Margin Trimmers - 79/80 {15-80-8-12}Mesial	RMT79/80M	R-10
3550	Excavator, Spoon Type- E1 [12-9-15]	REXCSTE1	R-11
3551	Excavator, Spoon Type E2 [15-9-15]	REXCSTE2	R-11
3552	Excavator, Spoon Type - E3 [25-9-15]	REXCSTE3	R-11
3553	Excavator, Spoon Type- 6	REXCST6	R-11
3554	Excavator, Spoon Type- 14	REXCST14	R-11

3555	Excavator, Spoon Type - 17	REXCST17	R-11
3556	Excavator, Spoon Type - 18	REXCST18	R-11
3557	Excavator, Spoon Type- 19	REXCST19	R-12
3558	Excavator, Spoon Type - 38/39 [11.5-7-14]	REXCST38/39	R-12
3559	Excavator, Oval Spoon Type- 17W	REXCOT17W	R-12
3560	Excavator, Oval Spoon Type- 18W	REXCOT18W	R-12
3561	Excavator, Oval Spoon Type - 19W	REXCOT19W	R-12
3562	Excavator, Oval Spoon Type- 220/221 Darby -Perry	REXCOT220/221	R-12
3563	Excavator, Blade Type - 36/37 [15-8-14]	REXCBT36/37	R-13
3564	Excavator, Blade Type - 40/41 [11-7-14]	REXCBT40/41	R-13
3565	Excavator, Blade Type - 44/45 [20-9-14]	REXCBT44/45	R-13
3566	Excavator, Blade Type - 61/62 [20-9-12]	REXCBT61/62	R-13
3567	Excavator, Blade Type - 63/64 [15-8-12]	REXCBT63/64	R-13
3568	Excavator, Blade Type - 65/66 [10-6-12]	REXCBT65/66	R-13
3569	Excavator, English Pattern - 125/126	REXCEP125/126	R-13
3570	Excavator, English Pattern - 127/128	REXCEP127/128	R-13
3571	Excavator, English Pattern - 129/130	REXCEP129/130	R-14
3572	Excavator, English Pattern - 131/132	REXCEP131/132	R-14
3573	Excavator, English Pattern - 153/154	REXC153/154	R-14
3574	Excavator, English Pattern - Goldstein Interporximal Excavator	REXCEPGIE	R-14
3575	Excavator, English Pattern - 243	REXCEP243	R-14
3576	Excavator, English Pattern - 244	REXCEP244	R-14
3577	Excavator, English Pattern - 245	REXCEP245	R-14
3578	Excavator, English Pattern - 246	REXCEP246	R-14
3630	Placement Instruments - Calcium Hydroxide Placer	RPICHP	R-15
3631	Placement Instruments - 6061 Mini Spatula/ Placer	RPIMSP6061	R-15
3632	Placement Instruments - 10 Novatech Placer	RPINTP10	R-15
3633	Placement Instruments - 11 Novatech Placer	RPINTP11	R-15
3634	Placement Instruments - 1 Composite Brush Handle	RPICBH1	R-15
3635	Placement Instruments - Miller Articulating Paper Foceps	RMPAF2	R-15
3636	Placement Instruments - Miller Articulating Paper Foceps	RMPAF4	R-15
3637	Miller Mod Articulating Paper Foceps	RMPAM	R-15
3665	Amalgam Carrier - Mini(1.5mm), Teflon Tip	ACMTT1.5	R-16
3666	Amalgam Carrier - Regular(2.0mm), Teflon Tip	ACRTT2.0	R-16
3667	Amalgam Carrier - Large(2.5mm), Teflon Tip	ACLTT2.5	R-16
3668	Amalgam Carrier - Jumbo(3.2mm), Teflon Tip	ACJTT3.2	R-16
3669	Amalgam Carrier - Mini(1.5mm), SS Tip	ACMSST1.5	R-16
3670	Amalgam Carrier - Regular(2.0mm), SS Tip	ACRSST2.0	R-16
3671	Amalgam Carrier - Large(2.5mm), SS Tip	ACLSST2.5	R-16
3672	Amalgam Carrier - Jumbo(3.2mm), SS Tip	ACJSST3.2	R-16
3673	Amalgam Carrier - Mini(1.5mm)/, Mini(1.5mm), Teflon Tip	ACMTT1.5/1.5	R-17
3674	Amalgam Carrier - Regular(2.0mm)/ Regular(2.0mm), Teflon Tip	ACRTT2.0/2.0	R-17
3675	Amalgam Carrier - Large(2.5mm)/ Large(2.5mm), Teflon Tip	ACLTT2.5/2.5	R-17
3676	Amalgam Carrier - Jumbo(3.2mm)/ Jumbo(3.2mm), Teflon Tip	ACJTT3.2/3.2	R-17
3677	Amalgam Carrier - Mini(1.5mm)/, Mini(1.5mm), SS Tip	ACMSST1.5/1.5	R-17
3678	Amalgam Carrier - Regular(2.0mm)/ Regular(2.0mm), SS Tip	ACRSST2.0/2.0	R-17
3679	Amalgam Carrier - Large(2.5mm)/ Large(2.5mm), SS Tip	ACLSST2.5/2.5	R-17

3680	Amalgam Carrier - Jumbo(3.2mm)/ Jumbo(3.2mm), SS Tip	ACJSST3.2/3.2	R-17
3681	Amalgam Carrier - Mini(1.5mm)/, Regular(2.0mm), Teflon Tip	ACM1.5R/2.0TT	R-18
3682	Amalgam Carrier - Large(2.5mm)/ Jumbo(3.2mm), Teflon Tip	ACL2.5/J3.2TT	R-18
3683	Amalgam Carrier - Mini(1.5mm)/, Regular(2.0mm), SS Tip	ACM1.5R2.0SST	R-18
3684	Amalgam Carrier - Large(2.5mm)/ Jumbo(3.2mm), SS Tip	ACL2.5J3.2SST	R-18
3685	Amalgam Carrier, (Chrome Plated) 2mm, SS Tip	ACCP2.0SST	R-18
3686	Amalgam Carrier, (Chrome Plated) 3.2mm, SS Tip	ACCP3.2SST	R-18
3687	Amalgam Carrier, (Chrome Plated) 2mm, Teflon Tip	ACCP2.0TT	R-18
3688	Amalgam Carrier, (Chrome Plated) 3.2mm, Teflon Tip	ACCP3.2TT	R-18
3695	Amalgam Well	AW	R-18
3750	Restorative Pluggers - ½ Andrews	RPLGA1/2	R-19
3751	Restorative Pluggers - ½ Black	RPLGB1/2	R-19
3752	Restorative Pluggers - ½ Black, Serrated	RPLGBS1/2	R-19
3753	Restorative Pluggers - H1 Hollenback	RPLGHBH1	R-19
3754	Restorative Pluggers - H2 Hollenback	RPLGHBH2	R-19
3755	Restorative Pluggers - H3 Hollenback	RPLGHBH3	R-19
3756	Restorative Pluggers - H4 Hollenback	RPLGHBH4	R-19
3757	Restorative Pluggers - 0/1 Marquette Condenser	RPLG0/1MC	R-19
3758	Restorative Pluggers - 0/1 Marquette	RPLG0/1M	R-19
3759	Restorative Pluggers - 2 Mortonson	RPLG2M	R-20
3760	Restorative Pluggers - 2 Smith Serrated	RPLG2SS	R-20
3761	Restorative Pluggers - SPO Serrated	RPLGSPOS	R-20
3762	Restorative Pluggers - 1M Markley	RPLG1MM	R-20
3763	Restorative Pluggers - 2M Markley	RPLG2MM	R-20
3764	Restorative Pluggers - 4M Markley	RPLG4MM	R-20
3765	Restorative Pluggers - 5M Markley	RPLG5MM	R-20
3766	Restorative Pluggers - 6M Markley	RPLG6MM	R-20
3767	Restorative Pluggers - 1T Tanner	RPLGT1TT	R-21
3768	Restorative Pluggers - 1 Oregon	RPLGOG1	R-21
3769	Restorative Pluggers - 2 Oregon	RPLGOG2	R-21
3770	Restorative Pluggers - 3 Oregon	RPLGOG3	R-21
3771	Restorative Pluggers - 4 Oregon	RPLGOG4	R-21
3772	Restorative Pluggers - 1 Ward	RPLGW1	R-21
3773	Restorative Pluggers - 2 Ward	RPLGW2	R-21
3774	Restorative Pluggers - 2T Tanner	RPLG2TT	R-22
3775	Restorative Pluggers - 3T Tanner	RPLG3TT	R-22
3776	Restorative Pluggers - 9/10 Back Action	RPLG9/10BA	R-22
3777	Restorative Pluggers - 11/12 Back Action	RPLG11/12BA	R-22
3778	Restorative Pluggers - 153 English Pattern Serrated	RPLG153EPS	R-22
3779	Restorative Pluggers - 154 English Pattern Serrated	RPLG154EPS	R-22
3830	CARVER - Inter - proximal	RCVIP	R-23
3831	CARVER - Inter - proximal Off -Angle	RCVIPOA	R-23
3832	CARVER - 1/2-3 Hollenbac k	RCVHB1/2-3	R-23
3833	CARVER - 3S Hollenback	RCVHB3S	R-23
3834	CARVERS - 3 Hollenback	RCVHB3	R-23
3835	CARVERS - 8 Wiland Minn. IPC1	RCVWIM8IPC	R-23
3836	CARVERS - 1 Andrew	RCVA1	R-24

3837	CARVERS - 2 Andrew	RCVA2	R-24
3838	CARVER - 3 Andrew	RCVA3	R-24
3839	CARVER- 4 Clappison	RCVCL4	R-24
3840	CARVER- 2/3 Frahm	RCVFR2/3	R-24
3841	CARVERS- 104	RCV104	R-24
3842	CARVER - 1 Hollenback 2 Modified Cleoid Small	RCVHB1/MCS2	R-24
3843	CARVER -1 Hollenback 3 Modified Cleoid Small	RCVHB1/MC3	R-24
3844	CARVERS - 1 Hollenback	RCVHB1	R-25
3845	CARVERS - 2 Hollenback	RCVHB2	R-25
3846	CARVERS - L7 Levy	RCVL7L	R-25
3847	CARVERS - SHO-A Shosnan	RCVSHOAS	R-25
3848	CARVERS - 18 Sprengel	RCV18S	R-25
3849	CARVERS - 1 Tharp	RCV1TH	R-25
3850	CARVER - 2 Tharp	RCV2TH	R-25
3851	Cleoids-Discoid(CD) - 3 Tanner	RCV3TCD	R-26
3852	Cleoids-Discoid(CD) - 4 Tanner	RCV4TCD	R-26
3853	Cleoids-Discoid(CD) - 5 Tanner	RCV5TCD	R-26
3854	Cleoids-Discoid(CD) - 6 Tanner	RCV6TCD	R-26
3855	Cleoids-Discoid(CD) - 1 Wall	RCVWL1CD	R-26
3856	Cleoids-Discoid(CD) - 3 Wall	RCVWL3CD	R-26
3857	Cleoids-Discoid(CD) - 1 Ward	RCVWD1CD	R-26
3858	Cleoids-Discoid(CD) - 2 Ward	RCVWD2CD	R-26
3859	Cleoids-Discoid(CD) - ½	RCVCD1/2	R-27
3860	Cleoids-Discoid(CD) - 4/5	RCVCD4/5	R-27
3861	Cleoids-Discoid(CD) - 3/6	RCVCD3/6	R-27
3862	Cleoids-Discoid(CD) - 89/92	RCVCD89/92	R-27
3863	Cleoids-Discoid(CD) - 90B	RCVCD90B	R-27
3864	Cleoids-Discoid(CD) - W-A-C-D	RCVCDWACD	R-27
3910	BURNISHERS-18	RB18	R-28
3911	BURNISHERS-26/27S	RB26/27S	R-28
3912	BURNISHERS-27/29	RB27/29	R-28
3913	BURNISHERS-2	RB2SE	R-28
3914	BURNISHERS-2 Double-End	RB2DE	R-28
3915	BURNISHERS-2/29	RB2/29	R-28
3916	Burnishers-21	RB21	R-29
3917	Burnishers-21B	RB21B	R-29
3918	Burnishers-2 Ladmore	RB2LM	R-29
3919	Burnishers-PKT-3R (Round Cone)	BPKT3R	R-29
3920	BURNISHERS-Romerowski	RBROM	R-29
3921	BURNISHERS-Ladmore	RBLM	R-29
3922	Burnisher-29/34	RB29/34	R-29
3970	GINGIVAL RETRACTORS (Aesthetics) - Kincheloe	RGRK	R-30
3971	GINGIVAL RETRACTORS (Aesthetics) - GP 10 Goldman-Fox	RGRGF10	R-30
3972	GINGIVAL RETRACTORS (Aesthetics) - 1 Meinershagen	RGRM1	R-30
3973	GINGIVAL RETRACTORS (Aesthetics) - 2 Meinershagen	RGRM2	R-30
3974	GINGIVAL RETRACTORS (Aesthetics) - 3 Meinershagen	RGRM3	R-30
3975	GINGIVAL RETRACTORS (Aesthetics) - 4 Meinershagen	RGRM4	R-30

4030	GINGIVAL CORD PACKERS (Aesthetics) - CSI-1 Serrated	RGCPCSI1S	R-31
4031	GINGIVAL CORD PACKERS (Aesthetics) - 113 Serrated	RGCP113S	R-31
4032	GINGIVAL CORD PACKERS (Aesthetics) - BN1	RGCPBN1	R-31
4033	GINGIVAL CORD PACKERS (Aesthetics) - S6	RGCP6	R-31
4034	GINGIVAL CORD PACKERS (Aesthetics) - Balshi	RGCPBAL	R-31
4035	GINGIVAL CORD PACKERS (Aesthetics) - 1 Yardley	RGCP1YAL	R-31
4036	GINGIVAL CORD PACKERS (Aesthetics) - 7 Guyer Serrated	RGCPG7S	R-31
4040	Crown Spreader Forceps	RCSF	R-32
4042	Dental Crown Gripper, TC	RDCGTC	R-32
4044	Morrel Crown Remover With 3 Tips and Key	RMCR	R-32
4046	Automatic Gun Type Crown Remover With 3 Tips and Key	RAGTCR	R-32
4048	Dental Crown Splitter	RDCS	R-32
4050	Christensen Crown Remover 90°	RCCR90	R-33
4052	Christensen Crown Remover Straight	RCCRS	R-33
4054	Crown Spreader	RCR	R-33
4056	Goldstein Crown Remover, Straight	RGSCRS	R-33
4058	Goldstein Crown Remove 45°	RGSCR45	R-33
4060	Goldstein Crown Remove 90°	RGSCR90	R-33
4062	Goldstein Crown Remove Occlusal Surface	GSCROCS	R-33
4090	SPATULAS - 24 Flexible 1 3/4 44mm	RSP24F	R-32
4091	SPATULAS - 324 Rigid 2" 51mm	RSP324R	R-32
4092	SPATULAS - A6 Rigid 1" 25mm	RSPA6R	R-32
4093	SPATULAS - 5 Novatech Long/Fluted	RSPNT5LF	R-32
4094	SPATULAS - 6 Novatech Long/Blade	RSPNT6LB	R-32
4095	SPATULAS - 7 Novatech Short /Blade	RSPNT7SB	R-32
4096	SPATULAS - 8 Novatech, Long	RSPNT8L	R-32
4097	SPATULAS - 9 Novatech, Short	RSPNT9S	R-32
4100	Impression Trays - Set of 10 EA (Upper & Lower)	RITS10	R-33
4101	Impression Trays - Set of 8 EA (Upper & Lower)	RITS8	R-33
4102	Impression Trays - Set of 6 EA (Upper & Lower)	RITS6	R-33
4103	Impression Trays, Perforated - Set of 10 EA (Upper & Lower)	RITPS10	R-33
4104	Impression Trays, Perforated - Set of 8 EA (Upper & Lower)	RITPS8	R-33
4105	Impressio Trays, Perforated - Set of 6 EA (Upper & Lower)	RITPS6	R-33
4110	Impression Trays, Perforated - Set of 8 EA (Upper & Lower)	RITS8	R-33
4112	Impression Trays, Perforated - Set of 8 EA (Upper & Lower)	RITPS8	R-33
4114	IMPRESSION TRAY - Upper 65x39 mm (each)	RITU	R-33
4116	IMPRESSION TRAY - Lower 64x39 mm (each)	RITL	R-33
4118	IMPRESSION TRAY - Upper 65x39 m (each)	RITPU	R-33
4120	IMPRESSION TRAY - Lower 64x39 mm (each)	RITPL	R-33
4122	IMPRESSION TRAY - Upper Right Lower Left (each)	RITPUR	R-33
4123	IMPRESSION TRAY - Upper Left Lower Right (each)	RITPLR	R-33
4124	IMPRESSION TRAY - Adjustable Perforated (each)	RITPA	R-33
4125	IMPRESSION TRAY - Adjustable Solid (each)	RITSA	R-33
4130	IMPRESSION TRAY - Bite Checa Tray (each)	RITBCT	R-33
4180	Spade Elevator, Straight	SELSS	S-1
4181	Spade Elevator, Mesially Angled	SELSMA	S-1
4182	Spade Elevator, Distally Angled	SELSDA	S-1

4183	Spade Elevator, Left	SESL	S-1
4184	Spade Elevator, Right	SELSR	S-1
4185	2.5mm, Small, Straight	SELSS2.5	S-1
4186	2.5mm, Small, Curved	SELSC2.5	S-1
4187	2.5mm, Small, Mesially Angled	SEL SMA2.5	S-1
4188	2.5mm, Small, Distally Angled	SELSDA2.5	S-1
4190	3.75mm, Large, Straight	SELLS3.75	S-2
4191	3.75mm, Large, Curved	SELLC3.75	S-2
4192	3.75mm, Large, Mesially Angled	SELLMA3.75	S-2
4193	3.75mm, Large, Distally Angled	SELLDA3.75	S-2
4194	4.5mm, Large, Straight	SELLS4.5	S-2
4195	4.5mm, Large, Curved	SELLC4.5	S-2
4196	4.5mm, Large, Mesially Angled	SELLMA4.5	S-2
4197	4.5mm, Large, Distally Angled	SELLDA4.5	S-2
4198	Small, Tight, Mesial Angle, 2.5mm	SELSTMA2.5	S-2
4199	Tight Mesial Spade, 2.5mm	SELTMS2.5	S-2
4204	Lindo Levin, Straight 3mm	LLS3	S-3
4205	Lindo Levin, Straight 4mm	LLS4	S-3
4206	Lindo Levin, Curved 3mm	LLC3	S-3
4207	Syndesmotomes Bernard Spade 11, Straight	BS11	S-3
4208	Syndesmotomes Bernard Spade 6, Left	BR6	S-3
4209	Syndesmotomes Bernard Spade 7, Right	BL7	S-3
4210	Syndesmotomes Bernard Spade, Offset IN	BSOIN	S-3
4211	Syndesmotomes Bernard Spade, Offset OUT	BSOOT	S-3
4213	ELEVATORS - 1 Heidbrik Straight	E1	S-4
4214	ELEVATORS - 2 Heidbrik Straight	E2	S-4
4215	ELEVATORS - 40 Heidbrik Straight	E40	S-4
4216	ELEVATORS - 41 Heidbrik Straight	E41	S-4
4217	Elevator 46	E46	S-4
4218	Elevator 46R	E46R	S-4
4219	Elevator 77	E77	S-4
4220	Elevator 77R	E77R	S-4
4221	Elevator 21 Cryer Mini	EM21	S-5
4222	Elevator 22 Cryer, Mini	EM22	S-5
4223	Elevator, 44 Cryer, Small	E44	S-5
4224	Elevator, 45 Cryer, Small	E45	S-5
4225	Elevator, 25 Cryer, Medium	E25	S-5
4226	Elevator, 26 Cryer, medium	E26	S-5
4227	Elevator, A30 Cryer	EA30	S-5
4228	Elevator, A31 Cryer	EA31	S-5
4229	Elevator, 31 Cryer, Large	EA31	S-5
4230	Elevator, 32 Cryer, Large	EA32	S-5
4231	ELEVATORS - Luxating 2mm, Straight	SLES2	S-6
4232	ELEVATORS - Luxating 2mm, Curved	SLEC2	S-6
4233	ELEVATORS - Luxator - Straight 2.5mm	SLES2.5	S-6
4234	ELEVATORS - Luxator - Curved 2.5mm	SLEC2.5	S-6
4235	ELEVATORS - Luxator - Straight 3mm	SLES3	S-6

4236	ELEVATORS - Luxator - Curved 3mm	SLEC3	S-6
4237	ELEVATORS - Luxator - Straight 4mm	SLES4	S-6
4238	ELEVATORS - Luxator - Curved 4mm	SLEC4	S-6
4239	ELEVATORS - Luxator - Straight 5mm	SLES5	S-6
4240	ELEVATORS - Luxator - Curved 5mm	SLEC5	S-6
4241	ELEVATORS - 34S Seldin	E34S	S-7
4242	ELEVATORS - 1L Seldin	E1L	S-7
4243	ELEVATORS - 1R Seldin	E1R	S-7
4244	ELEVATORS - 4L Seldin	E4L	S-7
4245	ELEVATORS - 4R Seldin	E4R	S-7
4246	ELEVATORS - 34 Seldin	E34	S-7
4247	ELEVATORS - 25 Seldin	E25	S-7
4248	ELEVATORS - 35 Seldin	E3S	S-7
4249	ELEVATORS - 301W Seldin	E301W	S-7
4250	ELEVATORS - 304W Seldin	E304W	S-7
4252	ELEVATORS - 71 Miller Apexo	E71	S-8
4253	ELEVATORS - 72 Miller Apexo	E72	S-8
4254	ELEVATORS - 73 Miller Apexo	E73	S-8
4255	ELEVATORS - 74 Miller Apexo	E74	S-8
4256	ELEVATORS - 4 Schmeckebier Apexo	E4	S-8
4257	ELEVATORS - 5 Schmeckebier Apexo	E5	S-8
4258	ELEVATORS - 301 Apexo	E301	S-8
4259	ELEVATORS - 302 Apexo	E302	S-8
4260	ELEVATORS - 303 Apexo	E303	S-8
4261	ELEVATORS - 8 Crane Pick	E8	S-8
4262	ELEVATORS - A Cogswell	EA	S-8
4263	ELEVATORS - B Cogswell	EB	S-8
4264	ELEVATORS - 1C Coupland Gouge	ECG1C	S-9
4265	ELEVATORS - 2C Coupland Gouge	ECG2C	S-9
4266	ELEVATORS - 3C Coupland Gouge	ECG3C	S-9
4267	ELEVATORS - A1C Coupland Gouge	ECGA1C	S-9
4268	ELEVATORS - A2C Coupland Gouge	ECGA2C	S-9
4269	ELEVATORS - A3C Coupland Gouge	ECGA3C	S-9
42611	ELEVATORS - 11M MacMillan Gouge	EMG11M	S-9
42621	ELEVATORS - 12M MacMillan Gouge	EMG12M	S-9
42631	ELEVATORS - 13M MacMillan Gouge	EMG13M	S-9
4270	ELEVATORS - 19	E19	S-10
4271	ELEVATORS - 20	E20	S-10
4272	ELEVATORS - 23	E23	S-10
4273	ELEVATORS - 38 Serrated	E38	S-10
4274	ELEVATORS - 39 Serrated	E38	S-10
4275	ELEVATORS, 559/1-Bein, Large	EBE559/1	S-10
4276	ELEVATORS,559/2-Bein, Medium	EBE559/2	S-10
4277	ELEVATORS, 559/3-Bein, Small	EBE559/3	S-10
4278	ELEVATORS - 31F Friedman [Craig] Gouge Serrated	E31F	S-10
4279	ELEVATORS - 32F Friedman [Craig] Gouge Serrated	E32F	S-10
4280	ELEVATORS - 51 Curtis	EC51	S-10

4281	ELEVATORS - 52 Curtis	EC52	S-10
4282	Elevator H5, Hourigan -Serrated	H5S	S-10
4283	Elevator H6, Hourigan -Serrated	H6S	S-10
4284	ELEVATORS - 6X Potts	E6X	S-11
4285	ELEVATORS - 7X Potts	E7X	S-11
4286	ELEVATORS - 15 Woodward	EW15	S-11
4287	ELEVATORS - 16 Woodward	EW16	S-11
4288	ELEVATORS - 11 Stout	E11	S-11
4289	ELEVATORS - 11A Stout	E11A	S-11
4290	ELEVATORS - 80	E80	S-11
4291	ELEVATORS - 81	E81	S-11
4292	ELEVATORS - 81A	E81A	S-11
4293	ELEVATORS - 81E	E81E	S-11
4294	ELEVATORS - 92 Serrated	E92	S-11
4295	ELEVATORS - 93 Serrated	E93	S-11
4296	ELEVATORS - 98	E98	S-11
4297	ELEVATORS - 190	E190	S-11
4298	ELEVATORS - 191	E191	S-11
4300	ELEVATORS - 6X Potts, T-Bar Handle	E6XTBH	S-12
4301	ELEVATORS - 7X Potts, T-Bar Handle	E7XTBH	S-12
4302	ELEVATORS - 3 Hourigan, T-Bar Handle	E3HTBH	S-12
4303	ELEVATORS - 4 Hourigan, T-Bar Handle	E4HTBH	S-12
4304	ELEVATORS - Warwick - James Straight, T-Bar Handle	WWJS	S-12
4305	ELEVATORS - Warwick - James Left, T-Bar Handle	WWJL	S-12
4306	ELEVATORS - Warwick-James Right, T-Bar Handle	WWJW	S-12
4307	ELEVATORS - 1LX, T-Bar Handle	E1LXCBH	S-12
4308	ELEVATORS - 1RX, T-Bar Handle	E1RXC BH	S-12
4309	ELEVATORS - 11LX, T-Bar Handle	E11LXCBH	S-12
4310	ELEVATORS - 11RX, T-Bar Handle	E11RXC BH	S-12
4311	ELEVATORS - 12LX, T-Bar Handle	E12LXCBH	S-12
4312	ELEVATORS - 12RX, T-Bar Handle	E12RXC BH	S-12
4313	ELEVATORS - 14LX, T-Bar Handle	E14LXCBH	S-12
4314	ELEVATORS - 14RX, T-Bar Handle	E14RXC BH	S-12
4316	ELEVATORS - ROOT TIP PICKS - 1 Heidbrink	EHB1	S-13
4317	ELEVATORS - ROOT TIP PICKS - 2 Heidbrink	EHB2	S-13
4318	ELEVATORS - ROOT TIP PICKS - 3 Heidbrink	EHB3	S-13
4319	ELEVATORS - ROOT TIP PICKS - 4 Heidbrink	EHB4	S-13
4320	ELEVATORS - ROOT TIP PICKS - 5 Heidbrink	EHB5	S-13
4321	ELEVATORS - ROOT TIP PICKS - 13/14 Heidbrink	EHB13/14	S-13
4322	ELEVATORS - ROOT TIP PICKS - 9 Apical	E9	S-13
4323	ELEVATORS - ROOT TIP PICKS - 9R	E9R	S-13
4324	ELEVATORS - ROOT TIP PICKS - 9L	E9L	S-13
4325	ELEVATORS - ROOT TIP PICKS - 1 West	EW1	S-14
4326	ELEVATORS - ROOT TIP PICKS - 2 West	EW2	S-14
4327	ELEVATORS - ROOT TIP PICKS - 3 West	EW3	S-14
4328	ELEVATORS - ROOT TIP PICKS - 4 West	EW4	S-14
4329	ELEVATORS - ROOT TIP PICKS - 5 West	EW4	S-14

4330	ELEVATORS - ROOT TIP PICKS - 6 West	EW6	S-14
4331	ELEVATORS - ROOT TIP PICKS - 1 Howard	EHO1	S-14
4332	ELEVATORS - ROOT TIP PICKS - 2 Howard	EHO2	S-14
4333	ELEVATORS - ROOT TIP PICKS - 3 Howard	EHO3	S-14
4334	ELEVATORS - ROOT TIP PICKS - 4 Howard	EHO4	S-14
4335	ELEVATORS - ROOT TIP PICKS - 5 Howard	EHO5	S-14
4336	ELEVATORS - ROOT TIP PICKS - 2/3 Howard	EHO2/3	S-14
4337	ELEVATORS - ROOT TIP PICKS - 1 Davis	ED1	S-15
4338	ELEVATORS - ROOT TIP PICKS - 2 Davis	ED2	S-15
4339	ELEVATORS - ROOT TIP PICKS - 3 Davis	ED3	S-15
4340	ELEVATORS - ROOT TIP PICKS - 4 Davis	ED4	S-15
4341	ELEVATORS - ROOT TIP PICKS - 8 Davis	ED8	S-15
4342	ELEVATORS - ROOT TIP PICKS - 9 Davis	ED9	S-15
4343	ELEVATORS - ROOT TIP PICKS - 10 Davis Double-End	ED10	S-15
4344	ELEVATORS - ROOT TIP PICKS - 11 Davis Double-End	ED11	S-15
4345	ELEVATORS - ROOT TIP PICKS - 94 Davis Single-edge	E94	S-15
4346	ELEVATORS - ROOT TIP PICKS - 95 Davis, Serrated Single-edge	E95	S-15
4347	Periotomes Elevator - 1 Posterior	SPT1	S-16
4348	Periotomes Elevator - 3 Posterior Right	SPTR3	S-16
4349	Periotomes Elevator - 4 Posterior Left	SPTL4	S-16
4350	Periotomes Elevator - 2 Anterior	SAT2	S-16
4351	Periotomes Elevator - 5 Anterior-Angled	SATA5	S-16
4352	Periotomes Elevator - 6 Anterior-Straight	SATS6	S-16
4353	Periotomes Elevator - Wise	SPW	S-16
4400	Periosteal Elevators - Bennett	SPSB	S-17
4401	Periosteal Elevators - P149	SPS149	S-17
4402	Periosteal Elevators - Hopkins -1	SPSH1	S-17
4403	Periosteal Elevators - Hopkins -2	SPSH2	S-17
4404	Periosteal Elevators - Mead-3	SPSM3	S-17
4405	Periosteal Elevators - Mead-3Mini	SPSM3M	S-17
4406	Periosteal Elevators - Molt -9	SPSM9	S-17
4408	Periosteal Elevators - GF-14	SPSGF14	S-17
4409	Periosteal Elevators - 3/9A	SPS3/9A	S-18
4410	Periosteal Elevators - PR-3	SPSPR3	S-18
4412	Periosteal Elevators - P16	SPSP16	S-18
4413	Periosteal Elevators - P15	SSPSP15	S-18
4414	Periosteal Elevators - Howard-9	SPSH9	S-18
4415	Periosteal Elevators - Howard-9 Mini	SPSH9M	S-18
4416	Periosteal Elevators - Allen-9	SPSA9	S-18
4417	Periosteal Elevators - 24G	SPS24G	S-19
4419	Periosteal Elevators - Benque	SPSBEN	S-19
4420	Periosteal Elevators - Buser	SPSBUS	S-19
4422	Periosteal Elevators - KN-1	SPSKN1	S-19
4423	Periosteal Elevators - KN-152	SPSKN152	S-19
4424	Periosteal Elevators - Hourigan-2	SPSHOURI	S-19
4425	Periosteal Elevators - HF-20	SPSHF20	S-20
4427	Periosteal Elevators - Molt -5	SPSM5	S-20

4428	Periosteal Elevators - Selden-23	SPSSEL23	S-20
4429	Periosteal Elevators - Freer	SPSFR	S-20
4430	Periosteal Elevators - OHL-10 Periosteal Elevator	SPSOHL10	S-20
4431	Periosteal Elevators - Woodson-1 Modify	SPSWSONM	S-20
4432	Periosteal Elevators - Allen Periosteal Elevator		
4480	SURGICAL CURETTES (SPOON) - 84 Lucas	SSCL84	S-21
4481	SURGICAL CURETTES (SPOON) - 85 Lucas	SSCL85	S-21
4482	SURGICAL CURETTES (SPOON) - 86 Lucas	SSCL86	S-21
4483	SURGICAL CURETTES (SPOON) - 87 Lucas	SSCL87	S-21
4484	SURGICAL CURETTES (SPOON) - 88 Lucas	SSCL88	S-21
4485	SURGICAL CURETTES (SPOON) - 8 Miller	SSC8M	S-22
4486	SURGICAL CURETTES (SPOON) - 9 Miller	SSC9M	S-22
4487	SURGICAL CURETTES (SPOON) - 10 Miller	SSC10M	S-22
4488	SURGICAL CURETTES (SPOON) - 11 Miller	SSC11M	S-22
4489	SURGICAL CURETTES (SPOON) - 12 Miller	SSC12M	S-22
4490	SURGICAL CURETTES (SPOON) - Molt 2/4	SSC2/4M	S-22
4491	SURGICAL CURETTES - Molt 2 (3.5mm)	SM2-3.5	S-22
4492	SURGICAL CURETTES - Molt 4 (7.0mm)	SM4-7.0	S-22
4540	BONE CHISELS - 1 Gardner 3mm	SBCSE1G3	S-23
4541	BONE CHISELS - 2 Gardner 4mm	SBCSE2G4	S-23
4542	BONE CHISELS - 5 Gardner 4.3mm	SBCSE5G4.3	S-23
4543	BONE CHISELS - 6 Gardner 4mm	SBCSE6G4	S-23
4544	BONE CHISELS - 9 Gardner 4mm	SBCSE9G4	S-23
4545	BONE CHISELS - 1 Chandler 5mm	SBCSE1CH5	S-23
4546	BONE CHISELS - 2 Chandler 5mm	SBCSE2CH5	S-23
4547	BONE CHISELS - 4 Chandler 5mm	SBCSE4CH5	S-23
4548	BONE CHISELS - 7 Chandler, Straight 4.7mm	SBCDB7CH4.7	S-23
4549	BONE CHISELS - 8 Chandler, BI-B, Straight 6.1mm	SBCDB8CH6.1	S-23
4550	BONE CHISELS - 9 Chandler, BI-B, Straight 7.1mm	SBCDB9CH7.1	S-23
4553	BONE CHISELS - 1W Wakefield 2.7mm	SBCDB18M4.2	S-23
4554	BONE CHISELS - 2W Wakefield 3.0mm	SBCDB1W2.7	S-23
4555	BONE CHISELS - 18 McFarland 4.2mm	SBCDB18M4.2	S-23
4600	BONE FILES - 1X Miller -Collburn Cross -Cut	SBF1X	S-24
4601	BONE FILES - 2X Miller -Collburn Cross -Cut	SBF2X	S-24
4602	BONE FILES - 3X Miller -Collburn Cross -Cut	SBF3X	S-24
4603	BONE FILES - 21 Miller	SBF21M	S-24
4604	BONE FILES - 45 Miller	SBF45M	S-24
4605	BONE FILES - 52 Miller	SBF52M	S-25
4606	BONE FILES - 64 Miller	SBF64M	S-25
4607	BONE FILES - 1 Wahl	SBF1W	S-25
4608	BONE FILES - 2 Wahl	SBF2W	S-25
4609	BONE FILES - 12 Howard	SBF12H	S-25
4610	BONE FILES - 67 Howard	SBF67H	S-25
4720	SCISSORS - 5 Wagner Straight 11.5cm/4 1/2 "	SS5WS	S-26
4721	SCISSORS - 5 Wagner Curved 11.5cm/4 1/2 "	SS6WC	S-26
4722	SCISSORS - 5 Wagner Angled 11.5cm/4 1/2 "	SS7WA	S-26
4723	SCISSORS - 1 Kelly Straight /Pointer 16cm/6 1/2 "	SS1KSP	S-26

4724	SCISSORS - 2 Kelly Straight /Pointer 16cm/6 1/2 "	SS2KSP	S-26
4725	SCISSORS - Kelly Curved 18cm/7"	SSKC	S-26
4726	SCISSORS - Kelly Straight 18cm/7"	SSKS	S-26
4727	SCISSORS - 8 Quinby Curved 12.5cm/5"	SS8QC	S-27
4728	SCISSORS - 9 Dean 16.5cm/6 1/2"	SS9D	S-27
4729	SCISSORS - 10 Double-Curved 13.5cm/5 1/2 "	SS10DC	S-27
4730	SCISSORS - 11 Locklin Curved Handles 16cm/6 1/1 "	SS11LCH	S-27
4731	SCISSORS - 12 Loc klin 16cm/6 1/ 2"	SS12L	S-27
4732	SCISSORS - 14 LaGrange Double-Curved 11.5cm/4 1/2 "	SS14LGDC	S-27
4733	SCISSORS - 15 Sullivan 14cm/5 1/2 "	SS15S	S-27
4734	SCISSORS - 16 Goldman-Fox 12.5cm/5"	SS16GF	S-27
4735	SCISSORS - Joseph Curved 14cm/5 1/ "	SSJC	S-28
4736	SCISSORS - Joseph Straight 14cm/5 1/ " 2	SSJS	S-28
4737	SCISSORS - 3 Mayo Curved/Blunt 14.5cm/5 3/	SS3MCB	S-28
4738	SCISSORS - 4 Mayo Curved/Blunt 14.5cm/5 3/	SS4MCB	S-28
4739	SCISSORS - Strabismus Straight Tenotomy 11.5cm/4 1/ "	SSSBMST	S-28
4740	SCISSORS - 31 Castroviejo Curved 10cm/4"	SS31CASC	S-28
4741	SCISSORS - 17 Iris Straight /Delicate 11.5cm/4 1/ "	SS17ISD	S-28
4742	SCISSORS - 18 Iris Straight /Delicate 11.5cm/4 1/ "	SS18ICD	S-28
4743	SCISSORS - Metzenbaum Baby 12.5cm/5"	SSMB	S-28
4744	SCISSORS - 21 Straight / Pointed 14.5cm/5 3/ "	SS21SP	S-29
4745	SCISSORS - 22 Straight / Blunt -Sharp 14.5cm/5 3/ "	SS22SBS	S-29
4746	SCISSORS - 23 Curved/ Pointed 14.5cm/5 3/ "	SS24CP	S-29
4747	SCISSORS - Lister Bandage 14.5cm/5 / "	SSLB	S-29
4748	SCISSORS - 13 Suture 15cm/6"	SS13S	S-29
4749	SCISSORS - 13S Suture3 15cm/4 3/4 "	SS13SS	S-29
4750	SCISSORS - Crown & Gold Straight 15cm/4 3/4 "	SSC&GS	S-29
4751	SCISSORS - Crown & Gold Curved 15cm/4 3/4 "	SSC&GC	S-29
4752	SCISSORS - Crown & Gold Curved Pedodontic 15cm/3.5"	SSC&GP	S-29
4800	TUNGUSTON CARBIDE TIPPED SCISSORS - 10 Double-Curved 13.5cm/5 1/ "	SS10DC	S-30
4801	TUNGUSTON CARBIDE TIPPED SCISSORS - 14 Double-Curved 11.5cm/5 1/ "	SS14DC	S-30
4802	TUNGUSTON CARBIDE TIPPED SCISSORS - 16 Goldman-Fox 12.5cm/5	SS16GFTC	S-30
4803	TUNGUSTON CARBIDE TIPPED SCISSORS - 17 Iris Straight /Delicate 11.5cm/4 1/ "	SS17ISDTC	S-30
4804	TUNGUSTON CARBIDE TIPPED SCISSORS - 18 Iris Curved/Delicate 11.5cm/4 1/ "	SS18ICDTC	S-30
4805	TUNGUSTON CARBIDE TIPPED SCISSORS - Metzenbaum Straight /Blunt 14.5cm/5 1/ "	SSMSBTC	S-30
4806	TC SCISSORS - Metzenbaum Curved/Blunt 14.5cm/5 1/2 "	SSMCBTC	S-30
4807	TUNGUSTON CARBIDE TIPPED SCISSORS - Metzenbaum Straight / Pointed 14.5cm/5 1/2 "	SSMSPTC	S-30
4808	TUNGUSTON CARBIDE TIPPED SCISSORS - Metzenbaum Curved/ Pointed 14.5cm/5 1/2 "	SSMCPTC	S-30
4809	TUNGUSTON CARBIDE TIPPED SCISSORS - Metzenbaum Delicate- Curved/ Blunt 14.5cm/5 1/2 "	SSMDCBTC	S-30
4810	SCISSORS - Goldman-Fox Straight 13cm/5	SSGFSTC	S-31
4811	SCISSORS - Goldman-Fox Curved 13cm/5"	SSGFCTC	S-31
4812	SCISSORS - Wire Cutting Notched 12cm/4 1/2 "	SSWCNTC	S-31
4813	SCISSORS - Iris (17) Straight 11.5cm/4 1/2 "	SS17ISTC	S-31
4814	SCISSORS - Iris (17) Curved 11.5cm/4 1/2 "	SS17ICTC	S-31
4815	SCISSORS - Dean 16.5cm/6 1/2 "	SSDTC	S-31
4816	SCISSORS - Locklin Curved Handle 16.5cm/6 1/2 "	SSLCHTC	S-31
4817	SCISSORS - Locklin Straight Handle 16.5cm/6 1/2 "	SSLSHTC	S-31

4818	SCISSORS - Kelly Curved 16cm/6 1/2	SSKCTC	S-31
4819	SCISSORS - Kelly Straight 16cm/6 1/2	SSKSTC	S-31
4870	Tissue Pliers, 31 Semkin-Taylor, Straight 12.5cm/5"	TP31STS	S-32
4871	TISSUE PLIERS 32 Semkin-Taylor Curved 12.5cm/5"	TP32STC	S-32
4872	TISSUE PLIERS AND FORCEPS - 33 Semk in-Taylor 1x 2, Straight 12.5cm/5"	TP33STS	S-32
4873	TISSUE PLIERS AND FORCEPS - 33 Semk in-Taylor 1x 2, Curved 12.5cm/5"	TP34STC	S-32
4874	TISSUE PLIERS AND FORCEPS - 16- 3x4, 15cm/6"	TP16	S-32
4875	TISSUE PLIERS AND FORCEPS - 50-L Fox Locking 3x3 15cm/6"	TP50L	S-32
4876	TISSUE PLIERS AND FORCEPS - Kramer -Nevins Tissue Pliers 1x 2 Curved	TPKN	S-32
4877	TISSUE PLIERS AND FORCEPS - 20 Corn Suture Pliers 15cm/6"	TP20C	S-32
4878	TISSUE PLIERS AND FORCEPS - 41 Adson Plain 12cm/4 3/4 "	TPAP41	S-32
4879	TISSUE PLIERS AND FORCEPS - 42 Adson1x2 12cm/4 3/4 "		S-32
4880	TISSUE PLIERS AND FORCEPS - 43 Adson-Brown 7x7 12cm/4 3/4 "	TPAB43	S-32
4881	TISSUE PLIERS AND FORCEPS - Russian 15cm/6"	TPRN15	S-32
4882	TISSUE PLIERS AND FORCEPS - Russian 20cm/8"	TPRN20	S-32
4883	TISSUE PLIERS AND FORCEPS - Gerald Plain 18cm/7" Straight /Serrated	TPGP18	S-32
4884	TISSUE PLIERS AND FORCEPS - Gerald 1x2, 15cm/7" Straight	TPG15	S-32
4885	TISSUE PLIERS AND FORCEPS - Allison Baby 4x5 12.5cm/5"	FALIB12.5	S-32
4886	TISSUE PLIERS AND FORCEPS - Allison Straight 5x6 15cm/6"	FALIS15	S-32
4887	TC - Adson3 12cm/4 / "	TPAD12TC	S-33
4888	TISSUE PLIERS/FORCEPS - Adson 1x2 12cm/4 3/4 "	TPAD1X2TC	S-33
4889	TISSUE PLIERS/FORCEPS - Adson Brown3 12cm/4 / "	TPADB12TC	S-33
4890	TISSUE PLIERS/FORCEPS - General 15cm/6"	TPG15TC	S-33
4891	TISSUE PLIERS/FORCEPS - General 1x2 17cm/6"	TPG1X2TC	S-33
4892	TISSUE PLIERS/FORCEPS - Gerald 18cm/7"	TPGRD18TC	S-33
4893	TISSUE PLIERS/FORCEPS - Gerald 1x2 18cm/7"	TPGRD1X2TC	S-33
4894	TISSUE PLIERS/FORCEPS - Cushing 18cm/7"	TPGUSHTC	S-33
4895	TISSUE PLIERS/FORCEPS - Cushing Angled 18cm/7"	SLB	S-33
4896	TISSUE PLIERS/FORCEPS - Cushing Bayone, 19cm/7 1/2 "	SLB	S-33
4897	TISSUE PLIERS/FORCEPS - Potts Smith 18cm/7" (Shown)	SLB	S-33
4898	TISSUE PLIERS/FORCEPS - Potts Smith 14.5cm/5/1/4" (Shown)	SLB	S-33
4899	TISSUE PLIERS/FORCEPS - De Bakey 15cm/6"(Shown)	SLB	S-33
4900	TISSUE PLIERS/FORCEPS - De Bakey 19cm/7 1/2 "	SLB	S-33
4950	NEEDLE HOLDERS - Hegar-Baumgartner (Gardner 5) 14cm/5 1/2	NHHB	S-34
4951	NEEDLE HOLDERS - Webster Smooth 13cm/5"	NHWS13	S-34
4952	NEEDLE HOLDERS - Mayo-Hegar (Gardner 6) Grooved 16cm/6 1/2 "	NHMH	S-34
4953	NEEDLE HOLDERS - Crile-Wood Grooved 15cm/6"	NHCW	S-34
4954	NEEDLE HOLDERS - Olsen-Hegar Needle Holder /Scissors 17cm/6 3/4 "	NHOHHS	S-34
4955	NEEDLE HOLDERS - Boynton 12cm/4 3/4 "	NHBOY	S-34
4956	NEEDLE HOLDERS - Castroviejo 14cm/5 1/2 "	NHCV	S-34
4957	NEEDLE HOLDERS - Collier 12.5cm/5"	NHC	S-34
4958	NEEDLE HOLDERS - Derf Grooved 12.5cm/5"	NHDG	S-34
4959	NEEDLE HOLDERS - Mathieu (Machatt) 14cm/5 1/2 "	NHMM	S-34
4960	TC-NEEDLE HOLDERS - Derf 11.5cm/4 1/2 "	NHDTC	S-35
4961	TC-NEEDLE HOLDERS - Hegar -Baumgartner 11.5cm/4 1/2 "	NHHBTC	S-35
4962	TC-NEEDLE HOLDERS - Halsey Smooth 13cm/5"	NHHSTC	S-35
4963	TC-NEEDLE HOLDERS - Halsey 13cm/5"	NHHTC	S-35

4970	TC-NEEDLE HOLDERS - Castroviejo, Straight, Flat Handle, 13cm/5.15"	NHCVSFH13TC	S-36
4971	TC-NEEDLE HOLDERS - Castroviejo, Curved, Flat Handle, 13cm/5.15	NHCVFHC13TC	S-36
4972	TC-NEEDLE HOLDERS - Castroviejo, Straight, Flat Handle, 14cm/5.5"	NHCVFHS14TC	S-36
4973	TC-NEEDLE HOLDERS - Castroviejo Curved, Flat Handle, 14cm/5.5"	NHCVFHC14TC	S-36
4974	TC-NEEDLE HOLDERS - Castroviejo, Straight, Flat Handle, 15cm/5.9"	NHCVFHS15TC	S-36
4975	TC-NEEDLE HOLDERS - Castroviejo Curved, Flat Handle, 15cm/5.9"	NHCVFHC15TC	S-36
4976	TC-NEEDLE HOLDERS - Castroviejo, Straight, Flat Handle, 16cm/6.3"	NHCVFHS16TC	S-36
4977	TC-NEEDLE HOLDERS - Castroviejo, Curved, Flat Handle, 16cm/6.3"	NHMCVFC16TC	S-36
4978	TC-NEEDLE HOLDERS - Castroviejo Straight, Flat Handle 18cm/7"	NHCVFHS18TC	S-36
4979	TC-NEEDLE HOLDERS - Castroviejo Curved, Flat Handle 18cm/7"	NHCVFHC18TC	S-36
4982	TC-NEEDLE HOLDERS - Castroviejo Straight, Round Handle, 13cm/5.15"	NHCVRHS13TC	S-37
4983	TC-NEEDLE HOLDERS - Castroviejo Curved, Round Handle, 13cm/5.15	NHCVRHC13TC	S-37
4984	TC-NEEDLE HOLDERS - Castroviejo, Straight, Round Handle, 14cm/5.5"	NHCVRHS14TC	S-37
4985	TC-NEEDLE HOLDERS - Castroviejo, Curved, Round Handle, 14cm/5.5	NHCVRHC14TC	S-37
4986	TC-NEEDLE HOLDERS - Castroviejo, Straight, Round Handle, 15cm/5.9"	NHCVRHS15TC	S-37
4987	TC-NEEDLE HOLDERS - Castroviejo, Curved, Round Handle, 15cm/5.9"	NHCVRHC15TC	S-37
4988	TC-NEEDLE HOLDERS - Castroviejo, Straight, Round Handle, 16cm/6.3"	NHCVRHS16TC	S-37
4989	TC-NEEDLE HOLDERS - Castroviejo, Curved, Round Handle, 16cm/6.3"	NHCVRHC16TC	S-37
4990	TC-NEEDLE HOLDERS - Castroviejo, Straight, Round Handle, 18cm/7"	NHCVRHS18TC	S-37
4991	TC-NEEDLE HOLDERS - Castroviejo, Curved, Round Handle, 18cm/7"	NHCVRHC18TC	S-37
4995	TC-NEEDLE HOLDERS - Arruga Serrated 16cm/6 1/2 "	NHARRUS16TC	S-38
4996	NEEDLE HOLDERS - Webster Smooth 11.5cm/4 1/2 "	NHWEBS11.5TC	S-38
4997	NEEDLE HOLDERS - Crile-Wood 15cm/6"	NHCW15TC	S-38
4998	NEEDLE HOLDERS - Olsen-Hegar Needle Holder 17cm/6.25 "	NHSOH17TC	S-38
4999	NEEDLE HOLDERS - Olsen-Hegar Needle Holder 14cm/5-1/4 "	NHSOH14TC	S-38
5000	NEEDLE HOLDERS - DeBakey 18cm/7" Delicate	NHDEB18TC	S-38
5001	Needle Holders-Mayo-Hegar, 16cm/6.25"	NHMH16TC	S-38
5002	NEEDLE HOLDERS - Mayo-Hegar 18cm/7"	NHMH18TC	S-38
5005	NEEDLE HOLDERS - Matheu-Kocher 14cm/5 1/2 "	NHMK14TC	S-39
5006	NEEDLE HOLDERS - Matheu 14cm/5 1/2 "	NHM14TC	S-39
5007	NEEDLE HOLDERS - HMT Thin Jaw 14cm/5 1/2 "	NHFSP15TC	S-39
5008	NEEDLE HOLDERS - Fine Swedish Pattern 15cm/6"	NHMVP15TC	S-39
5009	NEEDLE HOLDERS - Micro-Vascular Pattern 15cm/6"	NHMCP18TC	S-39
5010	NEEDLE HOLDERS - Micro-Vascular Pattern 18cm/7"	NHMRP15TC	S-39
5011	NEEDLE HOLDERS - Mini-Ryder Pattern 15cm/6"	NHMRP18TC	S-39
5012	NEEDLE HOLDERS - Mini-Ryder Pattern 18cm/7"		
5030	HAEMOSTATIC FORCEPS - Hartman Straight 10cm/4"	HHS10	S-40
5031	HAEMOSTATIC FORCEPS - Hartman Curved 10cm/4"	HHC10	S-40
5032	HAEMOSTATIC FORCEPS - Adson Straight 19cm/7"	HAS19	S-40
5033	HAEMOSTATIC FORCEPS - Adson Curved 19cm/7"	HAC19	S-40
5034	HAEMOSTATIC FORCEPS - No.5 -1x2 Teeth Straight 12.5cm/5"	H5S1X2	S-40
5035	HAEMOSTATIC FORCEPS - Rochester-pean Curved 16cm/6"	HRPC16	S-40
5036	HAEMOSTATIC FORCEPS - 3 Halsted-Mosquito Curved 12cm/4"	H3HMC12	S-41
5037	HAEMOSTATIC FORCEPS - 3 Halsted-Mosquito Straight 12cm/4"	H3HMS12	S-41
5038	HAEMOSTATIC FORCEPS - Kocher 1x2 Teeth Curved 14cm/5 1/2 "	HK1X2C14	S-41
5039	HAEMOSTATIC FORCEPS - 1 Kelly Curved 14cm/5"	H1KS14	S-41
5040	HAEMOSTATIC FORCEPS - 2 Kelly Curved 14cm/5"	H2KC14	S-41

5041	HAEMOSTATIC FORCEPS - Kelly-Rankin Straight 16cm/6"	HKRS16	S-41
5042	HAEMOSTATIC FORCEPS - Kelly-Rankin Curved 16cm/6"	HKRC16	S-41
5050	Middeldorf Cheek Retractor 17x14mm		
5051	Middeldorf Cheek Retractor 22x20mm		
5052	9 Ramus Retractors 16x80mm		
5053	8 Surgical Retractor 11x70mm		
5059	1 Surgical Retractor 7x25mm		
5060	1 Surgical Retractor10x42		
5061	1 Surgical Retractor 12x55mm		
5066	1 Austin Retractor		
5067	2 Austin Retractor		
5068	1 Henahan Retractor		
5069	2 Henahan Retractor		
5070	3 Henahan Retractor		
5100	RETRACTORS - University of Minnesota	RUOM	S-42
5101	RETRACTORS - Lindo Serrated	RLS	S-42
5102	RETRACTORS - Front Round	RFR	S-42
5103	RETRACTORS - Front Cut	RFC	S-42
5104	RETRACTORS - Cawood Minnesota	RCM	S-42
5110	RETRACTORS - Columbia 14cm	RC14cm	S-42
5111	RETRACTORS - Orringer Hands Free	ROHF	S-42
5112	RETRACTORS - Bishop	RBISH	S-42
5115	RETRACTORS - Mouth Gag (adult size) 14cm	RMGA14cm	S-42
5116	RETRACTORS - Mouth Gag (Child size) 12cm	RMGC12cm	S-42
5117	RETRACTORS - Mouth Gag (Pedo size) 10cm	RMGP10cm	S-42
5120	RETRACTORS - ALM 7cm	RALM7cm	S-42
5121	RETRACTORS - Shuman adult	RSHUA	S-42
5122	RETRACTORS - Shuman Child	RSHUC	S-42
5125	RETRACTORS - Retractor Large	RRL	S-43
5126	RETRACTORS - Lip Retractor	RLIP	S-43
5127	RETRACTORS - Retractors	RHR	S-43
5128	Retractor Small		
5130	RETRACTORS - Weider Tongue 30mm	RWTS30	S-43
5131	RETRACTORS - Weider Tongue Large 40mm	RWTL40	S-43
5132	RETRACTORS - Weider Cheek & Tongue Retractor Large 40mm with suction tube (2 in 1)	RSW	S-43
5153	RETRACTORS - Desmarres Chalazion Surgical Retractors 10X10 mm	DCSR10X10	S-43
5154	RETRACTORS - Desmarres Chalazion Surgical Retractors 11X16 mm	DCSR11X16	S-43
5155	RETRACTORS - Desmarres Chalazion Surgical Retractors 11X21 mm	DCSR11X21	S-43
5159	RETRACTORS, Farabeuf 12.5cm	RFARA12.5	S-43
5160	RETRACTORS, Farabeuf 15cm	RFARA15	S-43
5161	RETRACTORS, Farabeuf 16.5cm	RROUX16.5	S-43
5200	SCALPEL HANDLES - 3, 12cm	SH3-12cm	S-44
5201	SCALPEL HANDLES - 3G, 12cm	SH3G-12cm	S-44
5202	SCALPEL HANDLES - 4, 13.5cm	SH4-13.5cm	S-44
5203	SCALPEL HANDLES - 5, 14.5cm	SH5-14.5cm	S-44
5204	SCALPEL HANDLES - 5A, 14.5cm	SH5A-14.5cm	S-44
5205	SCALPEL HANDLES - Double Head 14.5cm	SHDH-14.5cm	S-44

5206	SCALPEL HANDLES - 7, 16.5cm	SH7-16.5cm	S-44
5207	SCALPEL HANDLES - 3L, 21cm	SH3L-21cm	S-44
5208	SCALPEL HANDLES - Scalpel 15.2cm	SHS15.2cm	S-44
5209	SCALPEL HANDLES - Scalpel 180° 16.2cm	SHS180	S-44
5249	Coupland Handle 25cm/9.9"	SASPCH	S-45
5250	ASPIRATORS - 1 Coupland Tip 1.5mm	SASPCT1	S-45
5251	ASPIRATORS - 2 Coupland Tip 2.0mm	SASPCT2	S-45
5252	ASPIRATORS - 3 Coupland Tip 2.5mm	SASPCT3	S-45
5253	ASPIRATORS - 4 Coupland Tip 3.0mm	SASPCT4	S-45
5254	ASPIRATORS - 1 Coupland Tip 1.5mm/with side hole	SASPCT1SH	S-45
5255	ASPIRATORS - 2 Coupland Tip 2.0mm/with side hole	SASPCT2SH	S-45
5256	ASPIRATORS - 3 Coupland Tip 2.5mm/with side hole	SASPCT3SH	S-45
5257	ASPIRATORS - 4 Coupland Tip 3.0mm/with side hole	SASPCT4SH	S-45
5259	ASPIRATORS - Byrd Self-Cleaning 20cm/8", Tip size 1.5mm	SASPB1	S-45
5260	ASPIRATORS - Byrd Self-Cleaning 20cm/8", Tip size 2.0mm	SASPB2	S-45
5261	ASPIRATORS - Byrd Self-Cleaning 20cm/8", Tip size 2.5mm	SASPB3	S-45
5262	ASPIRATORS - Cogswell 3 Suction Tube, 23cm/9"	SASPCG3	S-45
5263	ASPIRATORS - Frazier Suction Tube, 6 Frazier 1.5mm	SASPF6	S-45
5264	ASPIRATORS - Frazier Suction Tube, 8 Frazier 2.0mm	SASPF8	S-45
5265	ASPIRATORS - Frazier Suction Tube, 10 Frazier 2.5mm	SASPF10	S-45
5266	ASPIRATORS - Frazier Suction Tube, 12 Frazier 3.0mm	SASPF12	S-45
5270	ASPIRATORS - Oral Surgery Suction tube 22cm/8.5" 1mm	SASP1	S-45
5271	ASPIRATORS - Oral Surgery Suction tube 22cm/8.5" 2mm	SASP2	S-45
5272	ASPIRATORS - Oral Surgery Suction tube 22cm/8.5" 3mm	SASP3	S-45
5273	ASPIRATORS - Oral Surgery Suction tube 22cm/8.5" 4mm	SASP4	S-45
5275	ASPIRATORS - Yankauer Oral Surgery 26.5cm/10.5"	SASPYD	S-45
5276	ASPIRATORS - Andrew-Pychon 23cm/9"	SASPAP	S-45
5280	Locking Pliers 16.5cm/6-1/2"	SLP	S-45
5320	Intraligmentary Syringe 1.8ml	SY11.8ml	S-46
5322	Syringe 1.8ml	SY1.8ml	S-46
5323	1 FINGER SYRING1.8ML	SY111.8ml	S-46
5325	Syringe 1.8ml	SYS1.8	S-46
5328	SYRINGES - Paroject, Stainless Steel	SYPARS	S-46
5329	SYRINGES - Paroject Plastic	SYPARP	S-46
5380	1 Iwanson Spring Caliper-Metal	CLP1	S-47
5381	2 Iwanson Spring Caliper-Wax	CLP2	S-47
5383	Marix Retainer Tofflemire Junior	TJ	S-47
5384	Marix Retainer Tofflemire Senior	TS	S-47
5385	Marix Retainer Tofflemire Universal	TU	S-47
5386	Stainless Steel Ruler 15cm/6"	SS6R	S-47
5387	Boley Gauge	BG	S-47
5430	ATRAUMATIC FORCEPS - Upper Universal	EF1	S-48
5431	ATRAUMATIC FORCEPS - Lower Universal	EF2	S-48
5432	ATRAUMATIC FORCEPS - Upper Anterior	EF3	S-48
5433	ATRAUMATIC FORCEPS - Lower Anterior	EF4N	S-48
5434	ATRAUMATIC FORCEPS - Upper 3rd Molar	EF6	S-48
5435	ATRAUMATIC FORCEPS - Lower Molar Modify Cow-Horn Bead	EF7C	S-48

5436	ATRAUMATIC FORCEPS - Lower Root Extra Long	E46XL	S-48
5437	ATRAUMATIC FORCEPS - Upper Root	E51XL	S-48
5438	ATRAUMATIC FORCEPS - Lower Wisdoms	E222	S-48
5439	ATRAUMATIC FORCEPS - Lower Roots	E301	S-48
5440	ATRAUMATIC FORCEPS - Upper Roots	E30	S-48
5441	EXTRACTION FORCEPS - 5 Lower Third Molars	E5	S-49
5442	EXTRACTION FORCEPS - 15 Lower Molars	E15	S-49
5443	EXTRACTION FORCEPS - 16S Lower Molars	E16S	S-49
5444	EXTRACTION FORCEPS - 17 Lower Molars	E17	S-49
5445	EXTRACTION FORCEPS - 16 Lower Molar s CowHorn Beak	E16	S-49
5446	EXTRACTION FORCEPS - 1AF Upper Anteriors	E1AE	S-49
5447	EXTRACTION FORCEPS - 18L-Lower Molars, Left Side	E18L	S-49
5448	EXTRACTION FORCEPS - 18R-Upper Roots, Right Side	E18R	S-49
5449	EXTRACTION FORCEPS - Splitting Forces	E6	S-49
5450	EXTRACTION FORCEPS - Upper Centrals Laterals & Cuspids	E1	S-49
5451	EXTRACTION FORCEPS - 99C UpperPremolars & Incisors	E99C	S-50
5452	EXTRACTION FORCEPS - 88L Upper Molars Left Side	E88L	S-50
5453	EXTRACTION FORCEPS - 88L Upper Molars Right Side	E88R	S-50
5454	EXTRACTION FORCEPS - 10S Upper Molars Straight -Handle	E10S	S-50
5455	EXTRACTION FORCEPS - 53L Upper Molars Left Side	E53L	S-50
5456	EXTRACTION FORCEPS - 32 Upper Molars & Premolars Strong Bead	E53R	S-50
5457	EXTRACTION FORCEPS - 23 Lower Molars Cow-Horn Beak	E23	S-50
5458	EXTRACTION FORCEPS - 18R Upper Molars Right Side	E10H	S-50
5459	EXTRACTION FORCEPS - 69 Fragments & Small Roots Upper & Lower	E69	S-50
5460	EXTRACTION FORCEPS - 65 Upper Incisors & Roots	E65	S-50
5461	EXTRACTION FORCEPS - Universal for Upper Teeth & Roots	EU	S-50
5462	EXTRACTION FORCEPS - 53R Upper Molars Right Side	E32	S-50
5463	EXTRACTION FORCEPS - 151A Lower Premolars Incisors & Roots	E151A	S-51
5464	EXTRACTION FORCEPS - 151XAS Children Lower Premolars Incisors & Roots Canines	E151XAS	S-51
5465	EXTRACTION FORCEPS - 151AS Split Lower Premolar Incisors & Roots Canines	E151AS	S-51
5466	EXTRACTION FORCEPS - 151AF Apical Beak Lower Premolars Incisors & Roots Canines	E151AE	S-51
5467	EXTRACTION FORCEPS - 101 Universal for Premolars	E101	S-51
5468	EXTRACTION FORCEPS - 150 Upper Premolars Incisors & Roots Canines	E150	S-51
5469	EXTRACTION FORCEPS - 150AF Apical Beak Upper Premolars Incisors & Roots Canines	150AF	S-51
5470	EXTRACTION FORCEPS - 150AS Split Beak Upper Premolars Incisors & Roots Canines	E150AS	S-51
5471	EXTRACTION FORCEPS - 150XAS Children Split Upper Premolars Incisors & Roots Canines	E150XAS	S-51
5472	EXTRACTION FORCEPS - 150A Upper Premolars Incisors & Roots Canines	E150A	S-51
5473	EXTRACTION FORCEPS - 150S Upper Premolars Incisors & Roots for Children	E150S	S-51
5474	EXTRACTION FORCEPS - 151 Upper Premolars Incisors & Roots	E151	S-51
5475	EXTRACTION FORCEPS - 151S Lower Premolars Incisors & Roots for Children	E151S	S-52
5476	EXTRACTION FORCEPS - 210H Upper Wisdoms with Hook Handle	E210H	S-52
5477	EXTRACTION FORCEPS - 210S Upper Wisdoms with Straight Handle	E210S	S-52
5478	EXTRACTION FORCEPS - 222 Lower Wisdoms		
5479	EXTRACTION FORCEPS - 222AF Apical Beak Lower Wisdoms	E222AF	S-52
5480	EXTRACTION FORCEPS - 286 Upper Roots Premolars & Incisors	E286	S-52
5481	EXTRACTION FORCEPS - 30 Upper Roots & Incisors	E30	S-53
5482	EXTRACTION FORCEPS - 22 Lower Molars Hawks Bil	E22	S-53

5483	EXTRACTION FORCEPS - 21 Lower Molars	E21	S-53
5484	EXTRACTION FORCEPS - 20 Lower Wisdoms	E20	S-53
5485	EXTRACTION FORCEPS - 8 Lower Premolars	E8	S-53
5486	EXTRACTION FORCEPS - 7 Lower Premolars	E7	S-53
5487	EXTRACTION FORCEPS - 3 Upper Incisors & Canines	E3	S-53
5488	EXTRACTION FORCEPS - 1 Upper Incisors & Canines	E1	S-53
5489	EXTRACTION FORCEPS - 19 Upper Wisdoms	E19	S-53
5490	EXTRACTION FORCEPS - 18A Upper Molars	E18A	S-53
5491	EXTRACTION FORCEPS - 18Upper Molar sLeft	E18	S-53
5492	EXTRACTION FORCEPS - 17 Upper Molars Right	E17	S-53
5493	EXTRACTION FORCEPS - 67 AF Apical Beak Upper Wisdoms	E67AF	S-54
5494	EXTRACTION FORCEPS - 68 Lower Roots	E68	S-54
5495	EXTRACTION FORCEPS - 73 Lower Molars (Hawks Bill)	E73	S-54
5496	EXTRACTION FORCEPS - 74 Lower Incisors Root	E74	S-54
5497	EXTRACTION FORCEPS - 31 Lower Roots	E31	S-54
5498	EXTRACTION FORCEPS - 32 Lower Molars	E32	S-54
5499	EXTRACTION FORCEPS - 33 Lower Roots	E33	S-54
5500	EXTRACTION FORCEPS - 33 Lower Roots 33M Lower Roots With Long Beak	E33M	S-54
5501	EXTRACTION FORCEPS - 67 Upper Wisdoms	E67	S-54
5502	EXTRACTION FORCEPS - 51A Upper Roots Narrow Beak	E51A	S-54
5503	EXTRACTION FORCEPS - 45 Lower Roots	E45	S-54
5504	EXTRACTION FORCEPS - 44 Upper Roots	E44	S-54
5505	EXTRACTION FORCEPS - 95 Upper Molars Left	E95	S-55
5506	EXTRACTION FORCEPS - 121 Upper Wisdoms	E121	S-55
5507	EXTRACTION FORCEPS - 167 Uni versal for Lower Teeth Molars & Premolars	E167	S-55
5508	EXTRACTION FORCEPS - 168 Upper Canines	E168	S-55
5509	EXTRACTION FORCEPS AP - 74N Lower Incisors & Roots Canines	E74N	S-55
5510	EXTRACTION FORCEPS - 74NF Lower Incisors & Roots	E74NF	S-55
5511	EXTRACTION FORCEPS AMERICAN PATTERN - 79 Lower Wisdoms	E79	S-55
5512	EXTRACTION FORCEPS AMERICAN PATTERN - 79A Lower Wisdom	E79A	S-55
5513	EXTRACTION FORCEPS AMERICAN PATTERN - 86C Lower Molar s	E86C	S-55
5514	EXTRACTION FORCEPS AMERICAN PATTERN - 89 Upper Molars Right	E89	S-55
5515	EXTRACTION FORCEPS AMERICAN PATTERN - 90 Upper Molars Left	E90	S-55
5516	EXTRACTION FORCEPS AMERICAN PATTERN - 94 Upper Molars Right	E94	S-55
5517	EXTRACTION FORCEPS AMERICAN PATTERN - F-1Upper Universal	EF-1	S-56
5518	EXTRACTION FORCEPS AMERICAN PATTERN - F-2 Lower Universal	EF-2	S-56
5519	EXTRACTION FORCEPS AMERICAN PATTERN - F-3Upper Anterior	EF-3	S-56
5520	EXTRACTION FORCEPS AMERICAN PATTERN - F-4N Lower Anterior	EF-4N	S-56
5521	EXTRACTION FORCEPS AMERICAN PATTERN - 169 Lower Cnines	E169	S-56
5522	EXTRACTION FORCEPS - 233 Lower Root s Parallel Beak s	E233	S-56
5523	EXTRACTION FORCEPS - 235 Upper Roots Parallel Beaks	E235	S-56
5524	AP - MD-1 Universal Upper Wisdom Premolars & Roots	EMD-1	S-56
5525	EXTRACTION FORCEPS AMERICAN PATTERN - MD-4 Lower Molar s	EMD-4	S-56
5526	EXTRACTION FORCEPS AMERICAN PATTERN - 22L Lower Wisdoms Left	E22L	S-56
5527	EXTRACTION FORCEPS AMERICAN PATTERN - 22R Lower Wisdoms Right	E22R	S-56
5528	EXTRACTION FORCEPS - F-7C Lower Molar Modi fy Cow-Horn Beak	EF-7C	S-57
5529	EXTRACTION FORCEPS AMERICAN PATTERN - F-8 Upper Universal	EF-8	S-57

5530	EXTRACTION FORCEPS AMERICAN PATTERN - F-9 Lower Universal	EF-9	S-57
5531	EXTRACTION FORCEPS AMERICAN PATTERN - F-10 Lower Molar	EF-10	S-57
5532	EXTRACTION FORCEPS AMERICAN PATTERN - F-11 Upper Molar Left	EF-11	S-57
5533	EXTRACTION FORCEPS AMERICAN PATTERN - F-12 Upper Molar Right	EF-12	S-57
5534	EXTRACTION FORCEPS AMERICAN PATTERN - F-13 Upper Molar Right	EF-13	S-57
5535	EXTRACTION FORCEPS AMERICAN PATTERN - F-14 Upper Molar Left	EF-14	S-57
5536	EXTRACTION FORCEPS AMERICAN PATTERN - F-4 Lower Anterior	EF-4	S-57
5537	EXTRACTION FORCEPS AMERICAN PATTERN - F-5 Lower Anterior	EF-5	S-57
5538	EXTRACTION FORCEPS AMERICAN PATTERN - F-6 Upper 3rd Molars	EF-6	S-57
5539	EXTRACTION FORCEPS AMERICAN PATTERN - F-7 Lower 3rd Molars	EF-7	S-57
5540	EXTRACTION FORCEPS - F-15 Upper Premolars Universal	EF-15	S-58
5541	EXTRACTION FORCEPS AMERICAN PATTERN - F-16 Upper Root	EF-16	S-58
5542	EXTRACTION FORCEPS AMERICAN PATTERN - F-17 Lower Root	EF-17	S-58
5543	EXTRACTION FORCEPS - F-18 Upper Roots Fragment Narrow	EF-18	S-58
5544	EXTRACTION FORCEPS - F-18TN Upper Roots Fragment Narrow	EF-18TN	S-58
5545	EXTRACTION FORCEPS - F-19 Lower Roots Fragment Narrow	EF-19	S-58
5546	AMERICAN PATTERN - F-19NW Lower Root s Fragment Narrow	EF-19NW	S-58
5547	EXTRACTION FORCEPS - 46XL Lower Root Extra Long	E46XL	S-58
5548	EXTRACTION FORCEPS - F-20 Upper Premolar 1ST & 2RD Molar s	EF-20	S-58
5549	EXTRACTION FORCEPS AMERICAN PATTERN - 49 Upper Roots	E49	S-58
5550	EXTRACTION FORCEPS AMERICAN PATTERN - 51XL Upper Roots	E51XL	S-58
5551	Pedodontic Extraction Forceps - F-1, Upeer Universal Pedo	EPF-1	S-59
5552	Pedodontic Extraction Forceps - F-2, Lower Universal Pedo	EPF-2	S-59
5553	Pedodontic Extraction Forceps - F-3, Apper Anterior Pedo	EPF-3	S-59
5554	Pedodontic Extraction Forceps - F-4, Lower Anterior Pedo	EPF-4	S-59
5555	Pedodontic Extraction Forceps - F-5, Lower Anterior Pedo	EPF-5	S-59
5556	Pedodontic Extraction Forceps - F-6, Upper 3rd Molar Pedo	EPFU-6	S-59
5557	Pedodontic Extraction Forceps - F-6, Lower 3rd Molar Pedo	EPFL-6	S-59
5558	Pedodontic Extraction Forceps - F-16 Lower Root Pedo	EPF-16	S-59
5559	Pedodontic Extraction Forceps - F-17, Lower Root pedo	EPF-17	S-59
5560	Pedodontic Extraction Forceps - B23, Lower Molar Pedo	EPB23	S-59
5561	Pedodontic Extraction Forceps - 150-SK, Upper Universal Pedo	EP150-SK	S-59
5562	Pedodontic Extraction Forceps - 151-SK, Lower Universal Pedo	EP151-SK	S-59
6000	Rongeurs - Friedman, 14cm/(5-1/2")	RFM	S-60
6001	Rongeurs - Mini Friedman, 14.5cm/(5-3/4")	RMFM	S-60
6002	Rongeurs - Mini Friedman 90°, 14.5cm/(5-3/4")	RMFM90	S-60
6003	Rongeurs - Blumental Standard 30°, 14.5cm/(5-3/4")	RBMS	S-60
6004	Rongeurs - Blumental Standard 45°, 14.5cm/(5-3/4")	RBM45	S-60
6005	Rongeurs - Blumental 90°, 14.5cm/(5-3/4")	RBM90	S-60
6006	Rongeurs - 1A Mead, 16.5cm/(6-1/2")	R1A-M	S-60
6007	Rongeurs - 4A Cleveland, 17cm/(6-3/4")	R4A-C	S-60
6008	Rongeurs - 3H Dean, 17cm/(6-3/4")	R3UD	S-60
6010	Rongeurs - 4 cleveland, 14cm/(5-1/2")	R5C	S-60
6012	Rongeurs - 5 Cleveland Side Cutting, 17cm (6-3/4")	R4C	S-60
6013	Rongeurs - 5S Cleveland Side Cutting, 15cm/(6")	R5SC	S-60
6014	Rongeurs - Kramer-Nevins, 14.5cm/(5-3/4")	RKN	S-60
6020	Rongeurs - Kerrison Rongeurs 3mm for widening the sinus window	RK3M	S-61

6022	Rongeurs - Kerrison Rongeurs 4mm for widening the sinus window	RK4M	S-61
6025	Rongeurs - 18M (Double Action) Side Cutting, 17cm/(6-3/4")	R18MDA17	S-61
6030	Rongeurs - Beyers (Double Action), 18cm/(7")	RBEYDA18	S-61
6035	Rongeurs - Boehler (Double Action), 15cm/(6")	RBOEDA15	S-61
6040	Rongeurs - Boelher (Double Action), 16cm/(6-1/4")	RBOEDA16	S-61
BL1000	Explorers 2	DEXD2	BL-1
BL1001	Explorers 3CH	DEXD3CH	BL-1
BL1002	Explorers 5	DEXD5	BL-1
BL1003	Explorers TU17/23	DEXDTU17/23	BL-1
BL1004	Explorers Old Dominion University	DEXD11/12ODU	BL-1
BL1005	Exprobes 23/CP-12	DEXD23/CP12	BL-1
BL1020	Endo Explorers DG16	EDG16	BL-2
BL1021	Endo Explorers DG16/17	EDG16/17	BL-2
BL1022	Endo Explorers DG16/23	EDG16/23	BL-2
BL1023	Excavators 31L	E31L	BL-2
BL1024	Excavators 31W	E31W	BL-2
BL1025	Excavators 32L	E32L	BL-2
BL1026	Excavators 31W	E33L	BL-2
BL1030	Spreaders D11TS, Thinnest 21mm	ERCSD11TS	BL-3
BL1031	Spreaders MA57, Heat Carrier 27mm	ERCSMA57	BL-3
BL1032	Spreaders D11T, Thinner 22m	ERCSD11TS	BL-3
BL1033	Spreaders 00, 25mm	ERCS00	BL-3
BL1034	Spreaders GP1, 20mm	ERCSGP1	BL-3
BL1035	Spreaders D11, 23mm	ERCSD11	BL-3
BL1036	Pluggers 1/3, 21mm (0.40/0.45mm)	ERCP1/3	BL-4
BL1037	Pluggers 5/7, 21mm (0.50/0.75mm)	ERCP5/7	BL-4
BL1038	Pluggers GL1, 18mm (0.45/GL1)	ERCPGL1	BL-4
BL1039	Pluggers 9/11, 21mm (1.00/1.15mm)	ERCP9/11	BL-4
BL1830	RUBBER DAM CLAMP - 8A	RDC8A	BL-5
BL1831	RUBBER DAM CLAMP - W8A	RDCW8A	BL-5
BL1832	RUBBER DAM CLAMP - 8AD	RDC8AD	BL-5
BL1833	RUBBER DAM CLAMP - 10	RDC10	BL-5
BL1834	RUBBER DAM CLAMP - 11	RDC11	BL-5
BL1835	RUBBER DAM CLAMP - 14	RDC14	BL-5
BL1836	RUBBER DAM CLAMP - 14A	RDC14A	BL-5
BL1837	RUBBER DAM CLAMP - W14A	RDCW14A	BL-5
BL1838	RUBBER DAM CLAMP - 18	RDC18	BL-5
BL1839	RUBBER DAM CLAMP - 24	RDC24	BL-5
BL1840	RUBBER DAM CLAMP - 25	RDC25	BL-5
BL1841	RUBBER DAM CLAMP - 26N	RDC26N	BL-5
BL1842	RUBBER DAM CLAMP - 56	RDC56	BL-5
BL1843	RUBBER DAM CLAMP - 56S	RDC56S	BL-5
BL1844	RUBBER DAM CLAMP - 4	RDC4	BL-5
BL1845	RUBBER DAM CLAMP - W4	RDCW4	BL-5
BL1846	RUBBER DAM CLAMP - 5	RDC5	BL-5
BL1847	RUBBER DAM CLAMP - 8	RDC8	BL-5
BL1848	RUBBER DAM CLAMP - 27	RDC27	BL-5

BL1850	RUBBER DAM CLAMP - 30	RDC30	BL-6
BL1851	RUBBER DAM CLAMP - 31	RDC31	BL-6
BL1852	RUBBER DAM CLAMP - 201	RDC201	BL-6
BL1853	RUBBER DAM CLAMP - 205	RDC205	BL-6
BL1854	RUBBER DAM CLAMP - 3	RDC3	BL-6
BL1855	RUBBER DAM CLAMP - W3	RDCW3	BL-6
BL1856	RUBBER DAM CLAMP - 7	RDC7	BL-6
BL1857	RUBBER DAM CLAMP - W7	RDCW7	BL-6
BL1858	RUBBER DAM CLAMP - 7A	RDC7A	BL-6
BL1859	RUBBER DAM CLAMP - 12A	RDC12A	BL-6
BL1860	RUBBER DAM CLAMP - 31A	RDC31A	BL-6
BL1861	RUBBER DAM CLAMP - 26	RDC26	BL-6
BL1862	RUBBER DAM CLAMP - 28	RDC28	BL-6
BL1863	RUBBER DAM CLAMP - 0	RDC0	BL-6
BL1864	RUBBER DAM CLAMP - 1	RDC1	BL-6
BL1865	RUBBER DAM CLAMP - 1A	RDC1A	BL-6
BL1866	RUBBER DAM CLAMP - 2	RDC2	BL-6
BL1870	RUBBER DAM CLAMP - W2	RDCW2	BL-7
BL1871	RUBBER DAM CLAMP - 2A	RDC2A	BL-7
BL1872	RUBBER DAM CLAMP - 2AS	RDC2AS	BL-7
BL1873	RUBBER DAM CLAMP - 27N	RDC27N	BL-7
BL1874	RUBBER DAM CLAMP - 29	RDC29	BL-7
BL1875	RUBBER DAM CLAMP - 6	RDC6	BL-7
BL1876	RUBBER DAM CLAMP - 9	RDC9	BL-7
BL1877	RUBBER DAM CLAMP - 9S	RDC9S	BL-7
BL1878	RUBBER DAM CLAMP - 212	RDC212	BL-7
BL1879	RUBBER DAM CLAMP - 212SA	RDC212SA	BL-7
BL1880	RUBBER DAM CLAMP - 00	RDC-00	BL-7
BL1900	Rompen Pluggers, Mesial/Distal	BLRPM/D	BL-8
BL1901	Rompen Pluggers, Buccal/Lingual	BLRPB/L	BL-8
BL1902	Jovanovic Pluggers	BLJP	BL-8
BL1903	Labanca Pluggers	BLLBP	BL-8
BL1904	Grisdale Bone Condensers	BLGBP	BL-8
BL1915	Sinus Curettes 1	ISC1	BL-9
BL1916	Sinus Curettes 2	ISC2	BL-9
BL1917	Kramer-Nevins, Obtuse	IKNO	BL-9
BL1918	Kramer-Nevins, Acute	IKNA	BL-9
BL1919	Sinus Lift 300	ISL300	BL-9
BL1930	1/2 Buser Modified Periodontal Chisel 4mm/5mm	I1/2BMPBC	BL-10
BL1931	2/3 Buser Modified Periodontal Chisel 5mm/6mm	I2/3BMPBC	BL-10
BL1932	1/3 Buser Modified Periodontal Chisel 4mm/6mm	I1/3BMPBC	BL-10
BL1933	P24G/13KL Palacci Modified Chisel	1PAL13KG	BL-10
BL1934	3/4 Solt Chisel	ISOLT3/4	BL-10
BL1940	15/16 Kirkland Knife	BLKK15/16	BL-11
BL1941	Allen End-Cutting Intrasulcular Knife	BLKAECI	BL-11
BL1942	1/2 Allen Orban Knife, Round	BLKAOR1/2	BL-11
BL1943	1/2 Allen Orban Knife	BLKAO1/2	BL-11

BL1944	11 Goldman-Fox Knife	BL11GFK	BL-11
BL1945	1/2 Orban Knife	BL1/2OK	BL-11
BL1950	Tunneling Knife 1	BLTK1	BL-12
BL1951	Tunneling Knife 2	BLTK2	BL-12
BL1952	Giles De Quincey Tunneling Knife 1	BLGDQTK1	BL-12
BL1953	Giles De Quincey Tunneling Knife 2	BLGDQTK2	BL-12
BL3230	COMPOSITE INSTRUMENTS - #1	RTNCIGFT1	BL-13
BL3231	COMPOSITE INSTRUMENTS - #2	RTNCIGFT2	BL-13
BL3232	COMPOSITE INSTRUMENTS - #3 Extra-Flexible	RTNCIGFT3	BL-13
BL3233	COMPOSITE INSTRUMENTS - #4 Extra Flexible	RTNCIGFT4	BL-13
BL3234	COMPOSITE INSTRUMENTS - #5 Flexi -Thin	RTNCIGF5	BL-13
BL3235	COMPOSITE INSTRUMENTS - #6 (Flexible Thin)	RTNCIGFT6	BL-13
BL3236	COMPOSITE INSTRUMENTS - Micro Mini	RTNCIPCS	BL-13
BL3237	COMPOSITE INSTRUMENTS - Mini 1	RTNCIGFTM1	BL-13
BL3238	COMPOSITE INSTRUMENTS - Mini 3 flexible-Thin	RTNCIGFTM3	BL-13
BL3239	Composite Instruments Mini 4 Extra-Thin	RTNCIGFTM4	BL-13
BL3240	COMPOSITE INSTRUMENTS - AB1 Boghosian	RTNPFIA1	BL-14
BL3241	COMPOSITE INSTRUMENTS - AB2 Boghosian	RTNPFIA2	BL-14
BL3242	COMPOSITE INSTRUMENTS - Inerproximal Carver	RTNCVIPC	BL-14
BL3243	COMPOSITE INSTRUMENTS - W3	RTNPFIW3	BL-14
BL3244	COMPOSITE INSTRUMENTS - A6	RTNPFIA6	BL-14
BL3245	COMPOSITE INSTRUMENTS - 4/5 Gregg	RTNPFIA8A	BL-14
BL3246	COMPOSITE INSTRUMENTS - Contouring Instrument	RTNCFIR/L	BL-14
BL3247	COMPOSITE INSTRUMENTS - 8A	RTNPFIA8A	BL-14
BL3250	COMPOSITE INSTRUMENTS - 21B	RTNBB21B	BL-15
BL3251	COMPOSITE INSTRUMENTS - 27/29	RTNBB2729	BL-15
BL3252	COMPOSITE INSTRUMENTS - 2 Ladmore	RTNBBL2	BL-15
BL3253	COMPOSITE INSTRUMENTS - 3 Ladmore	RTNBBL3	BL-15
BL3254	COMPOSITE INSTRUMENTS - Small / Medium Ball Burnisher	RTNBBS/M	BL-15
BL3255	COMPOSITE INSTRUMENTS - Freedman "Duckhead"	RTNPCCI	BL-15
BL3256	COMPOSITE INSTRUMENTS - Freedman Small Contact forming	RTNFCIS	BL-15
BL3257	COMPOSITE INSTRUMENTS - Freedman Large Contact Forming	RTNFCIL	BL-15
BL3260	COMPOSITE INSTRUMENTS - A Cosmetic Contouring	RTNCCIA	BL-16
BL3261	COMPOSITE INSTRUMENTS - B cosmetic Contouring	RTNCCIB	BL-16
BL3262	COMPOSITE INSTRUMENTS - F cosmetic Contouring	RTNCCIF	BL-16
BL3263	COMPOSITE INSTRUMENTS - C cosmetic Contouring	RTNCCIC	BL-16
BL3264	COMPOSITE INSTRUMENTS - G Marginal Ridge & Embrasure Shaping Instrument	RTNCCIG	BL-16
BL3265	COMPOSITE INSTRUMENTS - D cosmetic Contouring	RTNCCID	BL-16
BL3266	COMPOSITE INSTRUMENTS - H Occlusal Anatomy Instrument	RTNCCIH	BL-16
BL3267	COMPOSITE INSTRUMENTS - E cosmetic Contouring	RTNCCIE	BL-16
BL3268	COMPOSITE INSTRUMENTS - I Cosmatic Packing Instrument	RTNCCII	BL-16
BL3270	Composite Instruments 3 Extra-Flex	RTNCIGFT3	BL-17
BL3271	COMPOSITE INSTRUMENTS - OT Tanner	RTNPLGOT	BL-17
BL3272	COMPOSITE INSTRUMENTS - Mini 1	RTNCIGFTM1	BL-17
BL3273	COMPOSITE INSTRUMENTS - 3 Hollenback	RTNPLGH3	BL-17
BL3274	COMPOSITE INSTRUMENTS - Micro-Mini	RTNCIPCS	BL-17
BL3275	COMPOSITE INSTRUMENTS - 5A	RTNPLG5A	BL-17

BL3276	COMPOSITE INSTRUMENTS - Medium Placing Condensing Instrument	RTNCIPCM	BL-17
BL3277	COMPOSITE INSTRUMENTS - Small/Medium Contact Forming Instrument	RTNCFISM	BL-17
BL3278	COMPOSITE INSTRUMENTS - Large Placing Condensing Instrument	RTNCIPCL	BL-17
BL3279	COMPOSITE INSTRUMENTS - Medium/Large Contact Forming Instrument	RTNCFIM/L	BL-17
BL3300	Composite Instruments - Composite Spatula, Small 1-1.5mm	CS-1/1.5	BL-18
BL3301	Composite Instruments - Composite Spatula, Large	CSL	BL-18
BL3302	Composite Instruments - Former and Spatula	CF&S	BL-18
BL3303	Composite Instruments - Composite Spatula 1.7mm (Special bend for Mesial and Distal areas)	CS1.7	BL-18
BL3304	Composite Instruments - Composite Spatula (Combo Flattened-Beavertail)	CSCFB	BL-18
BL3305	Composite Instruments - Composite Spatula 1.6mm	CS1.6	BL-18
BL3306	Composite Instruments - Plastic Filling Instruments-Composite Fissure Black	CICFFB	BL-18
BL3310	Composite Instruments - Spatula 2.2mm/Burnisher 1.5mm	S2.2B1.5	BL-19
BL3311	Composite Instruments - Composite Pluggers 1	CP1	BL-19
BL3312	Composite Instruments - Composite Pluggers 2	CP2	BL-19
BL3313	Composite Instruments - Spatula 2.4mm/Burnisher 1.9mm	S2.4/B1.9	BL-19
BL3314	Composite Instruments - Composite Pluggers 3	CP3	BL-19
BL3315	Composite Instruments - Plastic Instrument/Composite Plugger	PICP	BL-19
BL3316	Composite Instruments - Plastic Instrument/Special Plugger	PISP	BL-19
BL3582	Explorer 23, 0.2mm	BLPT23	BL-20
BL3583	Heidemann Spatula	BLHST2	BL-20
BL3584	Explorer T11, 0.2mm	BLPT11	BL-20
BL3585	Pluggers T11 (Wiland Shape), 0.4mm	BLPWST11	BL-20
BL3586	Spatula T1	BLST1	BL-20
BL3587	Spatula T2, Cobra Shape, 0.4mm	BLSCST2	BL-20
BL3588	Spatula T4	BLST4	BL-20
BL3600	Pluggers (Straight)0.5/0.8mm	RPS	BL-21
BL3601	Pluggers T1, Bent Shape 0.5/0.8mm	RPBST1	BL-21
BL3602	Pluggers T2, Bent Shape, 1.0/1.3mm	RPBST2	BL-21
BL3603	Pluggers T3, Bent Shape, 0.5/0.8mm	RPBST3	BL-21
BL3604	Pluggers T4, Bent Shape, 1.0/1.3mm	RPBST4	BL-21
BL3610	Pluggers T5, Columbia Shape, 0.6mm	RPCST5	BL-21
BL3611	Pluggers T6, Columbia Shape, 0.8mm	RPCST6	BL-21
BL3612	Pluggers T12, Wiland Shape, 0.6mm	RPWST12	BL-21
BL3613	Pluggers T13, Wiland Shape, 0.8mm	RPWST13	BL-21
BL3614	Pluggers T14, Wiland Shape, 1.0mm	RPWST14	BL-21
BL3620	Excavator T1, Wiland Shape, 0.8mm	REWST1	BL-22
BL3621	Excavator T2, Wiland Shape, 1.0mm	REWST2	BL-22
BL3622	Excavator T3, Wiland Shape, 1.3mm	REWST3	BL-22
BL3623	Excavator T4, Columbia Shape, 0.8mm	RECST4	BL-22
BL3624	Excavator T5, Columbia Shape, 1.0mm	RECST5	BL-22
BL3625	Excavator T6, Columbia Shape, 1.3mm	RECST6	BL-22
BL3626	Excavator T6 1/2, Straight -1.0mm	REST6-1/2	BL-22
BL3627	Excavator T6 3/4, Straight, 1.3mm	REST6-3/4	BL-22
BL3628	Curette 13/14, Columbia Shape	RCCS-13/14	BL-22
BL3629	Scaler 13/14, Columbia Shape	RSCS-13/14	BL-22
BL3650	Buser Periosteal Elevator	BLBPE	BL-23
BL3651	Allen Periosteal elevator	BLAPE	BL-23

BL3652	Allen Periosteal elevator, Anterior	BLAPEA	BL-23
BL3653	9 Molt Periosteal Elevator	BL9MPE	BL-23
BL3654	3 Prichard Periosteal Elevator	BL3PPE	BL-23
BL3655	2/4 molt curette	BLMC2/4	BL-23
BL3670	Swiss Peio, Modified P24GSP Periosteal 1	BLSPMP24GSPP	BL-24
BL3671	Swiss Perio, Modified P24GSIC Periosteal	BLSPMP24GSICP	BL-24
BL3672	swiss Periosteal, Modified P24GSP Periosteal 2	BLSPMP24GSP2	BL-24
BL3673	Swiss Periosteal, Modified Schlee Periosteal	BLSPMSP	BL-24
BL3674	Swiss Periosteal Modified PH2 Perio	BLSPMPH2	BL-24
BL3675	Modified Labanca Perio	BLMLP	BL-24
BL3680	Urban Mini Me periosteal	BLUMMP	BL-25
BL3681	Buser/Modified PH2MBH, Small	BLBMPH2MBHS	BL-25
BL3682	Buser/Modified PH2MBH, Large	BLBMPH2MBHL	BL-25
BL3683	152 Kramer-Nevins Periosteal	BLKN152P	BL-25
BL3684	24G Periosteal	BL24GP	BL-25
BL3685	14 Goldman-Fox Periosteal	BL14GFP	BL-25
BL3690	Surgical Curette 86 Lucas-2.8mm	BLLSC86	BL-26
BL3691	Surgical Curette 84 Lucas-2.0mm	BLLSC84	BL-26
BL3692	Surgical Curette 85 Lucas-2.5mm	BLLSC85	BL-26
BL3693	Surgical Curette 88 Lucas-4.7mm	BLLSC88	BL-26
BL3694	Surgical Curette 87 Lucas-3.5mm	BLLSC87	BL-26
BL4180	Spade Elevator, Straight	BLEZLSS	BL-27
BL4181	Spade Elevator, Mesially Angled	BLEZLSMA	BL-27
BL4182	Spade Elevator, Distally Angled	BLEZLSDA	BL-27
BL4183	Spade Elevator, Left	BLEZLSL	BL-27
BL4184	Spade Elevator, Right	BLEZLSR	BL-27
BL4185	2.5mm, Small, Straight	BLEZLSS2.5	BL-27
BL4186	2.5mm, Small, Curved	BLEZLSC2.5	BL-27
BL4187	2.5mm, Small, Mesially Angled	BLEZLSMA2.5	BL-27
BL4188	2.5mm, Small, Distally Angled	BLEZLSDA2.5	BL-27
BL4190	3.75mm, Large, Straight	BLEZLLS3.75	BL-28
BL4191	3.75mm, Large, Curved	BLEZLLC3.75	BL-28
BL4192	3.75mm, Large, Mesially Angled	BLEZLLMA3.75	BL-28
BL4193	3.75mm, Large, Distally Angled	BLEZLLDA3.75	BL-28
BL4194	4.5mm, Large, Straight	BLEZLLS4.5	BL-28
BL4195	4.5mm, Large, Curved	BLEZLLC4.5	BL-28
BL4196	4.5mm, Large, Mesially Angled	BLEZLLMA4.5	BL-28
BL4197	4.5mm, Large, Distally Angled	BLEZLLDA4.5	BL-28
BL4198	Small, Tight, Mesial Angle, 2.5mm	BLEZLSTMA2.5	BL-28
BL4199	Tight Mesial Spade, 2.5mm	BLEZLTMS2.5	BL-28
BL4205	Retractors, University of Minnesota	BLRUOM	BL-29
BL4206	Retractors, Lindo-Serrated	BLRLS	BL-29
BL4207	Retractors, Front Round	BLRFR	BL-29
BL4208	Retractors, Front cut	BLRFC	BL-29
BL4209	Retractors, Cawood Minnesota	BLRCWM	BL-29
BL4210	Retractors Weider Tongue, Small (30mm)	BLWTRS	BL-29
BL4211	Retractors Weider Tongue, Large (35mm)	BLWTRL	BL-29

BL4212	Retractors Bishop	BLRBISH	BL-29
BL4220	Elevator Luxator, Straight-2mm	BLLXS2	BL-30
BL4221	Elevator Luxator, Curved-2mm	BLLXC2	BL-30
BL4222	Elevator Luxator, Straight-2.5mm	BLLXS2.5	BL-30
BL4223	Elevator Luxator, Curved-2.5mm	BLLXC2.5	BL-30
BL4224	Elevator Luxator, Straight-3mm	BLLXS3	BL-30
BL4225	Elevator Luxator, Curved-3mm	BLLXC3	BL-30
BL4226	Elevator Luxator, Straight-4mm	BLLXS4	BL-30
BL4227	Elevator Luxator, Curved-4mm	BLLXC4	BL-30
BL4228	Elevator Luxator, Straight-5mm	BLLXS5	BL-30
BL4229	Elevator Luxator, Curved-5mm	BLLXC5	BL-30
BL4250	Elevator Slim Luxating Long, Straight 2mm	SLELS2	BL-31
BL4251	Elevator Slim Luxating Long, Curved 2mm	SLELC2	BL-31
BL4252	Elevator Slim Luxating Long, Straight 3mm	SLELS3	BL-31
BL5253	Elevator Slim Luxating Long, Curved 3mm	SLELC3	BL-31
BL4255	Elevator Slim Luxating Long, Straight 4mm	SLELS4	BL-31
BL4256	Elevator Slim Luxating Long, Curved 4mm	SLELC4	BL-31
BL4257	Elevator Slim Luxating Long, Straight 5mm	SLELS5	BL-31
BL4258	Elevator Slim Luxating Long, Curved 5mm	SLELC5	BL-31
BL4316	1 Heidbrink Root Tip Pick	BLEHB1	BL-32
BL4317	2 Heidbrink Root Tip Pick	BLEHB2	BL-32
BL4318	3 Heidbrink Root Tip Pick	BLEHB3	BL-32
BL4321	13/14 Heidbrink Root Tip Pick	BLEHB13/14	BL-32
BL4323	9R Apical Root Tip Pick	BLE9R	BL-32
BL4324	9L Apical Root Tip pick	BLE9L	BL-32
BL4343	10 Davis Double-end Root Tip Pick	BLED10	BL-33
BL4344	11 Davis Double ended Root Tip Pick	BLEDD11	BL-33
BL4345	2 Howard Root Tip Pick	BLEHO2	BL-33
BL4346	3 Howard Root Tip Pick	BLEHO3	BL-33
BL4347	Periotomes, Posterior (DE)	BLPPDE	BL-34
BL4348	Periotomes, Posterior, Right (SE)	BLPPRSE	BL-34
BL4349	Periotomes, Posterior, Left (SE)	BLPPLSE	BL-34
BL4350	Periotomes, Anterior (DE)	BLPADE	BL-34
BL4351	Periotomes, Anterior, Angled (SE)	BLPAASE	BL-34
BL4352	Periotomes, Anterior, Straight (SE)	BLPASDE	BL-34
BL4353	Periotomes, Wise (DE)	BLPW	BL-34